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KeyGuard™ 3D PRINTING RESINS | FOR PROFESSIONAL USE ONLY (EN)

Intended Purpose KeyPrint® KeyGuard™ is a biocompatible photopolymer resin for the fabrication of custom-made athletic mouth guards.

Product Description KeyGuard™ is designed for additive manufacturing in vat polymerization 3D printers utilizing wavelengths between 385nm-405nm. KeyGuard™ is for the fabrication of custom-made athletic mouth guards. KeyGuard™ is intended to be used within a computer-aided design and manufacturing (CAD/CAM) system that includes a 3D scanner, design software, 3D printer, and post-cure unit.

.....

Warnings & Precautions

KeyGuard™ is contraindicated for patients with known hypersensitivity or severe allergic reactions to acrylate monomers and oligomers.

Review the product Safety Data Sheet (SDS) prior to use.

KeyPrint® materials contain components that may cause skin dryness or irritation, sensitization (allergic contact dermatitis), or other allergic reactions in susceptible persons. Wear protective gloves, eye protection, and protective clothing when handling liquid resin.

Skin Contact: Wash thoroughly with soap and water. If skin sensitization occurs, discontinue use. If dermatitis or other symptoms persist, seek medical assistance. **Eye Contact:** Flush eyes promptly with copious amounts of water for 15 minutes and consult a physician. Wash skin with soap and water. **Inhalation:** Move subject to fresh air. Give oxygen or artificial respiration as required. **Ingestion:** Contact your regional Poison Control Center immediately.

Any serious incident related to the product should be reported to the manufacturer and the competent authority according to local and national regulations.

In case of Emergency: INFOTRAC: 800-535-5053 Outside the US: (001) 800-535-5053 Customer # 76987

Step-by-Step Instructions

KeyGuard™ is compatible with certified KeyPrint Compatible printers, wash stations, and cure units. Visit *https://keyprint.keystoneindustries.com/compatible/* for a current list of certified KeyPrint Compatible hardware.

Caution: Using non-certified equipment with KeyGuard™ may result in athletic mouth guard that is not safe an effective. Keystone Industries is not liable for product failures resulting from usage of non-certified equipment.Follow applicable product labeling and manufacturer Instructions for Use for all CAD/CAM components that are used in conjunction with the KeyGuard™.

.....

A. Design

KeyGuard™ is compatible with all Computer-Aided-Design (CAD) software that support the fabrication of dental devices. Keystone recommends the use of occlusal splint design workflows for digital design of athletic mouth guards.

Design Recommendations for KeyGuard™

Blockout: For maximum engagement tooth engagement, Keystone recommends not using undercut removal. Minimum thickness: 2 mm. Keystone recommends a facial and occlusal thickness of at least 3.5-4 mm. Radius Compensation: Off

B. Print

Caution: Wear proper personal protective equipment when handling uncured KeyPrint® resin. Do not mix different batches of KeyGuard™. Use a dedicated resin vat/cartridge for each lot of KeyGuard™ for optimal results.

Note: Use only certified KeyPrint Compatible printer settings. Visit https://keyprint.keystoneindustries.com/compatible/ for a current list of certified KeyPrint Compatible 3D Printers. For certified printer settings, reach out to your printer manufacturer.

- Thoroughly shake the bottle 1 hour prior to use for a minimum of 5 minutes, to ensure all pigments are dispersed. Allow the bottle to sit for 5 minutes to allow dissipation of bubbles.
- Leave the bottle at ambient temperature (20-25°C/68-77°F) prior to printing.
- In the unlikely event of a print failure, filter the liquid resin through a mesh screen with pore size < 200 microns. It is good practice to filter the resin in a vat periodically to prevent print failures.

Nesting Recommendations for KeyGuard™

KeyGuard™ devices should be printed at a minimum nesting angle of 75 degrees vertical for optimal accuracy and print quality. Supports should have a minimum of 0.9 mm contact width and minimum overall contact width of 1.5 mm.

C. Wash

Caution: Do not use denatured alcohol or ethanol for washing. Printed mouth guards should not exceed 3 minutes of IPA exposure as the parts may start to degrade.

Note: Although not required, Keystone Industries recommends the use of a post-wash agitation. Use of a small brush can be used for additional cleaning. Printed mouth guards should be washed with Isopropanol alcohol (> 97% purity) within approximately 8 hours from the completion of the print.

Post-Wash Recommendations for KeyGuard™

Two consecutive 1 -minute washes in fresh Isopropyl alcohol (>97% concentration). Dry KeyGuard™ mouth guards with compressed air dry in between wash 1 and wash 2 to help remove residual resin. The use of a soft bristled brush can also be used to remove residual resin.

D. Dry

Allow 15-20 minutes of air drying before post-cure.

E. Cure

Caution: For minimal surface defects or warping, allow parts to cool down before handling.

Note: Use only certified KeyPrint Compatible post cure settings. Visit https://keyprint.keystoneindustries.com/compatible/ for a current list of certified KeyPrint Compatible post cure equipment and curing times. For certified post cure settings, reach out to your equipment manufacturer.

G. Disposal

Unused surplus and non-recyclable liquid resin materials may be considered hazardous waste, dispose of in accordance with all local, regional, national, and international regulations.

Storage Store in a cool, dry place 15°C-30°C (59°F-86°F) and away from light. Ensure that the bottle is capped while not in use.

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