



DENTAL & MEDICAL SUPPLIES SOURCE

Summer 2024

Community Health & Safety Net Clinic



Rely On Us
for Business, Clinical,
Technology, and
Supply Chain Solutions.



Septodont
Septocaine® 4%
and Epinephrine
(228-8210)
Page 4



vVardis
International
Curodont™ Repair
Fluoride Plus
(687-0057)
Page 27



Crosstex®
Ultra® Sensitive
Face Masks,
Level 3
(107-7191)
Page 37



Abbott
Afinion™ HbA1c
Test Kit
(136-4596)
Page 53



INDEX

Abbott	51-52
ACCLEAN by Henry Schein.....	26
ADHA.....	21
AED.....	68
Anesthetics.....	4
Armor Dental	21
B. Braun.....	72
Baxter	64
BD.....	53-54, 69
Bionix.....	45
Butterfly.....	68
CareQuest.....	24
Cements & Liners.....	8
Centrix	25
Cepheid.....	55
Clorox Healthcare	48
Cold Storage	85
Cosmetic & Restorative Dentistry	15
Cranberry	33
CRITERION by Henry Schein	31
Crown & Bridge.....	16
Dental Top Sellers.....	10-11
Diagnostic Equipment.....	65-66
Exam Room Disposables	41-42
Exam Room Equipment	60, 62-63
Face Masks.....	37-38
GC America.....	3
Gloves.....	32
Halyard	34
Health & Beauty Aids	83
HealthFirst	30
Henry Schein Brand Products.....	18, 22, 26, 70, 73, 84
Hygiene Handpieces	20
Hypodermics.....	71
Infection Prevention.....	47, 49-50
Instruments	23
Integra Miltex.....	45
Kerr.....	6
Kuraray	5
Matrix Material	25
Metrex.....	46
Midmark.....	61
NSK.....	17
Office Supplies.....	85
Over-the-Counter (OTC)	82
Patient Monitoring	67
Personal Protection Equipment.....	39-40
Pharmaceuticals.....	78-81
Practice Green	35-36
Practice Marketing	87
Preventives.....	28
ProRepair	19
Pulpent	7
QuidelOrtho.....	56-57
Richmond.....	43
Sekisui.....	58-59
Septodont	5
Sherman Specialty.....	86
Shofu	14
Solvantum (Formerly 3M Healthcare).....	9, 12
Spring Health Products	21
Surgical Specialties	43
Ultrasound.....	68
Vaccines.....	74-77
Viome	27
Voco	Back Cover
vVardis.....	27
Waste Management	84
Wound Care	43-44
X-ray	29

WHAT'S INSIDE?

Check Out vVardis Curodont Repair Fluoride Plus!

See page 27

vVARDIS

SWITZERLAND

BIOMIMETIC DENTAL SCIENCE



(687-0057)

Just watch and wait to maximize patient outcomes with this unique brush-on system!

Look for CaviWipes from Metrex!

See pages 46 & 47

Metrex
INFECTION PREVENTION



(640-8482)



(640-0122)



(640-0024)



(140-3081)




(640-2805)

When it comes to keeping communities healthy, every decision matters. Make the easy choice with CaviWipes towelettes.




**LOOK FOR THIS ICON
THROUGHOUT THE FLYER FOR
SPECIAL LIMITED-TIME SAVINGS!**

THE HEALTH OF YOUR MOUTH AND YOUR HEART ARE CONNECTED



KNOW YOUR RISK

2 in 5 adults have gum or dental disease



IF YOU HAVE GUM DISEASE, YOU ARE

- ↔ **Twice as likely** to die from heart attack
- ↔ **Three times as likely** to die from stroke

And your risk increases over the age of 65

Teeth, Gums & Heart Health Connection


Gum disease and cavities can lead to serious problems for your heart health.



- ↔ Germs from the mouth enter the bloodstream and may add to plaque in the neck and heart arteries—these blockages can cause heart attack and stroke.
- ↔ Germs in the mouth can trigger inflammation that swells blood vessels, reduces blood flow, and increases risk of clots.

High Blood Pressure & the Health of Your Mouth

High blood pressure not only affects your heart health but can also impact your teeth and gums.



- ↔ Changes in blood flow can affect the health of gums.
- ↔ If you have high blood pressure, discuss your personal risks with your provider and get screened for other health conditions.

Prevention is Power is a new, multi-year public health awareness campaign to drive patient adoption of a preventive, integrated health care model to improve overall health outcomes. We will increase health literacy; foster an understanding of the connection between oral health and overall health for providers and patients; and promote the importance of maintaining routine preventive care to reduce the burden of noncommunicable diseases (NCDs), including heart disease, diabetes, and other conditions.

The campaign will include the Prevention is Power Health Care Screening Program, a special partnership with the National Association of Community Health Centers (NACHC) to promote utilization of integrated preventive health care, maintenance, and early detection among patients with chronic conditions through patient education and support. Selected CHCs receive cash, in-kind grants, and educational materials to support community-based screening events, campaigns, and programs related to NCDs and oral health.

The campaign will also include pilots in select geographies in partnership with providers to bring new patients into regular primary care, including oral care, as well as provide tools for clinicians in private practice and larger health systems.

Research is clear:

Regular access to preventive health care = Improved health for individuals, communities, and society = Higher productivity and economic growth

Overcoming Challenges

- 64%** of Americans only go to a physician or dentist when something feels “**extremely wrong**”
- 53%** say it is difficult to access health care
- 35%** say they would benefit from more information and education about the importance of screenings and routine care



FOR MORE INFORMATION



Learn more about Prevention is Power



Health Care Resources




Articaine 4% with Epinephrine

HENRY SCHEIN®

Indicated for local, infiltrative, or conductive anesthesia in both simple and complex dental procedures.

Articaine 4%.....50 per Box
Specify:
1:100,000, Gold **(112-5510)** 1:200,000, Silver **(112-5511)**



Lidocaine HCl 2% with Epinephrine Injection, USP

HENRY SCHEIN®

1.7-mL Anesthetic Cartridges

- Validated siliconization process: ensures smooth injection
- Meets all FDA requirements: allows total flexibility in brand selection

1:100,000, Red
(465-1205)..... 50 per Box
Contains: 5 blister packs of 10 cartridges per package.
1:50,000, Green
(465-0015)..... 50 per Box
Contains: 5 blister packs of 10 cartridges per package.



Septocaine® 4% and Epinephrine SEPTODONT

Articaine hydrochloride 4% with epinephrine.

Cartridges.....50 per Box
Specify:
1:100,000 Injection **(228-8210)**
1:200,000 Injection **(228-0944)**

Aspirating Syringe

HENRY SCHEIN®

- Stainless steel construction: provides strength
- Accepts 1.8cc carpules
- Accepts metal and plastic hub disposable needles
- Autoclavable



Astra Type
(100-4612)..... Ea.



C-W Type-Petite
(570-0410)..... Ea.



Premium Dental Needles-Plastic Hub

HENRY SCHEIN®

Features a softer and less aggressive hub and sheath. Rounded hub provides gentler contact against cheek and gum. No ribs on hub where blood can coagulate. Shorter butt end and larger bore for easier aspiration. Indicating mark on hub for simple identification of triple-bevel position. Round cutting edge protects gum. Thinner needle wall increases flexibility. Double-silicone coating reduces pain and tissue trauma. Seals and hubs are color coded. Plastic hub screws on both metric and imperial syringes. Sterilized.

Disposable..... 100 per Box
Specify:
30 Ga. X-Short-Blue ¾" **(900-1888)**
30 Ga. Short-Blue 1" **(900-1887)**
27 Ga. Short-Yellow 1" **(900-1886)**
27 Ga. Long-Yellow 1¾" **(900-1885)**
25 Ga. Short-Red 1" **(900-1884)**
25 Ga. Long-Red 1¾" **(900-1883)**



Benzo-Jel™ Topical Anesthetic Gel

HENRY SCHEIN®

20% benzocaine gel that provides effective relief of discomfort from local anesthetic injections, periodontal curettage, impression taking, and intraoral radiographs. Quick onset with no bitter aftertaste prevents excessive salivation with minimal wash-away. Gluten-free.

1 oz. Jar..... Ea.

Specify:

Banana **(570-0366)** Mint **(570-0363)**
Bubble Gum **(570-0362)** Piña Colada **(570-0364)**
Cherry **(570-0361)** Raspberry **(570-0365)**
Grape **(570-0344)** Strawberry **(570-0360)**



Lidocaine HCl 2% with Epinephrine COOK-WAITE

1.8 mL Injectable Anesthetic Cartridges Injection (lidocaine hydrochloride and epinephrine injection, USP)

Lidocaine.....50 per Box
Specify:
1:50,000 **(132-9848)**
1:100,000 **(132-9849)**

Oraqix® DENTSPLY SIRONA

(Lidocaine and Prilocaine Periodontal Gel 2.5%/2.5%)

Patient-friendly, needle-free application that is ideal for those requiring local anesthesia in periodontal pockets during scaling and/or root planing. Applied to the gingival margin and then to the periodontal pocket. Dispenses as a liquid and sets as a gel. Fast 30-second onset provides 20 minutes of anesthesia. Do not inject. Do not use on children.

Oraqix® Gel
(555-0626)..... Ea.
Contains: 20-17 Gm. cartridges & 20 blunt-tipped applicator tips.
Autoclavable Dispenser (Sold Separately)
(555-4520)..... Ea.

kuraray

Noritake



FAST, STRONG, SEAMLESS
THE BOND EVOLUTION

SIMPLIFIED UNIVERSAL RESTORATIVE SYSTEM FOR BONDING

CLEARFIL™ Universal Bond Quick 2

CLEARFIL™ Universal Bond Quick 2 provides new chemistry for next generation single-bottle adhesive dentistry.

Learn More



Doing Well by Doing Good



Visit www.HenrySchein.com/HSCares
to find out more about the Calendar of Caring Programs.



Thank you for helping us raise more than \$3 million
to help health happen!

Advanced Biotechnology for Easy Obturation



BIOROOT® FLOW BIOCERAMIC SEALER

- **Adaptable** for your technique warm or cold
- Hydrophilic for **reaching auxiliary canals**
- Resin-free bioceramic **provides a hermetic seal** without microleakage
- Highly pure nanotechnology for **outstanding biocompatibility**
- Creates hydroxyapatite aiding in healing with **no post-op sensitivity**
- **Easy to use** pre-mixed syringe



CASE STUDY
Dr. Lauren Kuhn
Nuth's newest
BioRoot Flow case

#228-6136:
1 (2g) syringe, 1 finger grip, 20 application tips



Questions about Septodont products? Visit septodontusa.com to learn more.



septodontusa.com
from the manufacturer of Septocaine®

NEW!

SimpliShadeTM

BULK FILL & BULK FILL FLOW

- ✓ One Shade
- ✓ Ease of Sculpting
- ✓ Seamless Color Blending
- ✓ Lower Shrinkage Stress
- ✓ High Adaptability
- ✓ Exceptional Esthetics



Ordering Information

SimpliShade Bulk Fill Flow

Item Code	Description
(145-4044)	SimpliShade Bulk Fill Flow Composite Syringe
(145-4043)	SimpliShade Bulk Fill Flow Composite Unidose

SimpliShade Bulk Fill

Item Code	Description
(145-7204)	SimpliShade Bulk Fill Composite Syringe
(145-7205)	SimpliShade Bulk Fill Composite Unidose

Order now through your Henry Schein Sales Representative

kerrendental.com

©2024 Kerr Corporation. All trademarks are property of Kerr Corporation.
All Rights Reserved. MKT-24-0009-HSD Rev 0

HENRY SCHEIN[®]
SOLUTIONS FOR HEALTH CARE PROFESSIONALS

Kerr[™]



Merging Restorative and Preventive Dentistry

Premier Partner in Caries Management



ACTIVA™ BioACTIVE-RESTORATIVE™

Dual cure ensures complete depth of cure - Injection technique eliminates voids
Proven fracture and wear resistance - Mineral apatite formation repairs teeth

ACTIVA™ Restorative/Kids Starter Kit - A1 (384-0135), A2 (384-0136), A3 (384-0161)

ACTIVA™ Restorative/Kids Value Refills - A1 (384-0149), A2 (384-0150), A3 (384-0151), A3.5 (384-0152)



ACTIVA™ BioACTIVE-CEMENT™

Continuously seals margins - Ideal for CAD/CAM, implants and all C&B

ACTIVA™ Cement Value Refills - Transparent (384-0168), A2 (384-0167)



ACTIVA™ BioACTIVE-BASE/LINER™

Seals dentin - Absorbs occlusal stress - Replenishes lost calcium and phosphate

ACTIVA™ Base/Liner Value Refills - VB2 (384-0157)

ACTIVA™ Presto™

Universal Stackable Composite • Light Cure

Mineral-enriched, ionic - Moisture friendly

ACTIVA™ Presto™: 2 x 1.2mL/2 gm syringes + 20 applicator tips

A1 (384-0194), A2 (384-0195), A3 (384-0196), A3.5 (384-0197),

A4 (384-0200), A6 (384-0202), B1 (384-0203), BW (384-0204)



POWERED BY



CRYSTA™
MCP TECHNOLOGY
SUPPORTS REMINERALIZATION

Lime-Lite™ ENHANCED

LIGHT CURE CAVITY LINER

Improved moisture-friendly formula:

Patented rubberized resin for shock absorption like Dentin

Greater release of Calcium, Phosphate and Fluoride

Lime Lite™ Enhanced - (384-0184)



embrace™ wetbond™ pit & fissure sealant

Moisture friendly ionic resin chemistry ensures exceptional marginal adaptation and effortless placement. Seals against microleakage and caries.

Embrace™ wetbond Pit & Fissure Sealant: 4 x 1.2mL syringes + tips - Natural (384-2864) or Off-White (384-3758)

Bulk Kit: 20 x 1.2mL syringes + tips - Natural (384-4671) or Off-White (384-6587)





Vitrebond[™] Plus 3M

Light-cure glass-ionomer liner/base. Packaged in the convenient Clicker[™] dispensing system for a faster, easier, and more consistent mix. Fluoride release of a true glass ionomer.

Single Pack #7571

(777-5630)..... Ea.
Contains: 1–10 Gm. Clicker[™] dispenser, 1 large mixing pad, 1 technique card & instructions.

Double Pack #7572

(777-6083)..... Ea.
Contains: 2–10 Gm. Clicker[™] dispensers, 1 large mixing pad, 1 technique card & instructions.



Gluma[®] Desensitizer KULZER

Fast and effective results with no mixing or curing. Can be used for treating cervical hypersensitivity; pretreating, priming, and bonding direct-composite restorations; pretreating prior to placing compomer or amalgam restorations; after crown & bridge preparation; and pretreating prior to crown & bridge cementing and bonding.

Standard Package

(665-8411)..... 5 mL

Clinic Package

(665-5804)..... Ea.

Contains: 3–5 mL bottles, 150 brushes & 1 tray.

Single-Dose Applicators—0.075 mL

(665-0090)..... 40 per Pkg.



TempBond[®] Clear[™] KERR

Dual-cure, temporary, resin-based cement with triclosan. Used for cementation of temporary and provisional restorations such as crowns, bridges, inlays, and onlays. Flexible and translucent when cured.

Handheld Automix Syringe

(123-4847)..... Ea.

Contains: 1–6 Gm. automix syringe & 10 mixing tips.

Automix Tips

(123-8703).....50 per Box



RelyX[™] Luting Plus Automix 3M

Resin-modified, glass-ionomer cement packaged in a handheld automix syringe. Tack light-cure option for faster removal of excess cement.

Trial Kit #3535TK

(777-0455)..... Ea.

Contains: 1–8.5 Gm. handheld automix syringe & 12 automix mix tips.

Value Package #3535

(777-0454)..... Ea.

Contains: 3–8.5 Gm. handheld automix syringes & 36 mixing tips.

Mixing Tips

(777-0456)..... 12 per Pkg.



GC Fuji PLUS[®] Capsules GC AMERICA

No-Etch Capsules

Starter Package

(333-3524)..... Ea.

Contains: 48 capsules & Capsule Applicator III.

Refill Package

(333-3522)..... 48 per Box

Capsule Applicator

(333-2891)..... Ea.

Capsule Applicator IV

(333-6531)..... Ea.



GC FujiCEM[®] 2 GC AMERICA

Resin-modified, glass-ionomer luting cement. Powered by F(2) Flex Fuse Technology, which improves indirect restorations with increased strength, high fluoride release, low film thickness, and excellent marginal integrity. Use on all types of metal, resin, alumina, and zirconia-based inlays, onlays, crowns, and bridges. Radiopaque.

Starter Kit

(333-3667)..... Ea.

Contains: 1–13.3 Gm. cartridge, 20 mixing tips & 1 FujiCEM[™] 2 dispenser.

Refill Package with Tips

(333-3664)..... Ea.

Contains: 2–13.3 Gm. cartridges & 42 mixing tips.

Refill Package without Tips

(333-3665)..... Ea.

Contains: 2–13.3 Gm. cartridges & 1 mixing pad.

Slide & Lock Mixing Tips

(333-3971).....20 per Pkg.

Intraoral Mixing Tips SL

(333-3969)..... 10 per Pkg.

Dispenser

(333-3270)..... Ea.



RelyX[™] Luting Plus 3M

Resin-modified, glass-ionomer cement. Packaged in the convenient Clicker[™] dispenser. Easy to dispense, mix, and load. Chemically bonds to tooth and restoration. Tack light-cure option for faster removal of excess cement.

Trial Kit #3525TK

(777-2737)..... Ea.

Contains: 1–11.5 Gm. Clicker dispenser of cement, technique guide & set of instructions.

Refill Package #3525

(777-6701)..... Ea.

Contains: 2–11.5 Gm. Clicker dispensers of cement & set of instructions.



Compare to Ultradent Ultra Blend[®] Plus

Lime-Lite[™] Enhanced PULPDENT

Light-cure, fluoride-releasing cavity liner/base material formulated for use with adhesives. Improved moisture-friendly formula and greater adhesion to dentin. Greater release of calcium and phosphate. Offers high compressive strength and wear resistance without any sensitivity. For use with adhesives, composites, and conventional restorative materials.

Complete Kit

(384-0184)..... Ea.

Contains: 4–1.2 mL syringes & 20 applicator tips.

Syringe Refill

(384-0183)..... Ea.

Contains: 1–3 mL syringe & 20 applicator tips.



GC Fuji I[®] Glass Ionomer Luting Cement GC AMERICA

For final cementing of crown & bridge restorations. Bonds to tooth structure and metal, providing excellent strength and marginal integrity. For securing inlays, onlays, posts, and orthodontic bands. Low film thickness, biocompatible, and high fluoride release.

1:1 Package

(333-5075)..... Ea.

Contains: 1 of each: 35 Gm. powder, 25 Gm. liquid, spatula & mixing pad.

Powder

(333-6051)..... 35 Gm.

Liquid

(333-2339).....25 Gm.

GC Fuji I[®] Capsules

(333-3518).....48 per Box

Capsule Applicator

(333-2891)..... Ea.

Capsule Applicator IV

(333-6531)..... Ea.

3M™ Clinpro™ Clear Fluoride Treatment

Clear and simple.

Coming soon: 3M™ Clinpro™ Clear Fluoride Treatment

The most remarkable thing about the next generation of fluoride treatment is what's not included. Featuring a clear, water-based formula and an easy, press-and-go delivery system, new Clinpro Clear Fluoride Treatment makes protecting smiles fast, comfortable, and effective:



Improved smell, taste and mouthfeel.



A smooth, clean appearance and feeling on teeth.



Minimum contact time of just 15 minutes.



The same fluoride uptake as traditional varnishes, but with a lower dosage.



Sign up to be the first to know about the launch of this exciting new product!
Scan QR code or visit: go.3M.com/4Tid

HenryScheinDental.com/3M

Ask your Henry Schein Dental Sales Consultant for complete ordering information.

Solventum and Clinpro are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. © Solventum 2024.

DENTAL TOP SELLERS



Filtek™ Supreme Flowable Restorative 3M

Flowable composite with unique nanotechnology that offers enhanced aesthetics and high wear resistance. Updated ergonomic syringe designed for improved comfort and easier injection of material. Use in anterior and posterior restorations, including liner/base and all Classes I–V. Flow-on-demand handling, excellent polish retention, and improved fluorescence.

Syringe Refills #6032

2 Gm..... 2 per Pkg.

Specify:

A1	(777-2466)	B2	(777-2472)
A2	(777-2467)	C2	(777-2473)
A3	(777-2468)	D2	(777-2475)
A3.5	(777-2469)	White	(777-2478)
A4	(777-2470)	Extra White	(777-2479)
B1	(777-2471)	Opaque A3	(777-2477)

Contains: 2–2 Gm. syringes & 20 dispensing tips.

Capsule Refills #6033

0.2 Gm..... 20 per Pkg.

Specify:

A1	(777-0313)	B2	(777-0320)
A2	(777-0314)	C2	(777-0321)
A3	(777-0315)	D2	(777-0326)
A3.5	(777-0316)	White	(777-0330)
A4	(777-0317)	Extra White	(777-0327)
B1	(777-0319)	Opaque A3	(777-0329)



QUATTROcare™ plus Spray KAVO

For cleaning and care of instruments/air motors in the QUATTROcare™ system. Designed specifically to withstand sterilization. Relubrication after sterilization is unnecessary, eliminating the risk of recontaminating sterilized handpieces.

QUATTROcare™ Plus Spray

(628-9277)..... 500 mL



TPH Spectra® ST DENTSPLY SIRONA

Universal composite restorative material. Filler technology, Sphere TEC™, has advanced granulated spherical filler that works together with the resin matrix to enhance handling and aesthetics. Simplified 5-shade system covers the entire VITA® shade range. Excellent cavity adaptation, sculptability, and stain resistance. Fast and easy to polish.

Compules® Tips

Introductory Kit..... Ea.

Specify:

LV–Low Viscosity(222-2278) HV–High Viscosity (222-2281)

Compules® Tip Refills

0.25 Gm..... 20 per Box

Specify:

Low Viscosity (LV)

A1	(222-2226)	A3.5	(222-2229)
A2	(222-2227)	A4	(222-2230)
A3	(222-2228)	BW	(222-2231)

High Viscosity (HV)

A1	(222-2232)	A3.5	(222-2239)
A2	(222-2233)	A4	(222-2240)
A3	(222-2237)	BW	(222-2241)

Syringe Refill..... 3 Gm.

Specify:

Low Viscosity (LV)

A1	(222-2243)	A3.5	(222-2247)
A2	(222-2244)	A4	(222-2249)
A3	(222-2246)		

High Viscosity (HV)

A1	(222-2251)	A3.5	(222-2254)
A2	(222-2252)	A4	(222-2255)
A3	(222-2253)		

Shade Guide

(222-2257)..... Ea.

Compule® Tips Gun

(222-4148)..... Ea.



Etch Gel

HENRY SCHEIN®

40% phosphoric dentin and enamel etch gel. Viscous solution stays in place once applied and does not run. Vibrant dark blue color, easy-to-use dispensing syringe, and precision syringe tip make it easy to place. Rinses clean with zero residue left behind.

40% Etch Gel

Standard Package

(570-2002)..... Ea.

Contains: 4–1.2 mL syringes & 20 tips.

Economy Package

(570-1998)..... Ea.

Contains: 12–1.2 mL syringes & 60 tips.

12 Gm. Bulk Syringe

(570-2003)..... Ea.

Contains: 1–12 Gm. syringe & 3 tips.

60 Gm. Jumbo Syringe Kit

(570-1999)..... Ea.

Contains: 1–60 Gm. syringe, 5–3cc empty refill syringes & 50 tips.

60 Gm. Jumbo Syringe Refill

(570-2000)..... Ea.



GUM® Soft-Picks®

SUNSTAR AMERICAS, INC.

Innovative, comfortable, and easy-to-use between-teeth cleaners for the on-the-go removal of dental plaque and food particles from between the teeth and around crowns, bridges, implants, dentures, and orthodontic appliances. Flexible and tapered synthetic rubber bristles fit easily between most teeth, removing dental plaque where gum disease often starts and toothbrushes alone cannot reach.

Tapered

5 per Pkg..... 72 Pkg. per Box

Specify:

Original (712-0176) Wider Spaces (712-0171)

DENTAL TOP SELLERS



Surface Barriers

HENRY SCHEIN®

Barrier Sleeves for Operating Lights, T-Style 4"W x 5¾"L

Allow a custom fit for most T-style light handles. Design keeps sleeves from falling off while in use, yet allows them to remove easily when finished. Clear.

(570-1316).....500 per Box

Barrier Sleeves for Air/Water Syringes

2½"W x 10"L

Provide protection for air/water syringes. Clear.

(570-1320).....500 per Box

Air/Water Syringe Covers

2½"W x 10"L

Provide protection for air/water syringes with a unique precut hole, allowing for easier application. Clear.

(104-7423).....500 per Box

Barriers for Dental Chairs

27½"W x 24"L

Feature a universal, tapered design that conforms neatly to the chair.

(570-1933).....225 per Box

Universal Adhesive Barrier Film

4"W x 6"L

Provides protection for all areas and surfaces. Unique backing adheres when placed and leaves no residue when removed. Perforated sheets.

Sheets.....1200 per Roll

Specify:
Clear (570-1315) Blue (570-1318)

Curing Light Sleeves

2"W x 9"L

Tapered from front to maximum width: ¾" x 1½"—2".

(112-7093).....100 per Box

Tray Barriers

10½"W x 14"L

Provide complete protection for flat or divided trays. Size B. Clear.

(570-1321).....500 per Box

Barrier Sleeves For X-Ray Tubes

15"W x 26"L

Universal size that provides complete protection for the entire X-ray tube head. Perforated sheets easily dispense from roll. Clear.

(570-1319).....250 per Roll



Impregum™ Penta™

3M

Precision polyether impression material for automatic mixing and dispensing. Accurate and hydrophilic, resulting in outstanding detail, even in a moist environment.

For Use in Pentamix 3 Mixing Unit

Refill Package—Medium Body #31642

(378-4898).....Ea.

Contains: 4 polybags (2 of each: 300 mL base & 60 mL catalyst).



TF Adaptive Nickel-Titanium Files

KERR

TF Adaptive Files: Small Canals

4 per Pkg.

Size	23mm	27mm
SM1	123-2212	123-2213
SM2	123-2220	123-2221
SM3	123-2218	123-2219

TF Adaptive Files: Medium/Large Canals

4 per Pkg.

Size	23mm	27mm
ML1	123-2226	123-2227
ML2	123-2222	123-2223
ML3	123-2229	123-2230

TF Adaptive Files: Procedure Pack

6 per Pkg.

Size	23mm	27mm
Small	123-2263	123-2264
Med/Large	123-2265	123-2266



Multilink® Automix Next Generation Universal Resin Cement

IVOCLAR

Dual-cure, self-etch, universal resin cement. Provides the security of high bond strength with a fast, simple application technique. Automix delivery provides accurate application with no mixing. Indicated for cementing metal, metal-free, and zirconia crowns as well as inlays, onlays, bridges, and posts.

NG System Package.....Ea.

Specify:

Transparent (945-0462) White (945-0460)
Yellow (945-0463) Opaque (945-0464)

Contains: 1–9 Gm. syringe, 1–3 Gm. bottle of each primer (A & B), 1–5 mL bottle of Monobond Plus, 1–2.5 Gm. liquid strip, 5 liquid strip application cannulas, 15 tapered mixing tips, 5 root-canal mixing tips, 50 microbrushes (25 small & 25 large), 1 mixing plate, 2 Optra Sticks & 1 mixing pad.

NG Refill Syringe.....Ea.

Specify:

Transparent (945-0466) White (945-0465)
Yellow (945-0467) Opaque (945-0468)

Contains: 1–9 Gm. syringe & 15 automix mixing tips.

Try-in Paste Syringe.....1.7 Gm.

Specify:

Transparent (945-0469) White (945-0472)
Yellow (945-0470) Opaque (945-0471)

Tapered Mixing Tips

(945-0473).....15 per Pkg.

Root Canal Tips

(945-0476).....5 per Pkg.

Primer Kit

(945-0474).....Ea.

Contains: 1–3 Gm. bottle of each primer (A & B).

Primer Refill.....3 Gm.

Specify:

A (945-0477) B (945-0479)

Monobond Plus

(945-0200).....5 Gm.

Applicator Refills

(945-7157)50 per Pkg.

**Look for more
Top Sellers in our
upcoming flyers!**

3M™ Filtek™ Easy Match Universal Restorative

4+1 Buy 4 Refills,
Get 1 Refill!*

1+1 Buy 1 Intro Kit,
Get 1 Refill!*

Offers Valid 7/1/24 – 9/30/24. Redeem by 10/10/24.
Redeem to OralCarePromos@solventum.com.

3 shades. 1 trusted technology.

With just 3 shades, 3M™ Filtek™ Easy Match Universal Restorative is designed to streamline single-shade restorations from the start. Just observe and select from three shades – Bright, Natural, and Warm – for an excellent match to almost any patient's tooth shade, anterior or posterior.

- ✓ Intuitive shade selection
- ✓ Naturally-adaptive opacity
- ✓ Reduced inventory



[HenryScheinDental.com/3M](https://henryschein.com/3M)

Ask your Henry Schein Dental Sales Consultant for complete ordering information.

*Bonus goods must be equal or lesser value of purchased product. To obtain bonus goods, you can redeem by email. Proof of purchase on one Henry Schein Dental invoice dated July 1, 2024 – September 30, 2024 is required (multiple invoices are accepted from a single purchase order when order cannot ship complete). Redeem by October 10, 2024. Redeem by email: OralCarePromos@solventum.com. Specify bonus product selections and include your business phone. Allow 3–5 weeks for delivery. Offers valid in the United States. Offer void if purchased product is returned. Solventum and Clinpro are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. © Solventum 2024.

Capsule refill includes:
20 – 0.2g capsules

(777-2552)	Bright
(777-2553)	Natural
(777-2555)	Warm

Capsule intro kit includes:
60 – 0.2g capsules (20 each shade)

(777-2560)	Assorted
------------	----------

Syringe refill includes:
1 – 4g syringe

(777-2556)	Bright
(777-2558)	Natural
(777-2559)	Warm

Syringe intro kit includes:
3 – 4g syringes (1 each shade)

(777-2562)	Assorted
------------	----------

GC FujiCEM® Evolve

Resin Modified Glass Ionomer Cement



Scan for a Lunch & Learn!

The solution is evolution

Building on the legacy of GC FujiCEM® 2, GC FujiCEM® Evolve provides more value for the same price!¹

- Superior bond strength to zirconia than a leading competitor²
- New tack-cure feature reduces cleanup time to seconds
- High radiopacity for easy visualization
- Ideal for cementation of zirconia, PFM and lithium disilicate restorations
- Rechargeable fluoride release and moisture tolerance

Try new GC FujiCEM® Evolve and join the evolution.

(333-5102) GC FujiCEM® Evolve Triple Automix Package:

Contains: Three 9.2g / 5.0mL syringes and 45 GC Push and Click Tips (Regular).

(333-5101) GC FujiCEM® Evolve Triple Handmix Package:

Contains: Three 9.2g / 5.0mL syringes and one mixing pad.

(333-5100) GC FujiCEM® Evolve Single Package:

Contains: One 9.2g / 5.0mL syringe and 15 GC Push and Click Tips (Regular).

(333-5103) GC Push and Click Tips (Regular x 15)

(333-5099) GC Push and Click Tips (Endo x 15)



www.gc.dental/america
www.gc.dental/america/training • © 2024 GC America Inc.

¹Same price per application for 'Automix' version and per gram for 'Handmix' version as GC FujiCEM® 2. ²Fujimoto et al, Bonding Durability of Novel Resin-Modified Glass-Ionomer Luting Cement to Zirconia, J Dent Res 98 (Spec Iss A): #1315, 2019

UNIVERSAL COMPOSITE

Beautiful II

Nano-hybrid universal composite with smooth, "non-sticky" consistency provides aesthetics and durability.



FLOWABLE COMPOSITE

Beautiful Flow Plus

Nano-hybrid, high strength flowable composite provides aesthetics and durability. Available in two viscosities, zero flow (F00) and low flow (F03).



ABRASIVE POLISHERS

OneGloss PS

One-step composite finisher and polisher comes individually packaged to prevent cross-contamination. Apply firm pressure to remove composite and soft pressure to polish to a "satin finish".



Super-Snap Singles

Award-winning disk polishing system for all composite restoratives comes individually packaged to prevent cross-contamination. Designed without a metal hub. Each kit contains one disk of each grit and a plastic mandrel. Available in Standard (12mm) and Mini (8mm).




SECONDARY DECAY IS THE MAIN CAUSE OF RESTORATIVE FAILURE

**GREATLY REDUCE
SECONDARY DECAY USING
SHOFU'S CLINICALLY
PROVEN BIOACTIVE
GIOMER RESTORATIVES**



Bioactive Giomer Technology³



Conventional Composites



Beautiful Flow Plus X & Super-Snap | Paiman Lalla, DDS, AEGD, FICOI, FAAIP

SINGLE USE SOLUTIONS



**Anti-Bacterial
Protection**



**Prevents Secondary
Decay²**



**Acid
Neutralization**



**Up to 30% Savings
Over Other Brands**

8-Year independent clinical study recorded 100% retention rate, no secondary caries, failures, or postoperative sensitivity.²

BUY 4, GET 1 (A2/F00 A2) AT NO CHARGE!*

*Restoratives only. Must be of like item. Mix & Match composite shades/viscosities. Email this page with a copy of your Henry Schein Dental invoice to Customer-Service@shofu.com or Fax to 760-736-3276. Bonus products will be shipped to you directly from Shofu Dental Corporation, San Marcos, CA. Allow 6-8 weeks for delivery.



**BIOACTIVE.
GIOMER TECHNOLOGY.
THE BETTER CHOICE.**

Visit shofu.com

Shofu Dental Corporation San Marcos CA 92078



¹Mjor IA, Moorhead JE, Dahl JE. Reasons for replacement of restorations in permanent teeth in general dental practice. Int Dent J. 2000 Dec;50(6):361-6. doi:10.1111/j.1875-595x.2000.tb00569.x. PMID: 11197194. ²Gordan VV, Mondragon E, Watson RE, et al. A clinical evaluation of a self-etching primer and a giomer restorative material: Results at eight years. JADA. 2007;138(5):621-627. ³Asahi University School of Dentistry.



Filtek[™] Supreme Ultra Universal Restorative

3M

Universal nanocomposite. Improved polish retention and enhanced fluorescence make it even easier to create natural-looking aesthetics. Deliver optimized aesthetic properties for both single-shade and multishade restorations. Use for anterior and posterior restorations, core buildups, splinting, and indirect restorations, including inlays/onlays and veneers.

Capsule Refills #6029

0.20 Gm.....20 per Box

Specify:

A2B (777-0136) A3B (777-0137)



Omnichroma[®]

TOKUYAMA DENTAL AMERICA

Single-shade universal composite. Spherical fillers, combined with Smart Chromatic Technology, create structural color matching of surrounding tooth color. Minimal wear with high compressive and flexural strength. Stain resistant to maintain cured color with glossy results. Opaque white before curing, blending naturally after curing. Blocker material used as a lingual layer in large Class III and IV cases.

PLT Refills—0.20 Gm.

(721-1022).....20 per Pkg.



GC Fuji Automix LC

GC AMERICA

Resin-reinforced, light-cured glass ionomer. Bioactive material with high fluoride release. High bond strength, even in the presence of saliva. Virtually no postoperative sensitivity. No etchant or adhesive required. Radiopaque.

Cavity Conditioner

(333-7058).....6 Gm.



OptiBond[®] Solo Plus[™]

KERR

15% filled, single-component, light-cure adhesive. Optimal particle size allows it to be used for direct and indirect applications, providing the ultimate in convenience, protection, and strength.

Bottle Refill

(123-5690).....5 mL

Unidose[®] Refills

(123-3029).....100 per Pkg.



TPH Spectra[®] ST

DENTSPLY SIRONA

Universal composite restorative material. Filler technology, Sphere TEC[™], has advanced granulated spherical filler that works together with the resin matrix to enhance handling and aesthetics. Simplified 5-shade system covers the entire VITA[®] shade range. Fast and easy to polish.

Compules[®] Tip Refills

0.25 Gm.....20 per Box

Specify:

A2 (222-2227)

A2 (222-2233)



EQUIA Forte[®]

GC AMERICA

Bulk-fill, self-adhesive posterior restorative system with new glass hybrid technology. Improved flexural strength and wear resistant. Ideal for Class I, II, and V restorations. Optimal marginal seal offers long-term resistance to microleakage and discoloration. High fluoride release at tooth restorative interface.

EQUIA Forte Fill

Capsule Refill

A2

(333-3687).....48 per Pkg.

Contains: 48–0.10 mL capsules.



Filtek[™] Supreme Flowable Restorative

3M

Flowable composite with unique nanotechnology that offers enhanced aesthetics and high wear resistance. Updated ergonomic syringe designed for improved comfort and easier injection of material. Use in anterior and posterior restorations, including liner/base and all Classes I–V. Flow-on-demand handling, excellent polish retention, and improved fluorescence.

Syringe Refills #6032

2 Gm.

A2

(777-2467).....2 per Pkg.

Contains: 2–2 Gm. syringes & 20 dispensing tips.



GC Fuji II LC[®] Capsule Refills

GC AMERICA

Restorative in premeasured, unit-dose capsules. Ultraconvenient. Fewer steps reduces chair time.

Capsule Refills

A2

(333-3464).....48 per Box



Scotchbond[™] Universal Adhesive

3M

A combined total-etch, self-etch, and selective-etch adhesive. Consistent bond strength to both moist and dry etched dentin. Single-bottle system for both direct and indirect restorations. Provides consistent, reliable bond strengths with virtually no postoperative sensitivity.

Refill Pack—#41258

(378-0349).....5 mL



MaxiTemp® HP



Automix, temporary crown-and-bridge material for creating temporary crowns, bridges, inlays, onlays, and veneers. High-performance, rounded cartridge prevents contamination at the tip and requires less force for easier dispensing of material. 4:1 mix ratio.

Complete Kit.....Ea.

Specify:

A1	(900-1284)	A3.5	(900-1282)
A2	(900-1281)	B1	(900-4090)
A3	(900-2126)	Bleach	(900-4091)

Contains: 1-50 mL cartridge & 6 mixing tips.

Cartridge Refills..... 2 per Pkg.

Specify:

A1	(900-1283)	A3	(900-2100)
A2	(900-1279)	A3.5	(900-1280)

Maxima® HP Mixing Tips—Light Blue

(112-6865)..... 50 per Pkg.

DS50 Dispensing Gun with 4:1/10:1 Slide

(112-6776)..... Ea.



FluoroCore® 2+

DENTSPLY SIRONA

Dual-cure, core build-up material plus fast cement for endodontic posts. High depth of cure and strong adhesion to dentin. Radiopaque.

Syringe Refill.....Ea.

Specify:

Tooth Colored	(222-0239)	Blue	(222-0241)
---------------	------------	------	------------

Contains: 4-4.75 Gm. syringes, 25 mixing tips, 10 endo intraoral tips & 25 intraoral tips.

Mixing Tips

(222-0242).....25 per Pkg.

Intraoral Tips

(222-0243).....25 per Pkg.

Endo Intraoral Tips

(222-0250).....30 per Pkg.

Root Canal Applicator Tips

(222-0253).....50 per Box



Protemp® Plus Temporization Material

3M

Automix bis-acrylic temporary material. Use on crowns, bridges, inlays, onlays, veneers, and long-term temporization. High fracture resistance and toughness. Smooth, glossy surface requires no polishing or glaze. Sophisticated filler technology provides strength and aesthetics. Easy plaque removal to support gingival health and patient comfort. 4:1/10:1 ratio.

Introductory Kit #46962

(777-0090)..... Ea.

Contains: 1-67 Gm. cartridge of shade A2, 16 blue Garant mixing tips, 1-16 Gm. catalyst of Relyx Temp NE, 1-36 Gm. base of Relyx Temp NE, 1 mixing pad, 1-2 Gm. syringe of Filtek Supreme Ultra A2 with tips, 1-4:1/10:1 Garant dispenser & technique cards.

Refill Package..... Ea.

Specify:

A1	(378-0081)	A3.5	(378-0281)
A2	(378-0086)	B3	(378-0091)
A3	(378-0090)	Bleach	(378-0080)

Contains: 1-67 Gm. cartridge & 16 blue Garant mixing tips.

Garant Mixing Tips—Blue

(378-1461)..... 50 per Pkg.

Garant Dispenser—4:1/10:1 Ratio

(378-7744)..... Ea.



Luxatemp® Ultra DMG AMERICA

Temporary provisional material that sets a new standard for flexural strength and beauty. Provides unmatched stability and long-term durability for temporary crowns, bridges, and more. Highly accurate with incredible marginal adaptation. Superior aesthetics. Can only be used with Smartmix® Tips.

Cartridge Refill—10:1 Ratio

Refill.....Ea.

Specify:

A1	(290-0042)	A3.5	(290-0045)
A2	(290-0043)	B1	(290-0046)
A3	(290-0044)	Bleach Light	(290-0047)

Contains: 1-76 Gm. cartridge & 15 tips.

Cartridge Mixing Tips—Fine

(290-6663).....25 per Pkg.

Dispensing Gun—Type 50—4:1/10:1

(290-9411)..... Ea.

Smartmix®

Handheld Syringe.....Ea.

Specify:

A2	(290-0048)	Bleach Light	(290-0050)
B1	(290-0049)		

Contains: 1-15 Gm. syringe & 10 Smartmix® mixing tips.

Smartmix® Tips—Green-Blue

(290-0114)..... 25 per Pkg.



Integrity® Temporary Crown and Bridge Material DENTSPLY SIRONA

Twin-port cartridge and dispensing system for fast, strong, aesthetic temporaries. Provides perfect mix ratio with no cross contamination. Minimal shrinkage, great consistency. Polishes to a natural luster. With fluorescence.

Cartridge Complete Package—10:1 Ratio

Shade A2

(222-7850).....Ea.

Contains: 1-76 Gm. cartridge shade A2, 20 mixing tips & 1 dispenser with slide.

Cartridge Refill Package..... Ea.

Specify:

A1	(222-6781)	B1	(222-8079)
A2	(222-9568)	Bleach White	(222-5263)
A3	(222-0774)	C2	(222-0780)
A3.5	(222-8736)		

Contains: 1-76 Gm. cartridge & 20 mixing tips.

Cartridge Mixing Tips

(222-4922).....20 per Box

Dispenser and Slide

(222-0701)..... Ea.

Handheld Automix Syringe

Refill Package..... Ea.

Specify:

A1	(222-9440)	B1	(222-3623)
A2	(222-4062)	Bleach White	(222-1479)
A3.5	(222-3846)		

Contains: 1-15 Gm. mini-syringe & 10 mini-syringe mixing tips.

Mini-Syringe Mixing Tips

(222-4206).....20 per Pkg.



Luxatemp® Automix Plus DMG AMERICA

Create durable and long-lasting provisional restorations with this superior bis-acryl material. Luxatemp is highly accurate with excellent marginal adaptation. Dispensed in the quick and easy-to-use automix system, preventing cross contamination. 10:1 ratio. Can only be used with Smartmix® Tips.

Cartridge Introductory Kit

(290-2523).....Ea.

Contains: 1-76 Gm. cartridge of Luxatemp Plus (shade A2), 15 fine Luxatemp Plus mixing tips & 1 Luxatemp Plus dispensing gun.

Cartridge Refill..... Ea.

Specify:

A1	(290-8937)	B1	(290-2279)
A2	(290-1628)	Bleach (Light)	(290-7176)
A3.5	(290-9180)		

Contains: 1-76 Gm. cartridge & 15 mixing tips.

Cartridge Bulk Package—A2

(290-0514)..... Ea.

Contains: 5-76 Gm. cartridges & 75 fine Luxatemp Plus mixing tips.

Cartridge Mixing Tips—Fine

(290-6663).....25 per Pkg.



Unprecedented Power


NEW Ti-Max Z Series Air-Driven Handpieces

Most Powerful at 44W*

* Compared with the same handpiece class offered by other major manufacturers, as of July 2023

DYNAMIC POWER SYSTEM

Power comparison with competitive products
(in-house research)**



**Air supply pressure specified by each company
*** Air pressure at 0.42 MPa

Ti-Max Z990L
L (786-0345)
KL (786-0347)
WL (786-0351)

Ti-Max Z890L
L (786-0346)
KL (786-0349)
WL (786-0350)

Ti-Max Z micro
L (786-0348)
KL (786-0352)



Z micro



NSK America Corporation • www.nskdental.com
1800 Global Parkway, Hoffman Estates, IL 60192 • Fax: +1 (800) 838-9328

HANDPIECES



**UNSURPASSED
QUALITY & SELECTION
AT YOUR FINGERTIPS.**

Shop our trusted, complete range of tools which offers savings that are good for the whole practice.

- High-Speed Air-Driven
- Low-Speed Air-Driven
- Electric
- Hygiene
- Replacement Parts



EXCLUSIVE BRANDS
& PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS



ProRepair The Leaders in Handpiece Repair NOW Digital

TRY Our Easy-to-Use Online Portal

See How Easy it is to Use our Online Portal

www.prorepaironline.com



FAST

turnaround time



USER-FRIENDLY



TRANSPARENT

*pricing upfront and tracking
for all repairs*

Easily Submit Repairs in 3 Easy Steps!

1

Log in

2

Choose the handpiece you
would like to repair

3

Add to cart checkout
and print UPS label



Scan for Instructions
to Submit Your
Handpiece

 **HENRY SCHEIN®**
PROREPAIR
Handpiece and Small Equipment Solutions



For more information, contact Henry Schein ProRepair:
email prorepair@henryschein.com or call
800-367-3674 Option 1, Option 1

©2024 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors.



ACCLEAN® Hygiene Handpiece 2



- Stainless steel construction for increased durability and longevity
- Improved torque for more responsive and smoother operation throughout operating speed range
- Designed to be both lightweight and powerful
- Simple-to-change color bands for easy identification and organization
- Quick-disconnect model for easy changing
- No-lube motor for simple, easy maintenance
- Autoclavable
- Compatible with all standard DPAs
- 1-year warranty

ACCLEAN® Hygiene Handpiece 2
(570-2861)..... Ea.

ACCLEAN® Swivel Hygiene Handpiece 2
(570-2855).....Ea.



Nupro RDH™ Hygiene Handpiece DENTSPLY SIRONA

Ergonomic shape was designed for a more natural hand positioning. 360° swivel and nonslip grip. Removable, autoclavable bands to personalize your handpiece. 2-year warranty.

Handpiece for Disposable Angles
Single Pack
(772-1322)..... Ea.

Three Pack
(772-3557)..... Ea.

Handpiece with Propy Right Angle
Single Pack
(772-6041)..... Ea.

Three Pack
(772-7371)..... Ea.



Prophy Star® 3 Handpiece STAR DENTAL

At 65 Gm., it is lighter than any other prophy handpiece, which reduces hand and wrist fatigue. Convenient speed adjustment: 0–5000 rpm in both forward and reverse directions. Steam sterilizable in temperatures up to 135°C. Air flow of 2.1 SCFM at maximal power of 40 psi. 360° nose cone swivel is longer, giving the user better swivel capabilities. Improved click-on and removal with no sticking. Includes 4 autoclavable, ergonomic rubber grips in assorted colors. Lube-free, no maintenance required. 2-year warranty.

Handpiece Less Swivel
(808-9698).....Ea.
Contains: 1 handpiece & 4 rubber replacement grips.



iProphy Air® Propylaxis Air Motor NSK AMERICA

This lightweight, ergonomic, and quiet hygiene handpiece features a slim, short design that minimizes hand fatigue and provides good balance. Swivels 360° to provide improved maneuverability during polishing procedures. Unique 2-part construction allows for nose cone and motor separation for easy maintenance and sterilization. Compatible with all standard Doriot-style disposable or metal angles, and connects directly to any 4-hole tubing. 1-year warranty.

Complete Unit
(786-0235)..... Ea.

Contains: 1 of each: motor, nose cone & set of 4 color rings.

Nosecone
(786-0237)..... Ea.



ProphyPal™ Hygiene Handpiece DENTICATOR

Specifically designed for disposable prophy angles. Features a 360° swivel and Sure-fit nose cone that allows for maximum control and easy transfer between buccal and lingual surfaces. Lightweight and balanced construction provides excellent ergonomic positions for all hand sizes. Motor operates at optimal speed and torque to provide tough polishing on heavy stains without stalling. Its “whisper quiet” air-powered operation reduces patient stress. Autoclavable up to 275°F/135°C. Made in the USA. 2-year warranty.

Complete Handpiece.....Ea.

Specify:
Pink (136-0007)

Silver (136-0006)

Wild Green (136-0008)

Contains: 1 of each: handpiece & lubricant.



PROPHYwiz Hygiene Handpiece KAVO

Ergonomic hygiene handpiece is the lightest, quietest, and smoothest-running handpiece available. Features 360° swivel that eliminates hand fatigue during polishing. Well-balanced to provide ultimate comfort for all-day use. Excellent torque with speed range up to 3,000 rpm for optimal polishing speed. No coupler needed, the PROPHYwiz adapts to any 4-hole standard tubing and can be used with both metal and disposable prophy angles. Easy to maintain and sterilize. 1-year warranty.

Single Pack
(628-6937)..... Ea.

3-Pack
(628-8901)..... 3 per Pkg.



Nupro Freedom® Hygiene Handpiece DENTSPLY SIRONA

Lightweight, cordless prophy hygiene handpiece system with SmartMode™ technology. Eliminates cord drag, giving the handpiece a lightweight feel. SmartMode™ technology provides the latest advance in handpiece speed control without the need for buttons or foot pedals. Handpiece provides 25% more speed and power and all-day battery life with rapid 15-minute recharge. Quieter than either air-driven hygiene or low-speed handpieces, providing a more pleasant polishing experience for both the patient and clinician. 1-year warranty.

System Package with Foot Pedal
(772-1111)..... Ea.

Contains: 1 inner module, 3 autoclavable metal outer sheaths, 1 autoclavable handpiece cradle, 1 charging base, 1 wireless foot pedal, 1 power supply, 1 box of 500 Disposa-Shield™ barriers & 200 Nupro Freedom™ disposable prophy angles.

*Wireless Foot Pedal sold separately.

ORAL HEALTH & ALL THAT JAZZ



Last Chance for New Orleans!

ADHA's Annual Conference is around the corner. Don't miss your chance to get inspired and unite with your dental hygiene community in the Big Easy.



Register now at
adha2024.org

adha24
ANNUAL CONFERENCE
JULY 26 - 28 / NEW ORLEANS

SPRING
HEALTH PRODUCTS



American Manufacturing & Ingenuity
DECADES OF EXPERIENCE!

The Cure
TC-3

TC-3 LED Cordless
(138-0780)



Zirconia Kit #Zirc 5
(138-0435)



Single & Multi-Use
FG Diamond Burs
REGULAR DIAMONDS:
25-Pack
**LARGE HEAD
DIAMONDS:**
25-Pack
TURBO DIAMONDS:
10-Pack



LOW PROFILE HEAD
FOR GREATER ACCESS

Giving Care to the Caregiver

Soft Retraction. Non-Metallic. Superior Control & Safety.

Developed and supported by IDD experts globally, Armor Dental launches **Comfort Soft**, the first soft, non-metallic retractor for special healthcare needs patients.

- Individuals with IDD need different instruments for a better dental experience.
- Properly managing dental care fears of IDD patients and their caregivers.
- Improve your chairside techniques with better visibility, access and patient COMFORT.
- Prior to an office visit, then used routinely for home oral care.

Two Comfort Soft Versions:


Comfort Soft Single Patient Multi-Use: Use chairside then given to caregiver for home care. (Small - age 6 and under 145-4143, and Large - age 7 plus 145-4142)

Comfort Soft Autoclavable - Steam autoclavable up to 30 cycles for in-office use. (Small - age 6 and under 145-5979, and Large - age 7 plus 145-5978)

Modernizing special oral care needs for providers, caregivers and individuals.

ARMORDENTAL

Manufactured by Armor Dental
www.armor-dental.com



RELY ON US

FOR HIGH QUALITY **INSTRUMENTS**



Trust our unbeatable selection of high quality, high value **instruments** so you can focus on taking care of your patients.



EXCLUSIVE BRANDS
& PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS



EVEREDGE® 2.0 Scalers HUFRIEDYGROUP

Using an enhanced finishing process, Hu-Friedy's professional artisans are creating edges that are over 72% sharper than the next closest competitor. Less force means easier, more comfortable scaling for clinicians, which enables them to deliver the best clinical outcome every time.

#9 Handle.....Ea.

Specify:
H6/H7 (601-1862)
204S (601-1864)
Nevi 4 Posterior (601-1866)
Gracey 13/14 (601-1868)
Gracey 11/12 (601-1870)
Barnhart 5/6 (601-1872)
Gracey 1/2 (601-1878)
H5/33 (601-1882)
1 Nevi/H5 1 Sickle (601-1893)



Montana Jack® Hygiene Instruments PDT, INC.

Our smart methods of 100% made-in-the-USA manufacturing enable us to significantly improve the Gracey design in a way that practitioner can easily see and experience for themselves.

Double End.....Ea.

Specify:
Regular (731-1681)
Rigid (731-1639)



Double-Sided Mirror HUFRIEDYGROUP

Two mirror faces allow the clinician to maintain an indirect vision or indirect illumination for enhanced visibility, in any position. Thin-profile, streamlined design for less bulk. Unique retractor shank design for increased patient comfort.

6-Pack.....Ea.

Specify:
#4 (600-8627) #5 (600-1967)



EverEdge™ 2.0 Harmony HUFRIEDYGROUP


Same great EverEdge 2.0 working ends with the new Harmony handle. Designed with TrueFit™ technology, an advanced, sensor-based technology system that measures pinch force in the fingers and pressure applied to the tooth when scaling. Features a 30% longer silicone grip, a 67% deeper double-helix grip design with curved grooves, and an organic, rounded handle that easily adapts to all hand sizes.

Sickle Scaler.....Ea.


Specify:
204S (601-3074)
H6/H7 (601-3073)
Nevi 4 Posterior (601-3083)

IMS® Signature Series® Cassette HUFRIEDYGROUP

Electropolished stainless steel is durable and won't rust or stain. Rail design holds instruments ½" off bottom. Stainless steel latch.



5½" x 8" x 1¼"
Purple
(600-0179).....Ea.



8" x 11" x 1¼"
Large—16-Instrument Capacity. Blue.
(600-7201).....Ea.



Cone Socket Mirror Handle HUFRIEDYGROUP

#6 Satin Steel®
(600-4442).....Ea.



Harmony™ Ergonomic Scalers and Curettes HUFRIEDYGROUP

Designed with TrueFit™ Technology to provide real ergonomic benefits driven by scientific evidence. This advanced, sensor-based testing system measures over 2.8M data points for pinch force in the finger and pressure applied to the tooth when scaling. Features EverEdge 2.0 working ends that are 72% sharper than the next leading competitor. This allows clinicians to efficiently remove calculus while applying less pressure to the handle.

Harmony.....Ea.

Specify:
Gracey Curette DE 1/2 (601-3062)
Gracey Curette DE 11/12 (601-3066)
Gracey Curette DE 13/14 (601-3068)
Universal Barnhart DE 5/6 (601-3056)



IMS® Signature Series® Amalgam Kit HU-FRIEDY

Contains a comprehensive set of instruments for amalgam procedures, all in a Signature Series® cassette in a convenient kit. Features a compilation of clinicians' most popular procedural setups. The 16-instrument Signature Series® cassette is 100% stainless steel.

(600-2472).....Ea.

Contains: 1 of each: #6 Satin Steel® mirror handle, 3-pack of #5 front surface mouth mirrors, #18 excavator, calcium-hydroxide placement instrument, regular/large CFI amalgam carrier, 0/1 Marquette condenser, 1/2 black condenser, 3/6 discoid-cleoid carver, 89/92 discoid-cleoid carver, 1/2-3 Hollenback carver, 21B burnisher, 27/29 burnisher, #2 college pliers, Miller articulating paper forceps, CW aspirating anesthetic syringe, amalgam well, bur cushion with short lid (holds 12) & A/W syringe tip clip.



Disposable Mirrors

HENRY SCHEIN®

Plane surface, high-quality patient mirrors made of white, nonautoclavable ABS plastic.

(100-6640).....72 per Box

The path ahead needs leadership. We bring together people and organizations who can change oral health care for the better.

We partner with public and private stakeholders across the health continuum to create better oral health for all. Because good overall health depends on good oral health.

We're on the path to transforming oral health care. We hope you join us.

Learn more at carequest.org

CareQuest 
Institute for Oral Health®



NEW FROM CENTRIX!

SilverSense SDF™

38% Silver Diamine Fluoride Sensation!

More
Treatment Options.
More
Patient Choices.
**Healthier
Teeth!**



- Stops bacterial formation
- Paints on fast & easy
- No anesthesia, injections or drilling
- The tooth restores beautifully
- Best Value – less than \$0.80/patient!
(Applicators included.)



TRY SILVERSENSE SDF TODAY!

Order now from Henry Schein or visit centrixdental.com/silversense-sdf

5 mL Bottles, Applicators Included

1-Bottle Standard Kit & 20 Applicators (163-0196)
3-Bottle Clinic Kit & 50 Applicators (163-0197)

© 2024 Centrix, Inc. | #HSCP-Q3

centrix
MAKING DENTISTRY EASIER

MAKING DENTISTRY EASIER | CENTRIXDENTAL.COM

matrix material



Palodent® Plus Sectional Matrix System DENTSPLY SIRONA

Utilizes optimized ring, band, and wedge technology. Constantly provides tight, accurate contacts. WedgeGuard completes the system, providing adjacent tooth protection and a tight gingival seal to minimize flash and time spent on finishing. Nickel-titanium rings can be retensioned over time providing better ring retention on tooth. Autoclavable.

Matrices Refills.....50 per Pkg.

Specify:
4.5mm (222-0381) 5.5mm (222-0382)

Matrices Refills–5.5mm
(222-0378)..... 100 per Pkg.

Ring Refills.....2 per Pkg.

Specify:
Universal (222-0385) Narrow (222-0386)

Wedge Refills..... 100 per Pkg.

Specify:
Small (222-0391) Medium (222-0392)

3D Fusion™ Sectional Matrix System GARRISON DENTAL SOLUTIONS

Same great features as the 3D XR™ but with the added benefit of 3D Fusion Bands and 3D Fusion Wedges. Superior tooth separation force through advanced engineering and material selection. Soft-Face™ silicone ring tips with enhanced marginal ridge anatomy hug the pre-contoured dead-soft matrices to the tooth for flash free, contoured restorations. Composi-Tight® 3D Fusion™ features easy ring stacking, over-the-wedge placement and Ultra-Grip™ anti-ring slip extensions.

Starter Kit with Rally® Polishers

(300-0384).....Ea.

Contains: 3 rings (short, tall & wide), 70–3D Fusion™ Bands (20 of each: 4.4mm, 5.6mm & 6.6mm; & 5 of each: 6.0mm with extension & 8.7mm with extension), 80–3D Fusion™ Wedges (20 of each size: X-small, small, medium & large), 4 FenderWedges™ (1 of each: 1.0mm, 1.4mm, 1.9mm & 2.3mm), 12 Rally® Polishers (2 of each: coarse point, coarse cup, medium point, medium cup, fine point & fine cup) & 1–3D Fusion™ Ring Placement Forceps.

3D Fusion™ Wedges

Refills..... 50 per Pkg.

Specify:
X-Small–Yellow (300-0332)
Small–Blue (300-0331)
Medium–Orange (300-0339)

3D Fusion™ Ring

Refills
Short
(300-0373) 2 per Pkg.



Compare to Ultradent's Omni-Matrix®

Pro-Matrix® YOUNG DENTAL

Single-use. Band designed to perform all dental restorations using both amalgam and composite filling materials. The sleek, narrow design allows for better visibility during procedures. Screw stays stationary during tightening for better patient comfort. Lever points to upper and lower gum when applying band. Ready to use, no assembly required, and individually wrapped.

Straight Bands..... 50 per Box

Specify:
Narrow 4.5mm (247-0124)
Wide 6mm (247-0125)

Curved Bands.....50 per Box

Specify:
Narrow 4.5mm (140-7818)
Wide 6mm (140-7819)

ACCLEAN®

by  HENRY SCHEIN®

HYGIENE APPOINTMENTS **MADE EASIER**



NEW ACCLEAN® ULTRASONIC INSERTS

Lasts over 60% longer* providing efficient and effective procedures. Compatible with most 25k and 30k Ultrasonic Scaling Units.

ACCLEAN® Easy Grip

Universal Slim #50 25K	(570-1878)
Universal Slim #50 30K	(570-1879)
Triple Bend Slim #1000 25K	(570-1929)
Triple Bend Slim #1000 30K	(570-2001)
Universal Extra Thin 30K	(570-2238)

ACCLEAN® Easy Grip XL

Universal Slim #50 25K	(570-2016)
Universal Slim #50 30K	(570-2039)
Triple Bend Slim #1000 25K	(570-2138)
Triple Bend Slim #1000 30K	(570-2237)
Universal Extra Thin 30K	(570-2239)

ACCLEAN® Classic

Universal P-10 25K	(570-1624)
Universal P-10 30K	(570-1729)
Universal Slim P-100 25K	(570-1873)
Universal Slim P-100 30K	(570-1874)



ACCLEAN® ULTRASONIC SCALER (570-2771)

A new, purposeful design to provide a better experience for you and your patients.

* Wear Resistance Test, 2023. Data on File.



Everyday Preventives for Healthy Patients

www.HenrySchein.com/ACCLEAN

V I O M E

Early Cancer Detection Saves Lives

Oral Health Pro™

with CancerDetect® Oral & Throat

Elevate your practice and personalize patient treatments with comprehensive oral microbiome health insights and early cancer detection.

✓ This test detects biomarkers for oral and throat cancers at their earliest stages with a specificity of 95% and sensitivity of 90%.

✓ Oral Health Scores reveal insights into oral microbiome health and connections to deeper, systemic health concerns.

- Periodontal health
- Caries
- Gum inflammation
- Other risks of oral and systemic diseases

CLIA
CERTIFIED
LAB



Learn more at
viomepro.com



V VARDIS
SWITZERLAND

BIOMIMETIC DENTAL SCIENCE

No More "Watch and Wait" with Curodont™

Discover Curodont™ Repair Fluoride Plus: the unique brush-on system for the treatment of watch areas through Guided Enamel Remineralization. Maximize patient outcomes with Curodont™ Protect



Curodont™ Repair Fluoride Plus

Guided Enamel Remineralization of early caries lesions

- For all surfaces including interproximal
- Clinically proven and safe
- Does not stain
- 1x chairside application
- 1 applicator can treat up to 3 lesions
- Non-invasive, easy, and fast application
- Suitable for all patients, including children

Dental
PRODUCT
SHOPPER
BEST
PRODUCT

Dental
PRODUCT
READER'S
CHOICE
2024



(687-0057)

Curodont™ Protect

Anticavity & Antisensitivity dental gel

- Chairside & at home complement to the Curodont™ Repair Fluoride Plus treatment
- Helps to prevent dental cavities
- Relieves tooth sensitivity
- Topical application on all tooth surfaces
- For all patients, including children age 6 and up



(687-2868)

www.vwardis.com



Curodont™ Repair Fluoride Plus **VARDIS**

Locally applied professional mouth rinse for the in-depth anticavity treatment of white spots. Dentist or hygienist applies this drill-free, chairside product directly to the affected site to reverse beginning lesions; applied to the white spot, it guides the remineralization. It has the highest concentration of the X-Formula Curolox®, and it strengthens and remineralizes the entire enamel layer. Single-use applicators. (687-0057)..... 10 per Pkg.

SilverSense SDF™ **CENTRIX**

38% silver diamine fluoride used for treating dental hypersensitivity. Solution helps deliver better patient care by alleviating hypersensitivity, discomfort, and pain and providing benefits that will allow patients to comfortably wait for a restorative treatment. Noninvasive procedure means no drilling and anesthetic. Stops harmful bacteria and hardens and conditions softened dentin. Does not stain healthy enamel or dentin. For use in adults over the age of 21.



Standard Kit
(163-0196)..... Ea.
Contains: 1–5 mL bottle & 20 Benda Micro applicators.

Clinic Kit
(163-0197)..... Ea.
Contains: 3–5 mL bottles & 50 Benda Micro applicators.

Benda® Micro Applicators—Fine, Yellow
(163-7534).....576 per Box



Colgate® Whole-Mouth Clean Bundle **COLGATE ORAL PHARMACEUTICALS**

Includes a toothbrush with a 25% recycled plastic handle that significantly reduces gingival bleeding; a toothpaste with a breakthrough formula that fights bacteria on teeth, tongue, cheeks, and gums for whole-mouth health; and a floss developed to slide easily between teeth without shredding.* The Colgate® Whole-mouth Clean Bundle is a great option for patients who desire improved gum health.

- Colgate® 360® toothbrush
- Colgate® Total[®] Clean Mint toothpaste (0.88 oz.)
- Colgate® Dental Floss
- Patient sample bags

(543-3743)..... 72 per Box
*Compared to an unwaxed floss



Riva Star Aqua **SDI**

Patented, water-based, silver-fluoride (AgF)-based solution that is an antimicrobial, a biofilm eliminator, and a desensitizer. It is a noninvasive, 2-step system used to desensitize tooth pain in a few drops. This aqueous AgF solution, with an ammonia-free base, offers immediate sensitivity relief, reduces staining, and provides a 2-year sensitivity effect. It also has no risk of tissue burn, requires no gingival barrier, and has an improved scent.

10-Patient Kit
(164-0350)..... Ea.
Contains: 20 single-dose capsules (10 of each: Riva Star Step 1 silver fluoride & Riva Star Step 2 potassium iodide, green bottle) & 20 color-coded application brushes.



2-Step Bottle Kit
(164-0351)..... Ea.
Contains: 1 of each bottle: 1.5-mL Riva Star Aqua Step 1 & 3-mL Riva Star Step 2.



Clinpro™ Sealant **3M**

The first sealant with smart color-change technology—it goes on pink for easy-to-see application, and cures to a natural white. Fluoride releasing. Low viscosity flows easily into pits and fissures for less waste and less need for occlusal adjustment. Available in bottle or fine tip syringe.

Introductory Bottle Kit
(777-3665)..... Ea.
Contains: 2–6 mL bottles sealant, 9 mL etchant, 2 brush handles, 60 brush tips, 1 mixing well & instructions.

Refill Bottle
(777-9284)..... 6 mL

Introductory Syringe Kit
(777-7884)..... Ea.
Contains: 2–1.2 mL syringes sealant, 20 sealant dispensing tips, 1–3 mL syringe etchant & 25 etchant dispensing tips.

Refill Syringe
(777-0371)..... Ea.
Contains: 1–1.2 mL syringe sealant, 10 dispensing tips & instructions.

Sealant Dispensing Tips
(777-0760)..... 10 per Box



ProFluorid® Varnish **VOCO**

5% sodium-fluoride varnish delivers with enhanced flow characteristics for hard-to-reach areas. Single-dose delivery system adheres to moist surface and sets quickly after contact with saliva. White, transparent varnish without yellow discoloration immediately releases fluoride to relieve hypersensitivity. Contains no tree nuts, peanuts, corn, shellfish, eggs, milk protein, soy, gluten, triclosan, petroleum, red dye/artificial coloring, saccharin, or aspartame.

Single Doses

0.25 mL—Kids—Melon
(840-9282).....50 per Box

0.4 mL—Assorted Flavors
(999-1366)..... 48 per Box
Contains: Bubble Gum, Caramel, Cherry, Cola Lime, Melon & Mint.

0.4 mL..... 50 per Box
Specify:
Bubble Gum (999-2292) Melon (840-9285)
Caramel (999-1362) Mint (999-1364)
Cherry (999-1363) Piña Colada (999-2819)
Cola Lime (999-2713)

0.4 mL..... 200 per Box
Specify:
Bubble Gum (999-2293) Melon (999-1370)
Caramel (999-1367) Mint (999-1368)
Cherry (999-1369) Piña Colada (999-2820)
Cola Lime (999-2710)


10 mL Tube..... Ea.
Specify:
Caramel (999-1669) Melon (114-3932)
Cherry (999-1670) Mint (999-1671)
Cola Lime (999-2712) Piña Colada (999-2821)



Colgate® Dental Floss **COLGATE ORAL PHARMACEUTICALS**

Polypropylene fiber with a microcrystalline coating that slides easily between teeth without shredding, compared to an unwaxed floss. Container made with 25% recycled plastic. Mint flavor to help freshen your mouth.

Trial Size—3 Yd.
(543-0252).....72 per Box



XCP-DS FIT® Universal Sensor Holder System DENTSPLY RINN

Universal design: the self-adjusting clip stretches to fit any sensor. One holder for both Size 1 and Size 2 sensors. Silicone-cushioned edges help provide patient comfort and prevent sensor cover perforations. Use with either XCP-DS or XCP-ORA arms and rings. Autoclavable—low cost per use.

XCP-DS FIT® Complete Kit #559909
(189-0484)..... Ea.
Contains: 1 XCP-ORA arm, 1 XCP-ORA ring, 1 XCP endo arm, 1 XCP endo ring, 10 XCP-DS FIT bite blocks (2 of each: anterior, posterior, horizontal bitewing, vertical bitewing & endo) & 50 sensor covers.

XCP-DS FIT® Hygiene Kit #559900
(189-0471)..... Ea.
Contains: 1 XCP-ORA arm, 1 XCP-ORA ring, 8 XCP-DS FIT bite blocks (2 of each: anterior, posterior, horizontal bitewing & vertical bitewing) & 50 sensor covers.

XCP-DS FIT® Hygiene Kit w/XCP #559910
(189-0715)..... Ea.
Contains: 4 XCP arms & rings, (1 anterior, 2 posterior & 1 bitewing), 8 XCP-DS FIT bite blocks (1 anterior, 1 anterior thin, 2 posterior, 2 horizontal bitewings & 2 vertical bitewings) & 50 sensor covers.

XCP-DS FIT® Endo Kit #559908
(189-0481)..... Ea.
Contains: 1 XCP-ORA arm, 1 XCP-ORA ring & 2 XCP-DS FIT endo bite blocks.

XCP-DS FIT® Biteblocks Refills.....2 per Pkg.
Specify:
Anterior—Blue #559901 (189-0472)
Anterior Thin—Blue #559906 (189-0708)
Posterior—Yellow #559902 (189-0473)
Horizontal Bitewing—Red #559903 (189-0474)
Vertical Bitewing—Red #559904 (189-0477)
Endo—Green #559905 (189-0480)
Size 0—Black #559907 (189-0713)

XCP-DS FIT® Biteblock Refills..... 10 per Box
Specify:
Anterior—Blue #559921 (189-0641)
Anterior Thin—Blue #559926 (189-0718)
Posterior—Yellow #559922 (189-0642)
Horizontal Bitewing—Red #559923 (189-0643)
Vertical Bitewing—Red #559924 (189-0644)

XCP-DS FIT® Replacement Elastic Bands..... 6 per Pkg.
Specify:
Short (Posterior, VBW, Endo, Size 0) #559913 (189-0565)
Long (Anterior, HBW) #559914 (189-0566)



Compare to Schick CDR® Holders Uni-Grip® 360 Universal Sensor Holder DENTSPLY RINN


Universal, disposable sensor holder with an easy rotating tab. Increased positioning range allows you to take all X-rays: anterior, posterior, and bitewings. Only one holder is needed to complete a full-mouth series. Four adhesive tabs prevent sensor from slipping, providing additional stability. Fits any #1- or #2-sized digital sensor, film, or phosphor plate. Bite-piece channels are color coded to ensure correct positioning.

Uni-Grip® 360 Combo Kit..... Ea.
Specify:
Uni-Grip® 360 Combo Kit (189-0464)
Digital Sensor (189-0116)
X-ray Film/Phosphor Plates
Contains: 1 universal arm, 1 universal ring, 100 holders & 100 covers.

Uni-Grip® 360 Kit.....50 per Box
Specify:
For Digital Sensor (189-0463)
For X-ray Film/Phosphor Plates (189-9156)
Contains: 50 holders & 50 covers.

Uni-Grip AR® Positioning Arm & Ring
(189-0119).....Ea.
Contains: 1 universal arm & 1 universal ring.

Universal Sensor Covers
2" x 10"—Clear
(189-0254)..... 500 per Box



DEXIS™ Platinum Sensor Holders DEXIS, LLC


Platinum Holder Kit
(114-8040)..... Ea.
Contains: 12 holders (2 posterior, 2 anterior, 4 bitewing & 4 handheld), 2 universal rings & 3 aiming bars.

Holders.....Ea.
Specify:
Anterior (114-0290) Horizontal Bitewing (114-1521)
Posterior (114-5308) Vertical Bitewing (114-8448)

Handheld Holders..... 2 per Pkg.
Specify:
Vertical (Anterior) (114-8305)
Horizontal (Posterior) (114-6170)

Endodontic Holder Set
(114-3011)..... Ea.
Contains: 4 holders, 1 ring & 1 endo aiming bar.

Sensor Sheaths—Universal
(234-0135)..... 100 per Pkg.



XCP®/BAI Replacement Parts DENTSPLY RINN

Rinn® XCP® Standard Arm..... Ea.
Specify:
#54-0857 Anterior—Blue Prongs (189-6909)
#54-0858 Posterior—Yellow Prongs (189-3611)
#54-0927 Bite Wing—Red Prongs (189-4247)
#55-0597 Endo—Green Prongs (189-0291)

Short Arm
Specially designed for use with handheld X-ray units.
Short Arm..... Ea.

Specify:
#54-0871 Anterior—Blue Prongs (189-0898)
#54-0872 Posterior—Yellow Prongs (189-0897)
#54-0873 Bite Wing—Red Prongs (189-0896)
#54-0874 Endo (189-0907)

Rinn® XCP® Aiming Ring..... Ea.
Specify:
#54-0865 Anterior—Blue (189-6357)
#54-0860 Posterior—Yellow (189-7712)
#54-0934 Bite Wing—Red (189-3918)
#55-0598 Endo—Green (189-8717)



XIOS® Intraoral Digital Sensor Holders DENTSPLY SIRONA

Disposable, self-adhesive "sticky-backed" sensor holders. More Control with the peel-and-stick adhesive as it firmly secures the sensor to the holder. Easy and accurate placement of sensors makes even the most difficult shots easy to obtain. Holders are color coded. Fits all types of sensors. Holders can be used with your existing rings and positioning arms. Tabs can be used with hemostats.

Sensor Holders..... 100 per Pkg.
Specify:
Anterior Holder—Blue (811-2762)
Posterior Holder—Yellow (811-4246)
Bitewing Holder—Red (811-5465)

Universal Bite Tab—Green
(811-3121)..... 100 per Pkg.
Endo Tab—Grey
(811-9012).....50 per Pkg.

Hygienic Sleeves
Sensor Sheaths..... 300 per Box
Specify:
Regular Plus
Size 1 (811-5847) Size 1 (811-0044)
Size 2 (811-8372) Size 2 (811-0047)



Preparedness, Prevention & Compliance



COMPLIANCE & TRAINING

- OSHA/HIPAA
- Workplace Safety
- Emergency Preparedness
- Cardiovascular, BBP & Hypertension



INFECTION PREVENTION

- Biological Spore Testing
- Water Treatment & Testing
- Air Purifiers

LEARN MORE

Contact your Henry Schein representative or visit the link below for product specs and additional information.

henryscheindental.com/healthfirst



ENVIRONMENTAL WASTE MANAGEMENT

- RX
- Sharps
- Amalgam
- Lead



MEDICAL EMERGENCY RESPONSE

- STAT KIT® Emergency Medical Kits
- Crash Cart Auto Replenishment Program
- Oxygen & Accessories

ROBUST PROTECTION IN STYLISH BLACK

Criterion™ Nitrile Gloves
Now Available in Black



- Chlorinated for **reduced friction** and **easy donning**
- Textured fingers provide **enhanced grip**
- Chemotherapy & fentanyl tested to ASTM D6978

CRITERION™ NITRILE EXAM N200 GLOVES—BLACK

X-Small	(570-3487)	Large	(570-3490)
Small	(570-3488)	X-Large	(570-3491)
Medium	(570-3489)		

CRITERION™ NITRILE EXAM N300 GLOVES—BLACK

X-Small	(570-3482)	Large	(570-3485)
Small	(570-3483)	X-Large	(570-3486)
Medium	(570-3484)		

MAKE IT A MATCHING SET

Criterion™ Earloop Masks Levels 1 & 3
are also available in stylish black

Level 1: (570-3272), Level 3: (570-3273)



ORDER TODAY AT [HENRYSCHEIN.COM/SM](https://www.henryschein.com/sm)

by **HENRY SCHEIN®**





CRITERION™ Nitrile N100



Formerly known as *Criterion® Nitrile*. Same Great

Product, New *Criterion®* Packaging

Nitrile Powder-Free Exam

Durable nitrile construction for reliable protection and tactility.

Chlorinated for reduced friction and easy donning. Textured fingers provide an enhanced grip. Chemotherapy tested to the ASTM D6978 standard. Beaded cuff. 3.15-mil thickness at the palm. Blue.

Gloves..... 100/Box

Specify:

XS (111-8533) L (111-8537)

S (111-8535) XL** (111-8539)

M (111-8536)

**90/Box



CRITERION™ Nitrile N200



Nitrile Powder-Free Exam

Slightly thinner than the N100, the durable nitrile construction provides reliable protection and tactility. Chlorinated for reduced friction and easy donning. Textured fingers provide an enhanced grip. Chemotherapy tested to the ASTM D6978 standard. Beaded cuff. 2.36-mil thickness at the palm. Available in both Periwinkle Blue and Sky Blue.

Periwinkle Blue.....200/Box

Specify:

XS (900-7437) L (900-7440)

S (900-7438) XL** (900-7441)

M (900-7439)

Sky Blue.....200/Box

Specify:

S (570-2098) L (570-2100)

M (570-2099) XL** (570-2101)

**180/Box



CRITERION™ Nitrile N300



Formerly known as *Criterion N300 Ultra Nitrile*.

Same Great Product, New *Criterion®* Packaging

Nitrile Powder-Free Exam

Ultrathin formulation for excellent tactile sensitivity and dexterity.

Chlorinated for reduced friction and easy donning. Textured fingers promote sure handling of instruments and materials. Chemotherapy tested to the ASTM D6978 standard. Ecofriendly packaging—with 300 gloves in a box—reduces waste and save space. Beaded cuff. 2.0-mil thickness at the palm. Ice Blue.

Ice Blue..... 300/Box

Specify:

XS (112-6829) L (112-6832)

S (112-6830) XL** (112-6833)

M (112-6831)

**250/Box



CRITERION™ Nitrile Low Derma



Formerly known as *Criterion Pure Freedom Nitrile*.

Same Great Product, New *Criterion®* Packaging

Nitrile Powder-Free Exam

Made from a soft, durable nitrile formulation without the use of chemical accelerators or sulfur, preventing Type-IV allergic reactions. Reduces the potential for dermatitis for glove wearers, yet still provides softness and comfort and maintains superior protection. Chlorinated for reduced friction and easy donning. Textured fingers promote sure handling of instruments and materials. Chemotherapy tested to the ASTM D6978 standard. Beaded cuff. 2.0-mil thickness at the palm. Slate Blue.

Gloves..... 200/Box

Specify:

XS (570-0630) L (570-0633)

S (570-0631) XL** (570-0634)

M (570-0632)

**180/Box



CRITERION™ Nitrile Aloe + Vitamin E



Formerly known as *Criterion Radiant Razz Nitrile*.

Same Great Product, New *Criterion®* Packaging

Nitrile Powder-Free Exam

Extra-soft stretch with an enhanced inner coating of aloe and vitamin E to soothe, moisturize, and rejuvenate hands. Fast and easy donning. Comfortable fit and feel to provide a superior range of motion. Textured fingers promote sure handling of instruments and materials. Beaded cuff. 1.97-mil thickness at the palm.

Magenta..... 200/Box

Specify:

XS (112-7143) M (112-7146)

S (112-7144) L (112-7147)

BEYOND PROTECTION

**LOW DERMATITIS
POTENTIAL**

**FENTANYL
RESISTANCE**

**CHEMOTHERAPY
DRUG TESTED**

**ACCELERATOR
FREE**



Transcend Nitrile Gloves 300/box

High tensile strength for improved barrier integrity and tear resistance
XS (565-0048) M (565-0050) *XL (565-0052)
S (565-0049) L (565-0051) *250/box



Carbon Air Nitrile Gloves 300/box

Distinct black color to minimize visible stains during use
XS (137-2693) M (137-2695) *XL (137-2697)
S (137-2694) L (137-2696) *250/box



Revo 200 Nitrile Gloves 200/box

Comfort formulation for unprecedented softness
XS (353-0013) M (353-0015) XL (353-0017)
S (353-0014) L (353-0016)



Contour Nitrile Gloves 100/box

Increased thickness for durability and less puncturing
XS (107-6664) M (107-2149) XL (107-2898)
S (107-0621) L (107-4112)



GUESS WHAT CAME IN WITH YOUR PATIENTS? YOU DESERVE THE RIGHT PROTECTION!



THE PERFECT PAIRING

PINK UNDERGUARD* gloves have a special smooth finish that makes it easy to slip on the **PURPLE NITRILE-XTRA*** glove. Used together, the benefits really add up:

- Double gloving is faster and easier
- Excellent permeation performance



52

Be USP <800> Compliant With CHEMO GEAR¹.
Both gloves now approved for **52 chemo drugs** per ASTM D6978.²

1 Meets the USP <800> guidelines for safe handling of hazardous drugs; refer to section 7.1 of USP <800> for complete Gloves guidelines
2 Standard Practice for Assessment of Medical Gloves to Permeation by Chemotherapy Drugs

For more information: Contact your local Henry Schein Medical Representative
or visit henryschein.com/medical.



HALYARD* is an Owens & Minor brand

*Registered Trademark or Trademarks of O&M Halyard or its affiliates. *PURPLE and PURPLE NITRILE are Registered Trademarks of O&M Halyard or its affiliates. The COLORS LAVENDER, and PURPLE are Registered Trademarks of O&M Halyard or its affiliates. The COLOR PINK is a Registered Trademark owned by Colur World, LLC. ©2005 Colur World LLC. ©2024. All rights reserved.



Together, we can build a **healthier world for all.**

Practice Green is designed to empower the health care community to positively impact the future of the planet by reducing the ecological footprint and promoting sustainability—both in daily life and while at work.



A product included in the Practice Green portfolio has been carefully selected by Henry Schein as “green” based on the product information provided by the supplier partner that it is either reusable, contains recycled materials, can be recycled, is biodegradable, or is from a sustainable source.



Scan the QR code to learn more, or visit
www.henryschein.com/PracticeGreen

© 2024 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. 24DS3152_M



Precision Eco Bio Face Masks

Level 3

PADM MEDICAL USA

Made from biopolymers that are produced from renewable crop resources.

- Hypoallergenic
- Absorbent, non-woven inner layer and hydrophobic outer layer
- Non-latex and graphene-free
- Bendable nose-clip wire for custom fit

Earloop Masks 50 per Box

Specify:

White (145-3980) Black (145-3981)



Soft White Nitrile

ANSELL

Nitrile Powder-Free Exam

Provides superior comfort, enhanced tactile sensitivity, and exceptional value. Textured fingertips and chlorinated. 2.6-mil thick. White.

Gloves 100/Box

Specify:

XS	(565-4510)
S	(565-8087)
M	(565-7431)
L	(565-9481)
XL	(565-1575)



Bio Nitrile®

CRANBERRY

Nitrile Powder-Free Exam

Features advanced biodegradable properties formulated to attract microbes in landfills to break down 80% of the glove in less than 1½ years. Offers strong puncture resistance while it maintains high dexterity. Compact and recycled packaging to reduce material waste, storage space, and shipping cost.

Gloves 300/Box

Specify:

XS	(143-8648)	L	(143-8651)
S	(143-8649)	XL**	(143-8652)
M	(143-8650)		

**250/Box



BeeSure® Naturals Face Masks

Level 3

ECOBEE

Designed with unique, natural elements.

- More than 60% earth-friendly cellulose components
- Quad-Fold® design offers 15% more breathing volume
- 100% non-latex and fiberglass-free filter layers

Earloop Masks 50 per Box

Specify:

Canyon Red (142-4176) Glacier Blue (142-4177)
Cloud White (142-4178) Forest Green (142-4180)



BeeSure® Floral Earloop Masks

Level 3

ECOBEE

- Quad-Fold® design offers 15% more breathing volume
- RoundEdge® corners and cushioned earloops for comfort
- Non-latex CoolSilk® and dye-free inner layer for a smooth, nonirritating feeling

Earloop Masks 50 per Box

Specify:

Hibiscus Pink	(142-4158)
Plumeria Blue	(142-4159)
Orchid Lavender	(142-4163)
Fern Green	(142-4165)
Daisy Yellow	(142-4164)



BeeSure® SuperSlim

ECOBEE

Nitrile Powder-Free Exam

Textured fingertips for enhanced grip. Innovative soft and stretchy formulation delivers superior tactile sensitivity and donning ease. 60% less packaging waste. 2.0-mil thick. Cobalt Blue.

Gloves 300/Box

Specify:

XS	(353-0052)	L	(353-0055)
S	(353-0053)	XL	(353-0056)
M	(353-0054)		



Ultra® Sensitive Face Masks Level 3

CROSSTEX

Same excellent features as the Ultra® High Efficiency Mask. White.

- Free of chemicals, inks, and dyes
- Fluid resistant
- White hypoallergenic inner cellulose layer
- Will not lint, tear, or shred
- Not made with natural rubber latex
- Made in the USA

Earloop Masks

(107-7191)..... 50 per Box



Ultra® Sensitive Face Masks with SecureFit® Technology Level 3

CROSSTEX

Same excellent features as the Ultra® Sensitive Face Masks. White.

- Fluid resistant
- White hypoallergenic inner cellulose layer
- Patent pending aluminum nose and chin piece for a custom fit
- Not made with natural rubber latex
- Made in the USA

Earloop Masks

(774-0131)..... 50 per Box



Ultra® Sensitive Fog-Free™ Face Masks Level 3

CROSSTEX

Same excellent features as the Ultra® Fog-Free™ Face Mask. White.

- White, hypoallergenic inner cellulose layer
- Will not lint, tear, or shred
- Not made with natural rubber latex
- Made in the USA
- ASTM Level 3

Earloop Masks

(111-6641)..... 40 per Box



Ultra® Face Masks with SecureFit® Technology Level 3

CROSSTEX

Same excellent features as the Ultra® High Efficiency Mask. Blue.

- Fluid resistant
- Patent pending aluminum nose and chin piece for a custom fit
- Not made with natural rubber latex
- Made in the USA

Earloop Masks

(774-0093)..... 50 per Box



Ultra® Sensitive Fog-Free™ Face Masks with SecureFit® Technology Level 3

CROSSTEX

Same excellent features as the Ultra® Sensitive Fog-Free™ Face Masks.

- Fluid resistant
- Hypoallergenic
- Nonwoven Fog-Free™ strip across inner layer of mask
- Patent pending aluminum nosepiece
- Dual chin pieces for a custom fit
- Not made with natural rubber latex
- Made in the USA

Earloop Masks

(774-0132)..... 40 per Box



CRITERION™ Sensitive Skin Earloop Masks Level 3

CRITERION

Same great features as the Criterion™ Sensitive Skin Earloop Masks Level 1, but with the highest barrier of protection.

- Fluid resistant
- Earloops and flexible nosepiece for a secure fit
- Hypoallergenic and free of chemicals, inks, and dyes
- Made in the USA
- Not made with natural rubber latex

Masks..... 50 per Box

Specify:

Blue (570-3276) Pink (570-3277)



CRITERION™ Earloop Masks Level 3

CRITERION

- Light and comfortable for extended wear
- Filter media is surrounded by soft, nonirritating inner and outer facing
- Full-width, flexible nosepiece provides fit and security
- Not made with natural rubber latex

Masks..... 50 per Box

Specify:

Blue (570-2044) Teal (570-2042)
Lavender (570-2047) White (570-2043)
Pink (570-2045) Yellow (570-2046)



CRITERION™ Earloop Masks Level 2

CRITERION

Formerly known as Henry Schein Procedure Masks Level 2. Same Great Product, New Criterion™ Packaging

- Fluid resistant
- Pleated design for enhanced breathability
- Earloops and a flexible nosepiece for a secure fit
- Made in the USA
- Not made with natural rubber latex

Masks..... 50 per Box

Specify:

Blue (104-7799) Pink (104-9200)
Lavender (104-8071)



CRITERION™ Surgical Tie-On Masks Level 3

CRITERION

4-ply surgical mask adds an extra layer of protection and comfort. Same great features as Criterion™ Surgical Tie-On Masks Level 1, but with the highest barrier of protection.

- Fluid resistant
- Flexible nosepiece
- Pleated design for enhanced breathability
- Not made with natural rubber latex

Blue

(570-3294)..... 50 per Box



CRITERION™ Earloop Masks Level 1

CRITERION

Formerly known as Henry Schein Earloop Level 1 Procedure Masks. Same Great Product, New Criterion™ Packaging

- Fluid resistant
- Pleated design for enhanced breathability
- Earloops and a flexible nosepiece for a secure fit
- Not made with natural rubber latex

Masks..... 50 per Box

Specify:

Blue (104-3809) Yellow (104-2849)
Lavender (104-8600) White (104-6611)
Pink (104-3730)



CRITERION™ Sensitive Skin Earloop Masks Level 1

CRITERION

Formerly known as Henry Schein UltraSoft Procedure Masks Level 1. Same Great Product, New Criterion™ Packaging

- Fluid resistant
- Earloops and a flexible nosepiece for a secure fit
- Hypoallergenic and free of chemicals, inks, and dyes
- Made in the USA
- Not made with natural rubber latex

Earloop Masks..... 50 per Box

Specify:

Blue (104-8072) Pink (104-8073)



Earloop/Eye Shield Procedure Masks Level 3

HENRY SCHEIN®

- Optically clear shield is treated on both sides with antifog material
- Not made with natural rubber latex
- Teal
- ASTM Level 3

(104-7321)..... 25 per Box



CRITERION™ Anti-Fog Earloop Masks Level 3

CRITERION

Lightweight and comfortable, these 3-ply procedural face masks offer ASTM Level 3 protection. Each mask features an integrated antifog strip to prevent eyewear from fogging and a soft, nonirritating inner and outer layer.

- Fluid resistant
- Pleated design for enhanced breathability
- Earloops and a flexible nosepiece for a secure fit
- Not made with natural rubber latex

Masks..... 50 per Box

Specify:

Blue (570-3281) White (570-3283)
Lavender (570-3282) Pink (570-3284)

CRITERION™ Extra Breathing Room Earloop Masks

CRITERION

Formerly known as Maxi-Gard™ Air Procedure Masks With Air Space™ Technology Level 1. Same Great Product, New Criterion™ Packaging

- Adjustable, full-width nosepiece
- Fiberglass-free
- Not made with natural rubber latex
- ASTM Level 1

Level 1..... 50 per Box

Specify:

Blue (570-0722) Lavender (570-0724)
White (570-0723)

Pleated Surgical Masks Level 1

HENRY SCHEIN®

- Glass-Free and Non-Latex
- Fiberglass-free
- Not made with natural rubber latex

Tie-On Masks

(104-6809)..... 50 per Box





Safety Glasses



- Special spherical lens curved in both horizontal and vertical positions, similar to precision prescription eyewear: superior field of vision; eliminates distortion; optically correct; less eye fatigue
- Designed to fit over most prescription street-wear spectacles up to 58mm eye size: protection for those who need to wear prescription glasses
- Integral-molded, solid side shield: total lateral protection from foreign contaminants such as spray, splash, or splatter
- Contoured, molded brow guard: total protection from top entry

#5148002
(100-7382).....Ea.



Maxi-Gard® Protective Eyewear



- Stylish, wraparound protective eyewear for greater comfort and protection
- Clear, shatter-resistant polycarbonate lens
- 10-base curve for complete coverage
- Includes a scratch-resistant, antifog coating and foam inserts for added comfort
- Meets ANSI Z87+ standards

Eyewear..... Ea.

Specify:

Black Frame (570-0434) Pink Frame (570-0431)
Blue Frame (570-0433) Purple Frame (570-0432)



Maxi-Gard® Protective Eyewear



Clear, shatter-resistant polycarbonate lens with scratch-resistant coating. Meets ANSI Z87+ standards. 9-base curve allows for great wraparound protection and coverage. Stop hinge prevents temple tip from scratching lens. Protective case included. Metal.

- Shatter-resistant polycarbonate lens that meets ANSI Z87+ standards
- Scratch-resistant coating
- 9" base curve for great wraparound protection and coverage
- Stop hinge prevents temple tip from scratching lens
- Includes protective case

811 Series
(900-4437).....Ea.

Protective Eyewear



- Lightweight, wraparound, unisex protective eyewear
- 10-base curve for complete coverage
- Excellent optical quality
- Ideal for patient and staff use
- Gray lenses and temples

(570-0435).....Ea.



Maxi-Gard® Protective Eyewear



- Clear, shatter-resistant polycarbonate lens with 9" (metal) and 10" (blue and black) base curve for complete coverage
- Scratch-resistant coating
- Meets ANSI Z87+ standards
- Sure-grip rubber temple insert for comfort

806 Series..... Ea.

Specify:

Black (900-4436) Pink (900-4779)
Blue (900-4488) Red (900-4781)



Disposable Face Shields



- Optically clear, see-through visor provides eye protection from splash and splatter
- "Soft-strap" has an impenetrable foam barrier providing comfortable, wide coverage without compromising safety
- Each shield is packaged individually in a plastic bag intended for single-patient use to prevent cross contamination
- Not made with natural rubber latex

7½"L x 13"W
(101-2254).....24 per Box



Compare to *ValuMax Extra-Safe Jackets and Lab Coats*
Maxi-Gard® SMS Jackets



Premium protective jacket available in 9 fashionable, attractive colors. Stylish hip-length jacket is made from high-quality SMS fabric that meets OSHA fluid-resistance requirements for protection. Lightweight, breathable material for maximum wearer comfort. Detailed finishing includes a 3-pocket design and snap-front closure. Not made with natural rubber latex.
Sizes S–XL.....10 per Pkg.

Specify:

Small

Black	(570-2288)
Blueberry	(570-2289)
Ceil Blue	(570-2290)
Light Pink	(570-2291)
Purple	(570-2292)
Raspberry	(570-2293)
Sky Blue	(570-2294)
Teal	(570-2295)
White	(570-2296)

Medium

Black	(570-2298)
Blueberry	(570-2299)
Ceil Blue	(570-2300)
Light Pink	(570-2301)
Purple	(570-2302)
Raspberry	(570-2303)
Sky Blue	(570-2304)
Teal	(570-2305)
White	(570-2306)

Large

Black	(570-2327)
Blueberry	(570-2328)
Ceil Blue	(570-2329)
Light Pink	(570-2297)
Purple	(570-2278)
Raspberry	(570-2330)
Sky Blue	(570-2331)
Teal	(570-2332)
White	(570-2333)

X-Large

Black	(570-2318)
Blueberry	(570-2319)
Ceil Blue	(570-2320)
Light Pink	(570-2321)
Purple	(570-2322)
Raspberry	(570-2323)
Sky Blue	(570-2324)
Teal	(570-2325)
White	(570-2326)

Size 2XL–Ceil Blue
(570-2307)..... 10 per Pkg.

Size 3XL–Ceil Blue
(570-2308)..... 10 per Pkg.



Maxi-Gard® SMS Lab Coats



Premium protective lab coat available in 9 fashionable, attractive colors. Stylish knee-length coat is made from high-quality SMS fabric for protection. Lightweight, breathable, non-latex material for optimal wearer comfort. Detailed finishing includes a 3-pocket design and snap-front closure.
Sizes S–XL.....10 per Pkg.

Specify:

Small

Black	(570-2334)
Blueberry	(570-2335)
Ceil Blue	(570-2336)
Light Pink	(570-2337)
Purple	(570-2338)
Raspberry	(570-2339)
Sky Blue	(570-2340)
Teal	(570-2341)
White	(570-2342)

Medium

Black	(570-2346)
Blueberry	(570-2347)
Ceil Blue	(570-2348)
Light Pink	(570-2349)
Purple	(570-2350)
Raspberry	(570-2351)
Sky Blue	(570-2352)
Teal	(570-2353)
White	(570-2354)

Large

Black	(570-2279)
Blueberry	(570-2280)
Ceil Blue	(570-2281)
Light Pink	(570-2282)
Purple	(570-2283)
Raspberry	(570-2284)
Sky Blue	(570-2285)
Teal	(570-2286)
White	(570-2287)

X-Large

Black	(570-2309)
Blueberry	(570-2310)
Ceil Blue	(570-2311)
Light Pink	(570-2312)
Purple	(570-2313)
Sky Blue	(570-2315)
Raspberry	(570-2314)
Teal	(570-2316)
White	(570-2317)

Size 2XL–Ceil Blue
(570-2343)..... 10 per Pkg.

Size 3XL–Ceil Blue
(570-2344)..... 10 per Pkg.

Table Paper



HENRY SCHEIN®

- White
- Quiet, absorbent barrier
- Designed and manufactured to provide strength, protection, and value

White, Smooth, 18" x 225 ft.

(900-4978)..... 12 rolls per Case

21" x 225 ft.

(900-4979)..... 12 rolls per Case



Premium Table Paper Rolls

HENRY SCHEIN®

- Designed and manufactured to provide strength, protection, and value
- White

Crepe, Standard, 125 ft.

18" (100-9865)..... 12 per Case

21" (100-1657)..... 25 per Case

Smooth, Standard, 225 ft.

14" (104-4501)..... 12 per Case

18" (100-9491)..... 12 per Case

21" (100-0177)..... 12 per Case

Smooth, Heavy Duty, 225 ft.

18" (104-6092)..... 12 per Case

21" (104-5273)..... 12 per Case



Poly-Perf® Exam Table Paper

HENRY SCHEIN®

- Absorbent top tissue layer with poly back for fluid barrier
- For exam table coverage; can also be used as scale and table liners
- Perforated every 18" or 21"
- 125-ft rolls
- White

18" (570-0728)..... 9 per Case

21" (570-0727)..... 9 per Case

Disposable Headrest Covers

HENRY SCHEIN®

Plastic

- Act as a barrier to protect patient and equipment from cross contamination
- Economical and convenient dispenser box
- Disposable poly barrier protects patient and dental chair headrest
- Convenient space-saver dispenser box

Regular Size

11½" x 10".....250 per Box

Specify:

Clear (112-4871) White (112-4873)

Jumbo Size

14½" x 10".....250 per Box

Specify:

Clear (112-4872) White (112-4874)



Headrest Covers

HENRY SCHEIN®

Tissue/Poly

- Tissue outer surface provides greater comfort for your patients
- Poly lined to help protect your equipment

10" x 10"..... 500 per Case

Specify:

White (101-7104) Mauve (101-3198)

Blue (101-4595) Green (102-8451)

Lavender (102-5543)

Regular Size (13" Open Side)

13" x 10"..... 500 per Case

Specify:

White (101-8627) Mauve (101-1051)

Blue (101-8955) Green (102-5542)

Lavender (102-9374)



Pillowcases

HENRY SCHEIN®

Tissue/Poly/Tissue

- Provide linen-like softness for patient comfort and an inner poly barrier for extra protection
- Strong and durable

White, 21" x 30"

(104-9206)..... 100 per Case



Premium FABRICEL® Pillowcase

HENRY SCHEIN®

1-Ply Tissue/Poly, White, 21" x 30"

- Disposable pillowcase is super soft and strong because it is made from Fabricel® material
- By bonding supple tissue to polyethylene film with a heat-lamination process, a unique material that is soft, strong, and tear resistant is created
- Deluxe material prevents soak-through and provides the most durable pillowcase available

21" x 30"–White

(101-2526)..... 100 per Case



Facial Tissues

HENRY SCHEIN®

- Soft, absorbent, elegantly embossed white facial tissue
- Patients and employees will notice the softness and quality

9" x 8"–White

(900-4214)..... 100 per Box



C-Fold Towels

HENRY SCHEIN®

- Individually folded towels
- Whitened without chemicals containing chlorine bleach
- 1-ply, strong, and absorbent
- Fits a wide range of universal dispensers
- 150 sheets/pack
- Made in the USA
- White

12.87" x 10.5"–White

(570-0704).....2400 per Case

Contains: 150 sheets per pack, 16 packs per case.



Multifold Towels

HENRY SCHEIN®

- White
- 9½" x 9¼"
- 250 sheets per pack

9½" x 9¼"–White

(570-0703)..... 4000 per Case

Contains: 250 sheets per pack, 16 packs per case.



Plastic Cups

HENRY SCHEIN®

Embossed Drinking Cups with Rolled Lips

- Heavy, sturdy plastic: leakproof; will not crumple in hands
- Releases easily from dispenser: easily accessible
- Rolled lip: easier to drink
- Embossed wall: nonslip
- High-gloss finish and assorted colors: decorative

5 oz..... 1000 per Case

Specify:

Beige (100-3193)

Blue (100-4361)

Dusty Rose (100-2104)

Gray (100-0473)

Green (101-2443)

Lavender (101-2156)

Peach (101-6812)

White (100-3902)

Yellow (100-1596)





Premium Exam Drapes



2-Ply, All-Virgin Tissue

- Designed to blanket the patient in lightweight comfort

White, 36" x 48"	(101-3658)	100/case
White, 40" x 48"	(100-1797)	100/case
Blue, 40" x 48"	(100-3360)	100/case
Mauve, 40" x 48"	(100-6428)	100/case
White, 40" x 60"	(101-8521)	100/case



Premium Exam Drapes



3-Ply, All-Virgin Tissue

- Soft, opaque, and heavyweight
- Embossed tissue drape with dependable wet strength is ideal for short to moderate exams
- Conforms to the body and offers the ultimate in absorbency

White, 40" x 48"	(101-7084)	100/case
------------------	------------	----------



Pink Ribbon Exam Drapes



2 Ply, 40" x 48"

- Provide patient draping for privacy

(900-4275)	100/case
------------	----------



Dri-Gard® Towels



2-Ply Tissue/1-Ply Poly

- Provide positive moisture barrier and quick absorbency
- Unique horizontal embossing pattern protects patient's clothing and increases the strength of the towel
- Made in the USA
- Not made with natural rubber latex

13" x 19"	500/case
-----------	----------

Specify:

Aqua	(100-9276)	Lavender	(101-0476)
Beige	(100-3501)	Peach	(100-7384)
Blue	(100-7265)	Rose	(100-6341)
Gray	(100-1313)	White	(100-5715)
Green	(100-1350)	Yellow	(100-9704)

Dri-Gard® Plus Towels



3-Ply Tissue/1-Ply Poly

- Provide a positive moisture barrier with maximum absorbency
- Unique horizontal embossing pattern protects patient's clothing and increases strength of the towel
- Made in the USA
- Not made with natural rubber latex


13" x 19"	500/case
-----------	----------

Specify:

Aqua	(100-9851)	Lavender	(101-6810)
Beige	(100-5015)	Peach	(100-7395)
Blue	(100-9825)	Rose	(100-2024)
Gray	(100-9545)	White	(100-7908)
Green	(100-4479)	Yellow	(100-3395)



Dri-Gard® Towels



Teddy Bear Design, 2-Ply Tissue/1-Ply Poly

- Add color to the operatory and bring smiles to younger patients
- Made in the USA
- Not made with natural rubber latex

13" x 19"	(100-9712)	500/case
-----------	------------	----------



Standard Professional Towels




3-Ply Paper

- Specially treated to increase the wet and dry tensile strength of towel. Provides the perfect blend of softness and absorbency. Made in the USA. Not made with natural rubber latex.

13" x 19"	500/case
-----------	----------

Specify:

White	(100-3161)	Blue	(100-0183)
-------	------------	------	------------



Contour Bibs



Plastic-Backed, 3-Ply Paper/1-Ply Poly

- Provide 30% more coverage than standard dental bibs
- 3-ply tissue offers extra absorption
- Poly backing provides durable moisture barrier
- Perfect for facial exam procedures


18" x 22"	400/case
-----------	----------

Specify:

Blue, Ultimate	(102-8707)
Green, Ultimate	(100-6332)
Lavender, Ultimate	(103-2185)
White, Ultimate	(100-8260)



Protection Plus Disposable Underpads



- Made of waterproof, blue poly backing
- Absorbent core of tissue layers
- White, hypoallergenic spunbond facing sheet
- Not made with natural rubber latex

17" x 24"	(112-5633)	300/case
23" x 24"	(112-5634)	200/case
23" x 36"	(570-1182)	150/case
30" x 30"	(570-1183)	150/case
30" x 36"	(570-1181)	100/case



Richmond

BRINGING THE BEST
TO DENTISTRY SINCE 1895




Cotton Rolls | Pellets | Sponges | Face Masks
Reflective Shields Plus® | Dispensers

RichmondDentalMedical.net

For a sample go to: www.bit.ly/RDMsample

7 out of 10 Americans use
social media.* **DO YOU?**



STAY INFORMED
on the most up-to-date
news in the dental industry.

DEVELOP YOUR SKILLS
at our numerous
educational events.

*Pew Research, 2021

Supplies | Equipment | Practice Management Systems
Digital Technology | Technical Service | Practice Services

#RELYONHENRYSCHEIN

HENRY SCHEIN®
DENTAL

Rely on Us™

© 2024 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors.

wound care



Premium Paper Surgical Tape

HENRY SCHEIN®

- 10-yd Rolls
- Excellent adhesion
- Pinked, nonfraying edges
- Strong woven backing
- Gentle removal
- Hypoallergenic
- Not made with natural rubber latex

10-yd. Rolls.....Box

Specify:

½", 24 Rolls per Box **(101-4305)**

1", 12 Rolls per Box **(101-9673)**

2", 6 Rolls per Box **(101-1332)**



Premium Hypoallergenic Transparent Surgical Tape

HENRY SCHEIN®

- 10-yd. Rolls
- Excellent adhesion
- Porous, with flexible backing
- Tears easily in both directions
- Not made with natural rubber latex

Tape.....Box

Specify:

½", 24 per Box **(101-6655)**

1", 12 per Box **(101-4645)**

2", 6 per Box **(101-9226)**



Premium Hypoallergenic Cloth Surgical Tape

HENRY SCHEIN®

- 10-yd Rolls
- Reliable adhesion
- Highly porous
- Thin, flexible backing
- Gentle removal
- Not made with natural rubber latex

10-yd. Rolls.....Box

Specify:

½", 24 per Box **(101-3025)**

1", 12 per Box **(101-2323)**

2", 6 per Box **(101-5979)**

All-Gauze Sponges



Soft, pliable, and absorbent. Available in 3 sizes in both 8-ply and 12-ply styles.

Non-sterile.

- Precision sizes and folding: consistency within functions
- Soft and pliable: nonabrasive touch
- Highly absorbent: soaks up unwanted fluid
- Can be used for general cleaning, debriding, and prepping

8 Ply, 2" x 2"	(570-2657).....	200/pkg
12 Ply, 2" x 2"	(570-2658).....	200/pkg
8 Ply, 3" x 3"	(570-2650).....	200/pkg
12 Ply, 3" x 3"	(570-2648).....	200/pkg
8 Ply, 4" x 4"	(570-2656).....	200/pkg
12 Ply, 4" x 4"	(570-2630).....	200/pkg



Premium All-Gauze Sponges



- Manufactured from 100% USP type-VII cotton gauze that is bleached and processed to established standards
- Soft and pliable
- Protects as a wound covering as it provides maximum absorbency

Non-Sterile

8 Ply, 2" x 2"	(570-2647).....	200/pkg
12 Ply, 2" x 2"	(570-2628).....	200/pkg
12 Ply, 3" x 3"	(570-2640).....	200/pkg
8 Ply, 4" x 4"	(570-2646).....	200/pkg
12 Ply, 4" x 4"	(570-2645).....	200/pkg
Sterile		
8 Ply, 2" x 2"	(570-2642).....	100/pkg
Contains: 2 sponges per envelope, 50 envelopes per package.		
12 Ply, 3" x 3"	(570-2631).....	80/pkg
8 Ply, 4" x 4"	(570-2644).....	50/pkg
12 Ply, 4" x 4"	(570-2637).....	50/pkg



Exodontia Sponges



100% cotton gauze sponges with maximal absorbency cotton inner liner. Superior and consistent quality. Available in 3 sizes. 8-ply.

Non-sterile—Filled

2" x 2"	(570-2636).....	5000/case
3" x 3"	(570-2602).....	4000/case
4" x 4"	(570-2651).....	2000/case
Sterile—Filled		
2" x 2"	(570-2653).....	200/pkg



Premium Exodontia Sponges



Manufactured from all USP 100% cotton. Soft, absorbent, and pure white. Highest quality. 8-ply.

- Manufactured from all USP 100% cotton
- Soft, absorbent, and pure white
- Highest quality

Non-Sterile—Filled

2" x 2"	(570-2635).....	5000/case
Sterile—Filled		
2" x 2"	(570-2652).....	200/pkg



Type VII Gauze




- Sterile
- Not made with natural rubber latex
- 4" x 4"

8 Ply	(570-2655).....	10/pkg
12 Ply	(570-2633).....	10/pkg
16 Ply	(570-2638).....	10/pkg
12 Ply, X-ray Detectable	(570-2632).....	10/pkg
16 Ply, X-ray Detectable	(570-2634).....	10/pkg



Gauze Pads



Highly absorbent. Can be used for general cleaning, debriding, and prepping. 100% woven cotton. Sterile 1's in peel-back packages. 12-ply.

2" x 2"	(570-2659).....	100/pkg
3" x 3"	(570-2629).....	100/pkg
4" x 4"	(570-2649).....	100/pkg



CRITERION™ Non-Woven Sponges



Sterile 2's, 4 Ply

- High quality
- Excellent construction
- Rayon/poly-blend that is soft and gentle to the touch
- Not made with natural rubber latex
- Thicker and more absorbent than standard sponges
- Multiuse sponge that can be used for all applications of wound-care treatments
- Easy-to-open, peel-back packaging
- Non-sterile

2" x 2"	(104-7897).....	50/pkg
Contains: 2 sponges per package, 25 packages per box.		
3" x 3"	(104-4132).....	50/pkg
Contains: 2 sponges per package, 25 packages per box.		
4" x 4"	(104-8665).....	50/pkg
Contains: 2 sponges per package, 25 packages per box.		



CRITERION™ Non-Woven Sponges



Non-sterile, 4 Ply

- High-quality, excellent construction, rayon/poly-blend sponge that is soft and gentle to touch and very comfortable against skin
- Thicker and more absorbent than standard sponges
- Multiuse sponge can be used for all applications of wound-care treatments
- Easy-to-open, peel-back packaging
- Not made with natural rubber latex

2" x 2"	(104-2930).....	200/pkg
3" x 3"	(104-2410).....	200/pkg
4" x 4"	(104-2411).....	200/pkg



Premium Non-Woven Sponges



Non-sterile, 4 Ply

- General-use, non-woven rayon/polyester-blend material
- Quarter-fold design
- Less linting: minimizes lint contamination
- Soft and pliable: prevents irritation to mucosa
- Greater absorbency: reduces number of sponges required
- Closed weave
- Not made with natural rubber latex

2" x 2"	(101-9278).....	200/pkg
3" x 3"	(101-8748).....	200/pkg
4" x 4"	(101-4336).....	200/pkg

Integra® Miltex® Sterile, Single Use Instruments



Limit uncertainty with quality and craftsmanship from the brand you trust

Integra® Miltex®

Disposables

Scissors

Forceps

Needle Holders

Scalpels


Blades

GYN Disposables

Biopsy Punches

Dermal Curettes

Ear Curettes



For more information or to place an order,
please contact:

 **HENRY SCHEIN®**
MEDICAL

www.henryschein.com/medical

Integra, the Integra logo and Miltex are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. ©2024 Integra LifeSciences Corporation. All rights reserved. T-0505705-1-EN 4557016-1-EN (2024-04)

Right for You and Your Patients – Bionix® OtoClear® Ear Irrigation Tips and Systems

The right tools for clinical excellence in cerumen removal. The soft, comfortable 3-streamed OtoClear Tip eases patient discomfort, improves efficiency, prevents the risk of over-insertion and reduces messy clean up. Choose from any of our three systems – your staff and patients will thank you!

- | | |
|------------|---|
| (124-4312) | OtoClear SprayWash Kit (5 Tips) |
| (111-9727) | OtoClear SprayWash Kit (20 Tips) |
| (136-7604) | OtoClear Aquabot® Kit (5 Tips) |
| (136-7603) | OtoClear Aquabot® Kit (20 Tips) |
| (611-7508) | OtoClear Electronic Ear Irrigation System (5 Tips) |
| (110-1386) | OtoClear Electronic Ear Irrigation System (40 Tips) |
| (131-8021) | OtoClear Tip (20 Tips) |
| (967-1733) | OtoClear Tip (40 Tips) |



For more information, contact your local Henry Schein Medical Representative or visit henryschein.com/medical.

 **HENRY SCHEIN®**

 **bionix®**

CaviWipes™

Easy as **1 2 3**

When it comes to keeping communities healthy,
every decision matters.
Make the easy choice with CaviWipes towelettes.



CaviWipes™ HP

Hydrogen Peroxide Wipe
for alcohol-sensitive
equipment and
fast turnarounds

Part #	Item Description
(640-0122)	Regular (6" x 6.75") 160 count canister
(142-5192)	XL (9" x 12") 65 count canister



CaviWipes™ 2.0

Low Alcohol/Quat Wipe
for everyday surfaces and
most medical devices

Part #	Item Description
(140-3081)	Regular (6" x 6.75") 160 count canister
(140-3082)	XL (9" x 12") 65 count canister



CaviWipes™ Bleach

Bleach Wipe
for high-risk areas and
outbreak control

Part #	Item Description
(129-6430)	Standard (6" x 10.5") 90 count canister



Scan QR Code to Learn More



CaviWipes™ Disinfecting Towelettes KERR/METREX

Kills both TB in 3 minutes and HBV in 2 minutes. These disposable towelettes are presaturated with CaviCide Surface Disinfectant Cleaner and are recommended for cleaning and disinfecting all nonporous surfaces. Meet OSHA and CDC guidelines. Fragrance-free. Intermediate level.

- Towelettes that kill both TB in 3 minutes and HBV in 2 minutes
- Disposable and presaturated with CaviCide Surface Disinfectant Cleaner and are recommended for cleaning and disinfecting all nonporous surfaces
- Meet OSHA and CDC guidelines
- Fragrance-free
- Intermediate level

Canister Pkg.

Large—6" x 6¾"

(640-2805)..... 160 per Canister

X-Large—9" x 12"

(640-8482)..... 65 per Canister



CaviWipes™ Disinfecting Towelettes

KERR/METREX

Non-woven, disposable wipe that offers a 1-minute kill time for viruses, bacteria, and fungi, including TB, MRSA, VRE, influenza A, HIV-1, H1N1, and HBV. This low-alcohol, multipurpose wipe can be used on hard, nonporous surfaces. Bleach-, phenol-, and fragrance-free. Intermediate level.

- One-minute contact time for virucidal, bactericidal (including TB), and fungicidal activity
- One-step cleaning and disinfecting*
- Compatible with common hard, nonporous surfaces found in the clinical environment

Canisters

Large—6" x 6¾"

(640-0024)..... 160 per Canister

X-Large—9" x 12"

(640-0023)..... 65 per Canister

FlatPack

7" x 9" convenient counterspace-saving packaging with patented adhesive backing for optimal placement on most smooth surfaces.

(640-0026)..... 45 per Pkg.

Individually Wrapped

X-Large—9" x 12"

(640-0027)..... 50 per Box



CaviCide® KERR/METREX

- Effective against TB in 3 minutes, and HIV-1, HCV, HBV, and MRSA in 2 minutes
- Disinfectant, cleaner, and decontaminant all in one
- May be used on most medical device materials
- For use on hard, nonporous surfaces

Spray Bottle

(640-2261)..... 24 oz.

Refill

(640-0012)..... Gal.

5 Gallon Refill

(640-5006)..... 5 Gal.

Contains: 2 x 2½ Gal.



CaviWipes™ HP KERR/METREX

- Alcohol-free, hydrogen-peroxide wipes that do not require PPE and are safe for use for both clinicians and patients
- Offer a 1-minute universal contact time against over 50 pathogens
- Qualified for EPA emerging viral pathogen claims for enveloped viruses (e.g., influenza A and B, HIV-1, and HBV) as well as small/large nonenveloped viruses (e.g., norovirus and rhinovirus types 1A and 37)

6" x 6¾"

(640-0122)..... 160 per Canister

9" x 12"

(142-5192)..... 65 per Canister



CaviWipes™ 2.0 Disinfectant Wipes KERR/METREX

- Effective against over 40 pathogens
- 2-minute efficacy against norovirus, influenza A and B, SARS-CoV-2, TB, HIV-1, and more
- Qualified for the EPA's Emerging Viral Pathogen claim for all virus types
- Compatible with hard, nonporous surfaces

Large: 6" x 6¾"

(140-3081)..... 160 per Canister

X-Large: 9" x 12"

(140-3082)..... 65 per Canister



PROTECTING WITH PURPOSE

Clorox Healthcare®
Fuzion® Cleaner Disinfectant
HS Item #: (126-5240)
9/32 fl. oz. Trigger spray

Clorox Healthcare®
Hydrogen Peroxide Wipes
HS Item #: (437-0022)
6.75" x 9", 6/95 ct.

Clorox Healthcare®
Hydrogen Peroxide Cleaner
HS Item #: (141-4662)
9/32 fl. oz. Trigger spray




© 2024 Clorox Professional Products Company,
1221 Broadway, Oakland, CA 94612.
NI-63086

infection prevention

Clorox® Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes CLOROX HEALTHCARE

- 30–60 second contact time for most bacteria and viruses (see product label for specific kill claims)
- Kill claims for 32 different microorganisms, including germs of high concern such as 9 antibiotic-resistant bacterial strains (including MRSA and VRE), norovirus, rhinovirus, rotavirus, RSV, influenza A, and E. coli
- High-quality tear-resistant wipes for greater durability
- No harsh chemical odor or fumes
- Noncorrosive, nonflammable

Regular–6¾" x 5¾"
(437-0023)..... 155 per Canister
Large–6¾" x 9"
(437-0022)..... 95 per Canister



Clorox Healthcare™ Bleach Germicidal Cleaners CLOROX HEALTHCARE

Fast-Acting Cleaner Disinfectant

- EPA-registered to kill *C. difficile* spores in 5 minutes
- Kills norovirus in 1 minute
- Formulated with 1:10 bleach solution strength, which is recommended by CDC guidelines for effective disinfecting
- Meets OSHA bloodborne pathogen standards
- Cleans and disinfects in 1 step

#68978, 1 Gallon
(437-0035)..... Ea.
Refill
(120-8623)..... 32 oz.



Clorox® Commercial Solutions Clorox® Disinfecting Spray CLOROX HEALTHCARE

- Bleach-free, phenol-free, broad-spectrum disinfectant and deodorizer that meets AOAC Germicidal Spray Product Test efficacy standards for hospital disinfectants
- Kills 99.9% of bacteria and viruses including HIV, TB, hepatitis A, poliovirus, E. coli, influenza, and more
- Fresh scent

19-oz Can
(250-1203)..... Ea.



Clorox® Healthcare Disinfecting Wipes CLOROX HEALTHCARE

Fresh-scented wipes that clean, sanitize, and disinfect nonporous surfaces throughout the dental office, including metal surfaces. Meets AOAC efficacy testing requirements. Bleach-free, phenol-free, and EPA-registered disinfectant. Low level.

Fresh Scent
(250-2369)..... 35 per Canister



**MaxiSpray® Plus****Surface Disinfectant Cleaner**

- Easy, ready-to-use formula: no mixing, measuring, or rinsing required
- Compatible with most medical/dental device material
- Safe to use in neonatal and infant-care areas
- Contains no phenols
- Nontoxic and nonirritating
- Satisfies OSHA's Bloodborne Pathogens Standard

Pump Spray

(101-4114)..... 24 oz.

Refill

(101-0285)..... Gal.

**MaxiWipe® Germicidal Cloths**

Disinfectant wipe that kills TB, RSV, HBV, HCV, HIV-1, rhinovirus, MRSA, VRE, influenza A (H1N1) virus, in 2 minutes.

- Bactericidal, tuberculocidal, and virucidal
- Compatible with a broad range of surfaces and equipment in healthcare
- Thick, absorbent, material helps remove soil from surfaces without tearing.
- Unique deep-well lid on canister seals securely to prevent moisture loss.
- Meets CDC and OSHA regulations
- Intermediate level disinfectant

Large—6" x 6¾"

(106-6794)..... 160 per Canister

X-Large—7½" x 15"

(112-5305)..... 65 per Canister

**MaxiWipe® Alcohol-Free
Germicidal Disposable Wipes**

Formulated with care: Alcohol- and fragrance-free wipes for staff and patients with respiratory sensitivities. Efficient, one-step cleaner and disinfectant. Effective bactericidal, tuberculocidal, and virucidal, including SARS-CoV-2 (the virus that causes COVID-19), in 3 minutes. EPA Registration No. 9480-9, now qualified for List Q*. Intermediate level.

Large—6" x 6¾"

(570-3056)..... 160 per Canister

X-Large—7½" x 15"

(570-3057)..... 75 per Canister

*MaxiWipe® Alcohol-Free has demonstrated effectiveness against emerging viral pathogens, such as the monkeypox virus, on hard, nonporous surfaces.

**PURELL® Advanced Aloe Instant Hand
Sanitizer****GOJO**

Enhanced with 4 different skin-conditioning agents that help maintain skin health and improve overall skin feel. Kills 99.99% of the most common germs on hands, not just in the lab. Works in as little as 15 seconds without water or towels. 70% ethyl alcohol.

Pump Bottle

(590-0066)..... 12 oz.

**PURELL® Advanced Instant Hand
Sanitizer****GOJO**

Kills 99.99% of most common germs that may cause illness. Works in as little as 15 seconds. No water or towels needed. Specially formulated with moisturizers to leave hands feeling soft and refreshed. Alcohol-based formula.

Squeeze Bottle

(313-1626)..... 4 oz.



Hand Sanitizer Gel



Waterless, alcohol-based gel hand sanitizer. Ideal for use where soap and water are not readily available. Rinse-free: no water required. Kills 99.9% of most common disease-causing germs spread by hands. Soothing emollients and a pleasant, fresh, clean scent. Leaves hands soft and clean. 66.5% ethyl alcohol. Made in the USA.

- Contains alcohol (active ingredient) to kill germs
- Soothing emollients leave hands feeling soft and clean
- Pleasant citrus fragrance

Pump Bottle
(900-4443)..... 8 oz.
(900-4442)..... 16 oz.



Antibacterial Liquid Hand Soap



- Provides a luxurious feel
- Leaves hands feeling soft and smooth
- Mild, gentle, and moisturizing
- Pleasant Cucumber and Aloe fragrance
- Active ingredient is 0.2% triclosan

Pump Bottle
(900-4439)..... 8 oz.
(900-4438)..... 16 oz.

Refill Bottle
(900-4440)..... Gal.



Foaming Antibacterial Hand Soap



- Provides a rich, dense lather
- Rinses easily using minimal water
- Leaves hands feeling soft and smooth
- Pleasant Spring Water Clean fragrance
- Contains 0.6% chloroxylenol (PCMX) for antimicrobial action, and aloe vera to add moisture with each use

Pump Bottle
(900-4441)..... 18.6 oz.



Maxicens® Instant Hand Sanitizer and Dispensers



- Kills 99.9% of germs in 15 seconds without the use of soap and water
- Contains aloe vera to add moisture with every use
- Henry Schein's hand sanitizer helps you meet OSHA, APIC, and CDC hand-washing recommendations
- Citrus scent

800-mL Bag for Manual Dispenser
(112-5988)..... Ea.
800 mL Bag for Automatic Dispenser
(112-5991)..... Ea.
Manual Dispenser
(112-5984)..... Ea.
Automatic Dispenser
(112-5986)..... Ea.
Stand for Dispenser
(112-5989)..... Ea.



Quick Wipes



- Soft, strong towelettes: allow easy and gentle cleansing of hands and face
- Effective and gentle formulation that contains aloe: formulated for fast drying and gentleness
- Convenient packaging: towelettes easily dispensed through orifice; each towelette is saturated with cleansing liquid

Pop-up Container
(104-4176)..... 150 per Cont.



Maxicens® Antimicrobial Hand Wash



- Contains 0.6% PCMX for antimicrobial action and aloe vera to add moisture with each use
- Aids in reducing cross contamination
- Helps meet APIC and OSHA recommendations for hand washing

800-mL Bag for Manual Dispenser
(112-5987)..... Ea.
800-mL Bag for Automatic Dispenser
(112-5990)..... Ea.
Manual Dispenser
(112-5983)..... Ea.
Automatic Dispenser
(112-5985)..... Ea.
Stand for Dispenser
(112-5989)..... Ea.



Maxicens™ Antiseptic Skin Cleanser



4.0% chlorhexidine gluconate formulation. Fast acting in just 15 seconds, while inhibiting bacterial regrowth for up to 6 hours. Skin emollients protect against irritation, keeping hands soft and soothed. Residue-free, nonsticky formulation with no added dyes. All-in-one solution is tough on bacteria, yet gentle on the skin. Pleasant floral scent.

Bottle
(570-3268)..... 32 oz.
Refill
(570-3269)..... Gal.



The right treatment begins with the right diagnosis.

Test your patient with both the **BD Vaginal Panel** and the **BD CTGCTV2** assay to diagnose vaginitis and the most prevalent non-viral STIs with one swab.¹⁻⁷

References: 1. BD Vaginal Panel Package Insert (P0258) and BD Vaginal Panel for use with the BD COR™ System (L012454). 2. BD CTGCTV2 for BD MAX™ System Package Insert (P0237). 3. Van Der Pol B et al. *Sex Transm Dis.* 2021;48(2):134–40. 4. Broache M et al. *Obstet Gynecol.* 2021;138(6):853–9. 5. Brown H and Drexler M. *Popul Health Manag.* 2020;23(S1): S3–S12. 6. Hillier SL et al. *Clin Infect Dis.* 2021;72(9): 1538–43. 7. Miller JM et al. *Clin Infect Dis.* 2018;67(6):e1–e94.

BD and the BD Logo are trademarks of Becton, Dickinson and Company or its affiliates. © 2024 BD. All rights reserved. (BD-62308 3310-US-0424)



point-of-care testing



Affirm™ VPIII Kit

BD

Easy-to-perform, moderately complex DNA probe technology. Conclusive identification of *Candida* species, *Gardnerella vaginalis*, and *Trichomonas vaginalis* from a single vaginal sample.

- Accuracy up to 30% higher than slide microscopy
- Less than 2 minutes hands-on time required
- Easy-to-read results
- Built-in procedural controls for each sample
- **Made in the USA**

#446252

(987-7212)24/box



MAKING AN IMPACT IN COMMUNITY HEALTH TOGETHER

FOCUSING ON THOSE IN NEED

Health centers improve the health and well-being of underserved communities providing accessible, high-quality primary health care services. As a global leader in rapid testing, Abbott's point-of-care technology provides fast, actionable results to support appropriate care. Let Abbott and Henry Schein help you do what you do best, care for your community.

- Cost-effective solutions
- Accurate and reliable test results in minutes
- Optimal care for patient satisfaction

MAKING AN IMPACT AT THE POINT OF CARE

STEP
1

CHOOSE WHICH PRODUCTS WILL BEST SERVE YOUR COMMUNITY



☐ ID NOW™

A platform for detection of Flu A & B, Strep A, RSV, and COVID-19 infections, providing the fastest rapid molecular results within 2 to 13 minutes¹



☐ Afnion™ 2

Compact, rapid analyzer that provides fast and reliable HbA1c and ACR test results in as little as 3 to 5 minutes at the point-of-care



☐ Cholestech LDX™ Analyzer

Rapid, accurate results in 5 minutes for a complete lipid profile and glucose to identify CVD risks



☐ Determine™ HIV-1/2 Ag/Ab Combo

First and only FDA approved 4th generation rapid test, detecting both acute and chronic HIV infection with results in 20 minutes



☐ BinaxNOW™ COVID Ag Card

Designed to detect active COVID-19 infections with or without symptoms (authorized under EUA)*



☐ iCassette Urine Test Dx Fentanyl Drug Screen**

Help keep your community safe with a CLIA waived point-of-care urine test cassette for detection of Fentanyl with 1ng/mL cutoff

STEP
2

SELECT THE OPTION THAT IS RIGHT FOR YOU

☐ Capital Option

- No contracts
- No volume commitments required

☐ Placement Option

- Aggregate annual spend commitment
- One agreement for all products**
- No upfront instrument capital needed

STEP
3

CONTACT YOUR ABBOTT OR HENRY SCHEIN REPRESENTATIVE TO LEARN MORE

VISIT [GLOBALPOINTOFCARE.ABBOTT](https://globalpointofcare.abbott) FOR MORE INFORMATION ON OUR TESTING SOLUTIONS

1. <https://abbott.mediaroom.com/2018-10-25-Abbott-Introduces-the-Next-Generation-of-Influenza-A-B-and-Strep-A-Assays-with-Fastest-Ever-Time-to-Molecular-Results>

* The BinaxNOW™ COVID-19 Ag Card has not been FDA cleared or approved. It has been authorized by the FDA under an emergency use authorization. It has been authorized only for the detection of proteins from SARS-CoV-2, not for any other viruses or pathogens, and is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. § 360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

**Clinical Toxicology products covered under separate agreement

© 2024 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. COL-24891-01 05/24



The right outcome starts with the right test.

One system with multiple assays delivering reliable results for SARS-CoV-2,* Flu A, Flu B, RSV, and Group A Strep within 15 minutes at the point-of-care.

See what's new at BDVeritor.com

* In the USA, this product has not been FDA cleared or approved but has been authorized by the FDA under an Emergency Use Authorization for use by authorized laboratories; use by laboratories certified under the CLIA, 42 U.S.C. §263a, that meet requirements to perform moderate, high, or waived complexity tests. The product is authorized for use at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation.

This product has been authorized only for the detection of proteins from SARS-CoV-2, influenza A and influenza B, not for any other viruses or pathogens; and, in the USA, the emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID 19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. §360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

BD Life Sciences, 7 Loveton Circle, Sparks, MD 21152-0999, U.S.

BD, the BD Logo and BD Veritor are trademarks of Becton, Dickinson and Company or its affiliates.
© 2024 BD. All rights reserved. (BD-89766 3308-US-0424 April 2024)



point-of-care testing

BD Veritor™ Plus System

BD

The wireless BD Veritor™ Plus System provides health-care providers and laboratorians in physician offices, clinics, hospitals, and integrated delivery networks (IDNs) with objective, lab-quality immunoassay test results within minutes. This fast and accurate solution streamlines the point-of-care (POC) diagnostic workflow and allows providers to quickly review patient results and determine the appropriate treatment in a single consultation. The new system can help detect infectious diseases such as influenza A and B, respiratory syncytial virus (RSV), and group A strep while the patient is still on site, which can enable faster decision-making for an earlier intervention opportunity.

- **Faster access to test results:** Help inform appropriate and targeted treatment decisions
- **Improved workflow:** The BD Veritor™ Plus Analyzer improves workflow with an Analyze Now or Walk Away mode, providing test results (with a print-enabled function) and aiding in decision-making while the patient is on site
- **Improved traceability:** The BD Veritor™ InfoScan module can scan 3 different identifiers (operator ID, specimen ID, and test kit/lot) that can be locally recorded, printed, or downloaded, and reduces hands-on time through automated documentation



- **Modularity:** Provides a unique solution to meet each testing site's preference
- **Connectivity:** BD Veritor™ InfoSync module provides the same benefits as BD Veritor™ InfoScan. Connects and documents tests results across multiple sites and can integrate with a hospital electronic medical record or laboratory information system
- **Advanced particle technology:** Enhances sensitivity by using a proprietary process to produce highly stable, modified colloidal metal particles.
- **Adaptive read technology:** Compensates for background. Nonspecific binding helps improve specificity to reduce false-positive results

#256066, Veritor Plus Analyzer
(126-5527)..... ea

BD Veritor™ System Flu A+B Test Kit

BD

The first CLIA-waived RSV test referenced against a higher sensitivity standard than culture. Test device is ready to insert into reader after sample is added. The same sample processed for RSV can also be used for flu A+B. Digital results displayed in 10 seconds.

- Sample type: NP swab sample
- Positive percent agreement (PPA): 81.6%
- Negative percent agreement (NPA): 99.1%
- 6 years and younger

#256038, RSV Test
(121-7553) 30/box
CPT® Code: 87807QW†



BD Veritor™ System RSV CLIA-Waived Kit

BD

The first CLIA-waived flu A+B test referenced against PCR, a higher sensitivity standard than culture. Device is ready to insert into reader 10 minutes after sample is added.

- Positive percent agreement (A/B): 82.1%/74.6
- Negative percent agreement (A/B): 98.6%/99.6%
- Positive agreement vs. PCR for current visual read rapid tests ranges from 10% to 70%
- Streamlined workflow: Provides a digital result in <11 minutes, with <50 seconds of hands-on time

#256045, Flu A+B Test Kit
(992-0002) 30/box
CPT® Codes: Influenza A: 87804QW† / Influenza B: 87804-59†; Using the Objective and Differentiated Reader System: 87999†





As many as

50%

of women **with vaginitis**
are misdiagnosed,
creating a frustrating
cycle of *ineffective*
treatment, repeat visits,
and an increased
*risk of recurrence.*¹

**Xpert Xpress
MVP Reagent**
HS# (146-1059)

Women's health is our health.
We all deserve to know **more.**

1 Brown H and Drexler M. Popul Health Manag. 2020;23(S1):S3-S12.
US-IVD. In Vitro Diagnostic Medical Device.



Introducing **Xpert® Xpress MVP**, a CLIA-waived **multiplex vaginal panel** designed to aid clinicians in accurately diagnosing three distinct vaginal infections: Bacterial Vaginosis, Vulvovaginal Candidiasis, and Trichomoniasis.

www.Cepheid.com





Sofia® 2 Platform



The Sofia 2 analyzer represents the next generation of rapid diagnostic testing

Proven lateral-flow technology, proprietary advanced fluorescence chemistry and assay development techniques are integrated into two small and portable benchtop analyzers. Our high-performance analyzers deliver point-of-care testing for better patient outcomes.



Simplify your workflow

Our Sofia benchtop systems and fluorescent immunoassays deliver highly accurate, objective and automated results in three to 15 minutes.



Efficient data storage and retrieval

View results remotely from your computer with Virena® System wireless connectivity. Patient results are automatically stored and results never expire.



Learn more about Sofia by scanning the QR code, visiting quidelortho.com, or reaching out to your Henry Schein representative.



Convenient and easy to use

The Sofia systems are small in form, have intuitive interfaces and only require about one minute of total hands-on time.

New QuidelOrtho branding may not be available in all markets, subject to country-specific regulatory approval. Please confirm with your local commercial team.

© 2024 QuidelOrtho Corporation PR-10572-NA-EN-US-v1

point-of-care testing

Sofia® 2 Analyzer QUIDEL

Offers an easy-to-use interface, customizable settings, flexible workflow, barcode technology, and Virena® connectivity. Results are reported automatically on the screen and can be stored in Sofia 2, optionally printed on an external printer, or sent to Virena. Sofia 2 has 2 modes that can be used depending on your workflow requirements: Walk Away mode may be most convenient for reading a single patient sample. Advance Result Technology (ART) can provide results in as few as 3 minutes. Read Now mode may be most convenient for reading several patient samples.

Sofia 2 features:

- Accurate, objective and automated results in as few as 3 minutes
- Excellent performance compared to molecular and viral culture methods



- Flexible, dual-mode testing and simple, CLIA-waived procedures for a variety of laboratory environments
- Secure, remote instrument management through the Virena ecosystem
- Integrated data management automatically stores test and user history
- User and patient ID captured with onboard barcode scanner

#20299, Sofia 2 Analyzer
(128-5292).....ea

Sofia® Influenza A+B FIA Test QUIDEL

This CLIA-waived test uses advanced immunofluorescence-based lateral-flow technology to detect influenza A and B viral nucleoprotein antigens in nasal and nasopharyngeal swab and nasopharyngeal aspirate/wash specimens taken directly from symptomatic patients.

- Enhanced assay sensitivity
- Results in 15 minutes
- Accurate detection with nasal swab, NP swab and NPA/wash samples
- Nasal swab sensitivity (A/B): 90%/89%; specificity (A/B): 95%/96%
- NP swab sensitivity (A/B): 97%/90%; specificity (A/B): 95%/97%
- NPA/wash sensitivity (A/B): 99%/88%; specificity (A/B): 96%/96%
- Room temperature storage
- 21-month shelf life from date of manufacture
- For use with Sofia Analyzer (sold separately)
- CLIA-waived for nasal and nasopharyngeal swabs



#20218, Sofia Influenza A+B Test
(122-4092) 25/box
CPT® Codes: Influenza A: 87804QW†; Influenza B: 87804-59†

Sofia® Rapid Strep A Test QUIDEL

Easy-to-use test with enhanced assay sensitivity and results in 5 minutes.

With results available in minutes, you can rapidly diagnose and treat, all in the same visit.

- Accurate detection with throat swab samples
- Sensitivity: 90.6 %
- Specificity: 96.1%
- Room temperature storage
- 24-month shelf life from date of manufacture
- For use with Sofia Analyzer (sold separately)

#20253–Sofia Strep A
(122-4089)..... 25/box
#20274, Sofia Strep A+
(122-7865)..... 25/box
CPT® Code: 87880†



Sofia® RSV FIA Test QUIDEL

A 15-minute test that has a walk-away automation feature that allows user

to multitask. Offers accurate detection with nasopharyngeal swabs, or aspirate/wash samples. CLIA-waived for children less than 7 years of age, and moderately complex for pediatric patients from age 7 to less than 19 years of age.

- Single specimen may be used to run both Sofia influenza and RSV FIA tests
- Room temperature storage
- 24-month shelf life from date of manufacture
- For use with Sofia Analyzer (sold separately)

#20260, Sofia RSV Kit
(122-4091) 25/box
CPT® Code: 87807†



The Savanna HSV 1+2/VZV assay leverages real-time PCR technology with a differentiated VZV target to provide superior patient outcomes faster.

For more information, visit quidelortho.com or contact your Henry Schein Representative.

* Data on file.
New QuidelOrtho branding may not be available in all markets, subject to country-specific regulatory approval. Please confirm with your local commercial team.

© 2024 QuidelOrtho Corporation
All rights reserved.

Savanna® HSV 1+2/VZV Assay

The power of certainty

Lesions for HSV-1, HSV-2, and VZV are misdiagnosed by visual examination in up to 20 percent of cases.

Make an impact for your patients with Savanna HSV 1+2/VZV assay, a unique rt-PCR assay, by providing results up to 230 percent faster than a typical send out.*



QuidelOrtho™

point-of-care testing



Savanna® Platform—Quidel **QUIDEL**

This compact benchtop instrument utilizes real-time PCR to deliver automated, customized molecular results in minutes. Molecular sensitivity with rapid turnaround can detect multiple pathogens in a single sample. Eliminates the need for dedicated lab space. Facilitates decentralization of testing and allows multiplex, syndromic panel testing at point of care.

- Sample-to-result technology
- Decreased hands-on time with no upfront sample or reagent preparation required
- Superior molecular accuracy and confidence in a multiplex, real-time PCR assay
- Single-use, self-contained cartridge with no sample or reagent prep prior to testing
- Fully automated and easy to use with touchscreen display results in 24 minutes

#20382, Real-Time PCR Instrument
(146-1056)..... ea
#20445, HSV1+2/VZV Test Kit
(146-0986)..... 12/kit
#M409-1, HSV1+2/VZV Assay Control Set
(146-0987).....ea
#M409-4, HSV1+2/VZV Assay Control Set
(146-0985).....4/kit

Resilience + Passion

We make
diagnostics that

matter



We understand it is crucial to immediately diagnose and treat vulvovaginal disorders which are linked to a broad range of associated health complications, from preterm birth to an increase in Sexually Transmitted Infections.¹

The OSOM® BVBLUE® and OSOM® Trichomonas Tests are easy-to-use, CLIA-waived and produce accurate and objective results for two of the most common causes of vulvovaginal disorders so you can get the answers fast and your patients back to doing what they love.

Like you, we understand there is a patient behind every answer—and that's what matters most.

For more information, contact your local Henry Schein Medical Representative or visit henryschein.com/medical.



The Perfect Pair

OSOM® TRICHOMONAS RAPID TEST

OSOM® BVBLUE® TEST

SEKISUI
DIAGNOSTICS

Because every result matters™

¹ Brown, H. Improving the Diagnosis of Vulvovaginitis. Population Health Management. Vol. 23, suppl 1, 2020

© 2024 SEKISUI Diagnostics, LLC. All rights reserved. OSOM® is a registered trademark of SEKISUI Diagnostics, LLC. Because every result matters™ is a trademark of SEKISUI Diagnostics, LLC. BVBLUE® is a registered trademark of Gryphus Diagnostics, LLC.

point-of-care testing



OSOM® BVBLUE® Test SEKISUI

For *Gardnerella vaginalis*

Diagnose bacterial vaginosis (BV) in just 10 minutes for faster treatment and education in the same office visit.

- Test for the 4 major pathogens that cause BV in just one test—reacts with sialidase produced by *Gardnerella vaginalis*, *Bacteroides* spp., *Prevotella*, and *Mobiluncus* spp.
- **Simple to use**—Get easy-to-read results in 4 easy steps and no need for additional equipment
- Decrease the chance of human error with minimal hands-on requirements
- **Quicker results**—Read results in just 10 minutes
- **Easy to read**—Instant color change for accurate results
- **Get more patient tests per kit**—Exclusive QC Inside® provides 2 extra tests at no additional charge to run external quality control testing
- **Made in the USA**

#183, Test
(414-4170) 25/box

#184, Control
(327-4117) ea

CPT® Code: 87905QW*



OSOM® Trichomonas Rapid Test SEKISUI

The only CLIA-waived lateral-flow rapid test for the detection of trichomoniasis while providing 95% agreement with culture and wet mount combined.

- **Simple-to-use**—Get easy-to-read results in just 10 minutes and 4 easy steps
- **Eliminate the need for wet mount**—Antigen-based rapid test does not require visualization of live organisms
- CDC-approved method of testing for the detection of *T. vaginalis*—unaffected by other common infections
- **Trust your results**—95% agreement against Composite Reference Standard (culture and wet mount)
- **Easy-to-use**—Requires minimal training and no additional equipment
- **QC Inside®**—2 extra tests included at no additional charge to run external quality-control testing for more confidence in results
- Made in the USA

#181, Test
(406-7631) 25/box

#182, Control
(611-1157) ea

CPT® Code: 87808QW*

A Broad Range of Respiratory Health Point-of-Care Tests

We are committed to providing **high quality, molecular and antigen point-of-care tests for detecting the most common causes of respiratory infections**, so you can get the answers fast and your patients back to doing what they love.

Like you, we understand there is a patient behind every answer—and that’s what matters most.

Learn more about our COVID-19, Flu A&B, and Strep A tests.

For more information, contact your local Henry Schein Medical Representative or visit henryschein.com/medical.

The Aptitude Metrix® COVID-19 Test is provided by SEKISUI Diagnostics.

The Metrix COVID-19 Test has not been FDA cleared or approved but has been authorized for emergency use by FDA under an EUA. This product has been authorized only for the detection of nucleic acid from SARS-CoV-2, not for any other viruses or pathogens. The emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID-19 under Section 564(b) (1) of the Federal Food, Drug, and Cosmetic Act, 21 U.S.C. § 360bbb 3(b)(1), unless the declaration is terminated, or authorization is revoked sooner.

© 2024 SEKISUI Diagnostics, LLC. All rights reserved. OSOM® is a registered trademark of SEKISUI Diagnostics, LLC. Because every result matters™ is a trademark of SEKISUI Diagnostics, LLC. Metrix® is a registered trademark of Aptitude Medical Systems Inc.



SEKISUI
DIAGNOSTICS
Because every result matters™

point-of-care testing



Acucy™ System—Sekisui SEKISUI

The Acucy™ Influenza A&B Test on the Acucy™ System provides clinicians flexibility in workflow and accurate, standardized results for improved patient care.

- CLIA waived
- “Read Now” or “Walk Away” modes
- Results within 15 minutes
- Nasal and nasopharyngeal swab types
- Includes 2 additional tests for quality-control testing

#1039, Acucy Reader System
(136-8725)..... ea
#1010, Acucy Influenza A & B Test
(136-8722)..... ea
CPT® Code: Influenza A: 87804QW*, Influenza B: 87804-59†



OSOM® Ultra Strep A Tests SEKISUI

Faster results for faster treatment—statistically equivalent to single-swab culture with results in just 5 minutes.

- **Simple to use**—Get easy-to-read results in 4 easy steps and no need for additional equipment
- **Greater sensitivity**—Antibody-labeled color particles coated at 2 separate locations allows for greater sensitivity than conventional single-label technology without sacrificing specificity
- **Easy-to-read** 2-color results
- **Simplify reagent addition**—single reagent delivery system in 25-test kit eliminates the need to handle 2 separate reagents
- **QC Inside®**—2 extra tests included in the 50-test kit at no additional charge to run external quality-control testing for more confidence in results
- Made in the USA

#147
(392-7389) 25/box
#149 (Includes 2 additional tests for external QC testing.)
(102-7617) 52/box
CPT® Code: 87880QW
Ask about RGS. Genzyme compares to rigorous Gold Standard.



**OPTIGEN® Allergy Tests
MINARIS**

- #85063, OPTIGEN US Inhalant 36
(131-1154)20/kit
Contains: 20 OPTIGEN pettes & reagents.
- #85065, OPTIGEN US Mix 36
(131-8137)20/kit
Contains: 20 OPTIGEN pettes & reagents.
- #80109, Control Reagent Kit for OPTIGEN
(131-1156)ea
Contains: 3-mL positive/3-mL negative.
CPT® Code: 86003†



**CLA-1™ Luminometer
MINARIS**

This compact luminometer automatically prints a complete report that lists the tested allergens and the severity of the patient's response to each. Easy to use, the system is the simplest, quickest, and most convenient allergy testing method available today. It can analyze up to 180 allergens in 8 minutes. Made in the USA.
#90015
(126-3464) ea

**830 Series Procedure Chair
MTI**

The highly versatile 830 Series chairs are ADA compliant in all base versions (standard, swivel, or mobile) using the transfer support, providing easy patient ingress/egress, and equal access for patients with disabilities and caregivers. All 830 chairs come equipped with SmartTech® technology to enhance safety, reduce energy consumption, and provide smoother movement, additional memory options, and more user preferences.

- 725-lb patient load capacity
- 4 programmable positions
- Programmable home and ADA-compliant home positions
- Power lift, back, tilt, and foot extension
- Hand and foot controls
- Oval articulating headrest
- 22" or 25" seat width
- 18¾"-39" seat height
- Leg rest features knee-break with 7" power foot extension
- Optional patient arms: removable, floating (compensating) arms, or seat rail or back rail mounted arms
- Removable upholstery: seamless or designer
- Tapered backrest for easier caregiver access
- Sturdy all-steel frame
- Multiple G2 accessory options

- #829-001-80, Standard Base (22" seat width)
(135-8121) ea
- #829-002-80, Swivel Base (22" seat width)
(139-4043) ea
- #829-003-22, Mobile Base (22" seat width)
(140-6171)ea
- #830-001-80, Standard Base (25" seat width)
(133-8866) ea
- #830-002-80, Swivel Base (25" seat width)
(136-9044)ea
- #830-003-22, Mobile Base (25" seat width)
(133-9127) ea



**463 and 464 Exam/Procedure Chairs
MTI**

The most feature-packed, ADA compliant* chair on the market. Standard features include ComfortSync® coordinated footrest, EasyView™ chair-back mounted switches, removable slide-back floating armrests, push-button 333° swivel base and oval articulating headrest.

- 725-lb patient load capacity
- 4 programmable positions
- Programmable home and ADA-compliant home positions
- Power lift, recline, and tilt (no tilt on 463)
- 21" seat width
- 18¾"-37" seat height (2" lower than any other ENT chair on the market)
- Removable upholstery: seamless or designer
- Tapered backrest for easier caregiver access
- Sturdy, all-steel frame

- #463-002-80, Swivel Base
(136-8741) ea
- #464-002-80, Swivel Base
(136-8744) ea

*ADA compliant with transfer support accessory for easy patient ingress/egress



Experience the value of automated vital signs capture.

The **Midmark® Digital Vital Signs Device** brings automation to the point of care. Simplify vital signs and weight measurements using one screen to capture and import data directly into the patient record from a single, centralized location. The device now comes equipped with an **Exergen® TemporalScanner®**. Temporal temperature is preferred for its accuracy, speed, hygiene and cost-effectiveness.¹

For more information, call 1.800.P.Schein
or visit: henryschein.com/medical



Exergen and TemporalScanner are registered trademarks of Exergen Corporation.
1 <https://www.exergen.com/wp-content/uploads/2023/02/Exergen-ClinicalStudies-03222023.pdf>

© 2024 Midmark Corporation, Versailles, Ohio USA



Ritter 224 Barrier-Free® Examination Chair

MIDMARK

The number of elderly and obese in the U.S. is steadily growing, and over one fifth of adults in America are disabled. Assisting patients with mobility limitations onto exam chairs has become a frequently performed task. Although caregivers working in health-care facilities are in one of the most at-risk professions for musculoskeletal disorders, little has been done to evaluate and reduce risk. The Ritter 224 is designed to offer low and high positioning heights so you can provide the best care to all patients, whether treating the elderly, expectant mothers, or patients with disabilities.

#224-002, Without Clean Assist Roller System (133-0159)..... Ea.

#224-003, With Clean Assist Roller System (133-0268)..... Ea.



Ritter 230

MIDMARK

Universal Power Procedures Table

Features eight-way power for effortless patient positioning. Whether it's power height, back, foot or tilt, it does the lifting for the physician—minimizing back stress, reducing workplace injury, and keeping patients comfortable.

- Power height, back, foot and tilt
- 450-lb patient weight capacity
- 22½" to 40" seat height adjustment
- "Knee-chest" position configurable
- Hand control standard
- Seamless upholstery top
- Articulating headrest
- Adjustable stirrups
- 18" paper roll holder
- Meets UL 60601-1 and CAN/CSA C22.2 No.60101-M90 safety standards
- Hospital grade electrical outlets and table rotation also available
- Optional foot control
- Available in 8 seamless upholstery colors

#230-001, Ritter 230 Base Unit (108-1319).....Ea.

#230-002, Ritter 230 Base Unit With Receptacle (238-7678)..... Ea.

#230-003, Ritter 230 Base Unit with Rotation (475-6182).....Ea.



Ritter 225 Barrier-Free® Examination Chair

MIDMARK

The number of elderly and obese in the U.S. is steadily growing, and over one fifth of adults in America are disabled. Assisting patients with mobility limitations onto exam chairs has become a frequently performed task. Although caregivers working in health-care facilities are in one of the most at-risk professions for musculoskeletal disorders, little has been done to evaluate and reduce risk. The Ritter 225 is designed to offer low and high positioning heights so you can provide the best care to all patients, whether treating the elderly, expectant mothers, or patients with disabilities.

#225-002, Without Clean Assist Roller System (133-0279).....Ea.

#225-003, With Clean Assist Roller System (133-0290)..... Ea.



425 Air Lift Physician Stool

MIDMARK

Designer Series

Designed for your comfort and support all day, every day. Choose the features you want. Start from the "bottom up" with the frame and seat, then add the back and armrests. Order everything you want at once, or add on later, as you prefer.

(151-4236)..... Ea.



626 Barrier-Free® Examination Chair

MIDMARK

Many health-care providers are not aware of current standards for proper blood-pressure measurement. If this isn't addressed, it leaves health systems at serious risk. 9.8% of the population is affected by overestimation of high blood-pressure errors. Multiply the cost of treatment for hypertension by the typical patient panel size and the cost of overtreatment could be as much as \$135,000 per each practicing physician. The impact of overmedication can also result in complications such as dizziness, fainting, exhaustion, or even patient falls. The Midmark 626 is the first exam chair of its kind designed to promote proper patient positioning for a more accurate blood-pressure measurement.

- Easily position patients for a better BP measurement
- Powered movement of the back section ensures the patient's back is supported
- A Barrier-Free® low chair-height allows the majority of patients, even female patients with heights in the 3rd percentile, to place their feet flat on the floor

#626-001, Base Model (133-0301)..... Ea.

#626-002, With Heated Upholstery (133-0312)..... Ea.

#626-003, With IQ Scale (133-0323).....Ea.

#626-004, With Heated Upholstery and Scale (133-0334).....Ea.



Midmark 272 Exam Stool

MIDMARK

Ritter stools by Midmark provide optimal comfort, support, and stability. The star-shaped base design allows for closer access to the patient, and stable, easy maneuverability.

- Black, plastic composite 5-leg caster base and hood
- Pneumatic seat-height adjustment
- Many upholstery color options available

Exam Stool.....Ea.

Specify:

#272-001-845, Lunar Gray (132-2295)
#272-001-814, Stone (132-2291)

Clarus 7™ LED Exam Light MTI

MTI's brightest exam light uses 7 LEDs and has a slightly larger 3½" light head, providing maximum illumination and operational ease. Can be used with either a caster base or floor stand. 2-year limited warranty.

Physical dimensions: head diameter: 3½"; caster base height: 73½"; floor stand height: 71½"; caster base diameter: 23¾", floor stand diameter: 15".

- LED life: 40,000 hours
- Color temperature (CCT): 4,300° K
- Light intensity @ 16": 33,400 Lux
- Light intensity @ 12": 48,500 Lux
- Working distance: 16"
- Spot size at 10": 5"
- Spot size at 12": 6"
- Spot size at 16": 8"
- Spot size at 24": 11"
- CRI: 75
- Control box dimmable on light head
- Electrical input: 100–240 V, 50–60 Hz, 20W, 12V @ 1.67A
- Power cord: 3-wire grounded cord set, 8' long
- Listings: UL 153, UL 8750, CAN/CSA: C22.2, No. 12

Exam Light with Caster Base

(121-8211)ea

Exam Light with Floor Stand

(121-8212) ea

Exam Lamp DRIVE MEDICAL DESIGNS

This versatile 16" mobile lamp has a flexible gooseneck design and a dome-style shade. Its adjustable arm offers full 360° flexibility, allowing the user to direct light in any direction. With a weighted steel base for stability, the height can be easily adjusted from 48" to 72" to accommodate individual preferences. The lamp has a 3-prong plug, is UL listed, and supports incandescent light bulbs up to 60 watts, ensuring both safety and convenience. Limited lifetime warranty.

#13408MB

(139-7964) ea

Welch Allyn Green Series™ 600 Minor-Procedure Lights BAXTER

Featuring 3 LEDs and a wider head, the GS 600 is designed specifically for demanding viewing conditions. Whiter, brighter, and greener: 3 white LEDs offer greater intensity. Other features include touchless power controls and an adjustable light output. It also has a disposable sheath available to help reduce the risk of infection due to cross contamination. Rugged, ergonomic design allows for improved performance and long-lasting durability. Low cost of ownership due to elimination of bulb replacement and lower energy consumption. 5-year limited warranty to include the LED light source and power supply only.

#44600, With Mobile Stand

(566-0088)ea

#44610, With Table/Wall Mount

(566-0089) ea

MI 1000 Ceiling Surgical Light, Dual SYMMETRY SURGICAL/ BOVIE

The MI 1000 LED light is the first in the cutting-edge MI LED Series offered by Bovie. The light emits very low heat with an output rating of 100,000 lux (10,000 foot candles) at 4300°K color temperature. Engineered with green technology, the MI 1000 uses significantly less electricity than standard halogen lights. The long-life LEDs are rated at 50,000 hours, virtually eliminating the need for inconvenient bulb changes. Five intensity outputs are available and are controlled using the 5-stage dimming and on/off switch located on the sterilizable handle. 100 volts.

#XLD-DC

(117-2977)ea

MI 750 LED Exam Light SYMMETRY SURGICAL/ BOVIE

The MI 750 provides physicians with an economical and energy-efficient LED light for exam and procedure rooms. The MI 750 produces bright, white light for clear and accurate visibility. Designed with unmatched lifetime spring-arm technology, 50,000-hour LEDs, and output rated at 75,000 lux with a Kelvin color temperature of 4300°K, the MI 750 offers a high-quality lighting system that easily aligns with any practice's budget. To maintain the sterile field when necessary, the on/off and dimming switches are located on the autoclavable handle. The MI 750 LED is backed by a 3-year manufacturer warranty.

#XLDP-SC

(123-8002)ea

Welch Allyn Green Series™ 300 General Exam Lights BAXTER

Utilizing a single white LED, the GS 300 provides a whiter, brighter light than halogen and consumes less energy. It's the perfect combination of intense light output and cool operation in a compact, rugged design. Features a touchless on/off, a choice of mobile stand or wall/table mount for maximizing floor space, low cost of ownership due to elimination of bulb replacement and lower energy consumption, and a 5-year limited warranty that includes both the LED light source and power supply.

With Mobile Stand

(114-2692) ea

With Table/Wall Mount

(114-5632) ea

Welch Allyn Green Series™ Exam Light IV BAXTER

Perfect for the OB/GYN office or for specialties that require very focused light in specific areas, the GS Exam Light IV provides you with 3 times the typical light output and features a highly flexible light pipe that is twice as long as conventional pipes. And with a compact, unobtrusive head, it ensures pinpoint accuracy and provides true tissue color rendition during an exam. Utilizing white LEDs, the GS Exam Light IV provides a whiter, brighter light than halogen and consumes less energy. 5-year limited warranty to include the LED light source and power supply only.

- One LED
- Flexible pipe
- Touchless on/off
- Easy to assemble

Light with Table/Wall Mount

(114-2547) ea

(114-4177) ea

Disposable Light Sheath

(566-0989) 25/box

DIAGNOSTIC TECHNOLOGY, SMART TOOLS, INNOVATIVE TECHNOLOGY

Designed for Clinicians in Primary Care

Heart disease is the leading cause of death in the United States, and stroke ranks among the top five.¹

Your patients are facing heart related challenges daily. Accurate measurement and diagnosis are vital. Partner with us today to leverage our solutions to help meet the needs of your patients.

Contact your Henry Schein representative to learn more.



Check out our current promotions in patient monitoring and connected cardiology today:
Vitals: hillrom.com/2024vitalsrebate **Cardiology:** hillrom.com/2024cardiorebate

Rx Only. For safe and proper use of the products mentioned herein, please refer to the appropriate Operator's Manual or Instructions for Use.

Reference:

1. Kochanek, M.A., K. D., Murphy, B.S., S., Xu, M.D., J., & Arias, Ph.D., E. [2023, September 22]. National Vital Statistics reports - Centers for Disease Control. <https://www.cdc.gov/nchs/data/nvsr/nvsr72/nvsr72-10.pdf>

Baxter, CP 150, Connex, Diagnostic Cardiology Suite, ELI, Hillrom, Spot Vital Signs, Welch Allyn are trademarks of Baxter International Inc.

US-FLC45-240008 [v1.0] 04/2024



Welch Allyn Green Series™ 777 Integrated Wall System BAXTER

#777-PM2WAS-US, With PanOptic Basic LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer
(140-1989) ea

#777-PM3WAS-US, With PanOptic Plus LED Ophthalmoscope, MacroView Plus LED Otoscope for iExaminer, BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer
(140-2495) ea

#777-SM2WAS, With Coaxial LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer
(140-1801) ea

#777-SM2WAX, With Coaxial LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, and Ear Specula Dispenser
(140-1799) ea



Welch Allyn Connex® Spot Monitor BAXTER

- Designed for the physician-office environment with an easy-to-use, vivid touch-screen display
- More accurate hypertension diagnosis without slowing down your workflow thanks to automated blood-pressure averaging
- Improved capabilities for assessing respiratory conditions such as COPD, asthma, and community-acquired pneumonia (CAP), thanks to integrated professional-grade pulse oximetry
- Increase staff efficiency and eliminate data-entry errors with wireless transfer of vital signs directly into an EMR

7100 Base Models with Wired Connectivity

#71XX-B, SureBP Non-invasive Blood Pressure
(566-0464) ea

#71XT-B, SureBP, SureTemp
(672-0176) ea

#71XE-B, SureBP, Braun Thermoscan PRO 6000 Ear Thermometer
(124-6103) ea

#71WX-B, SureBP, Nonin SpO₂
(672-0177) ea

#71WT-B, SureBP, Nonin SpO₂, SureTemp
(566-0466) ea

#71WE-B, SureBP, Nonin SpO₂, Braun Pro6000
(124-5328) ea

7300 Bluetooth Connectivity Physician Office Models

#73XT-B, SureBP, SureTemp
(566-0465) ea

#73XE-B, SureBP, Braun PRO 6000
(124-8200) ea

#73WX-B, SureBP, Nonin SpO₂
(672-0178) ea

#73WT-B, SureBP, Nonin SpO₂, SureTemp
(672-0179) ea

#73WE-B, SureBP, Nonin SpO₂, Braun PRO 6000
(124-4801) ea

#73CT-B, SureBP, Covidien SpO₂, SureTemp
(672-0175) ea

#73CE-B, SureBP, Covidien SpO₂, Braun PRO 6000
(125-8900) ea

#73MT-B, SureBP, Masimo SpO₂, SureTemp
(566-0467) ea

#73ME-B, SureBP, Masimo SpO₂, Braun PRO 6000
(124-7063) ea

7400 Upgradable to WiFi Connectivity Hospital Models

#74CX-B, SureBP, Covidien SpO₂
(566-0469) ea

#74CT-B, SureBP, Covidien SpO₂, SureTemp
(566-0470) ea

#74CE-B, SureBP, Covidien SpO₂, Braun PRO 6000
(127-5208) ea

#74MX-B, SureBP, Masimo SpO₂
(566-0471) ea

#74MT-B, SureBP, Masimo SpO₂, SureTemp
(566-0478) ea

#74ME-B, SureBP, Masimo SpO₂, Braun PRO 6000
(127-6033) ea



Welch Allyn Spot™ Vision Screener

BAXTER

A handheld, portable device designed to help users quickly and easily detect vision issues on patients from 6 months of age through adult. Spot screens both eyes at once from a nonthreatening 3-ft distance. The touch-screen display allows for one-touch activation, simple management of patient data entry, and easy configuration for both vision screening and autorefractive applications.

Features:

- Supports AAP screening guidelines for early detection of amblyopic risk factors
- Automated screening provides thorough, objective, and easy-to-understand results
- Easy to use, easy to implement, with minimal user training required
- Adequately powered for mass screening events
- Captures readings 97% of the time, helping to screen otherwise difficult patients
- Wireless printing and WiFi-enabled for easy export of data
- Manual exports available via USB port
- Uses lights and sounds to help engage children
- Convenient carrying case

#VS100S-B, Vision Screener Set

(497-0004) ea



Welch Allyn Connex® Cardio ECG BAXTER

This ECG provides results that can be trusted in an easy-to-understand format designed with workflow in mind.

- Captures best 10 seconds of ECG data from memory, reducing need for repeat ECGs
- Complies with ACC, AHA, and HRS recommendations for adult and pediatric patients
- Streamlines workflow by capturing and saving an ECG in as few as 2 clicks

#CC-RXX-WAXX, PC-Based ECG with Wireless Acquisition Module (WAM)

(129-6992) ea

#CC-RXX-AAXX, PC-Based ECG with AM12™ (Wired) Cable

(129-6993) ea

#CC-RXX-AADX, ECG with AM12™ (Wired) Cable and DICOM Connectivity

(129-6995) ea

#CC-RXX-WADX, ECG with Wireless Acquisition Module (WAM) and DICOM Connectivity

(129-6994) ea



Green Series™ 777 Integrated Wall System



- All models include BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer.
- #777-PM2WAS-HS, With PanOptic Basic LED Ophthalmoscope and MacroView Basic LED Otoscope (570-2666)..... Ea.
- #777-PM3WAS-HS, With PanOptic Plus LED Ophthalmoscope and MacroView Plus LED Otoscope for iExaminer (570-2665)..... Ea.
- #777-SM2WAS-HS, With Coaxial LED Ophthalmoscope and MacroView Basic LED Otoscope (570-2667)..... Ea.
- #777-SM2WAX-HS, With Coaxial LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, and Ear Specula Dispenser (Does not include SureTemp Plus Thermometer.) (570-2668)..... Ea.



Pro Plus Series Stethoscopes



- Combination adult-size chest piece is precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability.
- Extra-large bell (1 3/8" diameter) for unsurpassed low-frequency response
 - Ultrasensitive diaphragm (1 1/8" diameter) for greater amplification and crisper high-frequency response
 - Color-coordinated nonchill bell and diaphragm retaining rim for patient comfort
 - Stainless steel binaurals with double-leaf internal spring fixed at 15° angle for greater comfort
 - Reinforcing yoke molded into flexible 22" PVC "Y" tubing
 - Deluxe PVC ear tips (small and large) for the ultimate in wearing comfort and acoustic seal
 - Spare diaphragm included
 - Available in 3 striking colors for departmental coding or personal identification
 - Not made with natural rubber latex
 - Scope ID tag included
 - Lifetime warranty
 - Weight: 6.5 oz.
 - Overall length: 31"
- Stethoscope..... Ea.
- Specify:
- Black (900-4815) Navy (900-4817)
- Burgundy (900-4816)



Otoscope Speculum Tips



- These disposable tips are carefully manufactured to be free of rough edges or burrs. They are designed for a snug fit on Welch Allyn®, Propper®, and Heine® otoscopes. They will not fit the Welch Allyn® Macroview®.
- Carefully manufactured to be free of rough edges or burrs
 - Designed for a snug fit on Welch Allyn®, Propper®, and Heine® otoscopes
 - Will not fit the Welch Allyn® Macroview®
- 2.5mm (112-5527)..... 1000 per Bag
- 4.0mm (112-5528)..... 1000 per Bag



The Essentials Series Pocket Aneroid



- Black enamel, 300mm-Hg, no-pin stop manometer
 - Navy blue nylon cuff with range markings to facilitate selection of correct cuff size, artery label, and gauge holder
 - Latex inflation bladder and bulb
 - Chrome-plated deflation valve
 - Complete with leatherette zippered carrying case and operating instructions
 - Inspected and packaged in the USA from foreign components
- Aneroid.....Ea.
- Specify:
- Large Adult (112-6076) Child (112-6080)
- Adult (112-6077)



Fingertip Pulse Oximeter



- Dual-color LED makes monitoring easy
 - Display modes include SpO2 (range 70%–100%) and pulse rate (range 30–235 bpm) along with pulse bar and plethysmogram
 - 2 AAA batteries provide 30 hours of continuous use
 - Visual and audible low-battery alarm
 - Suitable for both adult and pediatric patients
 - Electromagnetic compatibility
- (900-6368).....Ea.

ePM 10A Configured Patient Monitor MINDRAY



Exceptional user interface and enhanced parameter functionality, delivering a monitoring solution to fit multiple care environments. Immediate and measurable benefits of the smartphone-like user interface that delivers 98% of all common monitoring functions with just 2 gestures, reducing time clinicians spend navigating monitor functions, and more time on what matters most: patient care. Advanced monitoring options include 3- and 5-lead ECG, SpO₂, NIBP, respiration, and temperature supporting multiple clinical workflows. ePM 10A has a compact, ergonomic design that is perfectly suited for healthcare facilities managing a steady flow of patients daily. Add optional CO₂ capability to aid clinicians to identify signs of patient deterioration due to respiratory distress. Manufactured with industry-leading plastics that stand up to the most stringent medical-grade cleaners. The concealed grab-and-go handle eliminates hard-to-clean crevices and eases clinical workflow while on the move. Flexible mounting options for bedside and transport. Dimensions: 10.6" x 9.9" x 6.3"

Weight: 8.8 Lb.

Warranty: 3 years

Integrated
Sidestream CO₂, RecorderEa.

Specify:
Masimo SpO₂ (139-2003)
Nellcor SpO₂ (139-2005)

Recorder..... Ea.

Specify:
Masimo SpO₂ (139-1920)
Nellcor SpO₂ (139-1927)

iM50 Patient Monitor EDAN

Due to its lightweight and compact design, the new-generation iM50 patient monitor provides the perfect solution for transport or compact environments. Full parameters (NIBP, SpO₂, ECG, temperature, respiration, optional EtCO₂) are clearly displayed on a 8½" color touch screen.

#iM50.P, Patient Monitor with Printer
(974-0176)Ea.

#iM50.CS.P, Patient Monitor with CO₂ and Printer
(974-0183)Ea.

#iM50.PMT, Patient Monitor with Printer and Rolling Stand
(974-0188)Ea.

#iM50.RT, Rolling-Stand
(127-1084)Ea.



iM60 Patient Monitor EDAN

The iM60 monitor spans accurate monitoring across different care environments, acuity levels, and clinical requirements. All data is clearly displayed on a 10½" color touch screen while providing real-time readings on NIBP, ECG, SpO₂, respiration, temperature, and optional EtCO₂.

#iM60.P, Patient Monitor with Printer
(974-0177)Ea.

#iM60.CS.P, Patient Monitor with CO₂ and Printer
(974-0180)Ea.

#iM60.PMT, Patient Monitor with Printer and Rolling Stand
(974-0185)Ea.

#iM60.CS.PMT, Patient Monitor with CO₂, Printer, and Rolling Stand
(974-0186)Ea.



ePM 12MA Patient Monitor MINDRAY



Exceptional user interface and enhanced parameter functionality, delivering a monitoring solution to fit multiple care environments. Immediate and measurable benefits of the smartphone-like user interface that delivers 98% of all common monitoring functions with just 2 gestures, reducing time clinicians spend navigating monitor functions, and more time on what matters most: patient care. Advanced monitoring options are integrated: 3- and 5-lead ECG, SpO₂, NIBP, respiration, and temperature-supporting multiple clinical workflows. The ePM 12MA features a semimodular design with a plug-and-play module structure that adds monitoring functions such as CO₂ and multigas, where and when needed. The 12MA eGM multigas module automatically identifies and measures the inspired and expired values of O₂, CO₂, N₂O and 5 anesthetic agents, providing a compact, yet comprehensive OR monitoring solution. Manufactured with industry-leading plastics that stand up to the most stringent medical-grade cleaners. Concealed grab-and-go handle eliminates hard-to-clean crevices and eases clinical workflow while on the move. Flexible mounting options for bedside and transport.

• Dimensions: 12.2" x 11.4" x 6.7"

• Weight: 10.6 Lb.

• Warranty: 3 years

Semi-Configured/ModularEa.

Specify:
Masimo SpO₂ (139-2126)
Nellcor SpO₂ (139-2128)

ePM AG Module
(139-3276)Ea.

CO₂ Module Adult/Pediatric
(137-8484)Ea.

iM70 Patient Monitor EDAN

The iM70 Patient Monitor is embedded with the latest advancements in design to provide accurate monitoring across different applications that range from general care to anesthesia monitoring. Its high-resolution, 12" touch screen clearly displays all patient data (NIBP, ECG, SpO₂, respiration, temperature, EtCO₂ optional, anesthetic agents optional).

#iM70.P, Patient Monitor with Printer
(974-0178)Ea.

#iM70.CS.P, Patient Monitor with CO₂ and Printer
(974-0181)Ea.

#iM70_OR, Patient Monitor with Multi-Gas Analyzer and Printer
(974-0131)Ea.

#PM.RT, Rolling Stand
(127-1083)Ea.



iM80 Patient Monitor EDAN

The iM80 patient monitor is equipped with a massive, 15" color touch screen to clearly display patient readings and waveforms.

Applications range from general cases to anesthesia monitoring. Parameters include: NIBP, ECG, SpO₂, respiration, temperature, EtCO₂ optional, anesthetic agents optional.

#iM80.P, Patient Monitor with Printer
(974-0179)Ea.

#iM80.CS.P, Patient Monitor with CO₂ and Printer
(974-0182)Ea.

#iM80_OR, Patient Monitor with Multi-Gas Analyzer and Printer
(974-0135)Ea.

#iM80.RT, Rolling Stand
(127-1139)Ea.





ZOLL AED Plus®
ZOLL

Bringing Lifesaving Defibrillation to Everyone

- Provides a one-piece electrode for accurate placement
- Uses consumer lithium batteries that are easy to replace, long-lasting, and inexpensive

Semi-Automatic Unit

(866-7559)..... ea
Contains: AED, electrodes, batteries, carry case, ECG, wall-mount bracket & operator's guide.

Fully Automatic Unit

(119-7910)..... ea
Contains: AED, electrodes, batteries, carry case, ECG, wall-mount bracket & operator's guide.

CPR-D padz®—Adult Defibrillator Pad

(253-2611)..... ea

Pedi padz® II—Pediatric Electrode Pad

(203-7054)..... ea

Type 123 Lithium Batteries

(610-4769)..... 10/pkg

Soft Carrying Case

(106-1937)..... ea

3D Wall Sign

(110-5384)..... ea

Wall Cabinet

(108-7520).....ea

AED Plus Simulator

(499-7385)..... ea



Powerheart G5 AED
ZOLL

The Powerheart® G5 AED with Intellisense™ CPR provides corrective feedback to ensure delivery of high-quality compressions. One-button dual-language functionality and comprehensive self-tests give rescuers the confidence to act.

- Intellisense™ CPR Feedback provides real-time voice and text corrective prompts for CPR quality
- Rescue Ready® conducts regular, automatic self-tests on pads, battery, and internal circuitry
- RescueCoach™ guides users through each critical step of the rescue with user-paced voice/text prompts
- Dual language change from primary to secondary language in the push of a button

#G5S80C-S, Semi-Automatic

(139-5101)..... ea

#G5A80C-S, Fully Automatic

(139-5102)..... ea



AED Pro® Automated
External Defibrillator

ZOLL

This external defibrillator provides advanced capabilities for BLS and ALS users with innovative Real CPR Help® technology for rate and depth of chest compressions. The See-Thru CPR® capability filters CPR chest compressions, displaying the underlying cardiac rhythm. The AED Pro is unmatched in its ruggedness and durability and is ideal for any environment.

- Available as semiautomatic only or semiautomatic with manual override for ALS professionals
- Built-in internal memory and external USB memory technology allows expanded data storage capacity and maximum flexibility for data management and transfer
- Electrodes and battery compatibility with ZOLL M Series®, E Series®, and R Series® professional defibrillators
- ECG monitoring with 3-lead cable and high-quality display
- Unrivaled battery options (long-life disposable or rechargeable) that meet specific needs

(869-6270)..... ea



Butterfly™

Meet the iQ3.™



Best-in-class
image clarity



Advanced
imaging tools



User-centric
ergonomics

With iQ3, diagnosing patients becomes a more intuitive process, allowing you to focus on what matters most: **delivering exceptional care.**

“I am a heavy user of POCUS, and have fully integrated it into my examinations and daily workflow. In this short time, I have concluded that the iQ3 is a game changer.”

*Steve Erickson, MD FAAFP RMSK, Family Medicine Physician**

*iQ3 provided by Butterfly.
980-24281-00 Rev A


Item (146-7168)
Butterfly iQ3 All-In Bundle

Supporting you with quality hypodermics and domestic manufacturing



Consistent quality. Reliable performance.

BD is the only major manufacturer of hypodermic syringes that manufactures primarily in the U.S. We're prepared to increase production at our domestic manufacturing facilities to help supply impacted providers—providing the quality, performance and reliable supply healthcare providers rely on to provide high quality patient care.



Manufacturing plants

Canaan, CT

- Conventional syringes
- Safety needles

Columbus, NE

- Conventional syringes
- Conventional needles
- Safety needles



Learn more about our syringe portfolio and our commitment to manufacturing quality.

BD, Franklin Lakes, NJ, 07417, U.S.

bd.com

BD and the BD Logo are property of Becton, Dickinson and Company or its affiliates.
© 2024 BD. All rights reserved. BD-113505 (01/24)



Safety & Conventional Hypodermic Needles



Features & Benefits:

- 3 single-handed safety activation methods: hard-surface, thumb, or finger
- Double-locking mechanism with an audible “click”
- Needle hub-compatible with all standard Luer lock and Luer slip-tip syringes
- ISO color coded—easy and quick identification of the needle gauge
- Bevel-up needle—allows for low-angle injections

ITEM#	DESCRIPTION	SIZE
NEEDLES		
(112-7106)	□ 16 ga x 1"	100/Box
(112-7116)	□ 16 ga x 1 ½"	100/Box
(900-4469)	■ 18 ga x 1"	100/Box
(900-4471)	■ 18 ga x 1 ½"	100/Box
(112-7109)	■ 20 ga x 1"	100/Box
(112-7110)	■ 20 ga x 1 ½"	100/Box
(112-7100)	■ 21 ga x 1"	100/Box
(112-7111)	■ 21 ga x 1 ½"	100/Box
(112-7107)	■ 22 ga x 1"	100/Box
(112-7101)	■ 22 ga x 1 ½"	100/Box

ITEM#	DESCRIPTION	SIZE
(900-4473)	■ 23 ga x 1"	100/Box
(900-4468)	■ 25 ga x ½"	100/Box
(112-7112)	■ 25 ga x ¾"	100/Box
(900-4470)	■ 25 ga x 1"	100/Box
(900-4472)	■ 25 ga x 1 ½"	100/Box
(112-7108)	■ 27 ga x ½"	100/Box
(112-7103)	■ 27 ga x 1/4"	100/Box
(112-6148)	■ 30 ga x ½"	100/Box
(112-6149)	■ 30 ga x 1"	100/Box

ITEM#	DESCRIPTION	SIZE
SAFETY NEEDLE		
(570-2693)	■ 18 ga x 1"	100/Box
(570-2694)	■ 18 ga x 1 ½"	100/Box
(570-0338)	■ 20 ga x 1"	100/Box
(570-2697)	■ 20 ga x 1 ½"	100/Box
(570-2696)	■ 21 ga x 1"	100/Box
(570-2698)	■ 21 ga x 1 ½"	100/Box
(570-2699)	■ 22 ga x 1"	100/Box
(570-2700)	■ 22 ga x 1 ½"	100/Box

ITEM#	DESCRIPTION	SIZE
(570-0341)	■ 23 ga x 5/8"	100/Box
(570-2701)	■ 23 ga x 1"	100/Box
(570-0340)	■ 23 ga x 1 ½"	100/Box
(570-2705)	■ 25 ga x 5/8"	100/Box
(570-2703)	■ 25 ga x 1"	100/Box
(570-2707)	■ 27 ga x 1/4"	100/Box
(570-2709)	■ 30 ga x ½"	100/Box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

To Order: 1-800-772-4346 8:00am—8:00pm (et) • www.henryschein.com/medical

Doing Well by Doing Good



Visit www.HenrySchein.com/HSCares
to find out more about the Calendar of Caring Programs.

**Thank you for helping us raise more than \$2.8 million
to help health happen!**



© 2023 Henry Schein, Inc. Not responsible for typographical errors.

EZ FLO Needles

- Standard bevel, thin wall from 18 ga to 25 ga
- Regular wall in 27 ga and 30 ga
- Sterile
- Individually packaged
- Stainless steel and a polypropylene hub

#NEZ30012, 30 ga x ½"	100 per Box
(146-1193)	
#NEZ30100, 30 ga x 1"	100 per Box
(146-1194)	
#NEZ27114, 27 ga x 1¼"	100 per Box
(146-1074)	
#NEZ27200, 27 ga x 2"	100 per Box
(146-1192)	
#NEZ25058, 25 ga x ⅝"	100 per Box
(146-1068)	
#NEZ25112, 25 ga x 1½"	100 per Box
(146-1070)	
#NEZ23034, 23 ga x ¾"	100 per Box
(146-1065)	
#NEZ23100, 23 ga x 1"	100 per Box
(146-1066)	
#NEZ22034, 22 ga x ¾"	100 per Box
(146-0930)	
#NEZ22100, 22 ga x 1"	100 per Box
(146-1063)	
#NEZ22112, 22 ga x 1½"	100 per Box
(146-1064)	



#NEZ21100, 21 ga x 1"	100 per Box	#NEZ19112, 19 ga x 1½"	100 per Box
(146-0927)		(146-0926)	
#NEZ21112, 21 ga x 1½"	100 per Box	#NEZ18100, 18 ga x 1"	100 per Box
(146-0928)		(146-0920)	
#NEZ20100, 20 ga x 1"	100 per Box	#NEZ18112, 18 ga x 1½"	100 per Box
(146-0924)		(146-0921)	
#NEZ20112, 20 ga x 1½"	100 per Box	#NEZ18200, 18 ga x 2"	100 per Box
(146-0925)		(146-0922)	
#NEZ19100, 19 ga x 1"	100 per Box		
(146-0923)			



Hypodermic Disposable Needles

EXEL

Sterile, Nontoxic, Nonpyrogenic

Needles and boxes are color coded for easy identification. Each needle is individually packaged, sterilized, and ready for use.

Needles..... 100 per Box

Specify:

#26400, Gray, 27 Ga. x ½"	(337-2518)
#26427, Gray, 27 Ga. x 1¼"	(337-5234)
#26426, Gray, 27 Ga. x 1½"	(337-5455)
#26402, Beige, 26 Ga. x ½"	(337-2831)
#26403, Dark Blue, 25 Ga. x ⅝"	(337-9855)
#26405, Dark Blue, 25 Ga. x 1"	(337-5760)
#26406, Dark Blue, 25 Ga. x 1½"	(337-2640)
#26407, Light Blue, 23 Ga. x ¾"	(337-1604)
#26408, Light Blue, 23 Ga. x 1"	(337-2432)
#26424, Brown, 19 Ga. x 1"	
(337-1186)	2000 per Case



#26411, Black, 22 Ga. x 1"	(337-1881)
#26412, Black, 22 Ga. x 1½"	(337-3130)
#26414, Green, 21 Ga. x 1"	(337-3044)
#26416, Green, 21 Ga. x 1½"	(337-2738)
#26431, Yellow, 20 Ga. x ½"	(337-6962)
#26417, Yellow, 20 Ga. x 1"	(337-9317)
#26418, Yellow, 20 Ga. x 1½"	(337-2180)
#26438, Brown, 19 Ga. x 1½"	(337-8419)
#26419, Pink, 18 Ga. x 1"	(337-1562)
#26420, Pink, 18 Ga. x 1½"	(337-2091)





IV Bags Should be Safe for All

Protecting Patients and Our Environment

More than four decades ago, B. Braun recognized the patient and environmental risks posed by products containing PVC and DEHP. Medical devices containing these chemicals pose potential health hazards to patients and to the environment during waste disposal. As a result, B. Braun was one of the first medical device manufacturers to remove these harmful substances from many of our products, and continue to lead the way in offering a full line of IV Solutions, Pre-Mixed Drugs, and Irrigation Solutions containers that are manufactured using biologically inert, non-toxic plastic materials.

Contact your Henry Schein Representative for further information

Rx Only. ©2024 B. Braun Medical Inc., Bethlehem, PA. All Rights Reserved. 24-0135



I.V. MEDICATION DELIVERY

From A Name You Can Trust.

Medication Delivery Sets and Connectors

Henry Schein provides market-leading infection control technologies to help minimize bacterial contamination. We carry a wide selection of injection and needle-free medication delivery ports to meet a variety of medication delivery needs.

Features & Benefits:

- Flat top needle-free valve with tight septum seal helps to support effective disinfection
- Metal free for MRI compatibility
- DEHP-free, non-latex and sterile

Item Code	Type	Product Description	Qty/Unit
(570-0881)	Primary I.V. Administration Set	I.V. Administration Set. 15 drops/mL. 73" Vented, Slide Clamp, Roller Clamp, Injection Y-site, Rotating Male Luer Lock	Each
(570-0885)	Primary I.V. Administration Set	I.V. Administration Set. 15 drops/mL. 96" Vented, Roller Clamp, Needle-free Y-site, Rotating Male Luer Lock	Each
(570-0887)	Primary I.V. Administration Set	I.V. Administration Set. 15 drops/mL. 105" Vented, Back Check Valve, Slide Clamp, (2) Needle-free Y-sites, Roller Clamp, Rotating Male Luer Lock	Each
(570-0883)	Secondary I.V. Administration Set	Secondary I.V. Administration Set. 15 drops/mL. 40" Vented, Roller Clamp, Hanger, Rotating Male Luer Lock	Each
(570-0884)	Microbore Extension set	Microbore extension set. 8". Removeable Needle-free Connector, Slide Clamp, Rotating Male Luer Lock	Each
(570-0882)	Needle-free Connector	Needle-free Connector	Each

SAFETY IV CATHETER

From A Name You Can Trust.

Item Number	Description	Quantity
(570-2586)	24G x 0.75"	50/box
(570-2591)	24G x 0.5"	50/box
(570-2587)	22G x 1"	50/box
(570-2588)	20G x 1.16"	50/box
(570-2589)	20G x 1"	50/box
(570-2590)	18G x 1.16"	50/box

Triple Facet Back Out Needle
Provides ease of insertion and reduced trauma to puncture site

Push Off Tab
Assists with catheter advancement

Push Button Mechanism
Retracts needle fully and safely into catheter chamber

Fully Encapsulating Needle Shield
Reduces exposure to blood borne pathogens via needle stick injuries

Gently Tapered Catheter Tip
Provides smooth transition from needle to catheter during insertion

Contoured Housing
Ergonomically designed grips allow for one-handed or two-handed techniques

Large-Volume Chamber
Allows clear visual confirmation through flash-back and easy handling

Needle
(Featured in sizes 20g - 24g)
The notched needle promotes fast confirmation of venipuncture

Catheter Material
PUR material designed to remain kink free and soften within the body



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!



ABRYVO™
PFIZER
 Respiratory Syncytial Virus Vaccine
 0.5-mL Prefilled Syringe
(145-4487).....ea
(145-4488).....5/box
 CPT® Code: 90678⁺



Adacel®
SANOFI
 Tetanus Toxoid, Reduced Diphtheria Toxoid, and
 Acellular Pertussis Vaccine Adsorbed
 0.5-mL SDV
(546-7789).....10/pkg
 0.5-mL Prefilled Syringe Without Needle
 Not made with natural rubber latex.
(132-5841).....5/pkg
 CPT® Code: 90715⁺



AREXVY
GSK
 Respiratory Syncytial Virus Vaccine, Adjuvanted
 0.5-mL SDV
(145-4809).....10/box
 CPT® Code: 90679⁺



BEXSERO®
GSK
 Meningococcal Group B Vaccine
 0.5-mL Prefilled Syringe
(127-2151).....10/pkg
 CPT® Code: 90620⁺



BOOSTRIX®
GSK
 Tetanus Toxoid, Reduced Diphtheria, and Acellular
 Pertussis Vaccine, Adsorbed
 0.5-mL SDV
(254-6475).....10/pkg
 0.5-mL Tip-Lok® Prefilled Syringe
(254-0031).....10/pkg
 CPT® Code: 90715⁺



DAPTACEL®
SANOFI
 Diphtheria and Tetanus Toxoids and Acellular
 Pertussis (DTaP) Vaccine (Adsorbed), Thimerosal-
 Free
 0.5-mL SDV
(546-1136).....10/pkg
 CPT® Code: 90700⁺



Engerix-B®–Pediatric/Adolescent
GSK
 Hepatitis B Recombinant Vaccine
 10 µg, 0.5-mL Syringe
(254-0026).....10/pkg
 CPT® Code: 90744⁺




ENERIX-B®–Adult
GSK
 Hepatitis B (Recombinant) Vaccine,
 Preservative-Free
 20 µg/mL, 1-mL SDV
(254-8254).....10/pkg
 20 µg/mL, 1-mL Tip Lok™ Syringe without Needle
(254-0029).....10/pkg
 CPT® Code: 90746⁺




GARDASIL® 9
MERCK
 Human Papillomavirus 9-Valent Vaccine,
 Recombinant
 0.5-mL Prefilled Luer-lock Syringe
(558-0044).....10/pkg
 CPT® Code: 90651⁺
 Gardasil® 9 is a registered trademark of Merck & Co., Inc.



HAVRIX®–Adult
GSK
 Hepatitis A (Inactivated) Vaccine, Thimerosal-Free
 1440 ELU/mL, 1-mL Tip-Lok® Prefilled Syringe
(124-0016).....10/pkg
 CPT® Code: 90633⁺



HEPLISAV-B™ Prefilled Syringe
DYNAVAX
 Hepatitis B Vaccine (Recombinant), Adjuvanted
 Solution for Intramuscular Injection
 0.5-mL Prefilled Syringe
 20 mcg/0.5 mL
(132-9264).....5/box
 CPT Code: 90739⁺



Hibrix[®]
GSK
 Haemophilus b Conjugate Vaccine (Tetanus Toxoid Conjugate)
 0.5-mL Prefilled Syringe
(146-1811) 10/pkg
 CPT[®] Code: 90648[†]



INFANRIX[®]
GSK
 Diphtheria and Tetanus Toxoids and Acellular Pertussis (DTaP) Vaccine (Adsorbed), Thimerosal-Free
 0.5-mL Prefilled Syringe
(124-0010) 10/pkg
 CPT[®] Code: 90700[†]



IPOL[®]
SANOFI
 Poliovirus (Salk) (Inactivated) Vaccine
 5-mL, 10-Dose Vial
(546-7189) ea
 CPT[®] Code: 90713[†]
 No credit will be given for return of this product.



KINRIX[™]
GSK
 Diphtheria and Tetanus Toxoids and Acellular Pertussis Adsorbed and Inactivated Poliovirus Vaccine
 0.5-mL Prefilled TipLok Syringe
(254-0025) 10/pkg
 CPT Code: 90696[†]



MenQuadfi[™]
SANOFI
 Meningococcal (Groups A, C, Y, W) Conjugate Vaccine Solution for Intramuscular Injection
 0.5-mL SDV
(140-3989) 5/carton
 CPT[®] Code: 90619[†]



MENVEO One-Vial
GSK
 Meningococcal (Groups A, C, Y, and W-135) Oligosaccharide Diphtheria CRM197 Conjugate Vaccine
 0.5-mL SDV
(145-4145) 10/pkg
 CPT[®] Code: 90734[†]




M-M-R[®] II
MERCK
 Measles, Mumps, and Rubella Vaccine
 0.5-mL SDV with Diluent
(558-0110) 10/pkg
 CPT[®] Code: 90707[†]
 M-M-R[®] II is a registered trademark of Merck & Co., Inc.



PEDIARIX[®]
GSK
 Diphtheria and Tetanus Toxoids and Acellular Pertussis (DTaP) (Adsorbed), Hepatitis B (Recombinant), and Inactivated Poliovirus Vaccine (Combined)
 0.5-mL Tip-Lok[®] Prefilled Syringe
(254-0028) 10/pkg
 CPT[®] Code: 90723[†]



PedvaxHIB[®]
MERCK
 Haemophilus b Conjugate (Meningococcal Protein Conjugate) Liquid Vaccine
 0.5-mL SDV
(558-3778) 10/pkg
 CPT[®] Code: 90647[†]
 PedvaxHIB[®] is a registered trademark of Merck & Co., Inc.



PENBRAYA[™]
PFIZER
 Meningococcal (Groups A, B, C, W, and Y) Vaccine.
 0.5-mL Prefilled Syringe
(146-0233) 1/pkg
(146-0247) 5/pkg
 CPT[®] Code: 90623[†]

**Pentacel®
SANOFI**

Diphtheria and Tetanus Toxoids and Acellular Pertussis Adsorbed, Inactivated Poliovirus and Haemophilus b Conjugate (Tetanus Toxoid Conjugate) Vaccine

0.5-mL SDV
(112-3585).....5/pkg
CPT® Code: 90698*



**PNEUMOVAX® 23
MERCK**

Pneumococcal Vaccines (Polyvalent)

0.5-mL Prefilled Syringe
(122-2087).....10/pkg

CPT® Code: 90732*
Pneumovax® 23 is a registered trademark of Merck & Co., Inc.



**PreHevbrio
VBI VACCINES**

Hepatitis B Vaccine (Recombinant)

10 mcg/mL, 1-mL SDV
(142-3464).....10/box
CPT® Code: 90759*



**PREVNAR 20™
PFIZER**


Pneumococcal 20-Valent Conjugate Vaccine

0.5-mL Prefilled Syringe
(140-8153).....ea
(140-8152).....10/box
CPT® Code: 90677*

**PRIORIX
GSK**

Measles, Mumps, and Rubella Vaccine, Live


0.5-mL SDV
(142-3773).....10/box
CPT® Code: 90707*



**ProQuad®
MERCK**

Measles, Mumps, Rubella, and Varicella Virus Vaccine, Live, 0.5-mL SDV

0.5-mL SDV
(558-0053).....10/pkg
CPT® Code 90701*
ProQuad® is a registered trademark of Merck & Co., Inc.



**Quadracel™
SANOFI**

Diphtheria and Tetanus Toxoids and Acellular Pertussis Adsorbed and Inactivated Poliovirus Vaccine

0.5-mL SDV
(141-7656).....10/pkg
0.5-mL Prefilled Syringe
(141-7655).....10/pkg
CPT® Code: 90696*

**RECOMBIVAX HB® –Adult
MERCK**

Hepatitis B (Recombinant) Vaccines, Thimerosal-Free

10 µg/mL, 1-mL SDV
(558-0737).....10/pkg
1-mL Prefilled Luer-Lock Syringe Without Safety Needle
(122-2086).....10/pkg

CPT® Code: 90746*
RECOMBIVAX HB® is a trademark of Merck & Co., Inc.



**RECOMBIVAX HB®
Pediatric/Adolescent
MERCK**

Hepatitis B (Recombinant) Vaccine, Thimerosal-Free

5 µg/mL, 0.5-mL SDV
(558-6254).....10/pkg
0.5-mL Prefilled Luer-Lock Syringe without Safety Needle
(122-2084).....10/pkg

CPT® Code: 90744*
RECOMBIVAX HB® is a registered trademark of Merck & Co., Inc.



**ROTARIX®
GSK**

Rotavirus Vaccine, Live, Oral

1.5-mL Pediatric Prefilled Syringe
(143-4023).....10/pkg
CPT® Code: 90681*


RotaTaq[®]
MERCK


Rotavirus Vaccine (Live), Oral, Pentavalent

2-mL Single-Dose Tube

(558-8763).....10/pkg

(558-0009).....25/pkg


 CPT[®] Code: 90680[†]

 RotaTaq[®] is a registered trademark of Merck & Co., Inc.

SHINGRIX[®]
GSK

Zoster Vaccine Recombinant, Adjuvanted

0.5-mL SDV with Diluent

(129-6729).....10/pkg

 CPT[®] Code: 90750[†]

TDVAX[™]
GRIFOLS

Tetanus and Diphtheria Toxoids, Adsorbed

7+ Years, Not made with Natural Rubber Latex and Preservative-Free

0.5 mL SDV

(125-9839).....

 CPT[®] Code: 90714[†]


Trumenba[®]
PFIZER

Meningococcal Group B Vaccine

0.5-mL Prefilled Syringe

(248-0649).....5/pkg


(248-0648).....10/pkg

 CPT[®] Code: 90621[†]

TWINRIX[®] Adult
GSK

Hepatitis A (Inactivated) and Hepatitis B (Recombinant) Vaccine

 1-mL Tip-Lok[®] Prefilled Syringe without Needle

(254-0030).....10/pkg

 CPT[®] Code: 90636[†]

VAQTA[®] Adult
MERCK

Hepatitis A (Inactivated) Vaccines, Thimerosal-Free

50 U/1 mL Prefilled Luer-lock Syringe Without Safety Needle

(122-2081).....10/pkg

 CPT[®] Code: 90632[†]

VAQTA is a Registered Trademark of Merck & Co., Inc.


VAQTA[®] –Pediatric/Adolescent
MERCK

Hepatitis A (Inactivated) Vaccine, Thimerosal-Free, 25 U, 0.5-mL SDV

0.5-mL Prefilled Luer-lock Syringe Without Safety Needle

(122-2080).....10/pkg

 CPT[®] Code: 90633[†]

 VAQTA[®] is a registered trademark of Merck & Co., Inc.


VARIVAX[®]
MERCK

Varicella (Chicken Pox) Virus (Live) Vaccine, Thimerosal-Free, Lyophilized

0.5-mL SDV

(558-1592).....10/box

 CPT[®] Code: 90716[†]

 Varivax[®] is a registered trademark of Merck & Co., Inc.

Vaxelis[™]
SANOFI


Diphtheria and Tetanus Toxoids and Acellular Pertussis, Inactivated Poliovirus, Haemophilus b Conjugate and Hepatitis B Vaccine

0.5-mL SDV

(140-3160).....10/pkg

0.5-mL PFS

(140-3161).....10/pkg

 CPT[®] Code: 90697[†]

VAXNEUVANCE[™]
MERCK

Pneumococcal 15-Valent Conjugate Vaccine

(141-1002).....10/box

 CPT[®] Code: 90671[†]

 Vaxneuvance[®] is a trademark of Merck & Co., Inc.



Adrenalin® Injection
PAR STERILE PRODUCTS

1 mg/mL, 1-mL SDV
(248-0883)..... ea
(890-0257).....25/pkg
1 mg/mL, 30-mL MDV
(129-4847).....ea



Compare to ProAir®
Albuterol Sulfate Inhaler HFA
LUPIN

90 mcg
(140-9653)..... ea




Aplisol® Injectables
PAR STERILE PRODUCTS

Tuberculin PPD Stabilized Solution, 5 TU/0.1 mL
(110-5199).....1-mL Vial
(110-5099).....5-mL Vial



Auvi-Q Epinephrine Auto-Injectors
KALEO

0.1 mg
(142-4090).....2/pkg
0.15 mg
(142-4091).....2/pkg
0.3 mg
(142-4092).....2/pkg



Bicillin® L-A Injectable*
Pfizer

Penicillin G benzathine IM Suspension

Pediatric, 21 ga x 1"
600,000 U/mL, 1 mL
(321-6704).....10/pkg
Adult, 21 ga x 1½"
1,200,000 U/mL, 2 mL
(321-7445).....10/pkg
Adult, 18 ga x 1½"
2,400,000 U/mL, 4 mL
(321-1663).....10/pkg

*Tubex injector is not required.



Compare to Rocephin
Ceftriaxone Sodium Injectable
APOTEX

500-mg SDV
(135-9645).....10/box
1-gram Vial
(133-6436).....10/box
2-gram Vial
(135-7800).....10/box



Celestone® Soluspan® Injectable Suspension
ORGANON

6 mg/mL, 5-mL MDV
(140-7856)..... ea



Compare to Rubramin PC®
Cyanocobalamin (Vitamin B-12)
Injectable

AMERICAN REGENT
1000 mcg/mL
1-mL MDV
(244-1902).....25/pkg



Depo-Provera® Injectable
Pfizer

Medroxyprogesterone Acetate Suspension
150 mg/mL, 1-mL Prefilled Syringe
(908-0028)..... ea



Durolane® Injectable
BIOVENTUS
 60 mg/3 mL Prefilled Syringe
 (137-4043).....ea
 Series of injections: 1
 J Code: J7318



Compare to Adrenalin®
Epinephrine Injectable Ampules
BPI LABS
 1:1000, 1 mL
 (248-0899).....ea
 1:1000, 1 mL
 (128-9991).....10/box



Ethyl Chloride Topical Anesthetic
GEBAUER
 Skin Refrigerant, 3½-oz Glass Bottle with Nozzle
 Fine-Stream Spray, 3.9 oz
 (138-9019).....ea
 Medium-Stream Spray, 3.9 oz
 (138-5777).....ea
 Mist Spray Can, 3.9 oz
 (138-6406).....ea



Euflexxa® Prefilled Syringe
FERRING PHARMACEUTICALS
 2 mL
 (137-4434).....3/pkg
 Series of injections: 3
 J Code: J7323



Evenity® Injection
AMGEN
 105 mg/1.17 mL
 (214-0051).....2/box



Compare to Fluresc®
Fluorescein Sodium and Benoxinate Hydrochloride Ophthalmic Solution, 0.25%/0.4%
ALTAIRE PHARMACEUTICALS
 5-mL Bottle with Sterile Dropper
 (131-3074).....ea



Gelsyn-3™ Injection
BIOVENTUS
 Prefilled Syringe
 (136-8849).....ea
 Series of injections: 3
 J Code: J7328



HyperHEP B™ Injectable
GRIFOLS
 5-mL SDV
 (271-7848).....ea



Kenalog®-10 Injectable
BRISTOL MYERS SQUIBB
 Triamcinolone Acetonide
 10 mg/mL, 5-mL MDV
 (196-0545).....ea



Kenalog®-40 Injectables
BRISTOL MYERS SQUIBB
 Triamcinolone Acetonide
 40 mg/mL
 1-mL SDV
 (196-8300).....ea
 5-mL MDV
 (196-1737).....ea
 10-mL MDV
 (196-9429).....ea



Ketorolac Injection
FOSUN PHARMA
 15 mg/mL
 1-mL SDV
 (146-1679).....25/box
 30 mg/mL
 1-mL SDV
 (143-6324).....25/box
 2-mL SDV
 (143-6293).....25/box



Kyleena® IUD System
BAYER
 19.5 mg
 (126-7908)..... ea



Mirena® IUD System
BAYER
 52 mg
 (125-0996)..... ea




Compare to Xylocaine®
Lidocaine HCl Injection
PFIZER
 50-mL Vial MDV
 1%
 (248-0874).....ea
 2%
 (248-0973).....ea
 (104-9843)..... 25/box



Medroxyprogesterone
MYLAN INSTITUTIONAL
Compare to Depo-Provera®
 150 mg/mL, 1-mL SDV
 (146-6392)..... ea
 (146-6382)25/box



Compare to Depo-Provera®
Medroxyprogesterone Injectable
EUGIE
 150 mg/mL Prefilled Syringe
 (146-4156)..... ea



Compare to Depo-Medrol®
Methylprednisolone Acetate Injection
AMNEAL
 40 mg/mL, 1-mL SDV
 (140-9021).....ea
 40 mg/mL, 1-mL SDV
 (140-9026)..... 25/box
 80 mg/mL, 1-mL SDV
 (140-9048)..... 25/box



Monovisc®
DEPUY SYNTHES
 4 mL, 22 mg/mL, Prefilled Syringe
 (729-0037).....ea
 Series of injections: 1
 J Code: J7327



Ondansetron HCl Injection
APOTEX
 2 mg/mL, 2-mL SDV
 (248-0902)..... ea
 (125-9100)..... 25/box



Orthovisc®
DEPUY SYNTHES
 2-mL Single-Use Syringe
 (729-6418)..... ea
 Series of injections: 3 to 4
 J Code: J7324




Pain Ease® Mist Spray
GEBAUER
Topical Anesthetic Skin Refrigerant
 3.9-oz Mist Spray
 (791-0634)..... ea
 1-oz Medium-Stream Spray
 (132-9467).....ea
 3.9-oz Medium-Stream Spray
 (761-0029)..... ea
 Can be used in facilities that restrict use of flammable components.



Prolia Injectable Prefilled Syringe
AMGEN
 60 mg/mL
 (214-0006).....ea



Compare to Lexiscan®
Regadenoson Injection
HIKMA
 0.4 mg/5 mL, 5-mL Prefilled Syringe
 (143-6225).....ea



RhoGAM® Ultra-Filtered Plus
KEDRION
 Rho(D) Immune Globulin (Human), 300-µg Single-Dose Safety-Shield Syringes
 (623-0238).....ea
 (623-0237).....5/pkg
 (623-0235).....25/pkg
 CPT® Code: 90742†



Supartz FX™ Injectable
BIOVENTUS
 2.5-mL Prefilled Syringe
 (129-3653).....ea
 Series of injections: 3 to 5
 J Code: J7321



Compare to Depo-Testosterone
Testosterone Cypionate Injection
HIKMA
 200 mg/mL, 10-mL MDV
 (137-7452).....ea



Compare to Kenalog®-40
Triamcinolone Acetonide Injection
AMNEAL
 40 mg/mL
 1-mL SDV (142-2167).....ea
 1-mL SDV (142-3667).....25/box
 5-mL MDV (142-3668).....ea
 10-mL MDV (142-2227).....ea



Xylocaine® Plain Injectable
FRESENIUS KABI
 1%, 5-mL SDV, MPF (150-0069).....25/box
 1%, 30-mL SDV (150-0072).....25/pkg
 1%, 50-mL MDV (248-0849).....ea
 2%, 5-mL SDV (150-0101).....25/pkg



Zilretta™ Injectable
PACIRA
 32-mg SDV
 Triamcinolone Acetonide Extended-Release
 (142-9357).....ea
 Series of injections: 1
 J Code: J3304



Zynrelef® Injection
HERON THERAPEUTICS
 200 mg/6 mg, 10-mL SDV (140-5769).....10/pkg
 400 mg/12 mg, 20-mL SDV (140-5768).....10/pkg

Compare to Tylenol®
Acetaminophen
Regular-Strength Tablets
Unit-Dose Package



Nonaspirin pain reliever and fever reducer.

- 2 tablets per pouch
- Convenient dispenser box
- Made in the USA

325 mg.
(900-4689).....500 per Box
 Contains: 250 packets of 2.



Compare to Extra Strength Tylenol®
Acetaminophen
Extra-Strength Tablets
Unit-Dose Package



Nonaspirin pain reliever and fever reducer.

- Manufactured and packaged in the USA
- Highest quality-control standards under the FDA
- Packaged in 250 sealed-in pouches of 2 tablets each
- Convenient dispenser box
- Color-coded inner and outer packaging

500 mg.
(900-4690)..... 500 per Box
 Contains: 250 packets of 2.



Compare to Advil® or Motrin®
Ibuprofen Tablets
Unit-Dose Package



Powerful relief from pain and inflammation.

200 mg
(570-0177).....500 per Box
 Contains: 250 packets of 2 tablets.



Bacitracin Zinc Ointment, USP



Helps prevent infection in minor cuts and burns.
 Active ingredient: bacitracin zinc, 500 units.

1 oz. Tube
(570-3129).....Ea.
 0.9-Gm. Foil Packs
(570-3125).....144 per Box



Isopropyl Alcohol, 70%
HYDROX LABORATORIES

4 oz
(128-0491).....Ea.



Alcohol



70% Isopropyl
Plastic Container

(112-7067)..... 16 oz.
(112-7011)..... 32 oz.
(112-7068)..... Gal. ♦

♦ **Hazardous material**



Lubricating Jelly
HR PHARMACEUTICALS

Sterile

Product features:

- Sterile
- Microflora- and sperm-friendly
- Safe with Pap smears
- 6 FDA registrations
- 80-year trustworthy market history
- Kosher certified
- Made in the USA
- 90-day hold period post manufacturing
- Doubles as ultrasound gel
- Great viscosity for better resolution
- Water soluble for quick, easy cleanup
- Fragrance-, salt-, and alcohol-free
- Nontoxic and nonirritating
- Dye-free and not made with natural rubber latex
- PEG- and chlorhexidine-free


3-g Packets
(118-2072)..... 144 per Box
 5-g Packets
(118-2073)..... 144 per Box
 1¼-oz Packets
(121-6827)..... 48 per Box
 2-oz Flip-Cap Tubes
(118-7249)..... 12 per Pkg.
 4-oz Flip-Cap Tubes
(118-1350)..... 12 per Pkg.



Lubricating Jelly
HENRY SCHEIN®

Sterile

3 Gm. Pack
(570-3105)..... 144 per Box
 4-oz Flip-Top Tube
(570-3106).....4 oz. Tube



Compare to Dacriose®, Collyrium Eye Lotion
Medi-Wash Eye & Skin Irrigant
MEDIQUE PHARMACEUTICALS

Sterile, isotonic, buffered solution flushes foreign particles. Single use.

Plastic Bottle
(118-2533)..... 1 oz.
(151-4399).....4 oz.
(151-0009)..... 16 oz.

Active ingredients: sodium chloride, potassium chloride, calcium chloride, magnesium chloride, sodium acetate, sodium citrate, and preserved with 0.013% benzalkonium chloride.



Compare to Halls®
Medikoff® Drops–Bulk Refill
MEDIQUE

Menthol/eucalyptus vapor-action cough drops.

Cherry
(151-4400)..... 600 per Box

Compare to Naprosyn®, Aleve®
Mediproxen Tablets
Unit-Dose Package
MEDIQUE

Has the same pain-relief and fever-reduction ingredient found in nationally advertised brands Naprosyn® and Aleve®. Unit-dose packets are safe, hygienic, and tamper-resistant. Naproxen sodium, 220 mg.

(993-0948).....50 per Box
 Contains: 50 packets of 1.
 Naproxen sodium, 220 mg.
(151-9200)..... 100 per Box
 Contains: 100 packets of 1.





Aquaphor® Healing Ointment BEIERSDORF

Original Formula

Treats dry, damaged skin and minor burns.
Hypoallergenic: no fragrances or preservatives
that can irritate sensitive skin. Great for
turf burns.

- Dermatologist recommended
- Advanced therapy for dry, cracked, or irritated skin
- Clinically proven to restore smooth, healthy skin

1.75 oz. Tube

(911-0837)..... Ea.

#63608, 14-oz Jar

(911-3536)..... 14 oz. Jar



Aquaphor® Healing Ointment Unit-Dose Package BEIERSDORF

0.35-oz Packets

(911-0011)..... 144 per Box



Boudreaux's Butt Paste® PRESTIGE CONSUMER HEALTHCARE

4-oz Tube

(131-8864)..... Ea.



Compare to Blistex®

Lip-Guard MEDI-FIRST

Temporarily protects and helps
soothe chapped or cracked lips.
Provides temporary relief of pain and itching
associated with minor lip irritations and cold
sores. Contains 4 medications.

1/57-oz. Packets

(151-5789)..... 20 per Box

Contains: 20 packets of 1/57 oz.



Nix® Lice Killing Spray

PRESTIGE CONSUMER HEALTHCARE

5-oz Bottle

(131-8954)..... Ea.



Nix® Lice Crème Rinse

PRESTIGE CONSUMER HEALTHCARE

2-oz Bottle

(131-8942)..... Ea.

(131-8943)..... 2 per Pkg.



Nix Ultra®

Lice Removal Kit

PRESTIGE CONSUMER HEALTHCARE

(136-9432)..... Ea.

Contains: 1 of each: Nix Ultra® Lice Removal Comb, 3.4-oz
Nix Ultra® solution & 5-oz Nix® Lice Killing Spray.



Nix Ultra® Lice Treatment

PRESTIGE CONSUMER HEALTHCARE

- Kills lice and their eggs
- Kills super lice
- Fast and effective
- Pesticide-free and non-toxic
- Easy application
- Tough on lice; gentle on skin

Lice Shampoo

(135-5184)..... 4 oz per Box

Lice Treatment Liquid

(135-5618)..... 3.4 oz per Box



Polident® Denture Essential Kit

HALEON

Designed for patients to get started right away
with proper denture care. Polident® microclean
formula removes tough stains and kills 99.9%
of odor-causing bacteria. Super Poligrip® is
clinically proven to seal out more food particles,
even in well-fitting dentures.

Starter Kit

(146-3885) 12 Kits per Case

Contains: 1–0.3 oz. tube of Super Poligrip® adhesive,
6 Polident® 3-minute cleansing tablets, 1 denture brush
& 1 denture bath.



Safetec Sunscreen Lotion

SAFETEC OF AMERICA INC.

SPF 30, 1/8-oz Packets

(115-5375)..... 25 per Box



Polident® Tablets

GLAXOSMITHKLINE CONSUMER HEALTHCARE

Designed to fight stains and kill 99.9% of
odor-causing bacteria. Clean dentures daily by
soaking and brushing. Offers a complete clean
feeling with triple mint freshness.

3-Minute

(583-1271)..... 40 per Box

(583-0396)..... 84 per Box

Overnight Whitening

(583-1293)..... 40 per Box

(583-0339)..... 84 per Box

HENRY SCHEIN®
BRAND PRODUCTS

Rely on Us™

Sharps Containers

From a Name You Can Trust.

- For use by professional health-care givers
- Meets all NIOSH and OSHA regulations
- FDA approved
- Puncture-resistant
- Transparent body and/or top allows for easy viewing
- High-impact and shatter-resistant polypropylene



#0310-1500-HS, 1 qt
(112-7149) ea

#0310-150T-HS, 1.5 qt
(570-0376) ea

#0322-150R-HS, 2.2 qt
(570-0213) ea

#0354-150M-HS, 5 qt
(112-7152) ea

#0354-150B-HS, 5.4 qt
(112-7151) ea

#0319-1500-HS, 1 gal, Slide Lid
(570-0378) ea

#0319-150R-HS, 1 gal, Rotary Lid
(570-0379) ea

#0320-1500-HS, 2 gal, Slide Lid
(112-7150) ea

#0320-150R-HS, 2 gal, Rotary Lid
(570-0377) ea

#0330-150R-HS, 3 gal
(570-0212) ea

#0380-150R-HS, 8 gal
(570-0380) ea



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

Sharps Container AP MEDICAL

- 4½"H x 12"W x 10¾"L
- Body: 1.8mm; Lid: 1.53mm
- Counterbalance lid
- Irreversible permanent closure
- Maximum filling level indicated
- Resistant to sharp and cutting objects
- Red, 5.4 Qt.

#242468
(139-8802) Ea.



Sharps Containers DYNAREX

- Durable and puncture-resistant
- Horizontal drop maximizes container volume
- Locks for final disposal
- Transparent lid helps reduce risk of overfilling
- Red

Slide Lid
#4623, 2 qt
(139-8161) Ea.

Mailbox Lid
#4624, 5 qt
(139-8162) Ea.



Sharpstainer® Disposable Container MEDICAL ACTION

Red, autoclavable, and puncture-resistant. Tamper-proof locking cap. One-way safety insert.

Small—3½" x 3½" x 5" (0.7 Qt. Capacity)
(153-3569) Ea.

Tall Tray—1.7 Qt.
(153-8085) Ea.

10 qt
(120-4364) Ea.



Phlebotomy Collector BD

One-step system eliminates recapping. Rigid, puncture-resistant polypropylene container. May be autoclaved and compacted or incinerated.

Tray Size—1.5 Qt.
(987-5812) Ea.

i.Series® iPR120-GX Pharmacy Refrigerator—Helmer Scientific HELMER SCIENTIFIC

20.2-cu-ft Capacity, Glass Door, 2 to 10°C

GX Solutions and the i.C3® Information Center set i.Series® models apart by providing a professional medical-grade refrigerator with advanced monitoring and control. This pharmaceutical-grade refrigerator is designed for storing medications, vaccines, and other pharmaceuticals. It is NSF/ANSI 456 vaccine storage-certified and provides superior stability, uniformity, and recovery. It also has the added benefits of quiet operation, low heat output, and energy savings, all thanks to the advanced OptiCool® technology.

- Provides continuous temperature monitoring, event logs, temperature graphs, and more
- Designed specifically for the central pharmacy, medication rooms, clean rooms, satellite locations, and outpatient clinics
- 50–65% more efficient than traditional technology
- Quiet operation (47 dB) to create fewer distractions and better working environment
- Audible and visual alarms include power failure, high- and low-temperature alarms, door ajar, condenser temperature, and low battery
- Includes Certificate of Calibration generated with an ISO 17025 reference thermometer
- Energy Star® certified

#5115120-1

(141-8672) Ea.



i.Series® iLF105-GX Laboratory Freezer—Helmer Scientific HELMER SCIENTIFIC

5.3-cu-ft Capacity, Solid Door, -15 to -30°C

A GX Solutions professional medical-grade, -30°C freezer.

Designed for the unique needs of critical health care applications, such as storing samples, reagents, controls, and frozen vaccines. It delivers optimized control in 3 important areas: temperature, noise, and energy management.

- Powered by OptiCool® technology, this freezer provides improved performance and user experience through energy efficiency, low noise levels, and stable storage environments to safeguard valuable patient samples, frozen vaccines, and more
- 30% more efficient than traditional technology
- Digital microprocessor temperature controller, with built-in alarms, includes minimum/maximum display and reset
- Energy Star® certified
- Powerful i.C3® Information Center is exclusive to i.Series® models and provides constant temperature monitoring, multiple information logs, and security features
- GX Solutions is the first professional medical-grade line of freezers for clinical laboratories, pharmacies, and other health care facilities
- NSF/ANSI 456 Vaccine Storage Standard certified

#5222105-1

(141-9439) Ea.



office supplies



Trash Liners

HENRY SCHEIN®

These liners are a strong, durable, and affordable option, and fit many different shapes and sizes of trash cans. Made in the USA.

Clear

7–10 Gal., 6 mil

(570-1340) 1000 per Case

15–20 Gal., 8 mil

(570-1346) 1000 per Case

20–30 Gal., 0.39 mil

(570-1336) 250 per Case

20–30 Gal., 10 mil

(570-1341) 250 per Case

33 Gal., 13 mil

(570-1342) 250 per Case

40–45 Gal., 12 mil

(570-1343) 250 per Case

Black

33 Gal., 0.51 mil

(570-1337) 250 per Case

40–45 Gal., 0.51 mil

(570-1338) 250 per Case

60 Gal., 1.4 mil

(570-1339) 100 per Case

White

13 Gal., 0.8 mil

(570-1344) 100 per Case

20–30 Gal., 0.85 mil

(570-1345) 250 per Case

Lysol® Disinfecting Wipes

Helps stop the spread of germs and are perfect for touch-up cleaning. These strong premoistened wipes kill germs, disinfect and deodorize. Convenient pop-up container.

Lemon & Lime Blossom®

(903-8719) 80 Wipes/PK



Lysol® Professional Disinfectant Sprays

Kills 99.9% of germs on hard, nonporous surfaces. Effective against over 100 microorganisms when used as directed. Eliminates odors and prevents the growth of mold and mildew. EPA registered. 19-oz aerosol spray.

Crisp Linen Scent

(902-9209) EA

Fresh Scent

(902-6347)*** EA



***This product is not available for customers in Alaska and Hawaii

ShermanSpecialty COMPANY, INC.

TAKE HOME BAGS • STICKERS • TOYS

*Ask Your Henry Schein Representative
For More Fun Patient Giveaways!*



STANDARD TIMERS

50 per pack
(366-5162) - 2 Minute
(366-6564) - 3 Minute



TOOTH-SHAPED TIMERS

50 per pack
(367-6826) - 2 Minute
(367-0491) - 3 Minute



CLEAR SCATTER BAGS

7½" x 10" - 100 per pack
(367-5241) - Blue Teeth
(367-2869) - Cute Tooth Guy



CLEAR SCATTER BAGS

7½" x 10" - 100 per pack
(368-0795) - Stars and Stripes
(368-4026) - Bones and Yarn



PAPER BAGS

7½" x 10" - 100 per pack
(368-4036) - Flag Heart
(368-3493) - Thank You



PAPER BAGS

7½" x 10" - 100 per pack
(368-3494) - Smile Supplies
(368-3496) - Great Check-Up



WHITE PLASTIC BAGS

9" x 13" - 100 per pack
(367-4693) - Toothbrush Flag
(367-4698) - Stars Smile



PAPER PHARMACY BAGS

2"D x 5"W x 12"H - 100 per pack
(368-4051) - Patriotic Stars
(368-3864) - Thank You



PAPER "LUNCH" BAGS

3½"D x 5½"W x 11"H - 100 per pack
(368-4030) - Thank You Tooth
(368-4031) - Tooth Swirl



VINYL ZIPPER POUCHES

Assorted colored zips - 144 per pack
(367-5358) - 10" x 4"
(367-7937) - 10" x 6"



DENTAL SQUISHY MIX

Assorted styles. 2"-2½" - 50 per pack
(368-3668)



ANIMAL POPPER TOYS

Assorted styles - 24 per pack
(368-3756)



SQUISHY FUN ASSORTMENT

Assorted Styles. 2" - 50 per pack
(368-4010)



STANDARD TOY CHEST


200 assorted toys
(366-3455) - Standard Toy Mix Treasure
(366-0918) - Standard Toy Mix Refill




**BATMAN
STICKERS**
(368-0666)



**JUSTICE LEAGUE
STICKERS**
(368-0673)



**DESPICABLE ME 4
STICKERS**
(368-4062)



BARBIE STICKERS
(368-3914)



2½" STICKERS

6 Assorted Designs
100 Stickers per roll

POUCHES • SUNGLASSES • TREASURE CHESTS • RECALL CARDS • DENTAL TIMERS

⚠ WARNING: To avoid danger of suffocation, keep this plastic bag away from babies and children. Do not use this bag in cribs, beds, carriages, or play pens. This bag is not a toy. **CHOKING HAZARD** - All items contain/are a small part. Not for children under 3 years.



Chap Ice® Mini Lip Balms—Oralabs®

- A nice gift for post-op dry, chapped lips with vitamin E
- Contains 40 Kiwi, 40 Cherry, and 20 Citrus Orange

Assorted
(367-1110) 100 per Pkg.
5 Pkg. @



Lip Rageous® Mini Lip Balms—Oralabs®

- Delivers mega-moisture and flavor
- Formulated with aloe and vitamin E
- Contains 35 Blue Raspberry, 35 Watermelon, and 30 Strawberry

Assorted
(367-1112) 100 per Pkg.
5 Pkg. @



Plexiglas® Hand Mirror

 HENRY SCHEIN®

One side plain/one side magnifying. 5" diam.
(101-0055) Ea.



Dr. John's® Classic Fruit Collections

INGREDIENTS: Dextrin (dietary fiber), Erythritol, Xylitol, Citric Acid, Ascorbic Acid, Natural Flavors, Fruit and Vegetable Juice for Color, Spirulina Extract (Color), Turmeric (Color), Stevia Leaf Extract. Assorted Flavors: Cherry, Grape, Blue Raspberry, Orange, Lemon, and Lime. Individually wrapped.

(368-3271) 150 per Bag



Crayon Bubbles

(367-6087) 24 per Pkg.

 Small Parts




Tooth-Shaped Brushing Timers

Brushing timer lengths are approximate. 3½".

2-Min. Timer
(367-6826) 50 per Pkg.
3-Min. Timer
(367-0491) 50 per Pkg.

 Small Parts



(367-2869)

(367-7852)

(368-2323)


(367-5241)

(367-7588)

Scatter Print Bags

- Clear plastic bags with 2-sided print
- Transparent design allows you to see the supplies inside
- Great price!
- Available in one size: 7"W x 10"H

Scatter Bags 100 per Pkg.
5 Pkg.

 **WARNING: CHOKING HAZARD**

Item contains a small ball, small part and/or marble. Not for children under three years.





24SM2618

HENRY SCHEIN®
Henry Schein, Inc. • 135 Duryea Road • Melville, NY 11747

PRSRT STD
U.S. POSTAGE

PAID

HENRY SCHEIN, INC.


henryschein.com/smTo Order: **1-800-851-0400** 8am-8pm (et)

California customers please use website when ordering
for Prop 65 information.

© 2024 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. **To guarantee special pricing/offers, where applicable, please include the three-digit promo code when ordering.** Note: Participating in a promotional discount program (e.g., points, discount redemptions or other special awards) is only permissible in accordance with discount program rules. By participation in such program, you agree that, to your knowledge, your practice complies with the program requirements. For HS promotions: This promotion offers a bundled discount. The cost of any "buy/get" (or "no charge") items is a discount applied to your other purchases on the same invoice. The reportable value (unit price) of the "no charge" item will be identified on the invoice. To calculate the applicable discount on the other items on the same invoice, divide the total value (unit price(s)) of any no-charge item(s) by the total value (extension price(s)) of all other items purchased on the same invoice, and apply the resulting percentage equally to each item to determine net prices. For manufacturer promotions: The prices stated in the invoice or statement may reflect a discount or be subject to a rebate. You must fully and accurately report this stated discount price, or if applicable, any net pricing, after giving effect to any rebates, to Medicare, Medicaid, Tricare and any other federal or state program upon request by any such program. Accordingly, you should retain these records. It is your responsibility to review any agreements or other documents applicable to these prices to determine if they are subject to a rebate.

Built to Last

PROVEN PERFORMANCE,
ESTHETICS & LONGEVITY
BUILT ON TOOTHLIKE
PHYSICAL PROPERTIES



GrandISO

Nano-Hybrid Universal Composite

89% fillers by weight:

- Outstanding wear resistance
- Enhanced color stability

30-50% less resin compared to classic hybrid and pseudo-nano composite

- Low shrinkage of only 1.6 vol. %
- Non-sticky
- Prolonged working time under ambient light

Designer Nano particles:

- Excellent polishability and polish retention
- Tooth-like modulus of elasticity
- Tooth-like thermal expansion and flexural strength

Data on File at www.vocoamerica.com

ONE UNIVERSAL
COMPOSITE
FOR ALL YOUR
RESTORATIVE
NEEDS

MOST POPULAR SHADES

(999-1460) – Caps Kit	80 × 0.25 Gm.	
(999-1461) – Syringe Kit	5 × 4.00 Gm.	
(999-1407) – Capsules	16 × 0.25 Gm.	A1
(999-1408) – Capsules	16 × 0.25 Gm.	A2
(999-1409) – Capsules	16 × 0.25 Gm.	A3
(999-1410) – Capsules	16 × 0.25 Gm.	A3.5
(999-1411) – Capsules	16 × 0.25 Gm.	A4
(999-1417) – Capsules	16 × 0.25 Gm.	C2

(GrandISO is available in Caps and Syringe Refills in the following shades: A1, A2, A3, A3.5, A4, B1, B2, B3, C2, D3, Incisal, 0A1, 0A2, 0A3.5, Bleach, GA3.25 and GA5)

NO MORE
SLUMPING
OUT OF
CLASS V'S!



GrandISO Heavy Flow

Nano-Hybrid Universal Composite

- High filler content and exceptional physical properties enable indications beyond those of conventional flowables, including Class I-V restorations
- No undesirable runoff, does not slump; increased stability for easy application in comparison to conventional flowables
- High fill rate (83 % w/w) and low resin content, results in lower shrinkage (2.99 % v/v)*
- High radiopacity for easy x-ray identification
- Excellent polish and polish retention

Data on File at www.vocoamerica.com

HIGH VISCOSITY,
LIGHT-CURED
UNIVERSAL
COMPOSITE
SUITABLE
FOR CLASS I-V
RESTORATIONS

MOST POPULAR SHADES

(999-1441) – Capsules	16 × 0.25 Gm.	A1
(999-1448) – Syringes	2 × 2 Gm.	A1
(999-1442) – Capsules	16 × 0.25 Gm.	A2
(999-1449) – Syringes	2 × 2 Gm.	A2
(999-2162) – Syringes	2 × 2 Gm.	Incisal

GrandISO Heavy Flow is available in Caps and Syringe Refills in the following shades: A1, A2, A3, A3.5, A4, B1, B2, C2, 0A2, Incisal, Bleach, White Opaaue, GA5)