



DENTAL & MEDICAL SUPPLIES SOURCE

Fall 2024

Community Health & Safety Net Clinic



Rely On Us
for Business, Clinical,
Technology, and
Supply Chain Solutions.



**QuidelOrtho
Sofia Rapid Strep
A Test**
Page 11



**GOJO
PURELL Advanced
Aloe Instant Hand
Sanitizer**
Page 54



**3M
Filtek Easy
Match Universal
Restorative**
Page 73



**Colgate
Junior
Toothbrush**
Page 91



Special Markets

1.800.851.0400 8:00am – 8:00pm (et)
henryschein.com/sm



Medical

1.800.772.4346 8:00am – 8:00pm (et)
henryschein.com/medical

California customers please use website when ordering for Prop 65 information.

YOUR SOLUTIONS FOR INSTRUMENT REPROCESSING

HuFriedyGroup offers the products you need
to improve efficiency, safety, and compliance as your instrument
reprocessing goes from chairside to cleaning to sterilization to storage.



Learn More

Hu-Friedy Mfg. Co., LLC, 1666 E. Touhy Ave., Des Plaines, IL 60018 | HuFriedyGroup.com
All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies,
unless otherwise noted. Marks not registered in all jurisdictions.
©2024 Hu-Friedy Mfg. Co., LLC. All rights reserved. 2410

HuFriedyGroup
The Best In Practice

INDEX

A.Titan	83
Abbott	8, 12
ACCLEAN by Henry Schein	86
Accucold	31
ADHA	69
AED	42
Anesthetics	65-66
Arkray	89
Armor Dental	84
Baxter	37
BD	16-17
BeeSure	60
CareQuest	95
Cements & Liners	66-69
Cold Storage	31-32
Corza Medical	42
Cosmetic & Restorative Dentistry	69-74
Crown & Bridge	75-76
Cranberry	56
Dental Top Sellers	78-79
Diagnostic Equipment	37-39
DJO	42
DMG	69
Endodontics	77
EndoEdge by Henry Schein	84
Exam Room Disposables	36
Exam Room Equipment	33-35
GC America	70
Gloves	56-59
Graham Medical	36
Handpieces	80-82
Health & Beauty Aids (HBA)	28
HealthFirst	90
HemoCue	14
Henry Schein Brand Products	6-7, 18, 29-31, 45-47, 80
HuFriedyGroup	2, 48
Hypodermics	17-18
Infection Prevention	51-55
Instruments	84-85
Instrument Reprocessing	48
IV Medication Delivery	30-31
Kerr	67
Kuraray	69
Matrix Materials	83
Metrex	51
Midmark	33
Office Supplies	92
Orthopedics	42
Over-the-Counter (OTC)	27
Patient Monitoring	40
Personal Protection Equipment	59-65
Pharmaceuticals	23-26
Point-of-Care Testing	6-16
Practice Green	58-59
Practice Marketing	93-94
Premier Dental	87
Preventives	86-91
ProRepair	81
QuidelOrtho	10-11
Roche	9
Septodont	66
Sherman Specialty	93
Shofu	71
Solvatum (formerly known as 3M)	43, 72-73
Surgical Products	49-50
Ultrasound	41
Vaccines	19-22
Viome	88
Visby Medical	15
VOCO	Back Cover
vVardis	88
Waste Management	29
Wound Care	43-44
Wound Closure	44-47
X-Ray	92
Zoll	42

WHAT'S INSIDE?

Through our new optimization program, PROSPER, we aim to streamline the distribution channels

See page 5



To learn more about PROSPER or to enroll at **no charge**, please contact your local Field Sales Consultant.

HENRY SCHEIN
PROSPER
YOUR COMMUNITY HEALTH ONE DISTRIBUTION SOLUTION

Get point-of-care STI results in under 30 minutes with the Visby Medical Sexual Health Test!

visby medical™

See page 15




(143-5431)

The power of PCR in your hands

Check out SillHa Oral Wellness System and – Leading the Way in Salivary Analytics

See page 89

(143-5120)



SillHa™
Oral Wellness System

The smart choice is **mySaliva** to manage all patients who are not high risk, as high risk patients may need a lab test.

THE HEALTH OF YOUR MOUTH AND YOUR HEART ARE CONNECTED




Henry Schein Cares
FOUNDATION

**Prevention
is Power™**
TAKE CHARGE OF YOUR HEALTH

KNOW YOUR RISK

**2 in 5 adults
have gum
disease**



IF YOU HAVE GUM DISEASE, YOU ARE

- ▶ **Twice as likely** to have a heart attack
- ▶ **Three times as likely** to have a stroke

And your risk increases over the age of 65


Teeth, Gums & Heart Health Connection

Gum disease and cavities can lead to serious problems for your heart health.



- ▶ Germs from the mouth enter the bloodstream and may add to plaque in the neck and heart arteries—these blockages can cause heart attack and stroke
- ▶ Germs in the mouth can trigger inflammation that swells blood vessels, reduces blood flow, and increases risk of clots

High Blood Pressure & the Health of Your Mouth



High blood pressure not only affects your heart health but can also impact your teeth and gums.

- ▶ Changes in blood flow can affect the health of gums
- ▶ If you have high blood pressure discuss your personal risks with your provider and get screened for other health conditions

TAKE ACTION FOR YOUR DENTAL HEALTH & HEART HEALTH

THINGS TO DO EVERY DAY

- Brush your teeth at least twice a day and clean between your teeth
- Eat a healthy and balanced diet
- Limit sugary food and drink which can harm your teeth
- Avoid smoking or vaping
- Engage in at least 30 minutes of moderate physical activity
- Take your medications as prescribed

VISIT YOUR MEDICAL PROVIDER

- For blood pressure screening at least once a year
- Ask about the connection between the health of your teeth, gums and heart
- Ask if any of your medications have side effects
- Ask about other recommended health screenings and what you can do at home to improve your health

VISIT YOUR DENTIST*

- For cleanings at least twice a year
- Talk about your gums, teeth, and breath. Tell them if you have noticed any changes, bleeding or discomfort.
- Share your medical history and current medications
- Ask if they recommend any specialized treatment plans for you

*you may be referred to a specialist for additional care

FOR MORE INFORMATION



Learn more about Prevention is Power



Health Care Resources

Prevention is Power Landing page and Resource pages are now live: www.henryschein.com/preventionispower

HENRY SCHEIN CARES FOUNDATION PREVENTION IS POWER PARTNERS



IMPLEMENTATION PARTNERS



PROSPER

YOUR COMMUNITY HEALTH ONE DISTRIBUTION SOLUTION



Through our new optimization program, PROSPER, we aim to streamline the distribution channels for your medical and dental supplies, optimize inventory management, and enhance the overall resilience and effectiveness of your Community Health Center operations.



DENTAL & MEDICAL SUPPLIES

Competitive pricing and access to the most enhanced supplier contracts and portfolio in the industry.



BUSINESS SOLUTIONS

Negotiation of optimal discounts on essential CHC-focused business solutions.



EQUIPMENT & TECHNOLOGY

Offers a diverse array of medical and dental equipment and technology.



SUPPLY CHAIN EXCELLENCE

Supplier relationships and strategically located distribution centers facilitate timely access to supplies and services and minimize the risk of disruptions.



TECHNICAL SERVICE & REPAIRS

In-office technical service at 44% off local labor rates and a 15% discount on off-site repairs.

To learn more about PROSPER or to enroll at **no charge, please contact your local Field Sales Consultant.**



Help patients stay ahead of their health goals

WITH POINT-OF-CARE TESTING



Lipid Panel, Glucose, and A1C Results in Minutes

- Improve communication and collaboration between providers and patients
- Inspire patients to take a more active role in managing their health



CardioChek® Plus

Lipid panel and glucose results
in as little as 90 seconds

CardioChek Plus Starter Box

HS Item #: 570-1036

60 tests. Includes analyzer, lipid panel test strips, eGLU test strips, multi-chemistry controls, and capillary tubes

CardioChek Plus analyzer

HS Item #: 570-0411

Lipid Panel+eGLU Smart Bundle

HS Item #: 570-0412

PTS Collect Capillary Tubes

HS Item #: 570-0413

PTS Panels Multi-Chemistry Controls

HS Item #: 261-8361



A1CNow®+

HbA1c results in 5 minutes

The A1CNow+ test system - 20 tests

HS Item #: 570-3428

includes analyzer, test cartridges, and shaker kits

OneStep+ Ultra Flu Test

From a Name You Can Trust.

Intended for the qualitative detection of influenza type A and type B nucleoprotein antigens directly from nasal and nasopharyngeal swab specimens from patients with signs and symptoms of respiratory infection.



- CLIA waived
- Results in 10 minutes
- Sensitivity (A/B): 90.3%/88%
- Specificity (A/B): 96.7%/99.2%
- For use with nasal and nasopharyngeal swab specimens
- Made in the USA

(570-1886) 25 per Box

Contains: 25 test sticks, 25 sterile nasal swabs, 27 extraction buffer vials, 1 influenza A control swab & 1 influenza B control swab.

CPT® Codes: Influenza A, 87804QW⁺/Influenza B, 87804QW-59⁺



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!



DETERMINE™ HIV-1/2 Ag/Ab COMBO

RAPID AND EARLY IDENTIFICATION OF HIV

INNOVATIVE

First and only CLIA-waived* point-of-care test that detects both acute and chronic HIV infections

EARLY DETECTION

Incorporates the p24 antigen for earlier diagnosis compared to antibody-only tests

EFFICIENT

Get results in just 20 minutes and improve same-visit notification

FLEXIBLE

Use capillary and venous whole blood, serum, or plasma samples



Contact your Henry Schein Representative for additional information.

Determine HIV-1/2 Ag/Ab Combo (25 tests) (512-0022) Determine HIV-1/2 Ag/Ab Controls (512-0024)

*Fingerstick Collection

© 2024 Abbott. All rights reserved. All trademarks referenced are trademarks of either the Abbott group of companies or their respective owners. Any photos displayed are for illustrative purposes only. COL-18769-01 01/23

Determine™ HIV-1/2 Ag/Ab Combo

ABBOTT

The first FDA-approved rapid point-of-care test that detects both HIV-1/2 antibodies and free HIV-1 p24 antigen. This 4th-generation test has the ability to identify HIV sooner than conventional rapid tests, which rely solely on the presence of HIV-1/2 antibodies. It enables health care providers to diagnose HIV infection earlier.

- **Innovative:** It's the first rapid point-of-care test that detects both HIV-1/2 antibodies and free HIV-1 p24 antigen
- **Empowering:** Enables health care providers to diagnose HIV infection earlier, allowing individuals to seek medical care sooner
- **Dependable:** Built-in quality controls let you know the test is working
- **Efficient:** Requires minimal training, test in 3 simple steps, and presents clear results in just 20 minutes
- **Flexible:** Multiple sample types; fingerstick and venous whole blood, serum, or plasma samples

#7D2648
(512-0022).....25 per Box



cobas®

Roche

cobas® liat total system

PCR made for the point of care

The cobas® liat total system is compact and uses PCR technology to differentiate pathogens for timely patient treatment to prevent disease progression and the spread of infection.

Henry Schein Item # Assay Menu:


Time to Result:

146-4384	SARS-CoV-2 & Influenza A/B	20 minutes
147-3232	SARS-CoV-2	20 minutes
127-7742	Influenza A/B & RSV	20 minutes
127-7744	Strep A	15 minutes
142-5842	C. difficile**	20 minutes



Learn more

**This product has been cleared by the FDA for use by authorized laboratories certified under OLIA to perform moderate complexity testing.



cobas® SARS-CoV-2 & Influenza A/B Nucleic Acid Test

ROCHE

This test, for use on the cobas® Liat® System, is a multiplex, real-time polymerase chain reaction (PCR) test that detects and differentiates SARS-CoV-2, influenza A, and influenza B in 20 minutes from a single nasal or nasopharyngeal sample in just one test.

- CLIA waived
- 1 sample, 1 test, 20 minutes: know if it's COVID-19 or flu
- Take the workload out of testing with simple, intuitive user handling
- Nasopharyngeal swab (NPS) and anterior nasal swab (ANS) specimens

#0921101190, Assay Kit

(146-4384)..... 20 per Box

#09211128190, Positive/Negative Control

(146-5329)..... 9 per Box

cobas® Liat® PCR System ROCHE

Brings molecular testing with real-time PCR (polymerase chain reaction) to the point of care, enabling confidence in making the right treatment decision when patients need it most.

#07341920190, Liat PCR Analyzer

(127-7745)..... Ea.

#07341890190, Liat Influenza A/B Test Kit

(127-7408)..... 20 per Box

#07402660190, Liat Influenza A/B Control Kit, Positive/Negative

(127-7409)..... 3 per Kit

#08160104190, Liat Influenza A/B & RSV Test Kit

(127-7742)..... 20 per Box

#07402686190, Liat Influenza A/B & RSV Control Kit

(127-7810)..... 6 per Pkg.

#07341911190, Liat Strep A Test Kit

(127-7744)..... 20 per Box

#07402678190, Liat Strep A Control Kit

(127-7813)..... 6 per Pkg.

#07402759190, Liat Cleaning Tool Kit

(127-7741)..... Ea.





The Savanna HSV 1+2/VZV assay leverages real-time PCR technology with a differentiated VZV target to provide superior patient outcomes faster.

For more information, visit quidelortho.com, or contact your Henry Schein representative.

* Data on file.

New QuidelOrtho branding may not be available in all markets, subject to country-specific regulatory approval. Please confirm with your local commercial team.

© 2024 QuidelOrtho Corporation. All rights reserved.

PR-105745-NA-EN-US-v1

Savanna® HSV 1+2/VZV Assay

The power of certainty

Lesions for HSV-1, HSV-2 and VZV are misdiagnosed by visual examination in up to 20 percent of cases.

Make an impact for your patients with Savanna HSV 1+2/VZV assay, a unique rt-PCR assay, by providing results up to 230 percent faster than a typical send out.*







Savanna® Platform
QUIDELORTHO

This compact benchtop instrument utilizes real-time PCR to deliver automated, customized molecular results in minutes. Molecular sensitivity with rapid turnaround can detect multiple pathogens in a single sample. Eliminates the need for dedicated lab space. Facilitates decentralization of testing and allows multiplex, syndromic panel testing at point of care.

- Sample-to-result technology
- Decreased hands-on time with no upfront sample or reagent preparation required
- Superior molecular accuracy and confidence in a multiplex, real-time PCR assay
- Single-use, self-contained cartridge with no sample or reagent prep prior to testing
- Fully automated and easy to use with touchscreen display results in 24 minutes

#20382, Real-Time PCR Instrument	
(146-1056)	Ea.
#20445, HSV1+2/VZV Test Kit	
(146-0986)	12 per Kit
#M409-1, HSV1+2/VZV Assay Control Set	
(146-0987)	Ea.
#M409-4, HSV1+2/VZV Assay Control Set	
(146-0985)	4 per Kit



Sofia® 2 Platform

The Sofia 2 analyzer represents the next generation of rapid diagnostic testing

Proven lateral-flow technology, proprietary advanced fluorescence chemistry and assay development techniques are integrated into two small and portable benchtop analyzers. Our high-performance analyzers deliver point-of-care testing for better patient outcomes.

QuidelOrtho™



Simplify your workflow

Our Sofia benchtop systems and fluorescent immunoassays deliver highly accurate, objective and automated results in three to 15 minutes.



Efficient data storage and retrieval

View results remotely from your computer with Virena® System wireless connectivity. Patient results are automatically stored and results never expire.



Learn more about Sofia by scanning the QR code, visiting quidelortho.com, or reaching out to your Henry Schein representative.



Convenient and easy to use

The Sofia systems are small in form, have intuitive interfaces and only require about one minute of total hands-on time.

New QuidelOrtho branding may not be available in all markets, subject to country-specific regulatory approval. Please confirm with your local commercial team.

© 2024 QuidelOrtho Corporation PR-10572-NA-EN-US-v1

Sofia® Influenza A+B FIA Test QUIDELORTHO

This CLIA-waived test uses advanced immunofluorescence-based lateral-flow technology to detect influenza A and B viral nucleoprotein antigens in nasal and nasopharyngeal swab and nasopharyngeal aspirate/wash specimens taken directly from symptomatic patients.

- Enhanced assay sensitivity
- Results in 15 minutes
- Accurate detection with nasal swab, NP swab and NPA/wash samples
- Nasal swab sensitivity (A/B): 90%/89%; specificity (A/B): 95%/96%
- NP swab sensitivity (A/B): 97%/90%; specificity (A/B): 95%/97%
- NPA/wash sensitivity (A/B): 99%/88%; specificity (A/B): 96%/96%
- Room temperature storage
- 21-month shelf life from date of manufacture
- For use with Sofia Analyzer (sold separately)
- CLIA-waived for nasal and nasopharyngeal swabs

#20299, Sofia 2 Analyzer
(128-5292)..... Ea.

#20218, Sofia Influenza A+B Test
(122-4092)..... 25 per Box

CPT® Codes: Influenza A: 87804QW†; Influenza B: 87804-59†



Sofia® RSV FIA Test QUIDELORTHO

A 15-minute test that has a walk-away automation feature that allows user to multitask. Offers accurate detection with nasopharyngeal swabs, or aspirate/wash samples. CLIA-waived for children less than 7 years of age, and moderately complex for pediatric patients from age 7 to less than 19 years of age.

- Single specimen may be used to run both Sofia influenza and RSV FIA tests
- Room temperature storage
- 24-month shelf life from date of manufacture
- For use with Sofia Analyzer (sold separately)

#20260, Sofia RSV Kit
(122-4091)..... 25 per Box
CPT® Code: 87807†



Sofia® Rapid Strep A Test QUIDELORTHO

Easy-to-use test with enhanced assay sensitivity and results in 5 minutes. With results available in minutes, you can rapidly diagnose and treat, all in the same visit.

- Accurate detection with throat swab samples
- Sensitivity: 90.6 %
- Specificity: 96.1%
- Room temperature storage
- 24-month shelf life from date of manufacture
- For use with Sofia Analyzer (sold separately)

#20253–Sofia Strep A
(122-4089)..... 25 per Box

#20274, Sofia Strep A+
(122-7865)..... 25 per Box

CPT® Code: 87880†





Abbott

LAB-ACCURATE RESULTS. ON SITE. IN MINUTES.

***i*-STAT® SYSTEM + PICCOLO® XPRESS CHEMISTRY ANALYZER**

FAST

Results in minutes enable clinical decisions to be made without the need for additional lab appointments.

ACCURATE

Rapid, lab-quality point-of-care testing results across a comprehensive menu of assays including CLIA-waived tests.

VERSATILE

Supports testing in a variety of care settings, including those without a laboratory on-site.

***i*-STAT AND PICCOLO'S COMPREHENSIVE TEST MENUS FEATURE A RANGE OF ASSAYS THAT SUPPORT TESTING IN A VARIETY OF CARE SETTINGS**



DATA EASILY INTEGRATED INTO EMR/LIS



INTUITIVE FUNCTION AND DESIGN



INCREASES OPERATIONAL EFFICIENCY



***i*-STAT®**

BROAD TEST MENU

The *i*-STAT system range of diagnostics tests include lactate, blood gases, chemistries and electrolytes, cardiac markers, coagulation, hematology, and endocrinology, including CLIA-waived Creatinine



piccolo^{xpress}
chemistry analyzer

CLIA-WAIVED TESTS

CMP, BMP, Lipid Panel, Lipid Panel Plus, General Chemistry 13, General Chemistry 6, Kidney Check, Liver Panel Plus, Renal Function Panel, MetLyte 8, and Electrolytes



TO LEARN MORE CONTACT YOUR ABBOTT POINT-OF-CARE REPRESENTATIVE, YOUR HENRY SCHEIN SALES REPRESENTATIVE, OR VISIT WWW.GLOBALPOINTOFCARE.ABBOTT

©Abbott Point of Care Inc. | *i*-STAT is a trademark of Abbott. | For *in vitro* diagnostic use only. | This material is intended for a U.S. audience only. | 400 College Road East, Princeton, NJ 08540 USA | (609) 454-9000 | (609) 419-9370 (fax) | www.globalpointofcare.abbott | Piccolo Xpress is a registered trademark of Abaxis, Inc. and distributed by Abbott Point of Care. | Abaxis 888-9685 Rev A Henry Schein Ad— Abbott POC *i*-STAT/Piccolo | 5366.Rev1.APOC.EN-US 05/24

Piccolo Xpress® Chemistry Analyzer— ABBOTT

The first broad-menu clinical chemistry analyzer designed for near-patient testing. *Basic Panel* contains BUN, Ca, Cl-, CRE, GLU, K+, Na+, and tCO₂. *Comprehensive Panel* contains ALB, ALP, ALT, AST, BUN, Ca, Cl-, CRE, GLU, K+, Na+, TBIL, tCO₂, and TP. *Lipid Panel* contains CHOL, CHOL/HDL (calculated), HDL, LDL (calculated), TRIG, and VLDL (calculated).

- **Speed:** Add the sample, insert the reagent disc, and in approximately 12 minutes receive a report suitable for the patient file or transferred to a computer, printer, or LIS/EMR.
- **Ease:** A fully automated system requiring no special operating skills.
- **Exclusivity:** Offers the only CLIA-waived comprehensive metabolic panel (CMP) blood test on the market.
- **Accuracy:** Self-calibrates with each run. The onboard continuous intelligent Quality Control (iQCT) monitors the analyzer, reagent reactions, and sample to ensure chemistry and instrument integrity.

#07P05-01, Piccolo Xpress® Chemistry Analyzer
(106-6368)Ea.

Contains: 01 cc minipipette, minipipette tips, power cord, power supply, printer results paper & owners manual.

#07P02-08, Comprehensive Metabolic Panel:
ALB, ALP, ALT, AST, BUN, Ca, Cl-, CRE, GLU, K+,
Na+, TBIL, tCO₂, TP
(174-6978)10 per Box



#07P02-04, Basic Metabolic Panel: BUN, Ca,
Cl-, CRE, GLU, K+, Na+, tCO₂, eGFR*
(137-3987)10 per Box

#07P02-05, Lipid Panel: CHOL, CHOL/HDL*,
nHDLc*, HDL, LDL*, TRIG, VLDL* (* Calculated)
(146-0189)10 per Box

#07P02-12, Lipid Panel Plus: ALT, AST, CHOL,
CHOL/HDL*, nHDLc*, GLU, HDL, LDL*, TRIG,
VLDL* (* Calculated)
(104-6269)10 per Box

#07P02-10, Liver Panel Plus: ALB, ALP, ALT, AMY,
AST, GGT, TBIL, TP
(104-3998)10 per Box

#07P02-01, General Chemistry 6: ALT, AST, BUN,
CRE, GGT, GLU, eGFR*
(108-5450)10 per Box

#07P02-11, General Chemistry 13: ALB, ALP,
ALT, AMY, AST, BUN, Ca, CRE, GGT, GLU, TBIL,
TP, UA, eGFR*
(982-0407)10 per Box

#07P02-02, Electrolyte Panel: Cl-, K+, Na+, tCO₂
(106-1669)10 per Box

#07P02-09, Kidney Check: BUN, CRE, eGFR*
(109-7255)10 per Box

#07P02-16, MetLac 12 Panel: ALB, BUN, Ca, Cl-,
CRE, GLU, LAC, K+, Mg, Na+, PHOS, tCO₂, eGFR*
(631-0050)10 per Box

#07P02-15, BioChemistry Panel Plus: ALB,
ALP, ALT, AMY, AST, BUN, Ca, CRE, CRP, GGT,
GLU, TP, UA
(116-9038)10 per Box

#07P02-14, MetLyte Plus CRP: BUN, CK, Cl-, CRE,
GLU, K+, Na+, tCO₂, CRP
(115-5701)10 per Box

#07P02-13, Basic Metabolic Panel Plus: BUN, Ca,
Cl-, CRE, GLU, K+, Na+, tCO₂, Mg, Lactate Dehy-
drogenase
(104-5577)10 per Box

#07P02-07, Renal Function Panel: ALB, BUN, Ca,
Cl-, CRE, GLU, K+, Na+, PHOS, tCO₂, eGFR*
(166-0886)10 per Box

#07P02-06, Hepatic Function Panel: ALB, ALP,
ALT, AST, DBIL, TBIL, TP
(195-7347)10 per Box

#07P02-03, MetLyte 8 Panel: BUN, CK, Cl-, CRE,
GLU, K+, Na+, tCO₂, eGFR*
(114-1284)10 per Box

i-STAT® Handheld System ABBOTT

Real-time, lab-quality results in minutes—not hours. The i-STAT® Handheld is an advanced, easy-to-use blood analyzer that provides a wide range of tests using just 2 or 3 drops of blood

- **Speed:** With components of a full-scale lab in each cartridge, results occur in minutes
- **Accuracy:** A unique combination of biosensors and calibrants deliver high-quality results
- **Reliability:** Single-use cartridges reduce many of the problems inherent in multi-use systems

#04J60-20, iSTAT Handheld System
(111-3324)Ea.

Contains: analyzer, printer kit, downloader/recharger, simulator, and rechargeable battery pack.

CLIA Waived Chemistry Tests—Results in 2 Minutes

#03P84-25, Creatine, Crea
(111-3394)25 per Box

Moderately Complex Chemistry Tests—Results in 2 Minutes

#09P31-26, CHEM 8+: Na, K, Cl, iCa, TCO₂, Glu,
BUN/Urea, Crea, Agap*, Hct, Hgb* (*Calculated)
(138-1904)25 per Box

#03P83-26, Glucose: Glu
(146-5249)25 per Box

Blood Gas Tests—Results in 2 Minutes

#03P85-51, CG4+: pH, PCO₂, PO₂, TCO₂*, HCO₃*,
BE_{ecf}*, sO₂*, Lactate (*Calculated)
(146-0531)25 per Box

#03P78-26, G3+: pH, PCO₂, PO₂, TCO₂*, HCO₃*,
BE_{ecf}*, sO₂*, (*Calculated)
(139-7158)25 per Box

Blood Gas, Electrolytes and Hematology Tests—Results in 2 Minutes

#03P77-25, EG6+: Na, K, pH, PCO₂, PO₂, TCO₂*,
HCO₃*, BE_{ecf}*, sO₂*, Hct, Hgb* (*Calculated)
(111-3387)25 per Box

#03P76-25, EG7+: Na, K, iCa, pH, PCO₂, PO₂,
TCO₂*, HCO₃*, BE_{ecf}*, sO₂*, Hct, Hgb* (*Calculated)
(111-3386)25 per Box

#03P88-25, CG8+: Na, K, iCa, Glu, pH, PCO₂, PO₂,
TCO₂*, HCO₃*, BE_{ecf}*, sO₂*, Hct, Hgb* (*Calculated)
(115-6662)25 per Box

Cardiac Marker Tests—Results in 10 Minutes (CK-MB Results in 5 Minutes)

#03P90-25, cTnI: Troponin I
(136-9907)25 per Box

#03P92-25, CK-MB
(111-3408)25 per Box

#03P93-25, BNP
(111-3158)25 per Box



Coagulation Tests—Results in 17 Minutes (PT/INR Results in 5 Minutes)

#03P86-25, Celite ACT
(111-3409)25 per Box

#03P87-25, Kaolin ACT
(111-3412)25 per Box

#03P89-50, PT/INR: Prothrombin Time
(146-0532)24 per Box

Pregnancy Tests—Results in 10 Minutes

#05P58-25, β-hCG
(982-0021)25 per Box

Make every blood drop count

Point-of-care testing can be more convenient for providers and families while being less invasive for pediatric patients. Tests are performed using a finger stick sample, which can help minimize discomfort for young patients.



Hb 801 System

leadcare II



Learn more about point-of-care hemoglobin and lead testing for children at hemocue.us



Glucose 201 Analyzer–Hemocue HEMOCUE

The HemoCue Glucose Systems are based on a glucose-dehydrogenase method and consist of a small, dedicated analyzer and a unique, disposable microcuvette. They will produce lab-accurate results for diagnosing or monitoring diabetes, as well as monitoring neonatal blood glucose levels.

- More accurate and 68% higher reimbursement than glucose meters
- FDA-compliant unit when screening and diagnosing
- Analyzer: 0–444 mg/dL (0–24.6 mmol/L)
- Only 5 µL of capillary, venous, or arterial blood is needed
- No inherent limitations related to hematocrit or known drug interferences
- No recalibration needed between cuvette lots
- 3-year warranty

#120706, Analyzer
(104-2290) ea

#110706, Microcuvettes
(439-9158) 100 per Box



Hb 801 Analyzer HEMOCUE

The HemoCue® Hb 801 system seamlessly integrates point-of-care testing into your workflow.

- Wireless connectivity through middleware
- Results in under 1 second
- Reagent-free cuvettes, improving open-vial stability
- Large display with intuitive on-screen symbols
- Internal quality control self-test
- USB adapter (included) or 3 AA batteries (not included) power supply

Hb 801 Analyzer
(135-3349)..... Ea.

Hb 801 Hemoglobin Microcuvettes
(135-3352)..... 200 per Box
Contains: 25 per bottle, 4 bottles per box.



The power of PCR in your hands

Get point-of-care STI results in under 30 minutes with the Visby Medical Sexual Health Test

For more information, contact your local Henry Schein Medical Representative or visit henryschein.com/medical

visby medical™



Sexual Health Test VISBY MEDICAL

A CLIA-waived, single-use (disposable), instrument-free PCR *in vitro* diagnostic test intended for use in point-of-care or clinical laboratory settings for the rapid detection and differentiation of DNA from *Chlamydia trachomatis*, *Neisseria gonorrhoeae*, and *Trichomonas vaginalis* in self-collected female vaginal swab.

- CLIA waived
- PCR test for chlamydia, gonorrhea, and trichomoniasis
- Result in under 30 minutes
- Vaginal self-collected swab

#PS-400372, Sexual Health Tests
(143-5431)..... 20 per Case

#PS-400515, Sexual Health Starter Kit
(145-7222)..... Ea.

#PS-000715, Vaginal Specimen Collection Kit
(141-4570)..... 50 per Kit

* For details, see instructions for use.



For more information visit [womens-health-solutions.bd.com](https://www.womens-health-solutions.bd.com)

Shape the future of women's health

with a comprehensive diagnostic portfolio from BD



Elevate your HPV testing with **BD Onclarity™ HPV Assay**, the only FDA-approved HPV test that offers extended genotyping.

Get complete vaginitis and STI testing with the **BD Vaginal Panel Assay** and **BD CTGCTV2 Assay** on the **BD MAX™ System**



Bring automated molecular point-of-care vaginitis testing into your lab with **BD Affirm™ VPIII Microbial Identification System**

BD, the BD Logo, BD Affirm, BD MAX and Onclarity are trademarks of Becton, Dickinson and Company or its affiliates. © 2024 BD. All rights reserved. (BD-115852 3220-US-0224 February 2024)



Discover the difference **BD Vacutainer® UltraTouch™** push button blood collection sets can make.

The **BD Vacutainer® UltraTouch™** blood collection set improves the overall phlebotomy experience for both patients and clinicians:

- Reduction in perceived overall patient pain¹
- Faster tube fill times²
- Ability to use a 25-gauge needle without compromising sample quality³
- Up to 88% reduction in needlestick injuries⁴

Visit bd.com/Vacutainer/UltraTouch to learn more.

* Tube fill time evaluation is based on a bench study. Bench study results may not necessarily be indicative of clinical performance.
1. Mouser A, Uettwiller-Geiger D, Plokhoy E, Berube J, Ahuja AJ, Stankovic AK. Evaluation of pain and specimen quality by use of a novel 25-gauge blood collection set with ultra-thin wall cannula and 5-bevel tip design. *J Appl Lab Med*. 2017;2(2):201-210. 2. Evaluation of Tube Fill Time of the BD Vacutainer® UltraTouch™ Push Button Blood Collection Set with PentaPoint™ Comfort Bevel and RightGauge™ Cannula as Compared with the Current BD Vacutainer® Push Button Blood Collection Set with Thin Wall 3-Bevel Cannula. White paper VS9249-2. Franklin Lakes, NJ: BD; 2016. 3. Veronika D, Merrill, Assessing Phlebotomy Device Preference and Specimen Quality in an Oncology Outpatient Clinic. *JALIM*. 2021;1-9. 4. Hotaling M. A retractable winged steel (butterfly) needle performance improvement project. *Jt Comm J Qual Patient Saf*. 2009;35(2):100-105.

bd.com

BD, the BD logo, UltraTouch and Vacutainer are trademarks of Becton, Dickinson and Company or its affiliates. © 2024 BD. All rights reserved. (BD-120156 3354-US-0424 April 2024)

Ordering information

Henry Schein ref #	Needle gauge	Needle length	Wing color	Tubing length	Packaging (box/case)
(987-0051)	BD Vacutainer® one-use holder				
(137-9099) [†]	21	0.75	Green	12"	50/200
(137-9098) [†]	23	0.75	Light Blue	12"	50/200
(137-9100) [†]	25	0.75	Royal Blue	12"	50/200
(124-9458)	21	0.75	Green	12"	50/200
(124-9457)	23	0.75	Light Blue	12"	50/200
(124-9456)	25	0.75	Royal Blue	12"	50/200

[†] Pre-attached holder included





Ready, set, vaccinate

with BD hypodermic solutions

It's simple: the more immunity you can give your patients, the better. That means having enough of the right vaccines on hand *and* the high-quality needles and syringes you need to deliver them. BD hypodermic solutions are designed to deliver consistent quality and reliable performance.

BD SafetyGlide™ Needles are designed to help keep clinicians safe from accidental needlestick injuries and feature Activation-Assist™ Technology—making safe injections easier and more efficient.



To purchase your
BD vaccine devices,
scan this QR code.



BD, Franklin Lakes, NJ, 07417, U.S.

bd.com

BD, the BD Logo, Activation-Assist and SafetyGlide are trademarks of Becton, Dickinson and Company or its affiliates.
© 2024 BD. All rights reserved. BD-120838 (04/24)



Safety & Conventional Hypodermic Needles

From a Name You Can Trust.



Features & Benefits:

- 3 single-handed safety activation methods: hard-surface, thumb, or finger
- Double-locking mechanism with an audible “click”
- Needle hub-compatible with all standard Luer lock and Luer slip-tip syringes
- ISO color coded—easy and quick identification of the needle gauge
- Bevel-up needle—allows for low-angle injections

ITEM#	DESCRIPTION	SIZE
SAFETY NEEDLE		
(570-2693)	18 ga x 1"	100 per Box
(570-2694)	18 ga x 1½"	100 per Box
(570-0338)	20 ga x 1"	100 per Box
(570-2697)	20 ga x 1½"	100 per Box
(570-2696)	21 ga x 1"	100 per Box
(570-2698)	21 ga x 1½"	100 per Box
(570-2699)	22 ga x 1"	100 per Box
(570-2700)	22 ga x 1½"	100 per Box
(570-0341)	23 ga x ⅝"	100 per Box
(570-2701)	23 ga x 1"	100 per Box
(570-0340)	23 ga x 1½"	100 per Box
(570-2705)	25 ga x ⅝"	100 per Box
(570-2703)	25 ga x 1"	100 per Box
(570-2707)	27 ga x 1¼"	100 per Box
(570-2709)	30 ga x ½"	100 per Box

ITEM#	DESCRIPTION	SIZE
NEEDLES		
(112-7106)	16 ga x 1"	100 per Box
(112-7116)	16 ga x 1½"	100 per Box
(900-4469)	18 ga x 1"	100 per Box
(900-4471)	18 ga x 1½"	100 per Box
(112-7109)	20 ga x 1"	100 per Box
(112-7110)	20 ga x 1½"	100 per Box
(112-7100)	21 ga x 1"	100 per Box
(112-7111)	21 ga x 1½"	100 per Box
(112-7107)	22 ga x 1"	100 per Box
(112-7101)	22 ga x 1½"	100 per Box
(900-4473)	23 ga x 1"	100 per Box
(900-4468)	25 ga x ⅝"	100 per Box
(112-7112)	25 ga x 1"	100 per Box
(900-4470)	25 ga x 1"	100 per Box
(900-4472)	25 ga x 1½"	100 per Box
(112-7108)	27 ga x ½"	100 per Box
(112-7103)	27 ga x 1¼"	100 per Box
(112-6148)	30 ga x ½"	100 per Box
(112-6149)	30 ga x 1"	100 per Box



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

www.henryschein.com/medical


ABRYVO™
PFIZER

Respiratory Syncytial Virus Vaccine
 0.5-mL Prefilled Syringe
(145-4487)..... Ea.
(145-4488)..... 5 per Box
 CPT® Code: 90678†


Adacel®
SANOFI

Tetanus Toxoid, Reduced Diphtheria Toxoid,
 and Acellular Pertussis Vaccine Adsorbed
 0.5-mL SDV
(546-7789).....10 per Pkg.
 0.5-mL Prefilled Syringe Without Needle
 Not made with natural rubber latex.
(132-5841)..... 5 per Pkg.
 CPT® Code: 90715†


AREXVY
GSK

Respiratory Syncytial Virus Vaccine, Adjuvanted
 0.5-mL SDV
(145-4809).....10 per Box
 CPT® Code: 90679†


BEXSERO®
GSK

Meningococcal Group B Vaccine
 0.5-mL Prefilled Syringe
(127-2151).....10 per Pkg.
 CPT® Code: 90620†


BOOSTRIX®
GSK

Tetanus Toxoid, Reduced Diphtheria, and Acellular
 Pertussis Vaccine, Adsorbed
 0.5-mL SDV
(254-6475).....10 per Pkg.
 0.5-mL Tip-Lok® Prefilled Syringe
(254-0031).....10 per Pkg.
 CPT® Code: 90715†


**Pneumococcal 21-valent
 Conjugate Vaccine**
MERCK


Pneumococcal Conjugate Vaccine
 0.5-mL Prefilled Syringe
(147-2899) 10/pkg
 CPT® Code: 90684†
 CAPVAXINE® is a registered trademark of Merck & Co., Inc.


ENGRIX-B®–Pediatric/Adolescent
GSK

Hepatitis B Recombinant Vaccine
 10 µg, 0.5-mL Syringe
(254-0026).....10 per Pkg.
 CPT® Code: 90744†


ENGRIX-B®–Adult
GSK

Hepatitis B (Recombinant) Vaccine,
 Preservative-Free
 20 µg/mL, 1-mL SDV
(254-8254).....10 per Pkg.
 20 µg/mL, 1-mL Tip Lok™ Syringe without Needle
(254-0029).....10 per Pkg.
 CPT® Code: 90746†


HAVRIX®–Pediatric/Adolescent
GSK

Hepatitis A (Inactivated) Vaccine,
 Thimerosal-Free, 720 ELU/0.5 mL
(254-0023).....10 per Pkg.
 CPT® Code: 90633†



**HAVRIX® –Adult
GSK**

Hepatitis A (Inactivated) Vaccine, Thimerosal-Free
1440 ELU/mL, 1-mL Tip-Lok® Prefilled Syringe
(124-0016).....10 per Pkg.

CPT® Code: 90632†




**HEPLISAV-B™
DYNAVAX**

Hepatitis B Vaccine (Recombinant), Adjuvanted
Solution for Intramuscular Injection

0.5-mL Prefilled Syringe
20 mcg/0.5 mL
(132-9264).....5 per Pkg.

CPT Code: 90739†




**INFANRIX®
GSK**

Diphtheria and Tetanus Toxoids and Acellular
Pertussis (DTaP) Vaccine (Adsorbed),
Thimerosal-Free

0.5-mL Prefilled Syringe
(124-0010).....10 per Pkg.

CPT® Code: 90700†



**IPOL®
SANOFI**

Poliovirus (Salk) (Inactivated) Vaccine

5-mL, 10-Dose Vial
(546-7189)..... Ea.
CPT® Code: 90713†




**IXCHIQ®
VALNEVA**

Chikungunya Vaccine

0.5-mL SDV
(146-1426)..... Ea.

CPT® Code: 90589†




**KINRIX™
GSK**

Diphtheria and Tetanus Toxoids and Acellular
Pertussis Adsorbed and Inactivated Poliovirus
Vaccine

0.5-mL Prefilled TipLok Syringe
(254-0025)..... 10 per Pkg.

CPT Code: 90696†




**MenQuadfi™
SANOFI**

Meningococcal (Groups A, C, Y, W) Conjugate
Vaccine Solution for Intramuscular Injection

0.5-mL SDV
(140-3989)..... 5 per Pkg.

CPT® Code: 90619†




**M-M-R® II
MERCK**

Measles, Mumps, and Rubella Vaccine

0.5-mL SDV with Diluent
(558-0110).....10 per Pkg.

CPT® Code: 90707†

M-M-R® II is a registered trademark of Merck & Co., Inc.




**MENVEO ONE-VIAL
GSK**

Meningococcal (Groups A, C, Y, and W-135)
Oligosaccharide Diphtheria CRM197 Conjugate
Vaccine

0.5-mL SDV
(145-4145).....10 per Pkg.

CPT® Code: 90734†



**PEDIARIX®
GSK**

Diphtheria and Tetanus Toxoids and Acellular
Pertussis (DTaP) (Adsorbed), Hepatitis B
(Recombinant), and Inactivated Poliovirus
Vaccine (Combined)

0.5-mL Tip-Lok® Prefilled Syringe
(254-0028).....10 per Pkg.

CPT® Code: 90723†




PedvaxHIB® MERCK

Haemophilus b Conjugate (Meningococcal Protein Conjugate) Liquid Vaccine

0.5-mL SDV
(558-3778).....10 per Pkg.

CPT® Code: 90647†

PedvaxHIB® is a registered trademark of Merck & Co., Inc.




PENTACEL® SANOFI

Diphtheria and Tetanus Toxoids and Acellular Pertussis Adsorbed, Inactivated Poliovirus and Haemophilus b Conjugate (Tetanus Toxoid Conjugate) Vaccine

0.5-mL SDV
(112-3585).....5 per Pkg.

CPT® Code: 90698†



PNEUMOVAX® 23 MERCK

Pneumococcal Vaccines (Polyvalent)

0.5-mL Prefilled Syringe
(122-2087).....10 per Pkg.

CPT® Code: 90732†

Pneumovax® 23 is a registered trademark of Merck & Co., Inc.



PREVNAR 20™ PFIZER

Pneumococcal 20-Valent Conjugate Vaccine

0.5-mL Prefilled Syringe
(140-8153).....Ea.

0.5-mL Prefilled Syringe
(140-8152).....10 per Pkg.

CPT® Code: 90677†



PRIORIX GSK

Measles, Mumps, and Rubella Vaccine, Live

0.5-mL SDV
(142-3773).....10 per Box

CPT® Code: 90707†



ProQuad® MERCK

Measles, Mumps, Rubella, and Varicella Virus Vaccine, Live, 0.5-mL SDV

0.5-mL SDV
(558-0053).....10 per Pkg.

CPT® Code: 90701†

ProQuad® is a registered trademark of Merck & Co., Inc.



RECOMBIVAX HB®—Adult MERCK

Hepatitis B (Recombinant) Vaccines, Thimerosal-Free

10 µg/mL, 1-mL SDV
(558-0737).....10 per Pkg.

1-mL Prefilled Luer-Lock Syringe without Safety Needle
(122-2086).....10 per Pkg.

CPT® Code: 90746†

RECOMBIVAX HB® is a trademark of Merck & Co., Inc.



RECOMBIVAX HB® Pediatric/Adolescent MERCK

Hepatitis B (Recombinant) Vaccine, Thimerosal-Free

5 µg/mL, 0.5-mL SDV
(558-6254).....10 per Pkg.

0.5-mL Prefilled Luer-Lock Syringe without Safety Needle
(122-2084).....10 per Pkg.

CPT® Code: 90744†

RECOMBIVAX HB® is a registered trademark of Merck & Co., Inc.



ROTARIX® GSK

Rotavirus Vaccine, Live, Oral 1.5-mL Pediatric Prefilled Syringe

(143-4023).....10 per Pkg.

CPT® Code: 90681†

RotaTaq[®] MERCK

Rotavirus Vaccine (Live), Oral, Pentavalent
2-mL Single-Dose Tube
(558-8763).....10 per Pkg.
(558-0009).....25 per Pkg.
CPT[®] Code: 90680[†]
RotaTaq[®] is a registered trademark of Merck & Co., Inc.



SHINGRIX[®] GSK

Zoster Vaccine Recombinant, Adjuvanted
0.5-mL SDV with Diluent
(129-6729).....10 per Pkg.
CPT[®] Code: 90750[†]



TENIVAC[®] SANOFI

Tetanus and Diphtheria Toxoids Adsorbed
(546-0016)10 per Pkg.
(546-0017)10 per Pkg.
CPT[®] Code: 90714[†]




Trumenba[®] PFIZER

Meningococcal Group B Vaccine
0.5-mL Prefilled Syringe
(248-0649).....5 per Pkg.
(248-0648).....10 per Pkg.
CPT[®] Code: 90621[†]



TWINRIX[®] Adult GSK

Hepatitis A (Inactivated) and Hepatitis B
(Recombinant) Vaccine
1-mL Tip-Lok[®] Prefilled Syringe without Needle
(254-0030).....10 per Pkg.
CPT[®] Code: 90636[†]



VAQTA[®] Adult MERCK

Hepatitis A (Inactivated) Vaccines, Thimerosal-Free
50 U/1 mL Prefilled Luer-lock Syringe Without
Safety Needle
(122-2081).....10 per Pkg.
CPT[®] Code: 90632[†]
VAQTA is a Registered Trademark of Merck & Co., Inc.



VAQTA[®]–Pediatric/Adolescent MERCK

Hepatitis A (Inactivated) Vaccine, Thimerosal-Free,
25 U, 0.5-mL SDV
0.5-mL Prefilled Luer-lock Syringe Without Safety
Needle
(122-2080).....10 per Pkg.
CPT[®] Code: 90633[†]
VAQTA[®] is a registered trademark of Merck & Co., Inc.


VARIVAX[®] MERCK

Varicella (Chicken Pox) Virus (Live) Vaccine,
Thimerosal-Free, Lyophilized
0.5-mL SDV
(558-1592).....10 per Box
CPT[®] Code: 90716[†]
Varivax[®] is a registered trademark of Merck & Co., Inc.




Vaxelis[™] SANOFI

Diphtheria and Tetanus Toxoids and Acellular
Pertussis, Inactivated Poliovirus, Haemophilus b
Conjugate and Hepatitis B Vaccine
0.5-mL SDV
(140-3160).....10 per Pkg.
0.5-mL PFS
(140-3161).....10 per Pkg.
CPT[®] Code: 90697[†]



VAXNEUVANCE[™] MERCK

Pneumococcal 15-Valent Conjugate Vaccine
(141-1002).....10 per Pkg.
CPT[®] Code: 90671[†]
Vaxneuvance[™] is a trademark of Merck & Co., Inc.





Adrenalin® Injection
PAR STERILE PRODUCTS

1 mg/mL, 1-mL SDV
(248-0883)..... Ea.
1 mg/mL, 1-mL SDV
(890-0257)..... 25 per Pkg.
1 mg/mL, 30-mL MDV
(129-4847)..... Ea.



Compare to ProAir®
Albuterol Sulfate Inhaler HFA
LUPIN

90 mcg
(140-9653)..... Ea.



Aplisol® Injectables
PAR STERILE PRODUCTS

Tuberculin PPD Stabilized Solution, 5 TU/0.1 mL
1-mL Vial
(110-5199)..... Ea.
5-mL Vial
(110-5099)..... Ea.



Auvi-Q Epinephrine
Auto-Injectors
KALEO

0.1 mg
(142-4090)..... 2 per Pkg.
0.15 mg
(142-4091)..... 2 per Pkg.
0.3 mg
(142-4092)..... 2 per Pkg.



Bicillin® L-A Injectable*
Prefilled Syringes
PFIZER

Penicillin G benzathine IM Suspension
Pediatric, 21 ga x 1"
600,000 U/mL, 1 mL
(321-6704)..... 10 per Pkg.
Adult, 21 ga x 1½"
1,200,000 U/mL, 2 mL
(321-7445)..... 10 per Pkg.
Adult, 18 ga x 1½"
2,400,000 U/mL, 4 mL
(321-1663)..... 10 per Pkg.

*Tubex injector is not required.



Compare to Rocephin®
Ceftriaxone Sodium Injectable
APOTEX

250-mg SDV
(133-6543)..... Ea.
250-mg SDV
(131-4702)..... 10 per Pkg.
500-mg SDV
(133-6542)..... Ea.
500-mg SDV
(135-9645)..... 10 per Box
1-gram Vial
(133-6436)..... 10 per Box
2-gram Vial
(135-7800)..... 10 per Box



Celestone® Soluspan® Injectable
Suspension
ORGANON

6 mg/mL, 5-mL MDV
(140-7856)..... Ea.



Compare to Rubramin PC®
Cyanocobalamin (Vitamin B-12)
Injectable
AMERICAN REGENT

1000 mcg/mL
1-mL MDV
(244-1902)..... 25 per Pkg.



Depo-Provera® Injectable
PFIZER
Medroxyprogesterone Acetate Suspension
150 mg/mL, 1-mL Prefilled Syringe
(908-0028)..... Ea.



Durolane® Injectable
BIOVENTUS
 60 mg/3 mL Prefilled Syringe
 (137-4043).....Ea.
 Series of injections: 1
 J Code: J7318

Compare to Adrenalin®

Epinephrine Injectable
Ampules
BPI LABS
 1:1000, 1 mL
 (248-0899).....Ea.
 1:1000, 1 mL
 (128-9991).....10 per Box



Ethyl Chloride Spray
GEBAUER
 Topical Anesthetic Skin Refrigerant
 Fine Spray, 3.9 oz.
 (138-9019).....Ea.
 Medium Spray, 3.9 oz.
 (138-5777).....Ea.
 Mist Spray Can, 3.9 oz.
 (138-6406).....Ea.



Euflexxa® Prefilled Syringe
FERRING PHARMACEUTICALS
 2 mL
 (137-4434).....3 per Pkg.
 Series of injections: 3
 J Code: J7323



Evenity® Injection
AMGEN
 105 mg/1.17 mL
 (214-0051).....2 per Box



Fluorescein Sodium and Benoxinate Hydrochloride Ophthalmic Solution, 0.25%/0.4%
ALTAIRE PHARMACEUTICALS
 5-mL Bottle with Sterile Dropper
 (131-3074).....Ea.



Gelsyn-3® Injection
BIOVENTUS
 Prefilled Syringe
 (136-8849).....Ea.
 Series of injections: 3
 J Code: J7328



HyperHEP B™ Injectable
GRIFOLS
 1-mL SDV
 (140-0805).....Ea.
 5-mL SDV
 (271-7848).....Ea.



Kenalog®-10 Injectable
BRISTOL MYERS SQUIBB
 Triamcinolone Acetonide
 10 mg/mL, 5-mL MDV
 (196-0545).....Ea.



Kenalog® Injectables
BRISTOL MYERS SQUIBB
 Triamcinolone Acetonide
 40 mg/mL
 1-mL SDV
 (196-8300).....Ea.
 5-mL MDV
 (196-1737).....Ea.
 10-mL MDV
 (196-9429).....Ea.



Compare to Toradol
Ketorolac Injection
FOSUN PHARMA
 15 mg/mL, 1-mL SDV
 (146-1679).....25 per Box
 30 mg/mL, 1-mL SDV
 (258-0457).....Ea.
 30 mg/mL, 1-mL SDV
 (143-6324).....25 per Box
 30 mg/mL, 2-mL SDV
 (258-0456).....Ea.
 30 mg/mL, 2-mL SDV
 (143-6293).....25 per Box



Kyleena® IUD System
BAYER

19.5 mg
(126-7908)..... Ea.



Mirena® IUD System
BAYER

52 mg
(125-0996)..... Ea.



Compare to Xylocaine®
Lidocaine HCl Injection
PFIZER

50-mL Vial MDV
1%
(248-0874)..... Ea.
2%
(248-0973)..... Ea.
2% 50 mL Vial
(104-9843).....25 Vials per Box

Compare to Depo-Provera®
Medroxyprogesterone
MYLAN INSTITUTIONAL

150 mg/mL, 1-mL SDV
(146-6392)..... Ea.
150 mg/mL, 1-mL SDV
(146-6382).....25 per Box



Compare to Depo-Provera®
Medroxyprogesterone Injectable
EUGIA

150 mg/mL Prefilled Syringe
(146-4156)..... Ea.



Monovisc®
DEPUY SYNTHES

4-mL Prefilled Syringe, 22 mg/mL
(729-0037)..... Ea.
Series of injections: 1
J Code: J7327




Compare to Depo-Medrol®
Methylprednisolone Acetate Injection
AMNEAL

40 mg/mL, 1-mL SDV
(140-9021)..... Ea.
40 mg/mL, 1-mL SDV
(140-9026).....25 per box
80 mg/mL, 1-mL SDV
(140-9031)..... Ea.
80 mg/mL, 1-mL SDV
(140-9048).....25 per Box



Ondansetron HCl Injection
APOTEX

2 mg/mL, 2-mL SDV
(248-0902)..... Ea.
2 mg/mL, 2-mL SDV
(125-9100).....25 per Box



Orthovisc®
DEPUY SYNTHES

2-mL Single-Use Syringe
(729-6418)..... Ea.
Series of injections: 3 to 4
J Code: J7324



Pain Ease®
GEBAUER

Mist
(791-0634).....3.9 oz. Can
1-oz Medium-Stream Spray
(132-9467)..... Ea.
Medium Stream
(761-0029).....3.9 oz. Can




Prolia Injectable Prefilled Syringe
AMGEN

60 mg/mL
(214-0006)..... Ea.



Compare to Lexiscan®
Regadenoson Injection
HIKMA

0.4 mg/5 mL, 5-mL Prefilled Syringe
(143-6225)..... Ea.



RhoGAM® Ultra-Filtered Plus
KEDRION BIOPHARMA

Rho(D) Immune Globulin (Human), 300-µg
Single-Dose Safety-Shield Syringes
(623-0238)..... Ea.
(623-0237)..... 5 per Pkg.
(623-0235)..... 25 per Pkg.
CPT® Code: 90742†



Supartz FX™ Injectable
BIOVENTUS

2.5-mL Prefilled Syringe
(129-3653)..... Ea.
Series of injections: 3 to 5
J Code: J7321



Compare to Depo-Testosterone®
Testosterone Cypionate Injection
HIKMA

200 mg/mL, 10-mL MDV
(137-7452)..... Ea.



Compare to Kenalog®-40
Triamcinolone Acetonide Injection
AMNEAL

40 mg/mL
1-mL SDV
(142-2167)..... Ea.
1-mL SDV
(142-3667)..... 25 per box
5-mL MDV
(142-3668)..... Ea.
10-mL MDV
(142-2227)..... Ea.



Xylocaine® Plain Injectable
FRESENIUS KABI

1%, 2-mL SDV, MPF
(150-0113)..... 25 per Pkg.
1%, 5-mL SDV, MPF
(150-0069)..... 25 per Pkg.
(150-6188)..... 5 amp per Pkg.
1%, 20-mL MDV
(150-0107)..... 25 per Pkg.
50-mL Vial
1% Plain
(248-0849)..... Ea.
(150-0119)..... 25 per Pkg.
2%, 5-mL SDV
(150-0101)..... 25 per Pkg.



Zilretta™ Injectable
PACIRA

32-mg SDV
Triamcinolone Acetonide Extended-Release
(142-9357)..... Ea.
Series of injections: 1
J Code: J3304



Zynrelef™ Injection
HERON THERAPEUTICS

200 mg/6 mg, 10-mL SDV
(140-5769)..... 10 per Pkg.
400 mg/12 mg, 20-mL SDV
(140-5768)..... 10 per Pkg.



Compare to Tylenol®

Acetaminophen Regular-Strength Tablets Unit-Dose Package



Nonaspirin pain reliever and fever reducer.

- 2 tablets per pouch
- Convenient dispenser box
- Made in the USA

325 mg.
(900-4689).....500 per Box
Contains: 250 packets of 2.



Compare to Extra Strength Tylenol®

Acetaminophen Extra-Strength Tablets Unit-Dose Package



Nonaspirin pain reliever and fever reducer.

- Manufactured and packaged in the USA
- Highest quality-control standards under the FDA
- Packaged in 250 sealed-in pouches of 2 tablets each
- Convenient dispenser box
- Color-coded inner and outer packaging

500 mg.
(900-4690)..... 500 per Box
Contains: 250 packets of 2.



Compare to Advil® or Motrin®

Ibuprofen Tablets Unit-Dose Package



Powerful relief from pain and inflammation.

200 mg
(570-0177).....500 per Box
Contains: 250 packets of 2 tablets.



Isopropyl Alcohol, 70% HYDROX LABORATORIES

4 oz
(128-0491).....Ea.



Alcohol



70% Isopropyl Plastic Container

(112-7067)..... 16 oz.
(112-7011)..... 32 oz.
(112-7068)..... Gal. ♦
♦ Hazardous material

TheraTears® Extra® Dry Eye Therapy Lubricant Eye Drops

PRESTIGE CONSUMER HEALTHCARE

0.5-oz Bottle
(141-2098)..... Ea.
0.5-oz Bottle
(141-2446)..... 2 per Pkg.
Preservative-free, 0.6-oz Single-Use Vials
(141-2445)..... 30 per Box



Compare to Refresh Plus®

TheraTears® PRESTIGE CONSUMER HEALTHCARE

Unlike ordinary artificial tears, TheraTears replicates human tears with its patented electrolyte formula that provides soothing, long-lasting, dry-eye relief. Moisturizes and protects sandy, scratchy, gritty, irritated eyes.

Single-Use Package
(141-0403)..... 30 per Box
Contains: 6 foil pouches of 5 single-use packages.

TheraTears® Dry Eye Therapy Lubricant Eye Drops PRESTIGE CONSUMER HEALTHCARE

½-oz Bottle
(141-0402)..... Ea.
1-oz Bottle
(141-2088)..... Ea.



Lubricating Jelly



Sterile

3 gm. Pack
(570-3105)..... 144 per Box
4-oz Flip-Top Tube
(570-3106)..... 4 oz. Tube



Compare to Dacriose®, Collyrium Eye Lotion Medi-Wash Eye & Skin Irrigant MEDIQUE

Sterile, isotonic, buffered solution flushes foreign particles. Single-use plastic bottle.

1-fl-oz Plastic Bottle
(118-2533)..... 1 oz.
4-fl-oz Plastic Bottle
(151-4399)..... 4 oz.
16-fl-oz Plastic Bottle
(151-0009)..... 16 oz.



Compare to Halls®

Medikoff® Drops-Bulk Refill MEDIQUE

Menthol/eucalyptus vapor-action cough drops.

Cherry
(151-4400)..... 600 per Box



Compare to Naprosyn®, Aleve® Mediproxen Tablets Unit-Dose Package MEDIQUE

Has the same pain-relief and fever-reduction ingredient found in nationally advertised brands Naprosyn® and Aleve®. Unit-dose packets are safe, hygienic, and tamper-resistant. Naproxen sodium, 220 mg.

(993-0948).....50 per Box
Contains: 50 packets of 1.
Naproxen sodium, 220 mg.
(151-9200)..... 100 per Box
Contains: 100 packets of 1.

Aquaphor® Healing Ointment Unit-Dose Package
BEIERSDORF

0.35-oz Packets
(911-0011)..... 144 per Box



Aquaphor® Healing Ointment
BEIERSDORF

Original Formula
Treatment of dry, damaged skin and minor burns. Hypoallergenic: no fragrances or preservatives to irritate sensitive skin. Great for turf burns.

1.75 oz. Tube
(911-0837)..... Ea.

#63608, 14-oz Jar
(911-3536)..... 14 oz. Jar



Boudreaux's Butt Paste®
PRESTIGE CONSUMER HEALTHCARE

4-oz Tube
(131-8864)..... Ea.



Compare to Blistex®
Lip-Guard
MEDI-FIRST

Temporarily protects and helps soothe chapped or cracked lips. Provides temporary relief of pain and itching associated with minor lip irritations and cold sores. Contains 4 medications.

1/57-oz. Packets
(151-5789)..... 20 per Box
Contains: 20 packets of 1/57 oz.



Nix® Lice Killing Spray
PRESTIGE CONSUMER HEALTHCARE

5-oz Bottle
(131-8954)..... Ea.



Nix® Lice Crème Rinse
PRESTIGE CONSUMER HEALTHCARE

2-oz Bottle
(131-8942)..... Ea.

2-oz Bottle
(131-8943)..... 2 per Pkg.



Nix Ultra® Lice Removal Kit
PRESTIGE CONSUMER HEALTHCARE

(136-9432)..... Ea.
Contains: 1 of each: Nix Ultra® Lice Removal Comb, 3.4-oz Nix Ultra® solution & 5-oz Nix® Lice Killing Spray.

Nix Ultra® Lice Treatment
PRESTIGE CONSUMER HEALTHCARE

- Kills lice and their eggs
- Kills super lice
- Fast and effective
- Pesticide-free and non-toxic
- Easy application
- Tough on lice; gentle on skin

Lice Shampoo
(135-5184)..... 4 oz per Box

Lice Treatment Liquid
(135-5618)..... 3.4 oz per Box



Polident® Denture Essential Kit
HALEON

Designed for patients to get started right away with proper denture care. Polident® microclean formula removes tough stains and kills 99.9% of odor-causing bacteria. Super Poligrip® is clinically proven to seal out more food particles, even in well-fitting dentures.

Starter Kit
(146-3885)..... 12 Kits per Case
Contains: 1-0.3 oz. tube of Super Poligrip® adhesive, 6 Polident® 3-minute cleansing tablets, 1 denture brush & 1 denture bath.



Safetec Sunscreen Lotion
SAFETEC OF AMERICA, INC.

SPF 30, 1/8-oz Packets
(115-5375)..... 25 per Box



Polident® Tablets
HALEON

Designed to fight stains and kill 99.9% of odor-causing bacteria. Clean dentures daily by soaking and brushing. Offers a complete clean feeling with triple mint freshness.

3-Minute
(583-1271)..... 40 per Box
(583-0396)..... 84 per Box

Overnight Whitening
(583-1293)..... 40 per Box
(583-0339)..... 84 per Box

HENRY SCHEIN®
BRAND PRODUCTS

Rely on Us™

Sharps Containers

From A Name You Can Trust.

- For use by professional health-care givers
- Meets all NIOSH and OSHA regulations
- FDA approved
- Puncture-resistant
- Transparent body and/or top allows for easy viewing
- High-impact and shatter-resistant polypropylene



#0310-1500-HS, 1 qt
(112-7149) ea
#0310-150T-HS, 1.5 qt
(570-0376) ea
#0322-150R-HS, 2.2 qt
(570-0213) ea
#0354-150M-HS, 5 qt
(112-7152) ea

#0354-150B-HS, 5.4 qt
(112-7151) ea
#0319-1500-HS, 1-gal Slide Lid
(570-0378) ea
#0319-150R-HS, 1-gal Rotary Lid
(570-0379) ea
#0320-1500-HS, 2 gal, Slide Lid
(112-7150) ea

#0320-150R-HS, 2-gal Rotary Lid
(570-0377) ea
#0330-150R-HS, 3-gal
(570-0212) ea
#0380-150R-HS, 8 gal
(570-0380) ea



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!



Sharpstainer® Disposable Container MEDICAL ACTION

Red, autoclavable, and puncture-resistant.
Tamper-proof locking cap. One-way safety insert.

Large—3½" x 3½" x 10" (1.7 Qt. Capacity)
(153-3888) Ea.

Phlebotomy Container, 1-Quart BEMIS HEALTHCARE, INC

- Fits phlebotomy trays
- Dual openings; vertical drop
- Dimensions: 6½"H x 4¼"L x 4¼" W
- Uses holder #400 to secure container on tray or tabletop

#100-030, Translucent Red
(122-4680) Ea.



Sharpstainer® Disposable Container MEDICAL ACTION

Red, autoclavable, and puncture-resistant.
Tamper-proof locking cap. One-way safety insert.

Small—3½" x 3½" x 5" (0.7 Qt. Capacity)
(153-3569) Ea.

Large—3½" x 3½" x 10" (1.7 Qt. Capacity)
(153-3888) Ea.

Medium—6.2 Qt.
(9½"H x 10"W x 6"D)
(153-2566) Ea.

10 qt
(120-4364) Ea.



Sharps Containers DYNAREX

- Durable and puncture-resistant
- Horizontal drop maximizes container volume
- Locks for final disposal
- Transparent lid helps reduce risk of overfilling
- Red

Slide Lid

#4622, 1 qt
(139-8160) Ea.

#4626, 1 gal
(139-8163) Ea.

#4629, 6 gal
(139-8170) Ea.

#4625, 5.4 qt
(139-8169) Ea.

#4628, 3 gal
(139-8171) Ea.

Safety IV Catheter

From A Name You Can Trust.

Triple Facet Back Cut Needle
Provides ease of insertion and reduced trauma to puncture site

Push-Off Tab
Assists with catheter advancement

Push-Button Mechanism
Retracts needle fully and safely into catheter chamber

Fully Encapsulating Needle Shield
Reduces exposure to bloodborne pathogens via needle stick injuries

Gently Tapered Catheter Tip
Provides smooth transition from needle to catheter during insertion

Contoured Housing
Ergonomically designed grips allows for one-handed or two-handed techniques

Large-Volume Chamber
Allows clear visual confirmation through flash-back and easy handling



Needle
(Featured in sizes 20g–24g)
The notched needle promotes fast confirmation of venipuncture



Catheter Material
PUR material designed to remain kink free and soften within the body

Item Code

(570-2586)
(570-2591)
(570-2587)
(570-2588)
(570-2589)
(570-2590)

Description

24G x 0.75"
24G x 0.5"
22G x 0.1"
20G x 1.16"
20G x 1"
18G x 1.16"

Quantity

50 per Box
50 per Box
50 per Box
50 per Box
50 per Box
50 per Box




RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

I.V. Medication Delivery

From A Name You Can Trust.

Medication Delivery Sets and Connectors

Henry Schein provides market-leading infection control technologies to help minimize bacterial contamination. We carry a wide selection of injection and needle-free medication delivery ports to meet a variety of medication delivery needs.



Features & Benefits:

- Flat top needle-free valve with tight septum seal helps to support effective disinfection
- Metal free for MRI compatibility
- DEHP-free, latex-free and sterile

Item Code	Type	Product Description	Qty/Unit
(570-0881)	Primary I.V. Administration Set	I.V. Administration Set. 15 drops/mL. 73" Vented, Slide Clamp, Roller Clamp, Injection Y-site, Rotating Male Luer Lock	Each
(570-0885)	Primary I.V. Administration Set	I.V. Administration Set. 15 drops/mL. 86" Vented, Roller Clamp, Needle-free Y-site, Rotating Male Luer Lock	Each
(570-0887)	Primary I.V. Administration Set	I.V. Administration Set. 15 drops/mL. 105" Vented, Back Check Valve, Slide Clamp, (2) Needle-free Y-sites, Roller Clamp, Rotating Male Luer Lock	Each
(570-0883)	Secondary I.V. Administration Set	Secondary I.V. Administration Set. 15 drops/mL. 40" Vented, Roller Clamp, Hanger, Rotating Male Luer Lock	Each
(570-0884)	Microbore Extension set	Microbore extension set. 8" Removeable Needle-free Connector, Slide Clamp, Rotating Male Luer Lock	Each
(570-0882)	Needle-free Connector	Needle-free Connector	Each



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

accucold®

Laboratory Cold Storage Solutions



Advanced Temperature Control



Extended Warranty



Multiple Alarms



Quick Ship



Eco-Friendly



QUICK SHIP
With 2-3 Day
Lead Time



Pharma-Vac Refrigerators

GLASS HSM	SOLID HSM	GLASS FSI	SOLID FSI	H x W x D	CAPACITY
136-3278	137-4190	ARG1PV	ARS1PV	21.5" x 17.5" x 19.5"	1 cu.ft.
146-0362	146-0394	ARG2PV	ARS2PV	22.25" x 18.5" x 23.25"	1.94 cu.ft.
139-7157	139-7157	ARG31PVBIADA	ARS32PVBIADA	32" x 19.75" x 22.63"	2.83 cu.ft.
136-3313	136-3287	ARG3PV	ARS3PV	34" x 18.5" x 19"	3 cu.ft.
139-8154	139-8155	ARG61PVBIADA	ARS62PVBIADA	31.75" x 23.38" x 25"	6 cu.ft.
136-3244	138-5945	ARG6PV	ARS6PV	32.5" x 23.38" x 24.38"	6 cu.ft.
136-3503	136-3502	ARG8PV	ARS8PV	50" x 23.38" x 24.38"	8 cu.ft.
136-3507	136-3506	ARG12PV	ARS12PV	61.75" x 23.38" x 24.38"	12 cu.ft.
136-3547	136-3550	ARG15PV	ARS15PV	72" x 23.38" x 24.38"	15 cu.ft.
141-1336	141-1186	ARG18PV	ARS18PV	78.5" x 27.5" x 26"	18 cu.ft.

Pharma-Vac Freezers

SOLID HSM	SOLID FSI	H x W x D	CAPACITY
146-0222	AFZ1PV	22.25" x 18.5" x 23.25"	1.4 cu.ft.
145-4557	AFZ2PVBIADA	32.0" x 19.75" x 22.63"	2.47 cu.ft.
145-4556	AFZ5PVBIADA	31.75" x 23.38" x 25"	3.88 cu.ft.


Pharma-Lab Refrigerators

HSM	FSI	H x W x D	CAPACITY	DOOR(S)	DOOR SWING
136-3813	ARS23ML	83.75" x 27.5" x 31"	23 cu.ft.	(1) Stainless	RHD
138-8165	ARS23MLLH	83.75" x 27.5" x 31"	23 cu.ft.	(1) Stainless	LHD
136-3815	ARG23ML	83.75" x 27.5" x 31"	23 cu.ft.	(1) Glass	RHD
138-3121	ARG23MLLH	83.75" x 27.5" x 31"	23 cu.ft.	(1) Glass	LHD
136-4008	ARS49ML	83.75" x 55.25" x 31"	49 cu.ft.	(2) Stainless	French
136-4116	ARG49ML	83.75" x 55.25" x 31"	49 cu.ft.	(2) Glass	French

Pharma-Lab Freezers

HSM	FSI	H x W x D	CAPACITY	DOOR(S)	DOOR SWING
136-3814	AFS23ML	83.75" x 27.5" x 31"	23 cu.ft.	(1) Stainless	RHD
139-3861	AFS23MLLH	83.75" x 27.5" x 31"	23 cu.ft.	(1) Stainless	LHD
136-4115	AFS49ML	83.75" x 55.25" x 31"	49 cu.ft.	(2) Stainless	French

accucold.com



High-Performance Undercounter Medical Freezer

PHC CORP OF NORTH AMERICA

Designed to maintain temperatures between 5°F and -13°F (-15°C and -25°C) with a tight uniformity of ±5°C for storing critical vaccines and biologics, this model is engineered to meet CDC guidelines in a small footprint. Precise temperature control and uniformity protect vaccine and biologic viability whether the unit is loaded or empty.

- Capacity: 5 cubic ft
- 3 coated steel wire shelves: 19.2" x 14.75"
- Temperature control range: 5°F to -13°F (-15°C to -25°C) in 0.1°C increments
- Alarms: high/low temperature
- Refrigerant: 6-600a
- Door swing: right hinge
- Exterior dimensions: 34.1"H x 23.8"W x 23.7"D
- Interior dimensions: 25.6"H x 19.9"W x 17.1"D

#PF-L5181W-PA
(141-5626)..... Ea.



Premier Pharmacy/Vaccine Freestanding Undercounter Refrigerator

AMERICAN BIOTECH SUPPLY

Features a microprocessor temperature controller and display, allowing for precise temperature control, verification, and recovery versus mechanical-style thermostats.

- Glass swing door with keyed lock
- 3 adjustable shelves with guard rail on back to prevent product from being pushed against the back wall
- Bottom-footed shelf to allow air circulation
- Back wall evaporator cover
- Forced draft circulation
- Warranty: 2-years on parts and labor; 3 years on compressor parts only
- UL/C & UL listed
- Energy Star certified
- Digital temperature display
- Alarms: audible and visual high/low
- 2 probes: 1 in air and 1 in sample bottle
- Capacity: 5.2 cubic ft
- Temperature range: 36°F to 46°F (2°C to 8°C)
- Exterior dimensions: 32½" H x 23¾" W x 24" D

#PH-ABT-HC-UCFS-0504G
(125-4635)..... Ea.



Pharma-Vac Performance Series Vaccine Refrigerators ACCUCOLD

This line of refrigerators supports meeting regulations and guidelines from the CDC, State Agencies, World Health Organization, and other health ministries on vaccine storage and handling. Designed with a hydrocarbon refrigerant, they provide an energy-efficient solution to vaccine storage. Units are Pyxis, Omnicell, and AcuDose RX compatible.

- Certified to the NSF/ANSI 456 Vaccine Storage Standard
- Adjustable-cycle defrost operation ensures the unit avoids temperature spikes and ice buildup
- Battery-powered alarm sounds if the power connection is interrupted
- Digital controls with a large LED readout in Celsius or Fahrenheit for intuitive temperature management
- Multiple shelf positions spaced at ½" to accommodate virtually any size
- Ecofriendly hydrocarbon-based refrigerant for a greener footprint that meets current CARB regulations
- Environmentally friendly design without ozone-damaging chemicals

#ARG1PV456, 1-cu-ft Compact
(143-5669)..... Ea.
#ARG15PV456, 15-cu-ft Upright
(143-5665)..... Ea.
#ARG18PV456, 18-cu-ft Upright
(143-5753)..... Ea.



**i.Series® iLF105-GX Laboratory Freezer
HELMER SCIENTIFIC**

5.3-cu-ft Capacity, Solid Door, -15 to -30°C

A GX Solutions professional medical-grade, -30°C freezer. Designed for the unique needs of critical health care applications, such as storing samples, reagents, controls, and frozen vaccines. It delivers optimized control in 3 important areas: temperature, noise, and energy management.

- Powered by OptiCool™ technology, this freezer provides improved performance and user experience through energy efficiency, low noise levels, and stable storage environments to safeguard valuable patient samples, frozen vaccines, and more
- 30% more efficient than traditional technology
- Digital microprocessor temperature controller, with built-in alarms, includes minimum/maximum display and reset
- Energy Star® certified
- Powerful i.C3® Information Center is exclusive to i.Series® models and provides constant temperature monitoring, multiple information logs, and security features
- GX Solutions is the first professional medical-grade line of freezers for clinical laboratories, pharmacies, and other health care facilities
- NSF/ANSI 456 Vaccine Storage Standard certified

#5222105-1
(141-9439)..... Ea.



01




The Next Generation of Midmark® Steam Sterilizers

Midmark Steam Sterilizers are designed inside and out to make instrument processing as easy and as automated as possible.

- 01** Navigate cycle setup and other processes with a clear, fingerprint-resistant 5-inch touchscreen you can utilize while wearing gloves.
- 02** Determine cycle status—cycle in progress, cycle complete and cycle error—from a distance with a progressing color-coded LED light bar, large cycle countdown clock and audible cues.
- 03** Stay audit-ready with device reminders, user authentication, unlimited storage of routine care events for the life of the sterilizer, and automated cycle recordkeeping.

For more information, visit:
henryschein.com/medical



02



03



© 2024 Midmark Corporation, Versailles, Ohio USA

Brewer Assist Power Procedure Table™

BREWER

This procedure table can be lowered down to 19", making wheelchair transfers effortless. The 450-lb weight capacity and the additional 5" leg extension accommodates larger patients more comfortably for safer clinical outcomes. As with all Brewer power tables, the Assist comes with our standard 3-year warranty and 7-day lead-time. Available in 12 standard seamless upholstery colors.

Product benefits:

- 450-lb weight capacity
- Standard 3-year warranty
- Technologically advanced
- Lowers to 19" for lateral patient transfer
- Exclusive 3-point pivot headrest
- Seamless, easy-to-clean upholstery
- Footrest extends an additional 5"
- Intuitive foot controls
- 12 standard colors with optional upgrades
- Wide assortment of accessories available

Procedure Table..... Ea.

Specify:

#7000-21, Clamshell	(124-1394)
#7000-22, Gunmetal	(124-1395)
#7000-23, Feather	(124-1396)
#7000-24, Black Satin	(124-1398)
#7000-25, Saddle	(124-1399)
#7000-26, Tapestry Red	(124-1400)
#7000-27, Cabernet	(124-1401)
#7000-28, Cocoa	(124-1402)
#7000-29, Deep Sea	(124-1403)
#7000-30, Ivy	(124-1404)
#7000-31, Deep Sapphire	(124-1405)
#7000-32, Blue Fog	(124-1406)



**830 Series Procedure Chair
MTI**

The highly versatile 830 Series chairs are ADA compliant in all base versions (standard, swivel, or mobile) using the transfer support, providing easy patient ingress/egress, and equal access for patients with disabilities and caregivers. All 830 chairs come equipped with SmartTech™ technology to enhance safety, reduce energy consumption, and provide smoother movement, additional memory options, and more user preferences.

- 725-lb patient load capacity
- 4 programmable positions
- Programmable home and ADA-compliant home positions
- Power lift, back, tilt, and foot extension
- Hand and foot controls
- Oval articulating headrest
- 22" or 25" seat width
- 18¾"-39" seat height
- Leg rest features knee-break with 7" power foot extension
- Optional patient arms: removable, floating (compensating) arms, or seat rail or back rail mounted arms
- Removable upholstery: seamless or designer
- Tapered backrest for easier caregiver access
- Sturdy all-steel frame
- Multiple G2 accessory options

#829-001-80, Standard Base (22" seat width)
(135-8121)..... Ea.

#829-002-80, Swivel Base (22" seat width)
(139-4043)..... Ea.

#829-003-22, Mobile Base (22" seat width)
(140-6171)..... Ea.

#830-001-80, Standard Base (25" seat width)
(133-8866)..... Ea.

#830-002-80, Swivel Base (25" seat width)
(136-9044)..... Ea.

#830-003-22, Mobile Base (25" seat width)
(133-9127)..... Ea.



**626 Barrier-Free®
Examination Chair
MIDMARK CORPORATION**

Many health-care providers are not aware of current standards for proper blood-pressure measurement. If this isn't addressed, it leaves health systems at serious risk. 9.8% of the population is affected by overestimation of high blood-pressure errors. Multiply the cost of treatment for hypertension by the typical patient panel size and the cost of overtreatment could be as much as \$135,000 per each practicing physician. The impact of overmedication can also result in complications such as dizziness, fainting, exhaustion, or even patient falls. The Midmark 626 is the first exam chair of its kind designed to promote proper patient positioning for a more accurate blood-pressure measurement.

- Easily position patients for a better BP measurement
- Powered movement of the back section ensures the patient's back is supported
- A Barrier-Free® low chair-height allows the majority of patients, even female patients with heights in the 3rd percentile, to place their feet flat on the floor

#626-001, Base Model
(133-0301)..... Ea.

#626-002, With Heated Upholstery
(133-0312)..... Ea.

#626-003, With IQ Scale
(133-0323)..... Ea.

#626-004, With Heated Upholstery and Scale
(133-0334)..... Ea.



**Ritter 225 Barrier-Free®
Examination Chair
MIDMARK**

The number of elderly and obese in the U.S. is steadily growing, and over one fifth of adults in America are disabled. Assisting patients with mobility limitations onto exam chairs has become a frequently performed task. Although caregivers working in health-care facilities are in one of the most at-risk professions for musculoskeletal disorders, little has been done to evaluate and reduce risk. The Ritter 225 is designed to offer low and high positioning heights so you can provide the best care to all patients, whether treating the elderly, expectant mothers, or patients with disabilities.

#225-002, Without Clean Assist Roller System
(133-0279)..... Ea.

#225-003, With Clean Assist Roller System
(133-0290)..... Ea.



Clarus 4™ LED Exam Light

MTI

MTI's most economic exam light uses 4 LEDs and is equipped with a high-quality 3" light head. The small light head and flexible neck allow easy movement to many efficient and comfortable positions. Can be used with either a caster base or floor stand. Physical dimensions: head diameter: 3"; caster base height: 73", floor stand height: 71"; caster base diameter: 23¾", floor stand diameter: 15". 2-year limited warranty.

Features:

- 4 white LEDs
- LED life: 40,000 hours
- Color temperature (CCT): 4,300° K
- Light intensity @ 16": 14,100 Lux
- Light intensity @ 12": 21,500 Lux
- Working distance: 12"
- Spot size at 10": 2"
- Spot size at 12": 3"
- Spot size at 16": 4"
- Spot size at 24": 5"
- CRI: 75
- Control box dimmable on light head
- Electrical input: 100–240 V, 50–60 Hz, 40W, 12V @ 3.3A
- Power cord: 3-wire grounded cord set, 8' long
- Listings: UL 153, UL 8750, CAN/CS: C22.2, No. 12

Exam Light with Caster Base

(121-8209).....Ea.

Exam Light with Floor Stand

(121-8210).....Ea.

Exam Lamp

DRIVE MEDICAL DESIGNS

This versatile 16" mobile lamp has a flexible gooseneck design and a dome-style shade. Its adjustable arm offers full 360° flexibility, allowing the user to direct light in any direction. With a weighted steel base for stability, the height can be easily adjusted from 48" to 72" to accommodate individual preferences. The lamp has a 3-prong plug, is UL listed, and supports incandescent light bulbs up to 60 watts, ensuring both safety and convenience. Limited lifetime warranty.

#13408MB

(139-7964).....Ea.



Solid State Portable Headlight
WELCH ALLYN

20 lumens of bright, white light for true tissue color. Variable spot size. Battery is contained in the soft terry headband. 1½ hours of battery life. 10,000 hours of solid-state lamp life.

- 20 lumens of bright, white light for true tissue color
- Variable spot size
- Battery contained in soft terry headband
- 1½ hours of battery life
- 10,000 hours of solid-state lamp life

#46070, Headlight

(566-8627).....Ea.

#49642, Headband for Headlight

(566-0531)..... 10 per Case

#72270, Rechargeable Battery

(566-8948).....Ea.



Ritter 250-003 LED Exam Light

ITTER BY MIDMARK

Features LED technology for longer life, brighter intensity and reduced energy consumption. Provides up to 25,000 hr of lighting—more than 12 times longer than a traditional halogen exam-room light. Light head size: 4½" diam. x 4" long. Height: 36" long, flex length: 22".

- *Light source:* 12V, 3 LEDs, MR16 type (with glass cover). Illumination @ 18": 1,800 fc. Color temperature: 4,000°K. Life: 25,000 hours

Electrical: 115V AC, 50–60 Hz, 0.12A. Power cord: 8 ft hospital-grade plug

- *Mechanical:* Light head size: 4"L x 4½" diameter
- *Flex arm:* 22"L x 0.585" diameter. Neck material: high carbon steel covered in vinyl
- *Base tube:* 36"L. Outer tube diameter: 1". Tube material: steel
- *Standard base housing:* 2¼"H x 13" diameter
- *Weight:* 12 lb
- *Color:* Pebble Grey

#250-003, LED Exam Light

(132-3839).....Ea.



Welch Allyn Green Series™ 600
Minor-Procedure Lights

BAXTER

Featuring 3 LEDs and a wider head, the GS 600 is designed specifically for demanding viewing conditions. Whiter, brighter, and greener: 3 white LEDs offer greater intensity. Other features include touchless power controls and an adjustable light reduce the risk of infection due to cross proved performance and long-lasting durability. Low cost of ownership due to elimination of bulb replacement and lower energy consumption. 5-year limited warranty to include the LED light source and power supply only.

#44600, With Mobile Stand

(566-0088).....Ea.

#44610, With Table/Wall Mount

(566-0089).....Ea.





Quality. Convenience. Comfort.

Graham Medical® is a leading provider of quality single-use medical products designed to meet the demands of medical professionals while providing comfort and protection for patients.

graham
medical®
The Promise of Protection®



WELCH ALLYN AND BAXTER ARE ONE TEAM

JOIN THE FIGHT TO HELP PRESERVE VISION IN YOUR COMMUNITY



WELCH ALLYN SPOT VISION SCREENER:

Detects amblyopic risk
factors in children and adults.

WELCH ALLYN RETINAVUE CARE DELIVERY MODEL:

Detects diabetic retinopathy
in primary care settings.



Scan the QR code to get details on vision screening & diagnostic solutions
for Community Health Centers or visit Hillrom.com/CommunityHealth.

Contact your Henry Schein representative to learn more.

RetinaVue 700 Imager: Rx Only - For safe and proper use of the products mentioned herein, please refer to the appropriate Operator's Manual or Instructions for Use.

Baxter, Hillrom, RetinaVue, Spot and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.

US-FLC158-240059 (v1.0) 07/2024



Pro Pocket Aneroid Sphygmomanometer



A value-priced aneroid that meets the needs of the most demanding professionals. Precision-crafted movement is rated to 30,000 cycles and backed by a 20-year calibration warranty. Extended graduation marks eliminate parallax error. White numbers on contrasting black scale with red indicator needle for easier viewing. Proprietary Size Guide Marking system to enhance measurement ease and prevent miscuffing. Chrome-plated brass deflation valve. Inspected, assembled, and packaged in the USA. Leatherette zippered carrying case. Not made with natural rubber latex.



Pro Plus Series Palm Aneroid



This revolutionary palm-style manometer (#804) has an integral non-latex bulb and trigger-style air release valve for one-hand convenience. Features:

- Precision-crafted, 300mm-Hg, CE-marked, no-pin stop movement tested to 50,000 cycles
- Movement housed in a zinc safety cage for unsurpassed durability
- Ambidextrous infinite control trigger-release air valve
- Nylon cuff with proprietary size guide marking system, to prevent miscuffing
- Dip-molded latex inflation bladder withstands heavy use
- Quick-release Luer-slip connector
- Complete with leatherette carrying case with self-repairing nylon zipper
- Assembled, inspected, and packaged in the USA from foreign components

- Value-priced aneroid that meets the needs of the most demanding professionals
- Precision-crafted movement is rated to 30,000 cycles and backed by a 20-year calibration warranty
- Extended graduation marks eliminate parallax error
- White numbers on contrasting black scale with red indicator needle for easier viewing
- Proprietary Size Guide Marking system enhances measurement ease and prevent miscuffing
- Chrome-plated brass deflation valve
- Inspected, assembled, and packaged in the USA
- Leatherette zippered carrying case
- Not made with natural rubber latex

Child	
(112-6073).....	Ea.
Small Adult	
(112-6071).....	Ea.
Adult	
(112-6072).....	Ea.
Large Adult	
(112-6075).....	Ea.
Thigh	
(112-6074).....	Ea.

Adult, Black, Non-Latex	
(112-6067).....	Ea.
Large Adult	
(112-6086).....	Ea.
Pro Plus Cuff & Bladder, 1 Tube	
(112-6090).....	Ea.
Child	
(112-6087).....	Ea.
Adult	
(112-6088).....	Ea.
Large Adult	
(112-6092).....	Ea.
Thigh	
(112-6091).....	Ea.
Child	
(112-6093).....	Ea.
Small Adult	
(112-6094).....	Ea.
Adult	
(112-6095).....	Ea.
Large Adult	
(112-6096).....	Ea.
Thigh	
(112-6097).....	Ea.



Pro Plus Convertible Cardiology Stethoscope



Stainless steel chest piece. Drum converts from pediatric diaphragm to extra-deep adult bell.

- Sensitive multifrequency diaphragms
- Bilumen stainless steel binaurals
- Ultrasoft, threaded ear tips
- Nonchill bell and diaphragm rim for patient comfort
- Binaurals fixed at 15° angle for comfort and acoustic seal
- Extra-deep bell for low frequency response
- 19" bilumen PVC tubing for soft flexibility

28" Black	
(900-7457).....	Ea.



Digital Ear Thermometer



- For professional or personal use
 - Instant reading results in 1 second
 - Memory feature recalls last reading
 - Fever alarm
 - Low battery indicator
 - Replaceable lithium battery
 - °F or °C measuring mode
 - Thermometer includes 15 disposable probe covers
- | | |
|-----------------|-------------|
| Thermometer | |
| (900-6438)..... | Ea. |
| Probe Covers | |
| (900-6439)..... | 150 per Box |

PRO Digital Thermometer



- Measures body temperature by oral, axillary, or rectal read
 - 60-second response
 - 5 probe covers
 - More than 200 hours of continuous operation
- | | |
|-----------------|-----|
| (570-2055)..... | Ea. |
|-----------------|-----|

Advantage™ 6021N **Automatic Digital Blood** **Pressure Monitor** **AMERICAN DIAGNOSTIC** **CORPORATION**

This fully automatic upper-arm digital blood-pressure monitor is designed for home use. It measures systolic and diastolic pressure as well as pulse. Equipped with Smart Logic Modulation (SLM) technology, it ensures a gentle inflation process for optimal comfort and accuracy. The monitor features a self-adjusting, extended-range non-latex D-bar cuff that fits arms from 8.7" to 16.5".

Key features:

- Automatic inflation and deflation control
- Irregular heartbeat detector
- Stores the last 99 readings with date and time stamp for 2 users
- Automatic shut-off function
- Made without natural rubber latex
- Includes 4 AA batteries
- 5-year warranty

(122-6355)..... Ea.



Welch Allyn Green **Series™ 777 Integrated** **Wall System** **BAXTER**

#777-PM2WAS-US, With PanOptic Basic LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer
(140-1989)..... Ea.

#777-PM3WAS-US, With PanOptic Plus LED Ophthalmoscope, MacroView Plus LED Oscope for iExaminer, BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer
(140-2495)..... Ea.

#777-SM2WAS, With Coaxial LED Ophthalmoscope, MacroView Basic LED Oscope, BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer
(140-1801).....Ea.

#777-SM2WAX, With Coaxial LED Ophthalmoscope, MacroView Basic LED Oscope, BP Aneroid, and Ear Specula Dispenser
(140-1799)..... Ea.



Welch Allyn **SureTemp® Plus 690** **Electronic Thermometer** **BAXTER**

Combines the speed of ear thermometry with the accuracy and reliability of traditional oral thermometry. Takes fast, 4–6-sec. oral temperatures, 10-sec. pediatric axillary and rectal temperatures, and 15-sec. adult axillary temperatures.



- Accurate to within +/- 0.2°F (0.1°C)
- Interchangeable, removable probe wells to reduce the risk of cross contamination
- Large LCD that can display temperatures in Fahrenheit or Celsius
- Last-temperature recall
- Mode selection and recall buttons
- Four different temperature modes
- Convenient storage for 25 probe covers
- Takes about 6,000 readings on 3 AA batteries (included)

Thermometer, with Oral Probe, Well, and 4 ft. Cord
(203-3856).....Ea.

Thermometer, with Wall Mount, Oral Probe, Well, and 4 ft. Cord
(178-8374)..... Ea.

Thermometer with Wall Mount, Oral Probe, and 9 ft. Cord
(454-2689).....Ea.

Probe Covers
(566-0460)..... 250 per Box

Wall Holder Mount for Thermometer
(206-8985).....Ea.

Oral/Axillary Probe, Well, and 9 ft. Cord
(695-8954).....Ea.

Rectal Probe, Well, and 4 ft. Cord
(741-8045).....Ea.

Proscope™ 660 **and 670 Series** **Stethoscopes**

AMERICAN **DIAGNOSTIC** **CORP.**

Nurse-scopes combine lightweight design and sensitive acoustics. Similar to the 660 series, the 670 series features a combination diaphragm/bell dual-head chest piece with nonchill ring. The 660 series features:

- Lightweight aluminum chest-piece construction
- Diaphragm-only, nurse-style chest piece
- Adjustable chrome binaurals
- Flexible, 22" 1-piece molded Y PVC tubing

660 Series.....Ea.

Specify:
#660N, Navy Blue **(130-2573)**
#660BK, Black **(499-0725)**

670 Series..... Ea.

Specify:
#670N, Navy Blue **(130-4588)**
#670BK, Black **(102-2463)**
#670G, Gray **(185-8499)**



Masimo Set® LNCS® Sensors and Cables **MASIMO**

LNCS, or low-noise cabled sensors, are compatible with Masimo SET™ pulse oximeters and Nellcor sensor-compatible instruments without OxiMax™.

#1863, LNCS DC-I Adult Reusable Sensor
(114-1951)..... Ea.

#2653, LNCS DBI, Adult, Soft, 3 ft
(499-9823)..... Ea.

#1864, LNCS DCIP Pediatric Reusable Sensor
(116-0986)..... Ea.

#1859, LNCS Amtx Adult Sensor
(114-1911)..... 20 per Box

#2317, LNCS Amtx3, Adult, 3 ft
(700-1279).....20 per Box

#1860, LNCS Pmtx Pediatric Sensor
(499-6226).....20 per Box

#2318, LNCS Pmtx3, Pediatric, 3 ft
(700-1280)..... 20 per Box

#2918, LNCS E1®, Adult, Ear, 3 ft
(700-1785).....10 per Box

LNCS® Disposable, **Adhesive SpO2** **Sensors with** **Replacement Tape**

#2328, LNCS, Infant, 18"
(116-0769)..... 20 per Box

#2319, LNCS, Infant3, 3 ft
(119-3255)..... 20 per Box

#2329, LNCS Neo-L Neonatal Sensor
(116-0770)..... 20 per Box

#2320, LNCS, Neonatal/Adult, 3 ft
(119-1696).....20 per Box

LNCS® RED 20-Pin Patient Cables for Use with Masimo rainbow® SET® Devices
#2055, Red LNC-04 LNCS Patient Cable, 4 ft
(114-5834)..... Ea.

LNC Patient Cables for Use with Masimo SET® Devices
#2364, LNC1, 1 ft
(125-7805)..... Ea.

#2017, LNC-4 LNCS Patient Cable
(700-1281)..... Ea.





ePM 10A Configured Patient Monitor

MINDRAY

Exceptional user interface and enhanced parameter functionality, delivering a monitoring solution to fit multiple care environments. Immediate and measurable benefits of the smartphone-like user interface that delivers 98% of all common monitoring functions with just 2 gestures, reducing time clinicians spend navigating monitor functions, and more time on what matters most: patient care. Advanced monitoring options include 3- and 5-lead ECG, SpO₂, NIBP, respiration, and temperature supporting multiple clinical workflows. ePM 10A has a compact, ergonomic design that is perfectly suited for healthcare facilities managing a steady flow of patients daily. Add optional CO₂ capability to aid clinicians to identify signs of patient deterioration due to respiratory distress. Manufactured with industry-leading plastics that stand up to the most stringent medical-grade cleaners. The concealed grab-and-go handle eliminates hard-to-clean crevices and eases clinical workflow while on the move. Flexible mounting options for bedside and transport.

Dimensions: 10.6" x 9.9" x 6.3"

Weight: 8.8 Lb.

Warranty: 3 years

Integrated
Sidestream CO₂, Recorder.....Ea.

Specify:
Masimo SpO₂ (139-2003) Nellcor SpO₂ (139-2005)

Recorder.....Ea.

Specify:
Masimo SpO₂ (139-1920) Nellcor SpO₂ (139-1927)

ePM 12MA Patient Monitor

MINDRAY

Exceptional user interface and enhanced parameter functionality, delivering a monitoring solution to fit multiple care environments. Immediate and measurable benefits of the smartphone-like user interface that delivers 98% of all common monitoring functions with just 2 gestures, reducing time clinicians spend navigating monitor functions, and more time on what matters most: patient care. Advanced monitoring options are integrated: 3- and 5-lead ECG, SpO₂, NIBP, respiration, and temperature-supporting multiple clinical workflows. The ePM 12MA features a semimodular design with a plug-and-play module structure that adds monitoring functions such as CO₂ and multigas, where and when needed. The 12MA eGM multigas module automatically identifies and measures the inspired and expired values of O₂, CO₂, N₂O and 5 anesthetic agents, providing a compact, yet comprehensive OR monitoring solution. Manufactured with industry-leading plastics that stand up to the most stringent medical-grade cleaners. Concealed grab-and-go handle eliminates hard-to-clean crevices and eases clinical workflow while on the move. Flexible mounting options for bedside and transport.

• Dimensions: 12.2" x 11.4" x 6.7"

• Weight: 10.6 Lb.

• Warranty: 3 years

Semi-Configured/Modular.....Ea.

Specify:
Masimo SpO₂ (139-2126) Nellcor SpO₂ (139-2128)

ePM AG Module
(139-3276).....Ea.

CO₂ Module Adult/Pediatric
(137-8484).....Ea.



iM60 Patient Monitor

EDAN

The iM60 monitor spans accurate monitoring across different care environments, acuity levels, and clinical requirements. All data is clearly displayed on a 10½" color touch screen while providing real-time readings on NIBP, ECG, SpO₂, respiration, temperature, and optional EtCO₂.

#iM60.P, Patient Monitor with Printer
(974-0177).....Ea.

#iM60.CS.P, Patient Monitor with CO₂ and Printer
(974-0180).....Ea.

#iM60.P.MT, Patient Monitor with Printer and Rolling Stand
(974-0185).....Ea.

#iM60.CS.P.MT, Patient Monitor with CO₂, Printer, and Rolling Stand
(974-0186).....Ea.



iM70 Patient Monitor

EDAN

The iM70 Patient Monitor is embedded with the latest advancements in design to provide accurate monitoring across different applications that range from general care to anesthesia monitoring. Its high-resolution, 12" touch screen clearly displays all patient data (NIBP, ECG, SpO₂, respiration, temperature, EtCO₂ optional, anesthetic agents optional).

#iM70.P, Patient Monitor with Printer
(974-0178).....Ea.

#iM70.CS.P, Patient Monitor with CO₂ and Printer
(974-0181).....Ea.

#iM70_OR, Patient Monitor with Multi-Gas Analyzer and Printer
(974-0131).....Ea.

#PM.RT, Rolling Stand
(127-1083).....Ea.



Connex® Spot Monitor

HILLROM

Features an easy-to-use, vivid touch-screen display and provides accurate vital signs measurement, including blood-pressure averaging, spot checking, interval monitoring and custom scoring across patient populations. Device connects wirelessly to an EMR and accurately sends vitals to the patient chart from the point of care to help improve efficiency.

- Designed for the physician-office environment with an easy-to-use, vivid touch-screen display
- More accurate hypertension diagnosis without slowing down your workflow thanks to automated blood-pressure averaging
- Improved capabilities for assessing respiratory conditions such as COPD, asthma, and community-acquired pneumonia (CAP), thanks to integrated professional-grade pulse oximetry
- Increase staff efficiency and eliminate data-entry errors with wireless transfer of vital signs directly into an EMR

#71XX-B, Blood Pressure
(566-0464).....Ea.

#71XT-B, BP, SureTemp
(672-0176).....Ea.

#71XE-B, BP, Braun PRO6000
(124-6103).....Ea.

#71WX-B, BP, Nonin SpO₂
(672-0177).....Ea.

#71WT-B, BP, Nonin SpO₂, SureTemp
(566-0466).....Ea.

#71WE-B, BP, SpO₂, Braun Pro 6000
(124-5328).....Ea.

Bluetooth Connectivity:

#73XT-B, SureBP, SureTemp
(566-0465).....Ea.

#73XE-B, SureBP, Braun PRO 6000
(124-8200).....Ea.

#73WX-B, BP, Nonin SpO₂
(672-0178).....Ea.

#73WT-B, PB, Nonin SpO₂, SureTemp
(672-0179).....Ea.

#73WE-B, SureBP, Nonin SpO₂, Braun PRO 6000
(124-4801).....Ea.

#73CT-B, BP, Covidien SpO₂, SureTemp
(672-0175).....Ea.

#73CE-B, SureBP, Covidien SpO₂, Braun PRO 6000
(125-8900).....Ea.

#73MT-B, SureBP, Masimo SpO₂, SureTemp
(566-0467).....Ea.

#73ME-B SureBP, Masimo SpO₂, Braun PRO 6000
(124-7063).....Ea.



Butterfly iQ+
Handheld Whole-
Body Ultrasound

BUTTERFLY NETWORK, INC.

A handheld, whole-body, point-of-care ultrasound system built to last through the toughest shifts. Powered by the Ultrasound-on-Chip™ technology, Butterfly iQ+ is a one-probe solution for use by qualified healthcare practitioners across the acute care and pre-hospital settings. With its 20+ clinical presets and a variety of AI tools, iQ+ is the ultimate system, providing access to point-of-care imaging. The probe can be used for diagnostic ultrasound imaging and measurement of anatomical structures and fluids of adult and pediatric patients across various clinical applications.

- #IQ+LBLK, Butterfly iQ+ Probe
With lightning connector and charging base with 4 regional adapters.
(143-7089).....Ea.
 - #IQ+CBLK, Butterfly iQ+ Probe
With USB-C connector and charging base with 4 regional adapters.
(143-7090).....Ea.
 - #SUB-00001-S, Ultrasound for You
Membership allows users to store and access studies from mobile/web. Includes single-user license, customized archives for unlimited storage, and anonymized study sharing.
(143-7083).....Ea.
 - #SUB-00016-T, Ultrasound for Teams
Membership allows users to store, access, and collaborate. Includes customized archives for unlimited storage and anonymized study sharing. Includes 5 users with the ability to add additional seats of 5 at a discount. Add using the "Team Membership - Additional User" product. Maximum seats for team is 10.
(143-7085).....Ea.
- Rx only (USA). For external use only by trained healthcare practitioners. To ensure safe and effective use, please refer to the product labeling for information on intended use, warnings, and contraindications. EMS practitioners should use iQ+ in accordance with relevant state licensing requirements.



Vscan Air™
GE

This wireless, handheld ultrasound system enables whole-body scanning and delivers crystal-clear images.

- The power of a high-performance ultrasound machine in a small, lightweight wireless solution
 - Optimized for the most common adult and pediatric physical exams to support a broad range of applications
 - The dual probe's 2-sided design allows for both shallow and deep exams
- #H8031VA, Vscan Air Base Package
(140-0608).....Ea.
 - #H45611AG, Vscan Air Protective Carrying Case
(140-0615).....Ea.



Ultrasound Gel



Multipurpose, hypoallergenic, medium-viscosity ultrasound gel that is bacteriostatic for patient satisfaction. Excellent conductive agent with superior lubricating and clinical features for the most common professional diagnostic and imaging procedures. Simple operation and disposal. Long shelf life.

8½-oz Bottle, Blue
(900-4094).....Ea.

5-liter Bottle, Blue
(900-4095).....Ea.



Clear Ultrasound Gels



Multipurpose, hypoallergenic, medium-viscosity gel is clear, dye-free, and bacteriostatic for patient satisfaction. It is also an excellent conductive agent with superior lubricating and clinical features, making it suitable for the most common professional diagnostic and imaging procedures. For use in therapeutic ultrasound, diagnostic imaging, and other external procedures that require a thick medium with excellent lubricating and glide characteristics. Will spread evenly without running or drying under standard clinical applications.

8½-oz Bottle, Clear
(900-4352).....Ea.

5-liter Bottle, Clear
(900-4353).....Ea.

ZOLL Medical's interactive AEDs guide you through the process of performing high-quality CPR with voice, text, and image prompts. Rigorous, automatic self-testing ensures they are always ready. Each AED has unique features to provide the best solution for your needs.

ZOLL®



AED Plus®

- * Real CPR Help® technology
- * Industry leading five-year battery and pad life
- * One-piece pad design

(866-7559) Semi-Automatic Unit

Package includes: AED, electrodes, batteries, carry case.



Powerheart® G5 AED

- * Intellisense® CPR Feedback (ICPR)
- * Non-polarized pads (can be placed in either position)
- * Dual language functionality

(139-5101) Semi-Automatic Unit

Package includes: AED, electrodes, battery.



ZOLL R Series®

- * Professional defibrillator/monitor
- * OneStep Complete™ electrodes with integrated three-lead and Real CPR Help® technology

(127-9979) Professional Defibrillator

Package includes: Monitor/defibrillator, battery, 2 sets of OneStep electrodes.

DONJOY® VELOCITY+

Ankle Brace - Moderate and Extra Support

Velocity® represents an evolution in DonJoy® ankle brace technology, utilizing several patented proprietary design features, resulting in enhanced ankle brace performance and providing the foundation for successful ankle injury prevention and treatment. The Velocity+ ES brace features a foot and ankle speed wrap to apply compression for acute ankle injuries and to control inflammation. (MS features ankle speed wrap only). The new proximal extension cuff provides added stability above the ankle and may be removed as needed during recovery process.

The Velocity+ brace is intended to prevent abnormal eversion and inversion.

EXTRA SUPPORT

HS ITEM CODE	DJO PART NUMBER-LEFT	HS ITEM CODE	DJO PART NUMBER-RIGHT	CALF CUFF	SIZE
(146-5857)	11-1517-2-06000	(146-5860)	11-1516-2-06000	Standard	S
(146-5858)	11-1517-3-06000	(146-5861)	11-1516-3-06000	Standard	M
(146-5859)	11-1517-4-06000	(146-5862)	11-1516-4-06000	Standard	L

MODERATE SUPPORT

HS ITEM CODE	DJO PART NUMBER-LEFT	HS ITEM CODE	DJO PART NUMBER-RIGHT	CALF CUFF	SIZE
(146-5869)	11-1513-2-06000	(146-5872)	11-1512-2-06000	Standard	S
(146-5870)	11-1513-3-06000	(146-5873)	11-1512-3-06000	Standard	M
(146-5871)	11-1513-4-06000	(146-5874)	11-1512-4-06000	Standard	L

For more information, please contact your Henry Schein Representative

enovis Parent company of **DJO®**

Individual results may vary. Neither DJO, LLC nor any of the Enovis companies dispense medical advice. The contents of this document do not constitute medical, consult your healthcare professional for information on the courses of treatment, if any, which may be appropriate for you.



3M Medical Securement Solutions

Drive efficiency in your clinical practice and reduce tape SKUs.

A choice that's critical to care. At Solventum, formerly 3M Health Care, we've made it easy to choose the right securement for each patient and application, with four classes of securement that can work harder and smarter.



Non-adhesive



3M™ Coban™ NL Non-Latex Self-Adherent Wrap†

Non-adhesive securement for clinical applications such as blood draws, dressings, immobilization, securement for difficult to dress areas (head, fingers, toes) and support, and mild compression for soft tissue injuries (i.e. strains, sprains).*



Multi-purpose



3M™ Micropore™ S Surgical Tape†

Multi-purpose securement for clinical applications such as blood draws, light-weight dressings, IV lines and tubing (secondary securement), and non-critical tubes.



Flexible



3M™ Medipore™ H Soft Cloth Surgical Tape

Flexible securement for clinical application such as dressing and added pressure, IV lines and tubing (secondary securement), when swelling or movement is anticipated, central venous catheters (secondary securement), chest tubes, and surgical drain tubes.



High-adhesion



3M™ Durapore™ Surgical Tape

High-strength securement for clinical applications such as urinary catheters, nasogastric tubes, orogastric tubes, and patient positioning.



Contact your Henry Schein representative or visit [3M.com/MedicalSecurement](https://www.3m.com/MedicalSecurement) to learn more.

†Skin performance: Designates products that deliver the securement power you need while minimizing damage to skin.

*May use as a component in a compression wrap system support only under supervision of a wound care specialist.

© Solventum 2024. All rights reserved. Solventum and the S logo are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. Other trademarks are the property of their respective owners.

3M™ Medipore™ H Soft Cloth Surgical Tape SOLVENTUM, FORMERLY 3M HEALTHCARE



- Conformability: bidirectional stretch accommodates swelling, distension, and movement
- Easy to tear at perforations: no scissors necessary
- Helps reduce risk: individually packaged, single-use length rolls can help protect patients by reducing the risk of cross contamination
- Not made with natural rubber latex

Standard Rolls, 10 yd

#2861, 1" Width (887-7957).....	24 per Box
#2862, 2" Width (139-3011).....	12 per Box
#2863, 3" Width (777-3265).....	12 per Box
#2864, 4" Width (471-1595).....	12 per Box
#2866, 6" Width (328-0293).....	12 per Box
#2868, 8" Width (110-8807).....	6 per Box

Single-Use Rolls, 2 yd

#2860S-1, 1" Width (137-6301).....	72 per Box
#2860S-2, 2" Width (137-6302).....	48 per Box
#2860S-4, 4" Width (137-6303).....	24 per Box
#2860S-6, 6" Width (137-6304).....	16 per Box



Transpore™ Hypoallergenic Transparent Tape

SOLVENTUM, FORMERLY 3M HEALTHCARE

10-yd Rolls

Hypoallergenic tape is not made with natural rubber latex, is perforated for bidirectional tear, and is easy to handle with gloves.

- Easy to handle with gloves
- Good adhesion to skin and tubing
- Transparent
- Porous
- Water-resistant
- Hypoallergenic: not made with natural rubber latex

Standard Rolls

10-yd. Rolls.....Box

Specify:

½" Width, 24 per Box	(777-9735)
1" Width, 12 per Box	(777-7305)
2" Width, 6 per Box	(777-9422)
3" Width, 4 per Box	(635-1724)

Single-Use Rolls

1½ yd.....Box

Specify:

#1527S-1, 1", 100/box	(777-2577)
#1527S-2, 2", 50/box	(123-8478)



Durapore™ Hypoallergenic Silk-Like Tape

SOLVENTUM, FORMERLY 3M HEALTHCARE

A silk-like cloth tape in the High-Strength Securement Class. For clinical jobs, including securing catheters/tubes and patient positioning.

- High strength for durability
- Adheres well to dry skin
- Not made with natural rubber latex and hypoallergenic for sensitive patients
- Convenient bidirectional tear for easy tailoring
- Conformable
- Versatile
- Water resistant

#1538-0, ½" Width

(777-6047)..... 24 per Box

#1538-1, 1" Width

(777-2892)..... 12 per Box

#1538-2, 2" Width, 6 per Box

(777-0941)..... 6 per Box

#1538-3, 3" Width, 4 per Box

(777-9487)..... 12 per Box

Single-Use Rolls, 1½ yd

#1538S-1, 1" Width

(777-0731)..... 100 per Box


#1538S-2, 2" Width

(777-2160)..... 50 per Box

Look SURGICAL SUTURE

Now with 302 Stainless Steel Needles.

A Legacy of Quality, Strength, and Durability



NEW LOOK,
BETTER PERFORMANCE.

corzamedical



TOPICAL SKIN ADHESIVE - OCTYL

High-Strength 2-Octyl Cyanoacrylate Formulation

Advantages

Rapid set time

Flexibility without a plasticizer


Ergonomic ampule designed for ease of handling and targeted application

Provides an effective microbial barrier and water-resistant protection

High tensile strength

Unique formulation will not solidify in ampule, allowing all 0.5 mL to be used

Gravity fed applicator delivery system—no need to squeeze ampule



Additional Information

Polymerization starts when adhesive comes in contact with tissue or fluid.

All 0.5 mL of adhesive can be used on the patient, as adhesive will not solidify in the ampule.

Topical Skin Adhesive will remain attached to the skin surface until natural cellular regeneration occurs (approx. 5–7 days). If needed, any standard adhesive remover product may be used to release the adhesive from the skin. Peel off the film, do not pull the skin apart.

Henry Schein® Topical Skin Adhesive — Octyl

0.5 mL ampules, 12/box.....(570-3554)

Order today at **HENRYSCHEIN.COM/MEDICAL**



EXCLUSIVE BRANDS
& PRODUCTS



QUALITY YOU
CAN RELY ON



SMART
SAVINGS



TOPICAL SKIN ADHESIVE - BUTYL

Strong Bonding and Reliable Protection

- Sterile, single-use plastic ampules of liquid skin adhesive consisting of n-butyl-2-cyanoacrylate
- Antibacterial barrier protects the wound
- Applicator tips provide precise, controlled application
- Will not clog or dry when opened
- Blue liquid color option available to enable easy control over quantity applied




EASY-TO-USE APPLICATION



Open the package and take out the tip.



Open the ampule by twisting off the ribbed tip.



Apply the tip to opened ampule.



Apply the glue in a thin layer to the wound edges and hold in position for 60 seconds.

Henry Schein® Topical Skin Adhesive

Flexible Ampule with Tip, Blue 0.5 mL, Ea.

(570-0351)

Ampule without Tip, Clear 0.5 mL, 10/box

(112-6106)

Flexible Ampule with Tip, Blue 0.5 mL, 10/box

(112-7158)

Ampule without Tip, Blue 0.5 mL, 10/box

(112-6111)

Order today at [HENRYSCHEIN.COM/MEDICAL](https://www.henryschein.com/medical)





SUPERIOR POINT

Needles are honed with a unique multiedged point to provide greater ease of penetration.

SUPERIOR ATTACHMENT

Needles are drilled for a seamless needle-suture juncture, providing a smoother transition with less tissue disruption.

SUPERIOR PACKAGING


Unique double-sterile barrier packaging ends costly waste.



ABSORBABLE SUTURES 12/Box

MATERIAL	NEEDLE	LENGTH	SIZE				
			6-0	5-0	4-0	3-0	2-0
Plain Gut	C-9 24mm 1/2 Circle	18"/45cm				612G* (100-6120)	
	C-3 13mm 3/8 Circle	18"/45cm		686G* (101-3192)			
	C-6 18mm 3/8 Circle	18"/45cm			(100-7771)		
	C-6 18mm 3/8 Circle	27"/70cm			H821H* (101-0668)	H822H* (101-2976)	
	DSCM-13 13mm 3/8 circle	18"/45cm	1916G* (900-7485)	1915G* (900-7471)			
Chromic Gut	C-3 13mm 3/8 Circle	18"/45cm		687G* (101-8824)			
	C-6 18mm 3/8 Circle	18"/45cm			1637G* (900-7637)	1638H* (900-7636)	
	C-6 18mm 3/8 Circle	27"/70cm			635H* (101-2787)	636H* (101-1582)	
	HR-26 26mm 1/2 Circle	30"/76cm				G122H* (570-0466)	
	DSM-13 13mm 3/8 Circle	18"/45cm			1654G* (570-0475)		
	DSM-19 19mm 3/8 Circle	18"/45cm			1637G* (570-0476)		

* - Ethicon Equivalent



Disposable Skin Stapler



- Specially designed for surgical skin suturing
- Novel structure and application security allow for quick speed when sewing up during operation, modest histological reaction, and good wound-matching and healing
- Assists in avoiding cross infection effectively
- No pain during removal
- Scarring is minimal

Skin Stapler, 35 Wide

(570-2956)..... Ea.


Skin Stapler, 35 Regular

(570-2957).....Ea.

Skin Staple Remover

(900-7011)..... Ea.

Wound Closure Strips



- Adhesive wound closures with an exceptionally porous, spun-bond fiber backing include bilateral nylon filaments for high tensile strength
- Each size is in a color-coded package for quick and easy size identification
- Hypoallergenic: Designed for patient comfort and longer-term exposure
- Porous: Allows breathing for advanced healing and moisture control
- Not made with natural rubber latex

1/8" x 3"

(900-4317)..... 250 per Box

Contains: 5 strips per envelope, 50 envelopes per box.

1/4" x 3"

(900-4315)..... 150 per Box

Contains: 3 strips per envelope, 50 envelopes per box.

1/4" x 4"

(900-4318)..... 500 per Box

Contains: 10 strips per envelope, 50 envelopes per box.

1/4" x 1 1/2"

(900-4316)..... 300 per Box

Contains: 6 strips per envelope, 50 envelopes per box.

1/2" x 4"

(900-4314)..... 300 per Box

Contains: 6 strips per envelope, 50 envelopes per box.

1" x 4"

(900-7506)..... 100 per Box

Contains: 4 strips per envelope, 25 envelopes per box.



**Cleaning Monitors
HUFRIEDYGROUP**

Assess the efficiency of the cleaning process in washer-disinfectors and ultrasonic cleaners using multiparameter testing. Colored test soil on the strip is designed to mimic the removal of blood and bioburden from an instrument surface. Reusable stainless steel holder secures the monitor in place.

Cleaning Monitors..... 50 per Box
Specify:
Ultrasonic (124-8460)
Washer-Disinfector (124-8461)
Cleaning Monitor Holder
(124-8469)..... Ea.



**Enzymax®
HUFRIEDYGROUP**

Specifically formulated for efficient ultrasonic cleaning or presoaking of cassettes and instruments. Contains Steelgard® agents to protect instruments from harmful minerals naturally found in water that contribute to buildup, spotting, and corrosion. Can also be used as an evacuation system cleaner.

Powder
In convenient bulk canister.
(600-4353)..... 800 Gm.
Liquid
Fast-acting and concentrated.
(600-0041)..... Qt.
Liquid, 1-gal Bottle Refill
(600-0548)..... Gal.
Pump Only for Gal. Refill
(600-8938)..... Ea.
1/3 oz. Liquid Packets
(600-0379)..... 40 per Box

**Enzymax® Spray Gel
HUFRIEDYGROUP**

Dual-enzyme gel with corrosion inhibitor acts as a ready-to-use instrument precleaner. Maintains moisture to help prevent debris from drying on instruments before cleaning and helps dissolve blood, tissue, and other debris. Can also be used as a linen spot remover.

Spray Bottle
(600-2236)..... 24 oz.



**Efferzyme®
Effervescent
Cleaning Tablets
CROSSTEX**

For Enzyme Ultrasonic Cleaners
Biodegradable, nonchlorine formulation contains fast-acting protease enzymes. Enhances removal of blood and other soils that may be resistant to chemical detergents. Natural pH formula. Two tablets make 1 gal. Pleasant Spearmint scent.

- Nontoxic and biodegradable
- Made in USA
- Pleasant fragrance
- Noncorrosive

(107-6395).....52 per Box



ConFirm™ Rapid Read 20-Minute System



When frequent sterilizer testing is needed and time is of the essence, release an instrument load into service in just 20 minutes.

- An enzyme-based detection system between the biological indicator and reader provides results in just 20 minutes
- Glass-free design helps ensure staff safety
- Simple Twist & Flick activation reduces contamination potential and eliminates need for a separate activator device
- Results are read automatically and indication given via lights at each well of the incubator

For more information, contact your local Henry Schein Medical Representative or visit henryschein.com/medical.

All company and product names are trademarks of Hu-Friedy Mfg. Co., LLC, its affiliates or related companies, unless otherwise noted.
2024 Hu-Friedy Mfg. Co., LLC. All rights reserved. HFL-908/2410





Alcohol Prep Pads



Sterile

- Soft, absorbent, non-woven pads impregnated with 70% isopropyl alcohol USP
- Primarily used as a topical antiseptic; disinfects prior to administering injections

Medium, 30mm x 30mm

(570-2219)..... 200 per Box

Large, 45mm x 45mm, Folded

(570-1307).....200 per Box



Fenestrated Drapes



Sterile, with 3" Round (Tissue/Poly/Tissue)

- Allows access to the operative site while it protects surrounding area
- 3" round fenestration option is die cut and heat sealed, leaving a smooth, even edge that prevents linting
- Generously sized
- Individually packed, 300 per case (6 dispensers of 50 each)
- Not made from natural rubber latex

18" x 26"—White

(900-4685).....50 per Box

O.R. Towels



- Prewashed, 100% high-quality cotton
- Highly absorbent and virtually lint-free
- All sterile towels are CSR-wrapped and have a special corner fold that allows for aseptic presentation into a sterile field
- Not made with natural rubber latex
- 17" x 27"
- Blue

Non-Sterile

(570-1113)..... 100 per Case



Povidone-Iodine Swabsticks and Prep Pads



Sterile

Topical antiseptic used for skin preparation helps to reduce bacteria that can cause infection. Saturated with 10% USP povidone-iodine.

Single Pack, 4"L

(570-1549)..... 25 Pkg. per Box

Contains: 1 swabstick per pkg.

Single Pouches, 1½" x 1½", Folded


(570-1543)..... 1000 per Case

Contains: 1 per pouch.

Triple Pack, 4"L

(570-1548)..... 25 Pkg. per Box

Contains: 3 swabsticks per pkg.



Nonfenestrated Drapes



Polylined, Sterile (Tissue/Poly/Tissue)

- Heat sealed along 2 outer edges to ensure a pinhole-free plastic layer that guarantees maximum protection against bacterial strike-through from below
- Features polyethylene lining between 2 layers of absorbent-strength tissue, virtually eliminating risk of bacterial migration and strike-through
- All sterile fields are available in convenient dispenser boxes of 50 each

18" x 26"—White

(900-4686)..... 50 per Case



Sterile

2s

(570-1117)..... 2 per Pkg.

4s

(570-1116)..... 4 per Pkg.

6s

(570-1114)..... 6 per Pkg.

BZK Towelettes



Premoistened antiseptic towelettes used for general cleansing for patients and staff.

(112-6735)..... 100 per Box



Ear Curettes



Disposable

- Color coded
- Individually wrapped
- Non-sterile
- Loop or spoon tip available

Curettes.....50 per Box

Specify:

Ring, Pink

(900-4669)

Loop, White

(900-4670)

Spoon, Gray

(900-4671)

Large Spoon, Peach

(900-4672)



Maxima® Safety Scalpels



Sterile, single use.

- Retractable safety scalpels
- Easily extend or retract the blades while avoiding contact with the sharp blade edges
- Color coded for easy identification of blade size

Safety Scalpels..... 10 per Box

Specify:

Size 10 (570-0243)

Size 12 (570-0244)

Size 11 (570-0245)

Size 15 (570-0246)



Disposable Scalpels



Individually peel-backed for convenience. Plastic handle with millimeter scale and durable stainless steel blade.

Sterile..... 10 per Box

Specify:

#10 (112-6196)

#15 (112-6187)

#11 (112-6189)

#15S (112-7296)

#12 (112-6194)

#20 (112-6193)



Disposable Vaginal Specula



Recommended for routine gynecological exams.

- Allows examination of vaginal mucosa at a 360° angle without rotation
- Provides wide range of adjustability (angular and elevation adjustments)
- One-handed manipulation for ease of operation
- Self-holding in various opening positions
- Decreases risk of contamination
- If needed, apply lubrication for better comfort during examination
- Made of hard, clear, nontoxic, high-strength plastic resin

Small	
(102-5396).....	10 per Pkg.
Medium	
(102-5397).....	10 per Pkg.
Large	
(102-9422).....	10 per Pkg.



Disposable Vaginal Specula – Light Source Compatible



- Disposable specula made of high-strength, clear, nontoxic plastic
- Allow examination of the vaginal mucosa at 360° angle without rotation
- Provide wide range of adjustability (angular and elevation adjustments), ease of operation and manipulation with one hand
- Self-holding in various opening positions
- Single use of speculum decreases the risk of cross contamination
- Compatible with corded light source of leading brands, eliminating the need for positioning a light source independently

Specula.....	25 per Case
Specify:	
Small (112-6154)	Large (112-6156)
Medium (112-6155)	



**Monsel's Thickened Solution
PREMIER**

May be applied to surgical site following excision of diseased tissue.

8-mL Bottle with Applicator
(100-9284)..... 12 per Box



Disposable Biopsy Punch



Seamless, stainless steel, sharpened blade for safe, clean, and easy sample extraction.

- Rigid, closed-end translucent guard protects blade from damage
- Ergonomic, ribbed-grip handle provides comfort, control, and a clear view of the biopsy site
- Handle is color coded for clear identification, and blade size is clearly marked
- Free from phthalates
- Sterile
- Single use
- Not made from natural rubber latex
- Available in a wide range of colors and blade sizes

Biopsy Punch Ea.

Specify:	
White, 1.5mm	(570-3435)
Ivory, 2mm	(570-3436)
Green, 2.5mm	(570-3437)
Grey, 3mm	(570-3438)
Lemon Yellow, 3.5mm	(570-3439)
Light Blue, 4mm	(570-3440)
Pink, 5mm	(570-3441)
Light Sea Green, 6mm	(570-3442)
Purple, 8mm	(570-3443)



**Disposable Uterine Sound
INTEGRA MILTEX**

Sterile

Uterine sounds are often used prior to many procedures to enable the clinician to determine the position, depth, and tilt of the uterus. Sounding is useful in identifying any anatomical abnormalities and to confirm the internal structure of the uterus without perforating the uterus.

- Sterile, single-use
- Flexible
- Graduated in centimeters
- Individually packaged
- 9¾" (24.7 cm)

(953-3457)..... 25 per Box




**KleenSpec® 590 Vaginal Specula
WELCH ALLYN**

These specula eliminate time and expense of reprocessing and reduce the chance of cross contamination. Each is designed and manufactured to deliver durability and reliability.

- Maneuverability: new wider handle provides improved ergonomics, better balance, and is much easier to manipulate during exams
- Superior illumination: newly designed light pipe provides a homogeneous spot and eliminates glare back; all of the light is projected forward right where it is needed
- Comfortable: smooth, molded acrylic will not pinch or bind tissue for greater patient comfort during exam
- Disposable: eliminates time and expense of reprocessing and reduces chance of cross contamination
- Safe and dependable: stronger grip, and the yoke is more flexible than low-cost substitutes
- Compatibility: compatible with current illumination system (#78810) as well as new illumination system (#79910)

Specula.....Box

Specify:	
#59000, Small, 24 per Box	(499-3427)
#59001, Medium, 24 per Box	(499-7957)
#59004, Large, 18 per Box	(566-0282)



**Dermablade
ACCUTEC**

Sterile

This flexible, one-piece shave biopsy instrument has easy-to-pick-up grips that provide stability and control of excision. MicroCoat® surgical blade has a polymer coating to provide sharp, smooth, resistance-free cutting. Individually packaged.

(143-7056).....50 per Box

CaviWipes™

Easy as **1 2 3**

When it comes to keeping communities healthy,
every decision matters.
Make the easy choice with CaviWipes towelettes.



CaviWipes™ HP

Hydrogen Peroxide Wipe
for alcohol-sensitive
equipment and
fast turnarounds

Part #	Item Description
(640-0122)	Regular (6" x 6.75") 160 count canister
(142-5192)	XL (9" x 12") 65 count canister



CaviWipes™ 2.0

Low Alcohol/Quat Wipe
for everyday surfaces and
most medical devices

Part #	Item Description
(140-3081)	Regular (6" x 6.75") 160 count canister
(140-3082)	XL (9" x 12") 65 count canister



CaviWipes™ Bleach

Bleach Wipe
for high-risk areas and
outbreak control

Part #	Item Description
(129-6430)	Standard (6" x 10.5") 90 count canister



Scan QR Code to Learn More



CaviWipes™ Disinfecting Towelettes

METREX

Kills both TB in 3 minutes and HBV in 2 minutes. These disposable towelettes are presaturated with CaviCide Surface Disinfectant Cleaner and are recommended for cleaning and disinfecting all nonporous surfaces. Meet OSHA and CDC guidelines. Fragrance-free. Intermediate level.

Canister Pkg.

Large—6" x 6¾"

(640-2805).....160 per Canister

X-Large—9" x 12"

(640-8482).....65 per Canister

FlatPack™—7" x 9"

Convenient counter-space-saving packaging with patented adhesive backing for optimal placement on most smooth surfaces.

(640-2566).....45 per Pkg.

Individually Wrapped

X-Large—9" x 12"

(640-2854).....50 per Box

Wall Mount

Saves counter space. Fits both Large and X-Large canisters.

(640-2856).....Ea.



CaviCide™

METREX

Ready-to-use, low-alcohol surface disinfectant that offers a 1-minute kill claim for TB and all other product labeled organisms. Nontoxic and bleach-, phenol-, and fragrance-free. Intermediate level.

Spray Bottles

(640-0020).....8 oz.

Spray Bottle

(640-0018).....24 oz.

Refill Bottle

(640-0019).....Gal.

(640-0017).....5 Gal.

Contains: 2 x 2½ Gal.



CaviWipes™ 2.0 Disinfectant Wipes

METREX

- Effective against over 40 pathogens
- 2-minute efficacy against norovirus, influenza A and B, SARS-CoV-2, TB, HIV-1, and more
- Qualified for the EPA's Emerging Viral Pathogen claim for all virus types
- Compatible with hard, nonporous surfaces

Large: 6" x 6¾"

(140-3081).....160 per Canister

X-Large: 9" x 12"

(140-3082).....65 per Canister



CaviWipes™ HP

METREX

- Alcohol-free, hydrogen-peroxide wipes that do not require PPE and are safe for use for both clinicians and patients
- Offer a 1-minute universal contact time against over 67 pathogens
- Qualified for EPA emerging viral pathogen claims for enveloped viruses (e.g., influenza A and B, HIV-1, and HBV) as well as small/large nonenveloped viruses (e.g., norovirus and rhinovirus types 1A and 37)

6" x 6¾"

(640-0122).....160 per Canister

9" x 12"

(142-5192).....65 per Canister



CaviWipes™ 1 Disinfecting Towelettes

METREX

Non-woven, disposable wipe that offers a 1-minute kill time for viruses, bacteria, and fungi, including TB, MRSA, VRE, influenza A, HIV-1, H1N1, and HBV. This low-alcohol, multipurpose wipe can be used on hard, nonporous surfaces. Bleach-, phenol-, and fragrance-free. Intermediate level.

Canisters

Large—6" x 6¾"

(640-0024).....160 per Canister

X-Large—9" x 12"

(640-0023).....65 per Canister

FlatPack

7" x 9" convenient counterspace-saving packaging with patented adhesive backing for optimal placement on most smooth surfaces.

(640-0026).....45 per Pkg.

Individually Wrapped

X-Large—9" x 12"

(640-0027).....50 per Box



CaviCide™

METREX

Convenient, ready-to-use, intermediate-level surface disinfectant. Safe for use on nonporous surfaces. When used as directed, it will also effectively clean and decontaminate critical and semicritical instrumentation. Broad spectrum of kill and nontoxic. Intermediate level.

Spray Bottle

(640-2261).....24 oz.

Refill

(640-0012).....Gal.

5 Gallon Refill

(640-5006).....5 Gal.

Contains: 2 x 2½ Gal.



Clorox® Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipes

HEALTHLINK

One-step cleaning and disinfecting wipes. Patented, activated hydrogen-peroxide formula is EPA registered to kill 13 antibiotic-resistant organisms. Intermediate level.

- 30–60 second contact time for most bacteria and viruses (see product label for specific kill claims)
- Kill claims for 32 different microorganisms, including germs of high concern such as 9 antibiotic-resistant bacterial strains (including MRSA and VRE), norovirus, rhinovirus, rotavirus, RSV, influenza A, and E. coli
- High-quality tear-resistant wipes for greater durability
- No harsh chemical odor or fumes
- Noncorrosive, nonflammable

Regular—6¾" x 5¾" 155 per Canister
(437-0023)

Large—6¾" x 9" 95 per Canister
(437-0022)



Clorox® Healthcare® Bleach Germicidal Wipes

HEALTHLINK

Fast-acting, ready-to-use bleach disinfectant wipes. EPA registered to kill a total of 58 microorganisms, including C. difficile spores and TB, in 3 min. or less. Proven efficacy against Influenza A and B viruses, measles virus, and norovirus. Engineered to leave less visible residue on surfaces. Intermediate level.

- Destroys tough bacteria, mold, and mildew at their source, eliminating the potential for transmission of various microorganisms, germs, viruses, and allergens
- Specifically formulated for cleaning hard, nonporous surfaces
- Offer expanded kill claims and a shorter contact kill time (30 sec. on most bacteria, 1 min. on most viruses)
- Complies with OSHA's Bloodborne Pathogen Standard and AOAC testing requirements for hospital disinfectant

8¾" x 9" 70 per Canister
(250-9245)

6" x 5" Ea.
(891-0015)

12" x 12", 110-Count Canister Ea.
(113-6270)

Refill for #30358, 110/pkg 2 per Case
(118-9223)



Clorox® Commercial Solutions Clorox® Disinfecting Spray

HEALTHLINK

Fresh scented spray that cleans, sanitizes, and disinfects nonporous surfaces throughout the dental office, including metal surfaces. Simply spray and allow the surface to air dry. Meets AOAC efficacy testing requirements. Bleach-free, phenol-free, EPA-registered disinfectant. Low level. Fresh scent.

- Bleach-free, phenol-free, broad-spectrum disinfectant and deodorizer that meets AOAC Germicidal Spray Product Test efficacy standards for hospital disinfectants
- Kills 99.9% of bacteria and viruses including HIV, TB, hepatitis A, poliovirus, E. coli, influenza, and more
- Fresh scent

19-oz Can Ea.
(250-1203)



Clorox® Healthcare Disinfecting Wipes

HEALTHLINK

Fresh-scented wipes that clean, sanitize, and disinfect nonporous surfaces throughout the dental office, including metal surfaces. Meets AOAC efficacy testing requirements. Bleach-free, phenol-free, and EPA-registered disinfectant. Low level.

Fresh Scent 35 per Canister
(250-2369)

Lemon Scent 35 per Canister
(142-1251)



Clorox® Healthcare® Bleach Germicidal Cleaner

HEALTHLINK

Ready-to-use bleach cleaner/disinfectant with broad-spectrum efficacy. EPA registered to kill a total of 50 microorganisms, including C. difficile spores, in 3 minutes or less. Meets CDC and OSHA guidelines for 1:10 bleach solutions for disinfecting health-care surfaces. Intermediate level.

Refill 32 oz.
(120-8623)



MaxiSpray® Plus



Surface Disinfectant Cleaner
Broad-spectrum, tuberculocidal, disinfectant/
decontaminant cleaner. Effective against viruses
(hydrophilic and lipophilic), bacteria (including
VRE), and fungi. Kills MRSA and hepatitis B and
C in 2 minutes and TB in 3 minutes. May be
used as an immersion solution for noncritical
instruments and as a surface disinfectant.
Can be used on plastics, stainless steel,
chrome, glass, Plexiglas®, and painted surfaces.
Intermediate level.

- Easy, ready-to-use formula: no mixing, measuring,
or rinsing required
- Compatible with most medical/dental
device material
- Safe to use in neonatal and infant-care areas
- Contains no phenols
- Nontoxic and nonirritating
- Satisfies OSHA's Bloodborne Pathogens Standard

Pump Spray
(101-4114).....24 oz.
Refill
(101-0285).....Gal.



PURELL® Advanced Aloe Instant
Hand Sanitizer
GOJO

Enhanced with 4 different skin-conditioning
agents that help maintain skin health and
improve overall skin feel. Kills 99.99% of the
most common germs on hands, not just in the
lab. Works in as little as 15 seconds without
water or towels. 70% ethyl alcohol.
Pump Bottle
(590-0066).....12 oz.
1000-mL NXT Space Saver Refills
(600-4239).....8 per Case



MAXIWIPE™ Germicidal Wipes



Disinfectant wipe that kills the COVID-19 virus in
1 minute on hard, nonporous surfaces. Effective
against 16 pathogens, including COVID-1, MRSA,
and H1N1. Provides fast disinfection, is proven
effective for use in health care settings, and
is compatible with a broad range of surfaces,
including stainless steel, plastic, Formica®,
and glass.

- Strong, real-cloth wipes presaturated with the
correct amount of solution
- No mixing, no measuring—just pull, wipe, and
throw away
- Eliminates hazards of airborne sprays
- Unique deep-well lid seals securely to prevent
moisture loss
- EPA registered

Large—6" x 6¾"
(106-6794).....160 per Canister
X-Large—7½" x 15"
(112-5305).....65 per Canister



PURELL® Advanced Instant
Hand Sanitizer
GOJO

Refreshing sanitizer that kills 99.9% of germs
on the hands without water or towels. Leaves
hands feeling soft. Contains 62% ethyl alcohol.
2-oz. Personal Squeeze Bottles with Flip Tops
(106-7973).....24 per Case
Squeeze Bottle
(313-1626).....4 oz.
Pump Bottle
(590-6751).....8 oz.



MAXIWIPE™ Alcohol-Free Germicidal
Disposable Wipes



Formulated with care: Alcohol- and fragrance-
free wipes for staff and patients with
respiratory sensitivities. Efficient, one-step
cleaner and disinfectant. Effective bactericidal,
tuberculocidal, and virucidal, including SARS-
CoV-2 (the virus that causes COVID-19), in 3
minutes. EPA Registration No. 9480-9, now
qualified for List Q®. Intermediate level.

Large—6" x 6¾"
(570-3056).....160 per Canister
X-Large—7½" x 15"
(570-3057).....75 per Canister

*MAXIWIPE® Alcohol-Free has demonstrated effectiveness against
emerging viral pathogens, such as the monkeypox virus, on hard,
nonporous surfaces.



PURELL® Advanced Hand Sanitizer
GOJO

Gentle and Free Foam
• Kills 99.99% of most common germs that may
cause illness in a health-care setting
• Advanced, fragrance-free hand sanitizer made
with naturally renewable ethanol
• Dye- and fragrance-free
• Clinically proven to maintain skin health
• Compatible with nitrile, vinyl, and latex gloves
• Hypoallergenic
#7751-02, 1200-mL Refill for Purell® ES8 Touch-Free
Hand Sanitizer Dispenser
(133-3627).....Ea.



Hand Sanitizer Gel



Waterless, alcohol-based gel hand sanitizer. Ideal for use where soap and water are not readily available. Rinse-free: no water required. Kills 99.9% of most common disease-causing germs spread by hands. Soothing emollients and a pleasant, fresh, clean scent. Leaves hands soft and clean. 66.5% ethyl alcohol. Made in the USA.

- Contains alcohol (active ingredient) to kill germs
- Soothing emollients leave hands feeling soft and clean
- Pleasant citrus fragrance

Pump Bottle
(900-4443)..... 8 oz.
(900-4442)..... 16 oz.



Maxicens™ Antiseptic Skin Cleanser



4.0% chlorhexidine gluconate formulation. Fast acting in just 15 seconds, while inhibiting bacterial regrowth for up to 6 hours. Skin emollients protect against irritation, keeping hands soft and soothed. Residue-free, nonsticky formulation with no added dyes. All-in-one solution is tough on bacteria, yet gentle on the skin. Pleasant floral scent.

Bottle
(570-3268)..... Qt.
Refill
(570-3269)..... Gal.



Foaming Antibacterial Hand Soap



- Provides a rich, dense lather
- Rinses easily using minimal water
- Leaves hands feeling soft and smooth
- Pleasant Spring Water Clean fragrance
- Contains 0.6% chloroxylenol (PCMX) for antimicrobial action, and aloe vera to add moisture with each use

Pump Bottle
(900-4441).....18.6 oz.



Quick Wipes



- Soft, strong towelettes: allow easy and gentle cleansing of hands and face
- Effective and gentle formulation that contains aloe: formulated for fast drying and gentleness
- Convenient packaging: towelettes easily dispensed through orifice; each towelette is saturated with cleansing liquid

Pop-up Container
(104-4176).....150 per Cont.

Maxicens® Instant Hand Sanitizer and Dispensers



- Kills 99.9% of germs in 15 seconds without the use of soap and water
- Contains aloe vera to add moisture with every use
- Henry Schein's hand sanitizer helps you meet OSHA, APIC, and CDC hand-washing recommendations
- Citrus scent

800-mL Bag for Manual Dispenser
(112-5988)..... Ea.
800 mL Bag for Automatic Dispenser
(112-5991)..... Ea.



Antibacterial Liquid Hand Soap



- Provides a luxurious feel
- Leaves hands feeling soft and smooth
- Mild, gentle, and moisturizing
- Pleasant Cucumber and Aloe fragrance
- Active ingredient is 0.2% triclosan

Pump Bottle
(900-4439)..... 8 oz.
(900-4438)..... 16 oz.
Refill Bottle
(900-4440)..... Gal.



Maxicens® Antimicrobial Hand Wash



- Contains 0.6% PCMX for antimicrobial action and aloe vera to add moisture with each use
- Aids in reducing cross contamination
- Helps meet APIC and OSHA recommendations for hand washing

800-mL Bag for Manual Dispenser
(112-5987)..... Ea.
800-mL Bag for Automatic Dispenser
(112-5990)..... Ea.
Manual Dispenser
(112-5983)..... Ea.
Automatic Dispenser
(112-5985)..... Ea.
Stand for Dispenser
(112-5989)..... Ea.



Manual Dispenser
(112-5984)..... Ea.
Automatic Dispenser
(112-5986)..... Ea.
Stand for Dispenser
(112-5989)..... Ea.

BEYOND PROTECTION

LOW DERMATITIS
POTENTIAL

FENTANYL
RESISTANCE

CHEMOTHERAPY
DRUG TESTED

ACCELERATOR
FREE



Transcend Nitrile Gloves 300/box

High tensile strength for improved barrier integrity and tear resistance
XS (565-0048) M (565-0050) *XL (565-0052)
S (565-0049) L (565-0051) *250/box



Carbon Air Nitrile Gloves 300/box

Distinct black color to minimize visible stains during use
XS (137-2693) M (137-2695) *XL (137-2697)
S (137-2694) L (137-2696) *250/box



Revo 200 Nitrile Gloves 200/box

Comfort formulation for unprecedented softness
XS (353-0013) M (353-0015) XL (353-0017)
S (353-0014) L (353-0016)



Xlim Nitrile Gloves 100/box

Engineered for optimal grip in wet conditions
XS (108-3835) M (108-4480) *XL (108-3559)
S (108-7932) L (108-9150) *250/box



© CRANBERRY.USA f CRANBERRYUSA ® CRANBERRYUSA.COM



CRITERION™ Nitrile N100

CRITERION

Durable nitrile construction for reliable protection and tactility. Chlorinated for reduced friction and easy donning. Textured fingers provide an enhanced grip. Chemotherapy tested to the ASTM D6978 standard. Beaded cuff. 3.15-mil thickness at the palm. Blue.

Gloves..... 100 per Box
Specify:

XS	(111-8533)
S	(111-8535)
M	(111-8536)
L	(111-8537)
XL**	(111-8539)

**90/Box



CRITERION™ Nitrile N200

CRITERION

Nitrile Powder-Free Exam

Slightly thinner than the N100, the durable nitrile construction provides reliable protection and tactility. Chlorinated for reduced friction and easy donning. Textured fingers provide an enhanced grip. Chemotherapy tested to the ASTM D6978 standard. Beaded cuff. 2.36-mil thickness at the palm. Available in both Periwinkle Blue and Sky Blue.

Periwinkle Blue..... 200 per Box

Specify:	
XS	(900-7437)
S	(900-7438)
M	(900-7439)
L	(900-7440)
XL**	(900-7441)

Sky Blue..... 200 per Box

Specify:	
XS	(570-2097)
S	(570-2098)
M	(570-2099)
L	(570-2100)
XL**	(570-2101)

**180/Box



CRITERION™ Nitrile Aloe + Vitamin E

CRITERION

Extra-soft stretch with an enhanced inner coating of aloe and vitamin E to soothe, moisturize, and rejuvenate hands. Fast and easy donning. Comfortable fit and feel to provide a superior range of motion. Textured fingers promote sure handling of instruments and materials. Beaded cuff. 1.97-mil thickness at the palm.

Magenta..... 200 per Box

Specify:		
XS	(112-7143)	M (112-7146)
S	(112-7144)	L (112-7147)



CRITERION™ Nitrile N300

CRITERION

Nitrile Powder-Free Exam

Ultrathin formulation for excellent tactile sensitivity and dexterity. Chlorinated for reduced friction and easy donning. Textured fingers promote sure handling of instruments and materials. Chemotherapy and fentanyl tested to the ASTM D6978 standard. Ecofriendly packaging—with 300 gloves in a box—reduces waste and saves space. Beaded cuff. 2.0-mil thickness at the palm. Black.

Gloves..... 300 per Box

Specify:	
XS	(570-3482)
S	(570-3483)
M	(570-3484)
L	(570-3485)
XL	(570-3486)



CRITERION™ Latex Chlorinated

CRITERION

Offers a comfortable fit and feel for a superior range of motion and a textured finish for an enhanced grip. Chlorinated for reduced friction and lower protein levels. Beaded cuff. 4.33-mil thickness at the palm. White.

Gloves..... 100 per Box

Specify:	
XS	(102-5354)
S	(102-8296)
M	(102-5355)
L	(102-7763)
XL**	(102-5357)

**90/Box



CRITERION™ Latex Aloe + Vitamin E

CRITERION

Latex Powder-Free Aloe Exam

Enhanced inner coating of aloe and vitamin E to soothe and moisturize. Fast and easy donning. Comfortable fit and feel provide superior range of motion. Textured finish for an enhanced grip. Beaded cuff. 3.15-mil thickness at the palm. White.

Gloves..... 100 per Box

Specify:	
XS	(112-5840)
S	(112-5841)
M	(112-5842)
L	(112-5843)
XL**	(112-5844)

**90/Box



Together, we can build a **healthier world for all.**

Practice Green is designed to empower the health care community to positively impact the future of the planet by reducing the ecological footprint and promoting sustainability—both in daily life and while at work.



A product included in the Practice Green portfolio has been carefully selected by Henry Schein as “green” based on the product information provided by the supplier partner that it is either reusable, contains recycled materials, can be recycled, is biodegradable, or is from a sustainable source.



Scan the QR code to learn more, or visit
www.henryschein.com/PracticeGreen

© 2024 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. 24DS3152

Bio Nitrile® CRANBERRY



Nitrile Powder-Free Exam

Features advanced biodegradable properties formulated to attract microbes in landfills to break down 80% of the glove in less than 1½ years. Offers strong puncture resistance while it maintains high dexterity. Compact and recycled packaging to reduce material waste, storage space, and shipping cost.

Gloves..... 300 per Box

Specify:

XS	(143-8648)	L	(143-8651)
S	(143-8649)	XL**	(143-8652)
M	(143-8650)		

**250 per Box



BeeSure® SuperSlim™ ECOBEE



Nitrile Powder-Free Exam

Textured fingertips for enhanced grip. Innovative soft and stretchy formulation delivers superior tactile sensitivity and donning ease. 60% less packaging waste. 2.0-mil thick. Cobalt Blue.

Gloves..... 300 per Box

Specify:

XS	(353-0052)	L	(353-0055)
S	(353-0053)	XL	(353-0056)
M	(353-0054)		



Microflex® Soft White Nitrile ANSELL



Nitrile Powder-Free Exam

Provides superior comfort, enhanced tactile sensitivity, and exceptional value. Textured fingertips and chlorinated. 2.6-mil thick. White.

Gloves..... 100 per Box

Specify:

XS	(565-4510)	L	(565-9481)
S	(565-8087)	XL	(565-1575)
M	(565-7431)		



BeeSure® Floral Earloop Masks Level 3 ECOBEE



- Packaging is recyclable
- Quad-Fold® design offers 15% more breathing volume
- RoundEdge® corners and cushioned earloops for comfort
- Non-latex CoolSilk® and dye-free inner layer for a smooth, nonirritating feeling

Earloop Masks..... 50 per Box

Specify:

Hibiscus Pink	(142-4158)
Plumeria Blue	(142-4159)
Orchid Lavender	(142-4163)
Fern Green	(142-4165)
Daisy Yellow	(142-4164)



BeeSure® Radiant Face Masks–Level 3 ECOBEE



Designed to maintain a noncollapsible space for maximum breathing room. Forms an ergonomic seal around nose, chin, and cheeks to minimize air leakage.

- Fluid resistant
- Not made with natural rubber latex
- Lint-free and dye-free
- Packaging is recyclable

Earloop Masks..... 40 per Box

Specify:

Aqua	(143-3260)	Orange	(143-3262)
Lime	(143-3261)	Violet	(143-3263)



BeeSure® Naturals Face Masks Level 3 ECOBEE



- Designed with unique, natural elements.
- More than 60% earth-friendly cellulose components
- Quad-Fold® design offers 15% more breathing volume
- 100% non-latex and fiberglass-free filter layers

Earloop Masks..... 50 per Box

Specify:

Canyon Red	(142-4176)	Cloud White	(142-4178)
Glacier Blue	(142-4177)	Forest Green	(142-4180)





Radiant Series Face Masks



40/box

3D Design:
Noncollapsible
space for maximum
breathing room

ASTM Level 3

Ergonomic Seal:

Unique shape to
seal and fit around
the face

**Lint Free
Inner Layer:**

Eliminate fluffy,
itchy feeling and
reduce skin
irritation

Radiant Colors:

Highlight your day
with 4 irresistible
colors



Aqua (143-3260)



Lime (143-3261)



Orange (143-3262)



Violet (143-3263)

Vibe Face Masks

ASTM Level 3 50/box



Cool Grey
(138-4278)



Pure White
(138-4277)



Luminous Blue
(353-0096)



Dazzling Green
(353-0097)





Glamorous Pink
(353-0099)

Anti-Fog
contouring nose flap
for uncompromised
visibility



Quad-Fold® Design
offers 15% more
breathing volume



DuoSeal nose and
chin adjustable strips
for contouring fit

**Soft & Stretchy
Ear-Loop** eliminates
pull & pressure on
ears

Floral Design Face Masks

ASTM Level 3 50/box

Quad-Fold design offers 15% more breathing volume, and Round Edge corners for added comfort.



Hibiscus Pink
(142-4158)



Plumeria Blue
(142-4159)



Orchid Lavender
(142-4163)



Daisy Yellow
(142-4164)



Fern Green
(142-4165)



CRITERION™ Earloop Masks Level 1

CRITERION™

Lightweight and comfortable for extended wear, these 3-ply procedural face masks offer ASTM Level 1 protection, safeguarding against low levels of moisture, spray, and splatter. Each mask features a soft, nonirritating inner and outer layer.

- Fluid resistant
- Pleated design for enhanced breathability
- Earloops and a flexible nosepiece for a secure fit
- Not made with natural rubber latex

Masks..... 50 per Box

Specify:

Black	(570-3272)	Pink	(104-3730)
Blue	(104-3809)	Yellow	(104-2849)
Lavender	(104-8600)	White	(104-6611)



CRITERION™ Earloop Masks Level 2

CRITERION™

Same great features as the Criterion™ Earloop Masks Level 1, but with ASTM Level 2 protection.

- Fluid resistant
- Pleated design for enhanced breathability
- Earloops and a flexible nosepiece for a secure fit
- Made in the USA
- Not made with natural rubber latex

Masks..... 50 per Box

Specify:

Blue	(104-7799)	Pink	(104-9200)
Lavender	(104-8071)		



CRITERION™ Earloop Masks Level 3

CRITERION™

- Light and comfortable for extended wear
- Filter media is surrounded by soft, nonirritating inner and outer facing
- Full-width, flexible nosepiece provides fit and security
- Not made with natural rubber latex

Masks..... 50 per Box

Specify:

Black	(570-3273)	Teal	(570-2042)
Blue	(570-2044)	White	(570-2043)
Lavender	(570-2047)	Yellow	(570-2046)
Pink	(570-2045)		

Blue—Made in the USA

(570-2448)..... 50 per Box



CRITERION™ Sensitive Skin Earloop Masks

CRITERION™

Lightweight, comfortable, 3-ply procedural face masks that offer ASTM Level 1 and Level 3 protection. Each mask features a hypoallergenic inner layer that will not lint, tear, or shred for users with sensitive skin.

- Fluid resistant
- Earloops and a flexible nosepiece for a secure fit
- Hypoallergenic and free of chemicals, inks, and dyes
- Made in the USA
- Not made with natural rubber latex

Level 1..... 50 per Box

Specify:

Blue	(104-8072)	Pink	(104-8073)
------	------------	------	------------

Level 3..... 50 per Box

Specify:

Blue	(570-3276)	Pink	(570-3277)
------	------------	------	------------



CRITERION™ Pediatric Masks

CRITERION™

Designed to fit smaller faces, these masks offer optimal respiratory protection for children ages 4–12. Each mask features a soft, non-woven material to provide enhanced breathability and comfort that children need.

- Fluid resistant
- Flexible nosepiece
- Pleated design for enhanced breathability
- Not made with natural rubber latex

Blue (570-2268)..... 50 per Box



CRITERION™ Extra Breathing Room Earloop Masks

CRITERION™

With a contoured design that enhances airflow and breathability, these 3-ply procedural face masks offer ASTM Level 1 and Level 3 protection. Each mask features a soft, nonirritating inner and outer layer.

- Fluid resistant
- Earloops and a flexible nosepiece for a secure fit
- Not made with natural rubber latex

Level 1..... 50 per Box

Specify:

Blue	(570-0722)	Lavender	(570-0724)
White	(570-0723)		

Level 3..... 50 per Box

Specify:

Blue	(570-3289)	Lavender	(570-3291)
White	(570-3290)		



Ultra® Sensitive Face Masks Level 3

CROSSTEX

Same excellent features as the Ultra® High Efficiency Mask. White.

- Free of chemicals, inks, and dyes
- Fluid resistant
- White hypoallergenic inner cellulose layer
- Will not lint, tear, or shred
- Not made with natural rubber latex
- Made in the USA

Earloop Masks
(107-7191)..... 50 per Box



Ultra® Face Masks with SecureFit® Technology Level 3

CROSSTEX

Same excellent features as the Ultra® High Efficiency Mask. Blue.

- Fluid resistant
- Patent pending aluminum nose and chin piece for a custom fit
- Not made with natural rubber latex
- Made in the USA

Earloop Masks
(774-0093)..... 50 per Box



Ultra® Sensitive Face Masks with SecureFit® Technology Level 3

CROSSTEX

Same excellent features as the Ultra® Sensitive Face Masks. White.

- Fluid resistant
- White hypoallergenic inner cellulose layer
- Patent pending aluminum nose and chin piece for a custom fit
- Not made with natural rubber latex
- Made in the USA

Earloop Masks
(774-0131)..... 50 per Box



Procedural Face Masks Level 2

CROSSTEX

Outer and inner spunbond layers with a white inner layer for sensitive skin.

- Fluid resistant
- Extra-long aluminum nosepiece
- Made in the USA
- Not made with natural rubber latex

Earloop Masks..... 50 per Box

Specify:

Blue (107-6443) Yellow (107-4286)
Pink (107-6070) Lavender (110-8821)

Blue
(933-0171)..... 100 per Box



Ultra® Sensitive Fog- Free™ Face Masks with SecureFit® Technology Level 3

CROSSTEX

Same excellent features as the Ultra® Sensitive Fog-Free™ Face Masks.

- Fluid resistant
- Hypoallergenic
- Nonwoven Fog-Free™ strip across inner layer of mask
- Patent pending aluminum nosepiece
- Dual chin pieces for a custom fit
- Not made with natural rubber latex
- Made in the USA

Earloop Masks
(774-0132)..... 40 per Box



Ultra® Sensitive Fog-Free™ Face Masks Level 3

CROSSTEX

Same excellent features as the Ultra® Fog-Free™ Face Mask. White.

- White, hypoallergenic inner cellulose layer
- Will not lint, tear, or shred
- Not made with natural rubber latex
- Made in the USA
- ASTM Level 3

Earloop Masks
(111-6641)..... 40 per Box



Safety Glasses



- Clear plastic lens with side shields for eye protection in the operatory or lab. Unisex.
- Special spherical lens curved in both horizontal and vertical positions, similar to precision prescription eyewear: superior field of vision; eliminates distortion; optically correct; less eye fatigue
 - Designed to fit over most prescription street-wear spectacles up to 58mm eye size: protection for those who need to wear prescription glasses
 - Integral-molded, solid side shield: total lateral protection from foreign contaminants such as spray, splash, or splatter
 - Contoured, molded brow guard: total protection from top entry

#5148002
(100-7382)..... Ea.



Maxi-Gard® Protective Eyewear



Clear, shatter-resistant polycarbonate lens with scratch-resistant coating. Meets ANSI Z87+ standards. 10-base curve allows for complete coverage. Sure-grip rubber temple insert provides extra comfort.

806 Series.....Ea.
Specify:
Black (900-4436) Pink (900-4779)
Blue (900-4488) Red (900-4781)



Maxi-Gard® Protective Eyewear



Stylish, wraparound protective eyewear for greater comfort and protection. Clear, shatter-resistant polycarbonate lens. 10-base curve for complete coverage. Includes a scratch-resistant, antifog coating and foam inserts for added comfort. Meets ANSI Z87+ standards.

Eyewear..... Ea.
Specify:
Black Frame (570-0434) Pink Frame (570-0431)
Blue Frame (570-0433) Purple Frame (570-0432)



Protective Eyewear



Lightweight, wraparound, unisex protective eyewear. 10-base curve for complete coverage. Excellent optical quality. Ideal for patient and staff use. Gray lens and temples.
(570-0435).....Ea.



Disposable Face Shields



Non-Latex
Optically clear, see-through, full-length visor provides increased protection against splash, spatter, and aerosols. Features a 1½" thick, 9" long foam headband that keeps the face shield off the face to provide maximum comfort and protection. Ultrasonically welded, 13" elastic headband for added strength and increased wearer comfort. Special fog-free coating on inside and outside of visor. Not made with natural rubber latex. Made in the USA.
7½"L x 13"W
(101-2254).....24 per Box



Maxi-Gard® Protective Eyewear



Clear, shatter-resistant polycarbonate lens with scratch-resistant coating. Meets ANSI Z87+ standards. 9-base curve allows for great wraparound protection and coverage. Stop hinge prevents temple tip from scratching lens. Protective case included. Metal.
811 Series
(900-4437).....Ea.



Compare to ValuMax Extra-Safe Jackets and Lab Coats

Maxi-Gard® SMS Jackets



Premium protective jacket available in 9 fashionable, attractive colors. Stylish hip-length jacket is made from high-quality SMS fabric for protection. Lightweight, breathable material for maximum wearer comfort. Detailed finishing includes a 3-pocket design and snap-front closure. Not made with natural rubber latex.

Sizes S–XL..... 10 per Pkg.

Specify:

Small		Large	
Black	(570-2288)	Black	(570-2327)
Blueberry	(570-2289)	Blueberry	(570-2328)
Ceil Blue	(570-2290)	Ceil Blue	(570-2329)
Light Pink	(570-2291)	Light Pink	(570-2297)
Purple	(570-2292)	Purple	(570-2278)
Raspberry	(570-2293)	Raspberry	(570-2330)
Sky Blue	(570-2294)	Sky Blue	(570-2331)
Teal	(570-2295)	Teal	(570-2332)
White	(570-2296)	White	(570-2333)
Medium		X-Large	
Black	(570-2298)	Black	(570-2318)
Blueberry	(570-2299)	Blueberry	(570-2319)
Ceil Blue	(570-2300)	Ceil Blue	(570-2320)
Light Pink	(570-2301)	Light Pink	(570-2321)
Purple	(570-2302)	Purple	(570-2322)
Raspberry	(570-2303)	Raspberry	(570-2323)
Sky Blue	(570-2304)	Sky Blue	(570-2324)
Teal	(570-2305)	Teal	(570-2325)
White	(570-2306)	White	(570-2326)

Size 2XL–Ceil Blue
(570-2307).....10 per Pkg.

Size 3XL–Ceil Blue
(570-2308)..... 10 per Pkg.



Maxi-Gard® SMS Lab Coats



Premium protective lab coat available in 9 fashionable, attractive colors. Stylish knee-length coat is made from high-quality SMS fabric for protection. Lightweight, breathable, non-latex material for optimal wearer comfort. Detailed finishing includes a 3-pocket design and snap-front closure.

Sizes S–XL..... 10 per Pkg.

Specify:

Small		Large	
Black	(570-2334)	Black	(570-2279)
Blueberry	(570-2335)	Blueberry	(570-2280)
Ceil Blue	(570-2336)	Ceil Blue	(570-2281)
Light Pink	(570-2337)	Light Pink	(570-2282)
Purple	(570-2338)	Purple	(570-2283)
Raspberry	(570-2339)	Raspberry	(570-2284)
Sky Blue	(570-2340)	Sky Blue	(570-2285)
Teal	(570-2341)	Teal	(570-2286)
White	(570-2342)	White	(570-2287)
Medium		X-Large	
Black	(570-2346)	Black	(570-2309)
Blueberry	(570-2347)	Blueberry	(570-2310)
Ceil Blue	(570-2348)	Ceil Blue	(570-2311)
Light Pink	(570-2349)	Light Pink	(570-2312)
Purple	(570-2350)	Purple	(570-2313)
Raspberry	(570-2351)	Sky Blue	(570-2315)
Sky Blue	(570-2352)	Raspberry	(570-2314)
Teal	(570-2353)	Teal	(570-2316)
White	(570-2354)	White	(570-2317)

Size 2XL–Ceil Blue
(570-2343)..... 10 per Pkg.

Size 3XL–Ceil Blue
(570-2344)..... 10 per Pkg.



Shoe Covers



Non-sterile
Nonconductive with a nonskid design. One size fits all. Not made with natural rubber latex. Blue.

One Size fits Most

(104-9452)..... 50 Pairs per Box

(104-6350)..... 150 Pairs per Case



Bouffant Caps



Full-sized cut provides total hair coverage and allows the elastic to be comfortable. Prevents hair loss from contaminating sensitive area or site. Elastic that is not made with natural rubber latex minimizes allergic reaction.

- Nonwoven polypropylene material: comfortable, lightweight, and fluid resistant
- Full-sized cut and elasticized edge provide total hair coverage and promote sterile environment
- Elastic is not made with natural rubber latex and will minimize allergic reaction
- Disposable: no laundering; no cross contamination

White

21"

(102-4151).....100 per Box

24"

(102-6614).....100 per Box

Blue

21"

(104-6032).....100 per Box

Blue

24"

(104-6058).....100 per Box

Evolution 4 Non-Reinforced Surgical Gowns

HALYARD

Sterile with Towel, SMS Fabric

Extremely comfortable with raglan sleeves and adjustable yellow neckline binding-hook and loop. Provides the right level of protection for shorter, less fluid-intensive procedures. Meets AORN standards for fluid-, flame- and lint-resistance. Blue.

- Raglan sleeves
- Adjustable neckline with yellow hook-and-loop-style neck binding



Large

(643-0809)..... Ea.

X-Large

(643-0871)..... Ea.



ORASOOOTHE®



ORASOOOTHE® "SOCKIT" GEL is an FDA cleared hydrogel oral wound dressing used to promote optimal healing and pain management after invasive dental procedures



ORASOOOTHE "Sockit" Gel Syringe (5pk): (228-6215)
ORASOOOTHE "Sockit" Gel Syringe (25pk): (228-6216)

EXTRACTIONS

An estimated **20 million teeth** are extracted every year – making extractions one of the most common dental procedures.^{1,2}

DRY SOCKETS

Up to 38% of third molar extractions can result in a dry socket, and approximately up to 5% of all extractions result in dry sockets.³ Dry sockets can cause the patient considerable pain and prolong the healing process.

DRY SOCKET PREVENTION

"Sockit" Gel actively **prevents dry sockets** by maintaining a healthy environment for the extraction site while providing immediate pain relief and antimicrobial defense.

1. Friedman JW. The prophylactic extraction of third molars: a public health hazard. Am J Public Health. 2007 Sep;97(9):1554-9. doi: 10.2105/AJPH.2006.100271. Epub 2007 Jul 31. PMID: 17666691; PMCID: PMC1963310.

2. 2023.1.9 • Colgate <https://www.colgate.com/en-us/oral-health/tooth-removal/when-surgical-extraction-of-teeth-is-necessary#>

3. Mamoun J. Dry Socket Etiology, Diagnosis, and Clinical Treatment Techniques. J Korean Assoc Oral Maxillofac Surg. 2018 Apr;44(2):52-58. doi: 10.5125/jkaoms.2018.44.2.52. Epub 2018 Apr 25. PMID: 29732309; PMCID: PMC5932271.



Orasoothe.com
from the manufacturer of Septocaine®



Plastic Hub Needles



Features a softer and less aggressive hub and sheath. Rounded hub provides gentler contact against cheek and gum. No ribs on hub where blood can coagulate. Shorter butt end. Indicating mark on hub for easy identification of bevel position. Color codes indicate needle diameter. Surgical stainless-steel tubing reduces the risk of breakage. Tamper-resistant packaging. Sterilized.

Disposable.....	100 per Box
Specify:	
30 Ga. Blue—X-Short 3/8"	(900-1893)
30 Ga. Blue—Short 1"	(900-1892)
27 Ga. Yellow—Short 1"	(900-1891)
27 Ga. Yellow—Long 1 3/8"	(900-1890)
25 Ga. Red—Short 1"	(900-1894)
25 Ga. Red—Long 1 3/8"	(900-1889)



ProTector® Needle Sheath Prop CERTOL

Use to safely uncap and recap needles using a one-handed method. Single-use, disposable, and adapts to most dental anesthetic sheaths.

(991-1785).....	100 per Box
(991-1756).....	500 per Box



LolliCaine® Unit-Dose LolliPack® CENTRIX

Single dose of 20% benzocaine topical gel. One clean swab and gel in peel-back packaging.

LolliPack.....	120 per Pkg.
Specify:	
Bubble Gum (163-0171)	Mint (163-0170)
Cherry (163-0169)	Piña Colada (163-0172)



Cetacaine® Topical Anesthetic Liquid CETYLITE

14% benzocaine, 2% tetracaine HCl, 2% butamben Rx anesthetic liquid controls pain for up to 60 minutes. Apply to mucous membrane or into periodontal pockets for dental and hygiene procedures. Includes Luer-lock cap. Rx only. Do not inject or use in the eye.

Clinical Kit (635-0327).....	Ea.
Contains: 1–24 Gm. bottle, 100 tips/syringes.	
Chairside Kit (635-0329).....	Ea.
Contains: 1–11 Gm. bottle, 20 tips/syringes.	
11 Gm. Bottle (635-0330).....	Ea.
24 Gm. Bottle (635-0328).....	Ea.



Mepivacaine HCl 3% Plain Injection, USP



1.7-mL Anesthetic Cartridges

- Validated siliconization process: ensures smooth injection
- Meets all FDA requirements: allows total flexibility in brand selection
- Color-coded stoppers and labeling: easy identification of each anesthetic formula
- Mylar® wrap: prevents breakage
- Blister packaging: no glass/glass contact

(467-3005).....	50 per Box
Contains: 5 blister packs of 10 cartridges per package.	

Approved for sale only to licensed practitioners in the USA.




Benzo-Jel™ Topical Anesthetic Gel



20% benzocaine gel that provides effective relief of discomfort from local anesthetic injections, periodontal curettage, impression taking, and intraoral radiographs. Great-tasting flavors increase patient acceptance. Quick onset with no bitter aftertaste prevents excessive salivation with minimal wash-away. Gluten-free.

1 oz. Jar.....	Ea.
Specify:	
Banana (570-0366)	Mint (570-0363)
Bubble Gum (570-0362)	Piña Colada (570-0364)
Cherry (570-0361)	Raspberry (570-0365)
Grape (570-0344)	Strawberry (570-0360)



Aspirating Syringe HUFRIEDYGROUP

Syringe.....	Ea.
Specify:	
A Type (600-3142)	C-W Type (600-2004)



Jenker™ Needle Stick Protector DENTSPLY SIRONA

Solid stainless steel receptacle for holding the needle sheath for safely recapping.

Autoclavable (117-9254).....	Ea.
---------------------------------	-----

New! SimpliCut™

Single-Patient Use Diamond Burs

Simply Convenient, Fast and Reliable

Pre-sterilized, single patient-use, and ready to go right out of the box, SimpliCut lets you skip the cleaning, sterilization, and processing steps so you can get to work right away.



A high level of sterilization.



Range of shapes to meet a variety of needs.



Color-coded grit sizes for easy identification.



Fast-cutting burs with shanks that minimize wear & tear.



Learn more at
go.kerrdental.com/SimpliCut

kerrdental.com

©2024 Kerr Corporation. All trademarks are property of Kerr Corporation.
All Rights Reserved. MKT-24-0908 Rev 0

Kerr™

Vitrebond™ Plus SOLVENTUM, FORMERLY 3M HEALTHCARE

Light-cure glass-ionomer liner/base. Packaged in the convenient Clicker™ dispensing system for a faster, easier, and more consistent mix. Fluoride release of a true glass ionomer.

Single Pack #7571
(777-5630)..... Ea.
Contains: 2–10 Gm. Clicker™ dispenser, 1 large mixing pad, 1 technique card & instructions.

Double Pack #7572
(777-6083)..... Ea.
Contains: 2–10 Gm. Clicker™ dispensers, 1 large mixing pad, 1 technique card & instructions.



Ketac™ Cem Maxicap™ SOLVENTUM, FORMERLY 3M HEALTHCARE

Capsulated glass-ionomer luting cement. Maxicaps are larger than Aplicaps and require the use of Maxicap Activator and Applier.

Refill Capsules #56021
(378-0529)..... 50 per Box

Maxicap™ Activator
(378-0842)..... Ea.

Maxicap™ Applier
(378-4890)..... Ea.



GC Fuji I® Glass Ionomer Luting Cement GC AMERICA

For final cementing of crown & bridge restorations. Bonds to tooth structure and metal, providing excellent strength and marginal integrity. For securing inlays, onlays, posts, and orthodontic bands. Low film thickness, biocompatible, and high fluoride release.

1:1 Package
(333-5075)..... Ea.
Contains: 1 of each: 35 Gm. powder, 25 Gm. liquid, spatula & mixing pad.

Powder
(333-6051)..... 35 Gm.

Liquid
(333-2339)..... 25 Gm.

GC Fuji I® Capsules
(333-3518)..... 48 per Box

Capsule Applier
(333-2891)..... Ea.

Capsule Applier IV
(333-6531)..... Ea.



3M™ RelyX™ Luting Plus SOLVENTUM, FORMERLY 3M HEALTHCARE

Resin-modified, glass-ionomer cement. Packaged in the convenient Clicker™ dispenser. Easy to dispense, mix, and load. Chemically bonds to tooth and restoration. Tack light-cure option for faster removal of excess cement.

Trial Kit #3525TK
(777-2737)..... Ea.

Contains: 1–11.5 Gm. Clicker dispenser of cement, technique guide & set of instructions.

Refill Package #3525
(777-6701)..... Ea.
Contains: 2–11.5 Gm. Clicker dispensers of cement & set of instructions.



GC Fuji PLUS® Capsules GC AMERICA

Starter Package
(333-3524)..... Ea.
Contains: 48 capsules & Capsule Applier III.

Refill Package
(333-3522)..... 48 per Box

Capsule Applier
(333-2891)..... Ea.

Capsule Applier IV
(333-6531)..... Ea.



Gluma® Desensitizer KULZER

Fast and effective results with no mixing or curing. Can be used for treating cervical hypersensitivity; pretreating, priming, and bonding direct-composite restorations; pretreating prior to placing compomer or amalgam restorations; after crown & bridge preparation; and pretreating prior to crown & bridge cementing and bonding.

Standard Package
(665-8141)..... 5 mL

Clinic Package
(665-5804)..... Ea.
Contains: 3–5 mL bottles, 150 brushes & 1 tray.

Single-Dose Applicators—0.075 mL
(665-0090)..... 40 per Pkg.



3M™ RelyX™ Luting Plus Automix SOLVENTUM, FORMERLY 3M HEALTHCARE

Resin-modified, glass-ionomer cement packaged in a handheld automix syringe. Tack light-cure option for faster removal of excess cement.

Trial Kit #3535TK
(777-0455)..... Ea.

Contains: 1–8.5 Gm. handheld automix syringe & 12 automix mix tips.

Value Package #3535
(777-0454)..... Ea.

Contains: 3–8.5 Gm. handheld automix syringes & 36 mixing tips.

Mixing Tips
(777-0456)..... 12 per Pkg.



GC FujiCEM® 2 GC AMERICA

Resin-modified, glass-ionomer luting cement. Powered by F(2) Flex Fuse Technology, which improves indirect restorations with increased strength, high fluoride release, low film thickness, and excellent marginal integrity. Use on all types of metal, resin, alumina, and zirconia-based inlays, onlays, crowns, and bridges. Radiopaque.

Starter Kit
(333-3667)..... Ea.
Contains: 1–13.3 Gm. cartridge, 20 mixing tips & 1 FujiCEM™ 2 dispenser.

Refill Package with Tips
(333-3664)..... Ea.
Contains: 2–13.3 Gm. cartridges & 42 mixing tips.

Refill Package without Tips
(333-3665)..... Ea.
Contains: 2–13.3 Gm. cartridges & 1 mixing pad.

Slide & Lock Mixing Tips
(333-3971)..... 20 per Pkg.

Intraoral Mixing Tips SL
(333-3969)..... 10 per Pkg.

Dispenser
(333-3270)..... Ea.



Compare to Ultradent Ultra Blend® Plus Lime-Lite™ Enhanced PULPDENT

Light-cure, fluoride-releasing cavity liner/base material formulated for use with adhesives. Improved moisture-friendly formula and greater adhesion to dentin. Greater release of calcium and phosphate. Offers high compressive strength and wear resistance without any sensitivity. For use with adhesives, composites, and conventional restorative materials.

Complete Kit
(384-0184)..... Ea.
Contains: 4–1.2 mL syringes & 20 applicator tips.

Syringe Refill
(384-0183)..... Ea.
Contains: 1–3 mL syringe & 20 applicator tips.



kuraray

Noritake



**FAST, STRONG, SEAMLESS
THE BOND EVOLUTION**

SIMPLIFIED UNIVERSAL ADHESIVE WITH NO WAITING TIME

CLEARFIL™ Universal Bond Quick 2

SCAN ME
TO LEARN
MORE






Initial situation



LuxaCrown provisionals

LuxaCrown:
**Maximum strength,
durable provisionals**

Rely on the chairside, long-lasting provisional crown and bridge solution with high fracture toughness (compared to ceramic) for everyday treatments and emergency cases.

Learn more about LuxaCrown:
www.dmg-america.com/luxacrown



Order Today!

Ordering Information
50ml Cartridge, 15 Automix Tips

Item#	Shade
(290-0136)	LuxaCrown - A1
(290-0137)	LuxaCrown - A2
(290-0138)	LuxaCrown - A3
(290-0139)	LuxaCrown - A3.5
(290-0140)	LuxaCrown - B1

A smile ahead together **DMG**

IT'S TIME TO GET CE SMART

**New dental hygiene CE subscription
just \$125 a year!**

Welcome to the CE Smart Subscription from ADHA, the nation's largest professional association dedicated to the dental hygiene profession. CE Smart offers unlimited access to FREE and deeply discounted on-demand courses and live webinars. Plus, the included CE Smart Tracker means you can manage ALL your CE, all in one place. That's more than smart, it's CE Smart!



Learn more + subscribe at
adha.org/cesmart

WHY SUBSCRIBE?

- Unlimited access to CE, both FREE and at the subscriber-discount rate
- CE Smart Tracker to organize and track the credits you earn (from all sources!)
- Expert-led courses and webinars on trending topics






A Family of composite restorative solutions



G-ænial Sculpt®

Light-Cured, Universal Nano Hybrid Compactable Composite

- Excellent polish and high gloss retention
- Self-polishing effect enables additional polish in mouth
- Extraordinary wear resistance due to uniform nano-filler dispersion
- Easy to shape, yet stays in place
- High radiopacity (>300% Al)
- Pairs perfectly with G-Premio BOND™, our universal, 8th generation bonding agent



G-ænial™ Universal Injectable

Universal Light-Cured Radiopaque Injectable Composite

- Injectable with total control - create anatomy as you inject with no slumping or stringing
- Superb handling with excellent adaptation to the cavity preparation, recommended for Class I, II, III, IV, and V restorations
- Strong and highly durable restorations, exhibiting outstanding flexural strength and wear-resistance
- Long-lasting esthetics - high gloss retention and resistance to discoloration



G-ænial™ BULK Injectable

Universal Light-Cured Radiopaque Injectable Composite for Bulk Fill Restorations

- Injectable composite with excellent viscosity and thixotropy for ideal bulk filling
- Can achieve highly durable restoration via high wear resistance and flexural strength (unlike competitive composite)
- Excellent polishability and gloss retention with self-polishing ability
- Guaranteed depth of cure of 4mm which allows for bulk filling in one layer



G-ænial™ A'CHORD

The advanced composite with unishade simplicity

- Achieves esthetics of 16 Classic VITA®† shades in just 5 shades
- Full-Coverage Silane Coating (FSC) and High-Performance Pulverized CERASMART® (HPC) fillers provide excellent physical and esthetic properties
- Truly universal material may be used in anterior or posterior, in restorations large or small, with no slumping
- Provides natural fluorescence and esthetics in ANY light
- Layering shades are available for more complex and challenging cases



Since 1921
100 years of Quality in Dental

www.gc.dental/america
www.gc.dental/america/training • © 2024 GC America Inc.

†VITA® is a registered trademark of VITA ZAHNFABRIK.

See the current Henry Schein Dental Supplies Catalog or visit www.henryschein.com/sm for a full product listing.

UNIVERSAL COMPOSITE

Beautifil II

Nano-hybrid universal composite with smooth, "non-sticky" consistency provides aesthetics and durability.



FLOWABLE COMPOSITE

Beautifil Flow Plus

Nano-hybrid, high strength flowable composite provides aesthetics and durability. Available in two viscosities, zero flow (F00) and low flow (F03).



ABRASIVE POLISHERS

OneGloss PS

One-step composite finisher and polisher comes individually packaged to prevent cross-contamination. Apply firm pressure to remove composite and soft pressure to polish to a "satin finish".



Super-Snap Singles

Award-winning disk polishing system for all composite restoratives comes individually packaged to prevent cross-contamination. Designed without a metal hub. Each kit contains one disk of each grit and a plastic mandrel. Available in Standard (12mm) and Mini (8mm).




SECONDARY DECAY IS THE MAIN CAUSE OF RESTORATIVE FAILURE

**GREATLY REDUCE
SECONDARY DECAY USING
SHOFU'S CLINICALLY
PROVEN BIOACTIVE
GIOMER RESTORATIVES**



Bioactive Giomer Technology¹



Conventional Composites



Beautifil Flow Plus X & Super-Snap | Paiman Lalla, DDS, AECD, FICOI, FAAIP



SINGLE USE SOLUTIONS



**Anti-Bacterial
Protection**



**Prevents Secondary
Decay²**



**Acid
Neutralization**



**Up to 30% Savings
Over Other Brands**

8-Year independent clinical study recorded 100% retention rate, no secondary caries, failures, or postoperative sensitivity.²

Buy 4 Beautifil II Tips, Get 1 (A2) at NO CHARGE!*

*Must be of like item. Mix & Match composite shades. Email this page with a copy of your Henry Schein Dental invoice to Customer-Service@shofu.com or Fax to 760-736-3276. Bonus products will be shipped to you directly from Shofu Dental Corporation, San Marcos, CA. Allow 6-8 weeks for delivery.



**BIOACTIVE.
GIOMER TECHNOLOGY.
THE BETTER CHOICE.**

Visit **shofu.com**



Shofu Dental Corporation San Marcos CA 92078

¹Mjor IA, Moorhead JE, Dahl JE. Reasons for replacement of restorations in permanent teeth in general dental practice. Int Dent J. 2000 Dec;50(6):361-6. doi: 10.1111/j.1875-595x.2000.tb00569.x. PMID: 11197194. ²Gordan VV, Mondragon E, Watson RE, et al. A clinical evaluation of a self-etching primer and a giomer restorative material: Results at eight years. JADA. 2007;138(5):621-627. Asahi University School of Dentistry.

3M Health
Care is now
Solventum

3M™ Clinpro™ Clear Fluoride Treatment

Fluoride. Simplified.

Introducing a fluoride treatment that makes it easy for patients to say “yes.”

3M™ Clinpro™ Clear Fluoride Treatment has a water-based formula that helps promote patient acceptance—which benefits your patients and your practice.



Fast and easy
L-Pop unit dose



Minimum contact time
of just 15 minutes



Effective fluoride uptake
with lower dosage



Suction-line
friendly formula

HenryScheinDental.com/3M

Ask your Henry Schein Dental Sales Consultant for complete ordering information.

3M Health Care is now Solventum. Learn how we enable better, smarter, safer healthcare to improve lives at [Solventum.com](https://www.solventum.com)



3M™ Clinpro™ Clear Fluoride Treatment

50 unit-dose packages
(0.50 mL L-Pop)

(777-2548) Flavorless
(777-2545) Mint Flavor
(777-2546) Watermelon Flavor

100 unit-dose packages
(0.50 mL L-Pop)

(777-2551) Flavorless
(777-2549) Mint Flavor
(777-2550) Watermelon Flavor

© Solventum 2024. Solventum, the S logo and Clinpro are trademarks of Solventum or its affiliates. 3M, the 3M logo are trademarks of 3M.

3M Health
Care is now
Solventum

3M™ Filtek™ Easy Match Universal Restorative

Experience shade matching made easy.

3M™ Filtek™ Easy Match Universal Restorative

The newest addition to the 3M™ Filtek™ Dental Restoratives line simplifies everyday restorations by allowing you to go from observation to shade selection – simply and intuitively. See how easy shade matching can be!

- ▶ Intuitive shade selection
- ▶ Naturally-adaptive opacity
- ▶ Reduced inventory

HenryScheinDental.com/3M

Ask your Henry Schein Dental Sales Consultant for complete ordering information.

3M Health Care is now Solventum. Learn how we enable better, smarter, safer healthcare to improve lives at [Solventum.com](https://www.solventum.com)

© Solventum 2024. Solventum, the S logo and Filtek are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M.



Capsule refill includes:
20 – 0.2g capsules

(777-2552) Bright

(777-2553) Natural

(777-2555) Warm

Capsule intro kit includes:
60 – 0.2g capsules (20 each shade)

(777-2560) Assorted

Syringe refill includes:
1 – 4g syringe

(777-2556) Bright

(777-2558) Natural

(777-2559) Warm

Syringe intro kit includes:
3 – 4g syringes (1 each shade)

(777-2562) Assorted



Filtek[™] Supreme Ultra Universal Restorative SOLVENTUM, FORMERLY 3M HEALTHCARE

Universal nanocomposite. Improved polish retention and enhanced fluorescence make it even easier to create natural-looking aesthetics. Deliver optimized aesthetic properties for both single-shade and multishade restorations. Use for anterior and posterior restorations, core buildups, splinting, and indirect restorations, including inlays/onlays and veneers.

Capsule Refills #6029
0.20 Gm.....20 per Box

Specify:

Body Shades

A1B	(777-0135)	A3B	(777-0137)
A2B	(777-0136)	A3.5B	(777-0138)



Futurabond[®] U VOCO

Dual-curing universal adhesive. Bonds without any extra primer or activator to all self-, dual- and light-cured composite materials in a self-, selective-, or total-etch procedure. Bonds to metal, zirconia, aluminum oxide, and silicate ceramic.

(999-1740).....200 per Pkg.



GC Fuji Automix LC GC AMERICA

Resin-reinforced, light-cured glass ionomer. Bioactive material with high fluoride release. High bond strength, even in the presence of saliva. Virtually no postoperative sensitivity. No etchant or adhesive required. Radiopaque.

Cavity Conditioner
(333-7058).....6 Gm.



GC Fuji II LC[®] Capsule Refills GC AMERICA

Restorative in premeasured, unit-dose capsules. Ultraconvenient. Fewer steps reduces chair time.

Capsule Refills
A2
(333-3464).....48 per Box



Omnichroma[®] TOKUYAMA DENTAL AMERICA

Single-shade universal composite. Spherical fillers, combined with Smart Chromatic Technology, create structural color matching of surrounding tooth color. Minimal wear with high compressive and flexural strength. Stain resistant to maintain cured color with glossy results. Opaque white before curing, blending naturally after curing. Blocker material used as a lingual layer in large Class III and IV cases.

PLT Refills—0.20 Gm.
(721-1022).....20 per Pkg.



Scotchbond[™] Universal Plus SOLVENTUM, FORMERLY 3M HEALTHCARE

Provides all of the benefits of the original 3M Scotchbond Universal Adhesive, but with even greater control and predictability. The first radiopaque universal adhesive, it bonds and seals caries-affected dentin with a formulation free from BPA derivatives such as BisGMA. It also increases bond strength to glass ceramics and provides full dual- and self-cure compatibility without the use of additional dual-cure activators. Features a streamlined vial with smooth surfaces and edges.

Vial Delivery System

Refill Pack #41294
(777-2497).....5 mL
Value Pack #41295
(777-2498).....Ea.
Contains: 3–5 mL vials of universal adhesive.



OptiBond[®] Solo Plus[™] KERR

15% filled, single-component, light-cure adhesive. Optimal particle size allows it to be used for direct and indirect applications, providing the ultimate in convenience, protection, and strength.

Bottle Refill
(123-5690).....5 mL
Unidose[®] Refills
(123-3029).....100 per Pkg.



TPH Spectra[®] ST DENTSPLY SIRONA

Universal composite restorative material. Filler technology, Sphere TEC[™], has advanced granulated spherical filler that works together with the resin matrix to enhance handling and aesthetics. Simplified 5-shade system covers the entire VITA[®] shade range. Fast and easy to polish.

Compules[®] Tip Refills
0.25 Gm.....20 per Box
Specify: Low Viscosity (LV)
A2 (222-2227) A3 (222-2228)
Specify: High Viscosity (HV)
A2 (222-2233) A3 (222-2237)



Filtek[™] Supreme Flowable Restorative SOLVENTUM, FORMERLY 3M HEALTHCARE

Flowable composite with unique nanotechnology that offers enhanced aesthetics and high wear resistance. Updated ergonomic syringe designed for improved comfort and easier injection of material. Use in anterior and posterior restorations, including liner/base and all Classes I–V. Flow-on-demand handling, excellent polish retention, and improved fluorescence.

Syringe Refills #6032
2 Gm.....2 per Pkg.
Specify:
A2 (777-2467) A3 (777-2468)
Contains: 2–2 Gm. syringes & 20 dispensing tips.

DIRECTA Luxator®



3+1*

MIX & MATCH

All that you need for a successful extraction!

Item Code	Description	Qty.
(201-0094)	Luxator 1mm Straight	Ea.
(201-3981)	Luxator 2mm Straight	Ea.
(888-0615)	Luxator 3mm Straight	Ea.
(888-9887)	Luxator 3mm Curved	Ea.
(201-2799)	Luxator 3mm Contra Angle	Ea.
(201-8321)	Luxator 3mm Inverted	Ea.
(201-0244)	Luxator 5mm Straight	Ea.
(201-6021)	Luxator 5mm Curved	Ea.



*Qualifying purchase must be on one invoice. Offer good for 3+1 equal or lesser value (mix & match) single instruments only. Date of purchase must be between 9/30/24 and 12/31/24. Pricing and offer valid during this period. Invoice must be received by 1/15/25. Please allow 3-6 weeks for delivery. To receive your no additional charge goods, email a copy of your Henry Schein invoice to freegoods@directadental.com.



Scan QR Code:

To unlock Henry Schein deals & learn more with a simple scan!

**ORDER FROM HENRY SCHEIN
CONTACT YOUR REP TODAY!**

Brought to you by

DIRECTA USA
A Division of DirectaDentalGroup | directausa.com

Doing Well by Doing Good



Visit www.HenrySchein.com/HSCares to find out more about the Calendar of Caring Programs.

**Thank you for helping us raise more than \$3 million
to help health happen!**



© 2024 Henry Schein, Inc. Not responsible for typographical errors.

MaxiTemp® HP



Automix temporary crown-and-bridge material for creating temporary crowns, bridges, inlays/onlays, and veneers. High-performance, rounded cartridge prevents contamination at the tip and requires less force for easier dispensing of material. Available in 1:1, 4:1, and 10:1 mix ratios.

1:1 Ratio Kit*..... Ea.

Specify:

A1	(570-3514)	A3.5	(570-3518)
A2	(570-3515)	B1	(570-3519)
A3	(570-3517)	Bleach	(570-3520)

4:1 Ratio Kit*.....Ea.

Specify:

A1	(900-1284)	A3.5	(900-1282)
A2	(900-1281)	B1	(900-4090)
A3	(900-2126)	Bleach	(900-4091)

4:1 Cartridge Refills.....2 per Pkg.

Specify:

A1	(900-1283)	A3	(900-2100)
A2	(900-1279)	A3.5	(900-1280)

10:1 Ratio Kit*.....Ea.

Specify:

A1	(570-3521)	A3.5	(570-3524)
A2	(570-3522)	B1	(570-3525)
A3	(570-3523)	Bleach	(570-3526)

Maxima® HP Mixing Tips for 4:1/10:1 Mix Ratio Cartridges 3.2mm—Light Blue
(112-6865).....50 per Pkg.

DS50 Dispensing Gun with 4:1/10:1 Slide
(112-6776).....Ea.

*Each kit contains: 1–50 mL cartridge & 6 mixing tips.



Protemp® Plus Temporization Material SOLVENTUM, FORMERLY 3M HEALTHCARE

Automix bis-acrylic temporary material. Use on crowns, bridges, inlays, onlays, veneers, and long-term temporization. High fracture resistance and toughness. Smooth, glossy surface requires no polishing or glaze. Sophisticated filler technology provides strength and aesthetics. Easy plaque removal to support gingival health and patient comfort. 4:1/10:1 ratio.

Introductory Kit #46962
(777-0090).....Ea.

Contains: 1–67 Gm. cartridge of shade A2, 16 blue Garant mixing tips, 1–16 Gm. catalyst of Relyx Temp NE, 1–36 Gm. base of Relyx Temp NE, 1 mixing pad, 1–2 Gm. syringe of Filtek Supreme Ultra A2 with tips, 1–4:1/10:1 Garant dispenser & technique cards.

Refill Package.....Ea.

Specify:

A1	(378-0081)	A3.5	(378-0281)
A2	(378-0086)	B3	(378-0091)
A3	(378-0090)	Bleach	(378-0080)

Contains: 1–67 Gm. cartridge & 16 blue Garant mixing tips.

Integrity® Temporary Crown and Bridge Material



DENTSPLY SIRONA

Twin-port cartridge and dispensing system for fast, strong, aesthetic temporaries. Provides perfect mix ratio with no cross contamination. Minimal shrinkage, great consistency. Polishes to a natural luster. With fluorescence.

Cartridge Complete Package—10:1 Ratio
Shade A2
(222-7850).....Ea.

Contains: 1–76 Gm. cartridge shade A2, 20 mixing tips & 1 dispenser with slide.

Cartridge Refill Package.....Ea.

Specify:

A1	(222-6781)	B1	(222-8079)
A2	(222-9568)	Bleach White	(222-5263)
A3	(222-0774)	C2	(222-0780)
A3.5	(222-8736)		

Contains: 1–76 Gm. cartridge & 20 mixing tips.

Cartridge Mixing Tips
(222-4922).....20 per Box

Dispenser and Slide
(222-0701).....Ea.

Handheld Automix Syringe

Refill Package.....Ea.

Specify:

A1	(222-9440)	B1	(222-3623)
A2	(222-4062)	Bleach White	(222-1479)
A3.5	(222-3846)		

Contains: 1–15 Gm. mini-syringe & 10 mini-syringe mixing tips.

Mini-Syringe Mixing Tips
(222-4206).....20 per Pkg.

Luxatemp® Ultra DMG AMERICA

Industry leading temporary provisional material that sets a new standard for flexural strength and beauty. Provides unmatched stability and long-term durability for temporary crowns, bridges, and more. Highly accurate with incredible marginal adaptation. Superior aesthetics. Can only be used with Smartmix® Tips.

Cartridge Refill—10:1 Ratio
Refill.....Ea.

Specify:

A1	(290-0042)	A3.5	(290-0045)
A2	(290-0043)	B1	(290-0046)
A3	(290-0044)	Bleach Light	(290-0047)

Contains: 1–76 Gm. cartridge & 15 tips.

Cartridge Mixing Tips—Fine
(290-6663).....25 per Pkg.

Dispensing Gun—Type 50—4:1/10:1
(290-9411).....Ea.

Smartmix®
Handheld Syringe.....Ea.

Specify:

A2	(290-0048)	Bleach Light	(290-0050)
B1	(290-0049)		

Contains: 1–15 Gm. syringe & 10 Smartmix® mixing tips.

Smartmix® Tips—Green-Blue
(290-0114).....25 per Pkg.



Luxatemp® Automix Plus DMG AMERICA

Create durable and long-lasting provisional restorations with this superior bis-acryl material. Luxatemp is highly accurate with excellent marginal adaptation. Dispensed in the quick and easy-to-use automix system, preventing cross contamination. 10:1 ratio. Can only be used with Smartmix® Tips.

Cartridge Introductory Kit
(290-2523).....Ea.

Contains: 1–76 Gm. cartridge of Luxatemp Plus (shade A2), 15 fine Luxatemp Plus mixing tips & 1 Luxatemp Plus dispensing gun.

Cartridge Refill.....Ea.

Specify:

A1	(290-8937)	B1	(290-2279)
A2	(290-1628)	Bleach (Light)	(290-7176)
A3.5	(290-9180)		

Contains: 1–76 Gm. cartridge & 15 mixing tips.

Cartridge Bulk Package—A2
(290-0514).....Ea.

Contains: 5–76 Gm. cartridges & 75 fine Luxatemp Plus mixing tips.



FluoroCore® 2+ DENTSPLY SIRONA

Dual-cure, core build-up material plus fast cement for endodontic posts. High depth of cure and strong adhesion to dentin. Radiopaque.

Syringe Kit
(222-2495).....Ea.

Contains: 4–4.75 Gm. syringes (2 of each shade; tooth colored & blue), 2–2 mL bottles (1 of each: Prime&Bond Elect® Adhesive & self-cure activator), 25 mixing tips, 10 endo intraoral tips, 25 intraoral tips, 1 dispensing well, 50 flocked applicator tips & 50 root—canal applicator tips.

Syringe Refill.....Ea.

Specify:

Tooth Colored	(222-0239)	Blue	(222-0241)
---------------	------------	------	------------

Contains: 4–4.75 Gm. syringes, 25 mixing tips, 10 endo intraoral tips & 25 intraoral tips.

Mixing Tips
(222-0242).....25 per Pkg.

Intraoral Tips
(222-0243).....25 per Pkg.

Endo Intraoral Tips
(222-0250).....30 per Pkg.

Root Canal Applicator Tips
(222-0253).....50 per Box





NeoMTA 2

AVALON BIOMED

Optimized for perforation repair, root-end filling (apexification) treatment and other microsurgies, as well as vital pulp therapy (pulp capping and pulpotomy). High versatility for use as an obturation material and sealer. Improved gel formula allows for easier mixing, saving workflow time. Nonstaining and stable when placed to ensure depth of seal and bioactive benefit. Forms hydroxyapatite on its surface to support healing of the pulp and root canal apices.

Professional Kit

(226-0056)..... Ea.

Contains: 1 of each: 2.5 Gm. powder, 4cc gel & scoop.



TF® Adaptive Files

KERR

Files..... 4 per Pkg.

Specify:

23mm		27mm	
SM1	(123-2212)	SM1	(123-2213)
SM3	(123-2218)	SM3	(123-2219)
SM2	(123-2220)	SM2	(123-2221)
ML2	(123-2222)	ML2	(123-2223)
ML1	(123-2226)	ML1	(123-2227)
ML3	(123-2229)	ML3	(123-2230)



Sodium Hypochlorite Solution

HENRY SCHEIN®

Used for endodontic irrigation, debridement, and cleansing of root canals during and after instrumentation. Contains antimicrobial properties and provides fast, effective pulp and tissue dissolution.

6%
(112-6516).....16 oz.

Sodium Hypochlorite 5%

(570-0770).....16 oz.



IRM® Ivory

DENTSPLY SIRONA

Polymer-reinforced, zinc oxide-eugenol restorative material. Provides sedative-like qualities on hypersensitive tooth pulp. Designed for intermediate restorations but can also be used as a base under nonresin cements.

Complete Package

(222-1135)..... Ea.

Contains: 38 Gm. powder & 14 mL liquid.

Powder

(222-5856)..... 38 Gm.

Liquid

(222-6509)..... 14 mL



RC-Prep®

PREMIER DENTAL

For chemo-mechanical preparation of root canals. Provides rapid cleansing and enlarging of canals due to EDTA, which softens walls and smooths calcium deposits.

Jar
(378-4499).....18 Gm.

Jar
(378-3842)..... 227 Gm.

Syringe Pkg.
(378-4342)..... Ea.

Contains: 2-9 Gm. syringes.

Pump
(378-9182).....18 Gm.

Syringe Kit-3cc
(378-7875)..... Ea.

Contains: 5-3cc syringes & 50 disposable tips.



Tempit® Moisture-Activated Temporary Filling and Sealing Material

CENTRIX

Ideal for endo-access openings and short-term temporaries. Moisture activated. Sets in 5 minutes or less and expands slightly when set. Seals out microleakage and bacterial ingress. Contains calcium sulfate and zinc oxide. Single-unit dose 0.35 Gm. syringe tips for easy placement.

0.35 Gm. Tips

(163-7705).....30 per Jar

Contains: 30-0.35 Gm. prefilled tips.



Hygenic Endo-Ice®

COLTENE

A fast and easy method to test pulp vitality. At -26.2°C, Endo-Ice is much colder than ice or ethyl chloride. Does not require refrigeration.

(547-2070).....6 oz. Can



Absorbent Points

HENRY SCHEIN®

Used in root-canal therapy. Removes moisture from root canal, introduces medicaments. Remains sterile to point of use in operatory. Sterile cell-strip package. Uniformly rolled to exacting tolerances.

Standard Sizes,

Style #504.....200 per Box

Specify:

XX-Fine (100-1372) Coarse (100-8170)

X-Fine (100-4373) X-Coarse (100-6045)

Fine (100-3941) Assorted (100-9271)

Medium (100-2563)

Contains: 200 points per box, 5 points per cell & 40 cells per box.

DENTAL TOP SELLERS



TPH Spectra® ST DENTSPLY SIRONA

Universal composite restorative material. Filler technology, Sphere TEC[®], has advanced granulated spherical filler that works together with the resin matrix to enhance handling and aesthetics. Simplified 5-shade system covers the entire VITA[®] shade range. Excellent cavity adaptation, sculptability, and stain resistance. Fast and easy to polish.

Compules[®] Tips

Introductory Kit.....Ea.

Specify:

LV—Low Viscosity (222-2278) HV—High Viscosity (222-2281)

Contains: 40–0.25 Gm. Compules[®] tips (5 of each shade: A1, A3.5, A4 & BW; 10 of each shade: A2 & A3), 1 Compules[®] tip gun, 1 shade guide & 1 i-shade label.

Compules[®] Tip Refills

0.25 Gm..... 20 per Box

Specify:

Low Viscosity (LV)

A1 (222-2226) A3.5 (222-2229)

A2 (222-2227) A4 (222-2230)

A3 (222-2228) BW (222-2231)

Specify:

High Viscosity (HV)

A1 (222-2232) A3.5 (222-2239)

A2 (222-2233) A4 (222-2240)

A3 (222-2237) BW (222-2241)

Syringe Refill..... 3 Gm.

Specify:

Low Viscosity (LV)

A1 (222-2243) A3.5 (222-2247)

A2 (222-2244) A4 (222-2249)

A3 (222-2246)

Specify:

High Viscosity (HV)

A1 (222-2251) A3.5 (222-2254)

A2 (222-2252) A4 (222-2255)

A3 (222-2253)



Articaine 4% with Epinephrine

HENRY SCHEIN[®]

Indicated for local, infiltrative, or conductive anesthesia in both simple and complex dental procedures. It is available in 2 formulations for you to decide which epinephrine concentration is right for your patient and the procedure.

Articaine 4%..... 50 per Box

Specify:

1:100,000, Gold (112-5510)

1:200,000, Silver (112-5511)



Colgate[®] Dental Floss

COLGATE ORAL PHARMACEUTICALS

Polypropylene fiber with a microcrystalline coating that slides easily between teeth without shredding, compared to an unwaxed floss. Container made with 25% recycled plastic. Mint flavor to help freshen your mouth.

Trial Size—3 Yd.

(543-0252)..... 72 per Box



GC Fuji II LC[®] Capsule Refills

GC AMERICA

Restorative in premeasured, unit-dose capsules. Ultraconvenient. Fewer steps reduces chair time.

Capsule Refills..... 48 per Box

Specify:

A1—Pedo (333-3463) B2 (333-3471)

A2 (333-3464) B3 (333-3473)

A3 (333-3466) B4 (333-3474)

A3.5 (333-3468) C2 (333-3477)

A4 (333-3470) C4—Geriatric (333-3479)

Assorted Refill Pkgs 24 per Box

Specify:

Light (333-3529) Dark (333-3530)



Lidocaine HCl 2% with Epinephrine Injection, USP

HENRY SCHEIN[®]

1.7-mL Anesthetic Cartridges

- Validated siliconization process: ensures smooth injection
- Meets all FDA requirements: allows total flexibility in brand selection
- Color-coded stoppers and labeling: easy identification of each anesthetic formula
- Mylar[®] wrap: prevents breakage
- Blister packaging: no glass/glass contact

Cartridges..... 50 per Box

Specify:

1:100,000, Red (465-1205)

1:50,000, Green (465-0015)

Contains: 5 blister packs of 10 cartridges per package.

Approved for sale only to licensed practitioners in the U.S.



Impregum[™] Penta[™] SOLVENTUM, FORMERLY 3M HEALTHCARE

Precision polyether impression material for automatic mixing and dispensing. Accurate and hydrophilic, resulting in outstanding detail, even in a moist environment.

For Use in Pentamix 3 Mixing Unit

Refill Package—Medium Body #31642

(378-4898)..... Ea.

Contains: 4 polybags
(2 of each: 300 mL base & 60 mL catalyst).



Hydrim[®] Washer Detergents SCICAN

Hydrim[®] Washer Detergent

For use in the Original Hydrim[®] instrument washer.

(138-4285)..... 2 per Case

Hydrim[®] C Washer Detergent

For use in the Hydrim[®] #C61WG4.

(138-0257)..... 8 per Case

Hydrim[®] L Washer Detergent

For use in the Hydrim[®] #L110W.

(138-0008)..... 2 per Case

DENTAL TOP SELLERS

Adhese® Universal IVOCLAR



Single-component, light-curing universal adhesive. For direct and indirect bonding procedures and compatible with all etching techniques: self-etch, selective-enamel etch, and total etch. High bond strength and virtually no postoperative sensitivity.

System Package

(945-2430).....Ea.

Contains: 1–2 mL VivaPen refill,
100 snap-on cannulas & 20 protective sleeves.

VivaPen Refill–2 mL

(945-0577).....Ea.

Bottle Refill–5 Gm.

(945-0580).....Ea.



CaviCide® KERR

Convenient, ready-to-use, intermediate-level surface disinfectant. Safe for use on nonporous surfaces. When used as directed, it will also effectively clean and decontaminate critical and semicritical instrumentation. Broad spectrum of kill and nontoxic. Intermediate level.

Spray Bottle

(640-2261).....24 oz.

Refill

(640-0012).....Gal.

5 Gallon Refill

(640-5006).....5 Gal.

Contains: 2 x 2½ Gal.

NXT Hg5® Amalgam Separator System Accessories SOLMETEX



NXT Hg5™ Collection Container with Recycle Kit

NXT-Hg5-002CR

Compatible with the NXT Hg5 Amalgam Separator.

(177-0030).....Ea.

Contains: 1 NXT Hg5 collection container, packaging
for recycling/shipping & online certificates for recycling/
compliance.



Oraqix® DENTSPLY SIRONA

Lidocaine/prilocaine periodontal gel (2.5%/2.5%). Patient-friendly, needle-free application, ideal for those requiring local anesthesia in periodontal pockets during scaling and/or root planing. Applied to the gingival margin and then to the periodontal pocket. Dispenses as a liquid and sets as a gel. Fast 30-second onset provides 20 minutes of anesthesia. Do not inject. Do not use on children.

Oraqix® Gel

(555-0626).....Ea.

Contains: 20–17 Gm. cartridges & 20 blunt-tipped
applicator tips.

Autoclavable Dispenser (Sold Separately)

(555-4520).....Ea.



DentaPure™ Cartridge HUFRIEDYGROUP

Simple installation for safe, compliant dental-unit water. Cartridges reduce patient and staff bacterial exposure by eluting nonallergic, elemental iodine into dental-unit water.

Independent Water Bottle System

Use with tap or distilled water. Quick installation onto existing dental-unit water bottle pick-up tube.

365-Day Cartridge

(778-0022).....Ea.

Municipal Water-System Cartridge

For municipally plumbed lines. Installed onto existing waterlines in the junction box or cabinet via included polyflo fittings.

365-Day Cartridge

(778-0023).....Ea.



RelyX™ Luting Plus Automix SOLVENTUM

Resin-modified, glass-ionomer cement packaged in a handheld automix syringe. Tack light-cure option for faster removal of excess cement.

Trial Kit #3535TK

(777-0455).....Ea.

Contains: 1–8.5 Gm. handheld automix syringe & 12
automix mix tips.

Value Package #3535

(777-0454).....Ea.

Contains: 3–8.5 Gm. handheld automix syringes & 36
mixing tips.



Integrity® Temporary Crown and Bridge Material DENTSPLY SIRONA

Twin-port cartridge and dispensing system for fast, strong, aesthetic temporaries. Provides perfect mix ratio with no cross contamination. Minimal shrinkage, great consistency. Polishes to a natural luster. With fluorescence.

Cartridge Complete Package–10:1 Ratio
Shade A2

(222-7850).....Ea.

Contains: 1–76 Gm. cartridge shade A2, 20 mixing tips & 1
dispenser with slide.

Cartridge Refill Package.....Ea.

Specify:

A1	(222-6781)	B1	(222-8079)
A2	(222-9568)	Bleach White	(222-5263)
A3	(222-0774)	C2	(222-0780)
A3.5	(222-8736)		

Contains: 1–76 Gm. cartridge & 20 mixing tips.

Handheld Automix Syringe

Refill Package.....Ea.

Specify:

A1	(222-9440)	B1	(222-3623)
A2	(222-4062)	Bleach White	(222-1479)
A3.5	(222-3846)		

Contains: 1–15 Gm. mini-syringe & 10 mini-syringe mixing
tips.

Reach Maxima® Momentum.

maxima®
by HENRY SCHEIN®

When your handpieces and small equipment deliver high performance, nothing can hold you back. Each Maxima piece of equipment keeps your practice operating at peak efficiency all day, with the reliability you demand. Plus, Maxima is backed by more than 90 years of expertise from Henry Schein.

Contact your Henry Schein Representative today.

HenrySchein.com/Maxima



Scan for more information



Performance Handpieces & Small Equipment



Rely On ProRepair...

10% off your next order when using the Portal.

Mention promo code **PORTAL10** for discount. Offer expires 12/31/24.



For quick support, follow the QR code or access us online at **www.prorepair.com**.

We are excited to introduce the ProRepair Online Portal, a **NEW** service available to expedite HP repairs.



- **ProRepair Online Portal**
- **Collect. Pack. Send.**
- Submit your handpiece repair online at **www.prorepaironline.com** any time.



Handpiece Repair

Our Handpiece Repair Service team has a combined experience of more than 150+ years, providing comprehensive solutions for all types of air and electric handpieces.



Trade-in Program

Our highly trained team offers trade-in programs and promotion solutions for replacement handpieces, ultrasonic inserts, and small equipment.

Instrument Sharpening Subscription Plan

Professionally sharpen any brand of dental instrument through our convenient sharpening program. Plus enjoy a complimentary Nordent® XDURA® Hygiene Kit when you sign up! *Keep your hand instruments sharp.*



Small Equipment

We're your one-stop source for all small equipment repair needs.



Ultrasonic Insert Replacement

Our exclusive program can replace your old insert to P-style or internal water flow inserts.



For more information, contact Henry Schein ProRepair:

1.800.367.3674 8am–7pm est,

email **prorepair@henryschein.com** or visit **www.prorepair.com**

©2024 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors.



ACCLEAN® Hygiene Handpiece 2

- ACCLEAN®
- Stainless steel construction for increased durability and longevity
- Improved torque for more responsive and smoother operation throughout operating speed range
- Designed to be both lightweight and powerful
- Simple-to-change color bands for easy identification and organization
- Quick-disconnect model for easy changing
- No-lube motor for simple, easy maintenance
- Autoclavable
- Compatible with all standard DPAs
- **1-year warranty**

ACCLEAN® Hygiene Handpiece 2 Fixed Back End (570-2861)..... Ea.

ACCLEAN® Swivel Model Hygiene Handpiece 2 (570-2855)..... Ea.



PROPHYwiz Hygiene Handpiece
KAVO

Ergonomic hygiene handpiece is the lightest, quietest, and smoothest-running handpiece available. Features 360° swivel that eliminates hand fatigue during polishing. Well-balanced to provide ultimate comfort for all-day use. Excellent torque with speed range up to 3,000 rpm for optimal polishing speed. No coupler needed, the PROPHYwiz adapts to any 4-hole standard tubing and can be used with both metal and disposable prophy angles. Easy to maintain and sterilize. 1-year warranty.

Single Pack (628-6937).....Ea.

3-Pack (628-8901).....3 per Pkg.



AeroPro Cordless Prophy Handpiece
PREMIER DENTAL

Lightweight, ergonomic, cordless, and pedal-less handpiece. Offers mobility, choice of prophy angles, and much-needed control. Innovative swiveling-angle design promotes hand and wrist health by providing maximum access with comfortable hand positioning. Quick-touch button offers fingertip speed control without pressure or excess heat. Quiet, yet powerful motor maintains a consistent speed during polishing and eliminates risk of over-spinning with a rheostat. ChargeSMART™ Battery Technology provides long-lasting power. Accepts most prophy angle brands. CDC-compliant with a convenient infection-prevention protocol. Patients get a great polish every time. Unsurpassed 3-year warranty on handpiece and battery.

Complete Package (378-0927).....Ea.

Contains: 1 inner module, 3 autoclavable outer sheaths, 1 charging station, 1 AC adapter, 500 AeroPro disposable barriers, 1 tray cradle & 1 replacement O-ring with tool.



Nupro Freedom® Hygiene Handpiece
DENTSPLY SIRONA

Lightweight, cordless prophy hygiene handpiece system with SmartMode™ technology. Eliminates cord drag, giving the handpiece a lightweight feel. SmartMode™ technology provides the latest advance in handpiece speed control without the need for buttons or foot pedals. Handpiece provides 25% more speed and power and all-day battery life with rapid 15-minute recharge. Quieter than either air-driven hygiene or low-speed handpieces, providing a more pleasant polishing experience for both the patient and clinician. 1-year warranty.

System Package with Foot Pedal (772-1111).....Ea.

Contains: 1 inner module, 3 autoclavable metal outer sheaths, 1 autoclavable handpiece cradle, 1 charging base, 1 wireless foot pedal, 1 power supply, 1 box of 500 Disposa-Shield® barriers & 200 Nupro Freedom® disposable prophy angles.

*Wireless Foot Pedal sold separately.



iProphy Air® Prophy Air Motor
NSK AMERICA

This lightweight, ergonomic, and quiet hygiene handpiece features a slim, short design that minimizes hand fatigue and provides good balance. Swivels 360° to provide improved maneuverability during polishing procedures. Unique 2-part construction allows for nose cone and motor separation for easy maintenance and sterilization. Compatible with all standard Doriot-style disposable or metal angles, and connects directly to any 4-hole tubing. 1-year warranty.

Complete Unit (786-0235).....Ea.

Contains: 1 of each: motor, nose cone & set of 4 color rings.

Nosecone (786-0237).....Ea.



Prophy Star® 3 Handpiece
STAR DENTAL

At 65 Gm., it is lighter than any other prophy handpiece, which reduces hand and wrist fatigue. Convenient speed adjustment: 0–5000 rpm in both forward and reverse directions. Steam sterilizable in temperatures up to 135°C. Air flow of 2.1 SCFM at maximal power of 40 psi. 360° nose cone swivel is longer, giving the user better swivel capabilities. Improved click-on and removal with no sticking. Includes 4 autoclavable, ergonomic rubber grips in assorted colors. Lube-free, no maintenance required. 2-year warranty.

Handpiece Less Swivel (808-9698).....Ea.

Contains: 1 handpiece & 4 rubber replacement grips.



Palodent® Plus Sectional Matrix System DENTSPLY SIRONA

Utilizes optimized ring, band, and wedge technology. Constantly provides tight, accurate contacts. WedgeGuard completes the system, providing adjacent tooth protection and a tight gingival seal to minimize flash and time spent on finishing. Nickel-titanium rings can be retensioned over time providing better ring retention on tooth. Autoclavable.

Matrices Refills..... 50 per Pkg.

Specify:
4.5mm (222-0381) 5.5mm (222-0382)

Matrices Refills—5.5mm

(222-0378)..... 100 per Pkg.

Ring Refills..... 2 per Pkg.

Specify:
Universal (222-0385) Narrow (222-0386)

Wedge Refills..... 100 per Pkg.

Specify:
Small (222-0391) Medium (222-0392)



Compare to Ultradent's Omni-Matrix®

Pro-Matrix® YOUNG DENTAL

Single-use, all-in-one matrix retainer and band designed to perform all dental restorations using both amalgam and composite filling materials. The sleek, narrow design allows for better visibility during procedures. Screw stays stationary during tightening for better patient comfort. Lever points to upper and lower gum when applying band. Ready to use, no assembly required, and individually wrapped.

Straight Bands..... 50 per Box

Specify:
Narrow 4.5mm (247-0124)
Wide 6mm (247-0125)

Curved Bands..... 50 per Box

Specify:
Narrow 4.5mm (140-7818)
Wide 6mm (140-7819)



Interdental Wood Wedges

HENRY SCHEIN®

Interproximal wood wedges provide tooth separation and perfect adaptation of the matrix. Offers optimal stability and compressibility. Color coded for easy size identification.

Single Sizes..... 100 per Bag

Specify:
Mini (Thin)—Green (572-2670)
Small—Pink (572-2671)
Medium—White (572-2672)
Large—Orange (572-2673)

ASSURED INJECTION ADVANTAGE

with the leading USA-made dental syringe

Unmatched Precision
with no removable tips. Accepts
1.7/1.8ml dental cartridges



Lightweight
to reduce hand stress
during injection

Made in the USA
with the highest
quality materials

Ergonomic Design
enhances control &
maneuverability

**Your Trusted
Instrument Partner**

✓ Surgical
✓ Operative

✓ Restorative
✓ Periodontal

✓ Sutures
✓ Endodontic



The ultimate choice for Dentists




Articulating Paper Forceps

 HENRY SCHEIN®

Stainless Steel—Cross Action

- Stainless steel: quality material
- Autoclavable: promotes sterility
- Special separations grip articulating paper firmly: paper will not slip

(100-5648)..... Ea.



Bite Blocks

A. TITAN INSTRUMENTS

Silicone.


Complete..... 2 per Pkg.

Specify:

Small (107-4819) Large (107-4521)

Medium (107-5858)

Soft Retraction. Non-Reflective. Superior Control.



Scan Mate Dental Soft Tissue Retractor

Scan Mate is a manual non-latex device used by a dental professional in a clinical environment as an isolation and retraction instrument. The device format is similar to that of a standard dental mirror, but with a flexible neck and protective tip to accommodate the intended tissue retraction needs. The device is intended to be used only in a clinical setting.

147-0224 | Autoclavable 12 per box
 Manufactured by Armor Dental. www.armor-dental.com  | 



OptraGate® 3D Lip & Cheek Retractors

IVOCLAR

Auxiliary aid that enables treatment field to be enlarged easily, effectively, and comfortably. Allows lips and cheeks to be retracted and ensures relative isolation. Comfortable to wear over long periods of time, as it is three-dimensionally flexible. Lips are completely covered and protected. Because it is not made with natural rubber latex, it can be used in patients with latex allergies.

Retractors

White..... 80 per Box

Specify:

Adult Regular (945-7975)

Adult Small (945-4785)

Junior (945-0074)

Assortment (945-4218)

Contains: 40 of each: Small & Regular.

ACCESS  **EDGEENDO®**

THE COMPLETE ASSORTMENT OF

now available through  HENRY SCHEIN®



Shop the full EdgeEndo portfolio

HenryScheinDental.com/EdgeEndo



EVEREDGE® 2.0 Scalers HUFRIEDYGROUP

Using an enhanced finishing process, Hu-Friedy's professional artisans are creating edges that are over 72% sharper than the next closest competitor. Less force means easier, more comfortable scaling for clinicians, which enables them to deliver the best clinical outcome every time. Made in the USA.

#9 Handle.....Ea.

Specify:	
H6/H7	(601-1862)
204S	(601-1864)
Nevi 4 Posterior	(601-1866)
Gracey 13/14	(601-1868)
Gracey 11/12	(601-1870)
Barnhart 5/6	(601-1872)
Columbia 13/14	(601-1874)
Gracey 1/2	(601-1878)
H5/33	(601-1882)
1 Nevi/H5 1 Sickle	(601-1893)




Cotton and Dressing Pliers



Premium quality. Satin finish with precise, nonslip tips and serrated handles. Rounded pin guide with stop. Made of surgical-grade stainless steel.

#317, Standard, Serrated Tips
(100-3313).....Ea.



Old Dominion University Explorer PDT, INC

Lightweight handle with improved knurling that extends to the steel working tip for better control. Color coded for easy identification.

Double End
11/12
(731-2977).....Ea.




Front-Surface Mirrors



Stainless steel. Sterilizable by autoclaving, dry heat, boiling, or cold-sterilizing solution.

Cone Socket
(104-9564).....12 per Box



Amalgam Carriers INTEGRA MILTEX

Lifetime guarantee against manufacturing defects and workmanship. Color-coded ring at the neck. Medium (White), Large (Yellow) & Jumbo (Red). Stainless steel.

Double-End Steel
Medium/Large
(953-8226).....Ea.

Safe-Lok® Tray Covers ZIRC

Clear cover for B-Lok trays with "lock-on" side grips. Durable, clear plastic. Meets OSHA standards for transporting blood-borne pathogens.

Locking
13 $\frac{3}{8}$ "W x 9 $\frac{7}{8}$ "D x $\frac{3}{4}$ "H
B-Size
(658-8835).....Ea.




IMS Signature® Series® Crown and Bridge Kit HUFRIEDYGROUP

Contains a comprehensive set of instruments for crown-and-bridge procedures, all in a Signature Series® cassette in a convenient kit. Features a compilation of clinicians' most popular procedural setups. The 16-instrument Signature Series® cassette is 100% stainless steel.

(600-2431).....Ea.

Contains: 1 of each: #6 Satin Steel® mirror handle, 3-pack of #5 front surface mouth mirrors, #5 explorer, CP 12 color-coded probe, #8 Wiland carver, #18 excavator, #113 serrated gingival cord packer, #2 Goldstein flexi-thin XTSTR composite instrument, #24 spatula, #2 college dressing plier, Miller articulating anesthetic syringe, #18 curved iris scissor, bur cushion with short lid (holds 12), 2-pack hinged instrument clips & A/W syringe tip clip.



IMS Signature Series® Collection™ Cassette HUFRIEDYGROUP

Combination of the most popular exam set-up cassette and instruments. Makes for an easy and efficient working environment. Instruments will not be lost or damaged during cleaning and sterilization. Stainless steel for long-lasting, low-maintenance instruments.

Exam Setup: 3"H x 8"W x 1 $\frac{1}{4}$ "D
(600-2467).....Ea.

Contains: 1 yellow Signature Series Exam 5 instrument cassette, 1 Satin Steel mirror handle, 1-#5 front-surface mouth mirror, 1-#5 explorer with Satin Steel handle, 1-#CP-12 color-coded probe with Satin Steel handle & 1-#U17 utility pick-up dressing pliers.



ACCLEAN®

DISPOSABLE PROPHY ANGLES


COVERED FROM ALL ANGLES.

New assortment provides ergonomic and design enhancements to improve clinical outcomes.




90° ANGLE

SNAP-ON CUP




90° ANGLE

SCREW-CUP



CONTRA ANGLE

SNAP-ON CUP



90° ANGLE

SNAP-ON CUP
EXTENDED REACH

90° ANGLE | SNAP-ON CUP

Soft, Short, 4 Web, Ribbed, with Outer Ridges	Gray	100/Box	(570-3192)
Soft, Short, 4 Web, Ribbed, with Outer Ridges	Gray	500/Box	(570-3170)
Soft, Short, 4 Web, Ribbed	Gray	100/Box	(570-3190)
Soft, Short, 4 Web, Ribbed	Gray	500/Box	(570-3169)
Firm, Short, 4 Web, Ribbed	White	100/Box	(570-3189)
Firm, Short, 4 Web, Ribbed	White	500/Box	(570-3168)
Tapered Brush	White	100/Box	(570-3191)

90° ANGLE | SCREW-CUP

Soft, Short, Webbed with Outer Ridges	Gray	100/Box	(570-3167)
Soft, Short, Webbed with Outer Ridges	Gray	500/Box	(570-3175)
Soft, Short, 4 Web, Ribbed	Gray	100/Box	(570-3166)
Soft, Short, 4 Web, Ribbed	Gray	500/Box	(570-3174)
Firm, Short, 4 Web, Ribbed	White	100/Box	(570-3165)
Firm, Short, 4 Web, Ribbed	White	500/Box	(570-3173)
Soft, Long, 4 Web, Ribbed	Gray	100/Box	(570-3164)
Soft, Long, 4 Web, Ribbed	Gray	500/Box	(570-3172)
Firm, Long, 4 Web, Ribbed	White	100/Box	(570-3163)
Firm, Long, 4 Web, Ribbed	White	500/Box	(570-3171)

ERGONOMIC CONTRA ANGLE | SNAP-ON CUP

Compact, Soft, Short, 4 Web, Ribbed, with Outer Ridges	Gray	100/Box	(570-3112)
Soft, Short, 4 Web	Gray	100/Box	(570-3101)
Firm, Short, 4 Web	White	100/Box	(570-3102)

EXTENDED REACH | ERGONOMIC LIGHTER WEIGHT 90° ANGLE | SNAP-ON CUP

Extended Reach, 90° Angle Soft, Short, 4 Web, M-Type	Gray	100/Box	(570-3110)
Extended Reach, 90° Angle Soft, Short, 4 Web, S-Type	Gray	100/Box	(570-3111)

Not made with rubber latex.

Solo Diamond®

Premium Diamond Burs
Superior Cut. Pre-Sterilized.
Now that's Efficiency!

- ✓ Sharp, high-quality diamonds
- ✓ Pre-sterilized = convenience + peace of mind
- ✓ Simplifies & eliminates unnecessary bur blocks
- ✓ Wide variety of clinical shapes & diamond grits
- ✓ 5-year pouch seal integrity validation

Infection
Control
Solutions



Available in Boxes of 25


See the current Henry Schein Dental Supplies Catalog
or visit www.henryschein.com/sm for a full product training.

ProFlare®

Articulating Disposable Propy Angle
Superior quality at an affordable price.

- ✓ Unique pivot & flare design*
– maximum surface adaptation
- ✓ Contours to each tooth – more efficient polishing
- ✓ Up to 10% greater flaring capacity¹
- ✓ Less wrist motion – better ergonomics
- ✓ Low profile – flares and flexes easily
- ✓ Non-Latex

NEW!



Innovative
articulating cup!



(378-1186) ProFlare® – Bag of 125 angles

AeroPro®

Cordless Propy Handpiece
New & Even Better!

- ✓ **New** Enhanced quick-touch button
– better ergonomics, control & protection
- ✓ **New** Sealing technology
– increased durability for safe reprocessing
- ✓ Lightweight, 100% mobile design
– promotes hand & wrist health
- ✓ Less aerosol & less splatter²
– precise finger-tip speed control



(378-0927) AeroPro® Complete Package

(3) Autoclavable Outer Sheaths, 1 each; Motor Component, Charging Station, AC Adapter, Stand, Box of 500 Disposable Barriers, and Accessories


(378-0928) AeroPro® Basic Kit

(1) Autoclavable Outer Sheath, 1 each; Motor Component, Charging Station, AC Adapter, Stand, Box of 500 Disposable Barriers, and Accessories

Enamel Pro®

Propy Paste
Removes Stains & Repairs Enamel.³

- ✓ Effective stain removal
– superior polish for brighter, whiter, more lustrous smiles⁴,⁵
- ✓ Delivers amorphous calcium phosphate (ACP)
– helps repair enamel through fluoride uptake
- ✓ Minimizes splatter & easy to rinse
– safely leaves teeth clean & smooth
- ✓ Fun swirls & great flavors
– for happy patients



Gluten Free

	Fine	Medium	Coarse
Mint**	(378-1675)	(378-7735)	(378-8187)
VanillaMint	(378-0305)	(378-0307)	(378-0308)
Bubblegum	(378-5477)	(378-5601)	(378-0038)
Cinnamon	(378-7588)	(378-5948)	(378-8688)
Strawberry	(378-2501)	(378-4974)	(378-3441)
Mixed Berry	(378-0699)	(378-0700)	(378-0701)
Grape	(378-3699)	(378-9480)	(378-8586)
RaspberryMint	(378-0432)	(378-0437)	(378-0440)
Mint (no fluoride)			(378-4426)

**Available in Mint Extra Coarse



*Patent pending. / 1. Data on file. / 2. Based on a University of California Irvine study, 2021. / 3. Dental Products Testing study on fluoride shows that Enamel Pro provides 31% greater fluoride uptake into the enamel than pastes without ACP. / 4. Data on file. / 5. Data on file. / 0824001

V I O M E

Early Cancer Detection Saves Lives

Oral Health Pro™

with CancerDetect® Oral & Throat

Elevate your practice and personalize patient treatments with comprehensive oral microbiome health insights and early cancer detection.



✓ This test detects biomarkers for oral and throat cancers at their earliest stages with a specificity of 95% and sensitivity of 90%.

✓ Oral Health Scores reveal insights into oral microbiome health and connections to deeper, systemic health concerns.

- Periodontal health
- Caries
- Gum inflammation
- Other risks of oral and systemic diseases



Learn more at
viomepro.com

(146-2957) Ea.

V VARDIS
SWITZERLAND
BIOMIMETIC DENTAL SCIENCE

No More "Watch and Wait" with Curodont™

Discover Curodont™ Repair Fluoride Plus: the unique brush-on system for the treatment of watch areas through Guided Enamel Remineralization
Maximize patient outcomes with Curodont™ Protect

Curodont™ Repair Fluoride Plus

Guided Enamel Remineralization of early caries lesions



- For all surfaces including interproximal
- Does not stain
- 1x chairside application
- 1 applicator can treat up to 3 lesions
- Non-invasive, easy, and fast application
- Suitable for all patients, including children

(687-0057)

Curodont™ Protect

Anticavity & Antisensitivity dental gel



- Chairside & at home complement to the Curodont™ Repair Fluoride Plus treatment
- Helps to prevent dental cavities
- Relieves tooth sensitivity
- Topical application on all tooth surfaces
- For all patients, including children age 6 and up


(687-2868)

www.vvardis.com

SillHa Oral Wellness System and mySaliva✓ — Leading the Way in Salivary Analytics

Add a powerful tool to your risk assessment
Combining mySaliva✓ results with your clinical assessment can help to pinpoint specific areas for improvement and support effective interventions.

SillHa is simple, immediate and offers a low cost option making it a valuable addition to oral health risk assessments to aid in comprehensive care and educating patients.




The smart choice is mySaliva✓ to manage all patients who are not high risk, as high risk patients may need a lab test.

arkray USA, inc.

Minneapolis, MN 55439 | sillhasaliva.com
©ARKRAY USA, Inc. All rights reserved.
SH-2641-v1 Rev. 08/24

Engage and Educate Patients



Measures cariogenic bacteria, acidity, buffer capacity, blood, leukocytes, protein and ammonia for a complete picture of oral health.

Item Code #	Mfr Part #	Description
(143-5120)	105522	SillHa Oral Wellness System and software
(143-5121)	105524	SillHa paper Test Strip Kit – 30 test strips, 30 vials rinse solution
(143-5122)	105526	Fine tip droppers 1.5mL


HealthFirst Helps Build Healthier Communities



Supporting Community Health With Our Medical And Dental Solutions

Automate Your Readiness Program Through OnTraq™

Simplify your facility's readiness with OnTraq, HealthFirst's smart automation platform. Ensure medications and devices are up-to-date, and save time with online reports, updates, and compliance management.



Featured OnTraq Solutions:



STAT KIT SM Series basic emergency medical kits can be a cost-effective way to be ready for a common medical emergency. Every kit includes access to OnTraq, providing insight into medication expiration dates, lot numbers, upcoming medication replenishment shipments, and more. **Item # 145-9634**



OSHA, HIPAA and Infection Control Compliance Solutions: Whether you have one or multiple locations, HealthFirst Compliance Solutions simplifies implementation, maintenance, and training for teams. Our digital solution proactively manages compliance risks, keeping practices current with federal and state guidelines. Easily onboard new employees and manage CE training and documentation via OnTraq's personalized compliance dashboard and benefit from expert-guided facility assessments and ongoing support.



11157 v001 07/2024

Scan the QR Code
for More
Information





Curodont® Repair Fluoride Plus VVARDIS

Locally applied professional mouth rinse for the in-depth anticavity treatment of white spots. Dentist or hygienist applies this drill-free, chairside product directly to the affected site to reverse beginning lesions; applied to the white spot, it guides the remineralization. It has the highest concentration of the X-Formula Curolox®, and it strengthens and remineralizes the entire enamel layer. Single-use applicators.

(687-0057)..... 10 per Pkg.



Riva Star Aqua SDI

Patented, water-based, silver-fluoride (AgF)-based solution that is an antimicrobial, a biofilm eliminator, and a desensitizer. It is a noninvasive, 2-step system used to desensitize tooth pain in a few drops. This aqueous AgF solution, with an ammonia-free base, offers immediate sensitivity relief, reduces staining, and provides a 2-year sensitivity effect. It also has no risk of tissue burn, requires no gingival barrier, and has an improved scent.

10-Patient Kit
(164-0350)..... Ea.

Contains: 20 single-dose capsules (10 of each: Riva Star Step 1 silver fluoride & Riva Star Step 2 potassium iodide, green bottle) & 20 color-coded application brushes.

2-Step Bottle Kit
(164-0351)..... Ea.

Contains: 1 of each bottle: 1.5-mL Riva Star Aqua Step 1 & 3-mL Riva Star Step 2.



REACH® Dental Floss LG H&H USA, INC.

Patient Trial Size
5 Yd.....144 per Box

Specify:
Waxed (281-0085) Mint-Waxed (281-0086)

Professional Refill
200 Yd..... Ea.

Specify:
Waxed (281-0083) Mint-Waxed (281-0082)



Colgate® Wave™ Toothbrush*

COLGATE ORAL PHARMACEUTICALS

- Ergonomically designed, contoured handle
- Soft, multilevel bristles
- Gently hugs tooth surfaces and gums to remove plaque and bacteria

Adult-Compact Head
(543-0227)..... 6 per Box

Adult-Ultra-Compact Head
(543-0286)..... 6 per Box

*Imprinting available.



Fuji TRIAGE® Capsules

GC AMERICA

Radiopaque, command-set glass ionomer with exceptional fluoride release used for pit-and-fissure sealing. Premeasured capsules provide ideal consistency and convenience. Can be used in a moist field, and is self-adhering: no bonding agent is required. Excellent treatment for toothbrush abrasion and root-surface sensitivity.

Starter Kit..... Ea.

Specify:
Pink (333-4608) White (333-3127)

Contains: 50 capsules, 1 Capsule Applier III, 1 Cavity Conditioner & 1 Fuji Coat LC.

Capsule Refills..... 50 per Box

Specify:
Pink (333-4763) White (333-7681)

Capsule Applier
(333-2891)..... Ea.

Mixing Pads #22
(333-1109)..... 5 per Box



Colgate® Junior Toothbrush*

COLGATE ORAL PHARMACEUTICALS

- Dr. Rabbit character graphics
- Small head
- For children with both baby and adult teeth

Kids (Ages 3+)—Extra Soft, Small Head
(543-0292)..... 6 per Box

*Imprinting available.



SilverSense SDF™

CENTRIX

38% silver diamine fluoride used for treating dental hypersensitivity. Delivers better patient care by alleviating hypersensitivity, discomfort, and pain. Noninvasive procedure means no drilling and anesthetic. Stops harmful bacteria and hardens and conditions softened dentin. Does not stain healthy enamel or dentin. For use in adults over the age of 21.

Standard Kit
(163-0196)..... Ea.
Contains: 1–5 mL bottle & 20 Benda Micro applicators.

Clinic Kit
(163-0197)..... Ea.
Contains: 3–5 mL bottles & 50 Benda Micro applicators.

Benda® Micro Applicators—Fine, Yellow
(163-7534)..... 576 per Box



Crest® Kids Sparkle Cavity Protection PROCTER & GAMBLE

Make brushing enjoyable and encourage the daily habit with fun graphic toothpaste. Plus, kids will be fighting cavities with clinically proven Fluoristat®, which is gentle on tooth enamel. Mint flavor.

Sparkle Fun—0.85 oz. Tubes
(109-3749)..... 72 per Case

Sparkle—4.6 oz. Tubes
(325-0308)..... 24 per Case



Colgate® My First™ Toothbrush

COLGATE ORAL PHARMACEUTICALS

- Soft bristles for baby's gums and first teeth
- Unique dosing dot to measure right amount of toothpaste
- Ergonomic handle for easy grip

Baby (Ages 0–2)—Extra Soft, Small Head
(543-0288)..... 12 per Case

Digital Sensor Sleeves



Disposable digital sensor barriers to help prevent cross contamination and protect sensor. Universal. Fit all size 1 and 2 digital sensors.

Sleeves..... 500 per Box

Specify:
Size 1 (112-5545)
Size 2 (112-5546)

Lead-Free Microfiber X-ray Apron KERR

Support breast cancer awareness with the charcoal lead-free X-ray apron. Xenolite™ alloy sheeting makes apron lightweight and environmentally safe while protecting patients from radiation exposure. Weighs up to 30% less than a conventional apron and can withstand most disinfectant liquids. Microfiber front and nonslip corduroy back prevent it from sliding off the patient. 100% recyclable.

With Collar
Adult
23"W x 26"L
Blue
(135-5891)..... Ea.



Schick® AimRight Adhesive-Positioning Systems DENTSPLY SIRONA

For optimal sensor placement to obtain the right shots easily and effectively. Adhesive backing allows position to be adjusted on the sensor to help ensure just the right image is captured, regardless of the patient's anatomy. Patient-friendly design adds no extra dimensions to the setup. Single use for optimal hygienic workflow. Rings and arms are autoclavable.

Adhesive Tab
Holders..... 50 per Pkg.
Specify:
Bite Tab—Red (811-0844)
Bitewing Tab—Blue (811-0846)
Endodontic Tab—Green (811-0849)
Anterior Periapical Tab—White (811-0843)
Posterior Periapical Tab—Yellow (811-0854)



office supplies



Trash Liners



These liners are a strong, durable, and affordable option, and fit many different shapes and sizes of trash cans. Made in the USA.

Clear

7–10 Gal., 6 mil
(570-1340)..... 1000 per Case
15–20 Gal., 8 mil
(570-1346)..... 1000 per Case
20–30 Gal., 0.39 mil
(570-1336)..... 250 per Case
20–30 Gal., 10 mil
(570-1341)..... 250 per Case
33 Gal., 13 mil
(570-1342)..... 250 per Case
40–45 Gal., 12 mil
(570-1343)..... 250 per Case

Black

33 Gal., 0.51 mil
(570-1337)..... 250 per Case
40–45 Gal., 0.51 mil
(570-1338)..... 250 per Case

60 Gal., 1.4 mil
(570-1339)..... 100 per Case

White

13 Gal., 0.8 mil
(570-1344)..... 100 per Case
20–30 Gal., 0.85 mil
(570-1345)..... 250 per Case

MetriMist—Kerr TotalCare METREX

Natural aromatic deodorizer that eliminates odors by cleaning the air, not masking the smell. May be sprayed directly in the air, on a surface, or added to laundry water to eliminate offensive odors. Nontoxic. Biodegradable and nonflammable.

(176-9931) 8 oz.



Lysol® Professional Disinfectant Sprays

Kills 99.9% of germs on hard, nonporous surfaces. Effective against over 100 microorganisms when used as directed. Eliminates odors and prevents the growth of mold and mildew. EPA registered. 19-oz aerosol spray.

Crisp Linen Scent
(902-9209).....Ea.

Fresh Scent
(902-6347)***Ea.

***This product is not available for customers in Alaska and Hawaii



STICKERS

**THESE STICKERS ARE
MADE IN THE USA,
FULLY RECYCLABLE,
BIODEGRADABLE,
AND NON-LATEX**



**ALL
NEW**



BLUE'S DENTAL
6 Assorted Designs. 2½"
(368-4077).....100 per Roll



BLUES CLUES AND YOU
6 Assorted Designs. 2½"
(368-4073).....100 per Roll



MICKEY MOUSE FUNHOUSE
6 Assorted Designs. 2½"
(368-4072).....100 per Roll



INSIDE OUT 2
6 Assorted Designs. 2½"
(368-4070).....100 per Roll



DISNEY JUNIOR'S ARIEL
6 Assorted Designs. 2½"
(368-4081).....100 per Roll



TRANSFORMERS ONE
6 Assorted Designs. 2½"
(368-4080).....100 per Roll



KUNG FU PANDA 4
6 Assorted Designs. 2½"
(368-4061).....100 per Roll



DESPICABLE ME 4
6 Assorted Designs. 2½"
(368-4062).....100 per Roll



FANTASY BARBIE
6 Assorted Designs. 2½"
(368-4074).....100 per Roll



ADOPT ME FROM ROBLOX
6 Assorted Designs. 2½"
(368-4100).....100 per Roll



BRATZ
6 Assorted Designs. 2½"
(368-4088).....100 per Roll



DISNEY PRINCESS ANIME
6 Assorted Designs. 2½"
(368-4075).....100 per Roll



BATMAN ANIME
6 Assorted Designs. 2½"
(368-4089).....100 per Roll



Dental Squish Mix

These individually wrapped dental squish toys are sure to bring a smile to your patient's face. Great for stress relief.

8 assorted designs.

(368-3668).....50 per Pkg.



Lip Rageous® Mini Lip Balms–Oralabs®

- Delivers mega-moisture and flavor
- Formulated with aloe and vitamin E
- Contains 35 Blue Raspberry, 35 Watermelon, and 30 Strawberry

Assorted

(367-1112)100 per Pkg.



Plexiglas® Hand Mirror

 **HENRY SCHEIN®**

One side plain/one side magnifying. 5" diam.

(101-0055).....Ea.



Dr. John's® Classic Fruit Sweet Originals Collection

INGREDIENTS: Isomalt, Dextrin (dietary fiber), Erythritol, Citric Acid, Natural Flavors, Fruit and Vegetable Juice for Color, Spirulina Extract (Color),Turmeric (Color), Stevia Leaf Extract. Assorted Flavors: Cherry, Grape, Blue Raspberry, Orange, Lemon, and Lime. Individually wrapped.

(368-3276).....150 per Bag



Tooth-Shaped Brushing Timers

Brushing timer lengths are approximate. 3½"

2-Minute Timer

(367-6826).....50 per Pkg.

3-Minute Timer

(367-0491).....50 per Pkg.



Sticky Hands–Sherman Specialty

5"

(778-5141) 72 per Pkg.



(367-2869)

(367-7852)

(368-2323)


(367-5241)

(367-7588)

Scatter Print Bags

- Clear plastic bags with 2-sided print
- Transparent design allows you to see the supplies inside
- Great price!
- Available in one size: 7"W x 10"H

Scatter Bags100 per Pkg.



WARNING: CHOKING HAZARD

Item contains a small ball, small part and/or marble. Not for children under three years.




The path ahead needs an oral health care system designed to work for everyone.

Together, we can stand up for those let down and left behind.

Millions of people, especially those in marginalized communities, are being left to suffer the medical, emotional, and economic consequences of poor oral health. It's a nationwide failure that we aim to address.

We're on the path to creating a system that is accessible, equitable, and integrated. We hope you join us.

Learn more at carequest.org



**Special Markets**

1.800.851.0400 8:00am – 8:00pm (et)
henryschein.com/sm

Medical

1.800.772.4346 8:00am – 8:00pm (et)
henryschein.com/medical

© 2024 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. To guarantee special pricing/offers, where applicable, please include the three-digit promo code when ordering. Note: Participating in a promotional discount program (e.g., points, discount redemptions or other special awards) is only permissible in accordance with discount program rules. By participation in such program, you agree that, to your knowledge, your practice complies with the program requirements. For HS promotions: This promotion offers a bundled discount. The cost of any "buy/get" (or "no charge") items is a discount applied to your other purchases on the same invoice. The reportable value (unit price) of the "no charge" item will be identified on the invoice. To calculate the applicable discount on the other items on the same invoice, divide the total value (unit price(s)) of any no-charge item(s) by the total value (extension price(s)) of all other items purchased on the same invoice, and apply the resulting percentage equally to each item to determine net prices. For manufacturer promotions: The prices stated in the invoice or statement may reflect a discount or be subject to a rebate. You must fully and accurately report this stated discount price, or if applicable, any net pricing, after giving effect to any rebates, to Medicare, Medicaid, Tricare and any other federal or state program upon request by any such program. Accordingly, you should retain these records. It is your responsibility to review any agreements or other documents applicable to these prices to determine if they are subject to a rebate.

HENRY SCHEIN®
Henry Schein, Inc. • 135 Duryea Road • Melville, NY 11747

PRSRT STD
U.S. POSTAGE

PAID
HENRY SCHEIN, INC.

SIMPLICITY IN 1 SHADE

+ X-TRA

- x-tra low shrinkage
- x-tra depth of cure
- x-tra high biocompatibility



Scan to
Request
Your Free
Sample



Admira Fusion x-tra

Single-shade **OMNI-CHROMATIC** restorative to simplify the procedure, save time and reduce inventory

- Only one omni-chromatic shade: no guess work and no more wasted shades
- 4 mm depth of cure – saves time
- Easy to polish, highly stain resistant with a high gloss
- Exceptional longevity with low initial shrinkage of only 1.25%
- 100% BPA-Free and no classic monomers for superior biocompatibility



(999-2208) – Syringe 1 × 3.0 Gm.
(999-2210) – Capsules 15 × 0.2 Gm.

What reaction do you want?



Scan to
Request
Your Free
Sample



VOCO Profluorid® Varnish

The Thin Transparent 5% Sodium Fluoride Varnish in a Non-messy SingleDose Delivery System

- High immediate fluoride release to relieve hypersensitivity
- Seven great tasting flavors without the unpleasant aftertaste of some other brands
- No yellow discoloration of teeth
- Adheres well to moist surfaces
- Sets quickly in seconds after contact with saliva



(999-2293) – SingleDose 200 × 0.40 mL bubble gum
(999-1367) – SingleDose 200 × 0.40 mL caramel
(999-1369) – SingleDose 200 × 0.40 mL cherry
(999-1370) – SingleDose 200 × 0.40 mL melon
(999-1368) – SingleDose 200 × 0.40 mL mint
(999-2710) – SingleDose 200 × 0.40 mL cola lime
(999-2820) – SingleDose 200 × 0.40 mL piña colada
(999-1366) – SingleDose 48 × 0.40 mL mixed
(840-9282) – SingleDose 50 × 0.25 mL melon
(999-2292) – SingleDose 50 × 0.40 mL bubble gum
(999-1362) – SingleDose 50 × 0.40 mL caramel
(999-1363) – SingleDose 50 × 0.40 mL cherry
(840-9285) – SingleDose 50 × 0.40 mL melon
(999-1364) – SingleDose 50 × 0.40 mL mint
(999-2713) – SingleDose 50 × 0.40 mL cola lime
(999-2819) – SingleDose 50 × 0.40 mL piña colada