

BACK-TO-SCHOOL

PRODUCT GUIDE 2024

Summer/Fall Edition



Because So Many
Rely on You...
You Can
Rely on Us.





AS THE BACK-TO-SCHOOL SEASON APPROACHES, ENSURE YOU'RE WELL-EQUIPPED WITH THE NECESSARY SUPPLIES FOR VACCINATIONS AND SCREENINGS.

Henry Schein offers a comprehensive selection of pediatric supplies and vaccines from trusted partners, designed to keep children healthy this fall.

We're your trusted partner in pediatric care, providing essential tools for diagnosis, and treatment, and ensuring patient comfort. Our medical supplies, equipment, and vaccines cater specifically to infants, children, and adolescents. We simplify your job, empowering you to focus on delivering exceptional healthcare. Let's make a positive impact together on the children who rely on you.



TABLE OF
■ CONTENTS:

Diagnostic Equipment.....3–6
Exam Room Equipment7–11
AEDs 12
Exam Room Disposables.....13–16
Exam Glove.....17–20
Practice Green.....18–19
PPE.....19–22
Infection Prevention 22–25
Surgical Products26–28

Exam Room Instruments.....28–30
Orthopedics 31
Wound Care32–37
Wound Closure38–40
Waste Management.....40–41
Hypodermics42–43
Pharmaceuticals45–49
Vaccines.....50–53
OTC54–55

HBC.....55–57
Varnish 58
Cold Storage..... 59–62
Allergy Testing..... 62
Point-of-Care Testing..... 63–77
Urinalysis78
Specimen Collection 79

Baxter

WELCH ALLYN AND BAXTER ARE NOW ONE TEAM

ZERO CONTACT, FAST, ACCURATE VISION SCREENINGS



WELCH ALLYN SPOT VISION SCREENER

Contact your Henry Schein representative
to learn more about performing safe vision screenings.

Hill-Rom reserves the right to make changes without notice in design, specifications and models.
The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

Baxter, Hillrom, Spot and Welch Allyn are trademarks of Baxter International Inc. or its subsidiaries.
ALL RIGHTS RESERVED. US-FLC158-230016 V1 04/23



Baxter

Welch Allyn Green Series™ 777 Integrated Wall System

#777-PM2WAS-US, With PanOptic Basic LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer
(140-1989).....ea

#777-PM3WAS-US, With PanOptic Plus LED Ophthalmoscope, MacroView Plus LED Otoscope for iExaminer, BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer
(140-2495).....ea

#777-SM2WAS, With Coaxial LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer
(140-1801).....ea

#777-SM2WAX, With Coaxial LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, and Ear Specula Dispenser
(140-1799).....ea



Baxter

Welch Allyn Connex® Spot Monitor

- Designed for the physician-office environment with an easy-to-use, vivid touch-screen display
- More accurate hypertension diagnosis without slowing down your workflow thanks to automated blood-pressure averaging
- Improved capabilities for assessing respiratory conditions such as COPD, asthma, and community-acquired pneumonia (CAP), thanks to integrated professional-grade pulse oximetry
- Increase staff efficiency and eliminate data-entry errors with wireless transfer of vital signs directly into an EMR

7100 Base Models with Wired Connectivity
#71XX-B, SureBP Non-invasive Blood Pressure
(566-0464).....ea

#71XT-B, SureBP, SureTemp
(672-0176).....ea

#71XE-B, SureBP, Braun Thermoscan PRO 6000 Ear Thermometer
(124-6103).....ea

#71WX-B, SureBP, Nonin SpO₂
(672-0177).....ea

#71WT-B, SureBP, Nonin SpO₂, SureTemp
(566-0466).....ea

#71WE-B, SureBP, Nonin SpO₂, Braun Pro6000
(124-5328).....ea

7300 Bluetooth Connectivity Physician Office Models
#73XT-B, SureBP, SureTemp
(566-0465).....ea

#73XE-B, SureBP, Braun PRO 6000
(124-8200).....ea

#73WX-B, SureBP, Nonin SpO₂
(672-0178).....ea

#73WT-B, SurePB, Nonin SpO₂, SureTemp
(672-0179).....ea

#73WE-B, SureBP, Nonin SpO₂, Braun PRO 6000
(124-4801).....ea

#73CT-B, SureBP, Covidien SpO₂, SureTemp
(672-0175).....ea

#73CE-B, SureBP, Covidien SpO₂, Braun PRO 6000
(125-8900).....ea

#73MT-B, SureBP, Masimo SpO₂, SureTemp
(566-0467).....ea

#73ME-B SureBP, Masimo SpO₂, Braun PRO 6000
(124-7063).....ea

7400 Upgradable to WiFi Connectivity Hospital Models
#74CX-B, SureBP, Covidien SpO₂
(566-0469).....ea

#74CT-B, SureBP, Covidien SpO₂ SureTemp
(566-0470).....ea

#74CE-B, SureBP, Covidien SpO₂, Braun PRO 6000
(127-5208).....ea

#74MX-B, SureBP, Masimo SpO₂
(566-0471).....ea

#74MT-B, SureBP, Masimo SpO₂, SureTemp
(566-0478).....ea

#74ME-B, SureBP, Masimo SpO₂, Braun PRO 6000
(127-6033).....ea



Baxter

Welch Allyn Spot™ Vision Screener

A handheld, portable device designed to help users quickly and easily detect vision issues on patients from 6 months of age through adult. Spot screens both eyes at once from a nonthreatening 3-ft distance. The touch-screen display allows for one-touch activation, simple management of patient data entry, and easy configuration for both vision screening and autorefraction applications. Features:

- Supports AAP screening guidelines for early detection of amblyopic risk factors
 - Automated screening provides thorough, objective, and easy-to-understand results
 - Easy to use, easy to implement, with minimal user training required
 - Adequately powered for mass screening events
 - Captures readings 97% of the time, helping to screen otherwise difficult patients
 - Wireless printing and WiFi-enabled for easy export of data
 - Manual exports available via USB port
 - Uses lights and sounds to help engage children
 - Convenient carrying case
- #VS100S-B, Vision Screener Set
(497-0004).....ea



Baxter

Welch Allyn Connex® Cardio ECG

This ECG provides results that can be trusted in an easy-to-understand format designed with workflow in mind.

- Captures best 10 seconds of ECG data from memory, reducing need for repeat ECGs
- Complies with ACC, AHA, and HRS recommendations for adult and pediatric patients
- Streamlines workflow by capturing and saving an ECG in as few as 2 clicks

#CC-RXX-WAXX, PC-Based ECG with Wireless Acquisition Module (WAM)
(129-6992).....ea

#CC-RXX-AAXX, PC-Based ECG with AM12™ (Wired) Cable
(129-6993).....ea

#CC-RXX-AADX, ECG with AM12™ (Wired) Cable and DICOM Connectivity
(129-6995).....ea

#CC-RXX-WADX, ECG with Wireless Acquisition Module (WAM) and DICOM Connectivity
(129-6994).....ea



You may never have to change an otoscope or ophthalmoscope bulb again.

HEINE® is the first company to replicate the color-rendering index properties of halogen bulbs in a cost-effective, durable LED solution. So, you can **see colors as they are in their natural state without the maintenance costs of halogen bulbs.**

For more information, call 1.800.P.Schein
or visit: henryschein.com/medical

California customers please use website when ordering for Prop 65 information.

HEINE is a registered trademark of HEINE Optotechnik GmbH + Co. KG.
© 2024 Midmark Corporation, Versailles, Ohio USA

seca

Precision for health



MEET OUR WORKHORSES

Back-to-school with our best physician scales

seca 777

HS Item # 139-2675

Analog measuring rod with
4" – 90" measuring range

- Capacity: 550 lbs / 250 kg
- Anti-tip column design
- Over-sized, low-profile full metal platform (no plastic)
- Rubber-coated castors for easy transport
- seca 477 handrail accessory for extra support (HS Item # 138-4356)
- Optional seca 459 Panda accessory (HS Item # 138-5120)
- Available in lbs/kg or kg-only

seca.com

seca 787

HS Item # 139-6690

Digital measuring rod with
4" – 90" measuring range

**seca**

Precision for health

seca 869 Digital Flat Scale

With Cable Remote Display

Suitable for stationary or mobile use. The display can be set up separately from the base and read from its position on the wall or desk due to its low-weight, compact size. Additional convenience with BMI calculation for nutritional conditions, HOLD to keep measurements displayed, and TARE for weighing small children held in an adult's arms.

- Cable remote display for wall-mounting or table installation
- Capacity: 550 lb/250 kg
- BMI function to determine nutritional status
- TARE, HOLD, and BMI functions
- Separate display and operating element with infinitely adjustable tilt angle
- Large, easy-to-read ¾"-high LCD figures

(116-0582).....ea

**seca**

Precision for health

seca 374 Digital Baby Scale

The seca 374 baby scale with extra-large weighing tray is perfect for babies and seated toddlers up to 44 lb. The adjustable damping function combined with the Auto-HOLD function ensures fast and precise weighing.

- Ergonomic tray with high side walls for increased safety
- Solid surfaces prevent liquid penetration, eliminating cross contamination
- Capacity: 44 lb
- Optional measuring rod (HS Item # 1165634) expands scale into a measuring station

(824-3615).....ea

seca

Precision for health

seca 222 Mechanical Telescopic Measuring Rod

With Large Measuring Range

Designed specifically for wall mounting, the seca 222 height measuring rod is aligned to work exactly to the millimeter from head to toe. The telescopic measuring rod ensures that the results can be easily read at eye level, even for the tallest patient. When not in use, the measuring slide can be folded down for safety.

- Measuring range: 2½"-90½"
- Includes heel positioner for ideal 3-point measurement
- Durable, abrasion-resistant aluminum measuring rod

(266-8199).....ea



HENRY SCHEIN®

Digital Electronic Eye-Height Scale

Compact, eye-height digital scale that features an integrated mechanical height rod, large backlit LCD display, and BMI calculator. Weight capacity up to 500 lb (225 kg). Low profile for easy access—platform is only 2½" high. Sta-Bilt stabilizing software ensures an accurate measure, even if the patient is moving. Autohold feature allows the weight to continue to be displayed after the patient has stepped off the scale. Integrated USB interface allows data to be transferred to a computer. Includes 2-year warranty.
22" x 15" x 53"
(112-7154).....ea



HENRY SCHEIN®

Mechanical Beam Scale with Height Rod


Standing mechanical beam scale with built-in height rod. Weight capacity of 450 lb (200 kg). 1-year warranty.
22¼" x 18¼" x 58¼"
(112-7157).....ea



HENRY SCHEIN®

Digital Floor Scale

Compact digital floor scale with a large backlit LCD display. Weighs accurately in both lb and kg, up to 400 lb (200 kg). Features antislip feet and a low profile for added safety. Automatically turns on. Includes 1-year warranty.
14" x 11½" x 1"
(112-7153).....ea



Medical Technologies

Easy on-PC Spirometry System

The Easy on-PC uses the power of your PC, laptop, or tablet, leveraging True Flow ultrasound technology for a complete spirometry solution. This PC spirometer uses pediatric incentives, displays inspiratory and expiratory real time curves, and offers challenge testing. Some of the essential features include:

- Flexible EMR capabilities
- Trending
- Selectable predicted values and interpretation
- Automated quality control
- Supports multiple languages

#2700-3, Easy on-PC Spirometry System
(114-6688).....ea
#2050-1HS, Spirette Mouthpiece
(570-0383)..... 50/case
#2050-5HS, Spirette Mouthpiece
(570-0384).....200/case



Medical Technologies

EasyOne Pro® LAB

This respiratory analysis system provides the most efficient and accurate way to perform complete lung function testing. With EasyOne Pro® LAB, the technician can fully focus on the patient. TrueCheck® ensures that each gas measurement is accurate without calibration and warm-up time.

- Both patient-friendly and easy to maintain—for MBW (multiple breath nitrogen washout), only an oxygen supply is needed: no expensive SF6 gas required
- Equipped with a modern, embedded PC solution that makes network and printer connections extremely easy and stable

#3100-00, EasyOne Pro LAB Only
(409-0003).....ea
#3050-2, Barriette Shields for EasyOne Pro
(119-7148).....100/case

Medical Technologies

EasyOne Pro® Portable Single-Breath DLCO System

Portable, single-breath DLCO device that allows physicians to provide their patients with prompt, accurate diagnosis and treatment. Performs PFTs in under 30 minutes, all in 1 square foot. Some of the many features include:

- Single-breath CO diffusion (DLCO)
- Spirometry
- Automatic calibration
- Flexible EMR capabilities
- No warm-up time
- Challenge testing
- No expensive service contracts
- High-resolution color touch-screen

#3000-00, Single-Breath DLCO System
(115-6219).....ea
Contains: 1 of each: single-breath DLCO system with flow sensor, gas supply tube, one-way valve, DLCO valve, pressure tube, temperature/humidity sensor & touch stylus; DLCO regulator & high-resolution color touch-screen.
#3000-10, Single-Breath DLCO System with Cart and Printer
(114-8519).....ea
#2050-1HS, Spirette Mouthpiece
(570-0383)..... 50/case
#2050-5HS, Spirette Mouthpiece
(570-0384).....200/case
#3050-1, Barriettes
(115-4279)..... 50/case
#3050-2, Barriette Shields for EasyOne Pro
(119-7148).....100/case



HENRY SCHEIN®

Green Series™ 777 Integrated Wall System

All models include BP Aneroid, Ear Specula Dispenser, and SureTemp Plus Thermometer.

#777-PM2WAS-HS, With PanOptic Basic LED Ophthalmoscope and MacroView Basic LED Otoscope (570-2666).....ea

#777-PM3WAS-HS, With PanOptic Plus LED Ophthalmoscope and MacroView Plus LED Otoscope for iExaminer (570-2665).....ea

#777-SM2WAS-HS, With Coaxial LED Ophthalmoscope and MacroView Basic LED Otoscope (570-2667).....ea

#777-SM2WAX-HS, With Coaxial LED Ophthalmoscope, MacroView Basic LED Otoscope, BP Aneroid, and Ear Specula Dispenser (Does not include SureTemp Plus Thermometer.) (570-2668).....ea



HENRY SCHEIN®

Otoscope Speculum Tips

These disposable tips are carefully manufactured to be free of rough edges or burrs. They are designed for a snug fit on Welch Allyn®, Propper®, and Heine® otoscopes. They will not fit the Welch Allyn® Macroview®.

2.5 mm (112-5527).....1000/bag
4.0 mm (112-5528).....1000/bag



HENRY SCHEIN®

PRO Digital Thermometer

• Measures body temperature by oral, axillary, or rectal read
• 60-second response
• 5 probe covers
• More than 200 hours of continuous operation (570-2055).....ea



HENRY SCHEIN®

Fingertip Pulse Oximeter

- Adult and pediatric use: ideal for all specialties
- Range for SpO₂: 70%–100%
- Range for pulse: 30–235 bpm
- Dual-color LED with display modes and brightness control
- Includes pulse bar and plethysmogram
- Visual and audible low-battery alarm
- 2 AAA alkaline batteries with 30 hours of continuous use
- Electromagnetic compatibility

(900-6368).....ea



HENRY SCHEIN®

Pro Plus Series Professional Stethoscope

Combination adult-size chest piece is precisely machined to exacting tolerances from surgical stainless steel for outstanding performance and rugged durability. Available in 3 striking colors for departmental coding or personal identification.

- Extra-large bell (1½" diameter) for unsurpassed low-frequency response
- Ultrasensitive diaphragm (1½" diameter) for greater amplification and crisper high-frequency response
- Color-coordinated nonchill bell and diaphragm retaining rim for patient comfort
- Stainless steel binaurals with double-leaf internal spring fixed at 15° angle for greater comfort
- Reinforcing yoke molded into flexible 22" PVC "Y" tubing
- Deluxe PVC ear tips (in small and large) for the ultimate in wearing comfort and acoustic seal
- Spare diaphragm included
- Available in 3 striking colors for departmental coding or personal identification
- Not made with natural rubber latex
- Scope ID tag included
- Lifetime warranty
- Weight: 6.5 oz
- Overall length: 31"

Stethoscope.....ea

Specify:

Black (900-4815) Navy (900-4817)
Burgundy (900-4816)



HENRY SCHEIN®

Digital Ear Thermometer

- For professional or personal use
- Instant reading results in 1 second
- Memory feature recalls last reading
- Fever alarm
- Low battery indicator
- Replaceable lithium battery
- °F or °C measuring mode
- Thermometer includes 15 disposable probe covers

Thermometer

(900-6438).....ea

Probe Covers

(900-6439).....150/box



HENRY SCHEIN®

The Essentials Series Pocket Aneroid

Navy, Latex

- Black enamel, 300mm-Hg, no-pin stop manometer
- Navy blue nylon cuff with range markings to facilitate selection of correct cuff size, artery label, and gauge holder
- Latex inflation bladder and bulb
- Chrome-plated deflation valve
- Complete with leatherette zippered carrying case and operating instructions
- Inspected and packaged in the USA from foreign components

Aneroid.....ea

Specify:

Child (112-6076) Large Adult (112-6080)
Adult (112-6077)

Brewer

Award Winning Versa Pediatric Exam Table

- 500 lb weight capacity
- 4.7 cubic ft of storage
- Converts within minutes
- 5 year warranty

NIGHTINGALE
AWARDS
2023



Gold Award



Brewer

Versa Exam Table

The Versa exam table provides the ability to readily and cost-effectively adapt to changing needs. When it comes to care, the Standard model provides a new standard of quality and efficiency for most practices. Easily convert to a pediatric top to safely and effectively care for children from newborns to adolescents.

- Improves the patient experience with a stable platform to support patients of all builds, and supports patient self-positioning
- Easy to connect and a removable upholstery top is available in standard adult and pediatric sizes to meet the changing needs of a facility
- Customizable to choose the options that fit within a budget
- Includes an ergonomically designed pass-through side drawer purposefully positioned at waist height to reduce bending and twisting for caregivers, while providing a generous 3.7 cubic feet of storage space



Standard..... ea

Specify:

#1000AX-XXDX-XXX-24, Black Satin (143-6121)

#1000AX-XXDX-XXX-31, Deep Sapphire (143-5875)

#1000AX-XXDX-XXX-32, Blue Fog (143-5879)

With Stirrups, Pelvic Tilt, Power Cord, Receptacle and Heater..... ea

Specify:

#1000AX-XSDT-RHX-31, Deep Sapphire (143-6136)

#1000AX-XSDT-RHX-32, Blue Fog (143-6172)

With Assist Handles, Stirrups, Pelvic Tilt, Power Cord..... ea

Specify:

#1000AX-HSDT-RXX-32, Blue Fog (143-9936)

#1000AX-HSDT-RXX-31, Deep Sapphire (143-9935)

#1000AX-HSDT-RXX-24, Black Satin (145-3150)

#1000PX-XXDX-XXX-32, Pediatric, Blue Fog (143-6065)..... ea

Brewer



Access™ High-Low Exam Table

This table's design provides an 18" profile in the seated position to allow level transfer from standard wheelchair heights, complete with upright seat-back and safety-grab bars to simplify reliable patient transfers. The system also offers abundant storage capacity with Brewer's signature pass-through drawer system and is backed by Brewer's standard 3-year warranty.

- 450-lb weight capacity
- Standard 3-year warranty
- Innovative chair design
- Chair-to-table design takes patient from consultation to examination

Exam Table.....ea

Specify:

#6501-31, Deep Sapphire (125-4768)

#6501-26, Tapestry Red (124-7582)

#6501-53, Mocha (550-0626)



Ritter® 224 Barrier-Free® Examination Chair

The number of elderly and obese in the U.S. is steadily growing, and over one fifth of adults in America are disabled. Assisting patients with mobility limitations onto exam chairs has become a frequently performed task. Although caregivers working in health care facilities are in one of the most at-risk professions for musculoskeletal disorders, little has been done to evaluate and reduce risk. The Ritter 224 is designed to offer low and high positioning heights so you can provide the best care to all patients, whether treating the elderly, expectant mothers, or patients with disabilities.

#224-002, Without Clean Assist Roller System (133-0159).....	ea
#224-003, With Clean Assist Roller System (133-0268).....	ea



Ritter® 230 Universal Power Procedures Table

Features eight-way power for effortless patient positioning. Whether it is power height, back, foot, or tilt, it does the lifting for the physician—minimizing back stress, reducing workplace injury, and keeping patients comfortable.

- Power height, back, foot, and tilt
- 450-lb patient weight capacity
- 22½" to 40" seat height adjustment
- "Knee-chest" position configurable
- Hand control standard
- Seamless upholstery top
- Articulating headrest
- Adjustable stirrups
- 18" paper roll holder
- Meets UL 60601-1 and CAN/CSA C22.2 No.60101-M90 safety standards
- Hospital grade electrical outlets and table rotation also available
- Optional foot control
- Available in 8 seamless upholstery colors

Tables	
#230-001, Ritter 230 Base Unit (108-1319).....	ea
#230-002, Ritter 230 Base Unit with Receptacle (238-7678).....	ea
#230-003, Ritter 230 Base Unit with Rotation (475-6182).....	ea



Ritter® 225 Barrier-Free® Examination Chair

The number of elderly and obese in the U.S. is steadily growing, and over one fifth of adults in America are disabled. Assisting patients with mobility limitations onto exam chairs has become a frequently performed task. Although caregivers working in health-care facilities are in one of the most at-risk professions for musculoskeletal disorders, little has been done to evaluate and reduce risk. The Ritter 225 is designed to offer low and high positioning heights so you can provide the best care to all patients, whether treating the elderly, expectant mothers, or patients with disabilities.

#225-002, Without Clean Assist Roller System (133-0279).....	ea
#225-003, With Clean Assist Roller System (133-0290).....	ea



626 Barrier-Free® Examination Chair

Many health care providers are not aware of current standards for proper blood-pressure measurement. If this isn't addressed, it leaves health systems at serious risk. 9.8% of the population is affected by overestimation of high blood-pressure errors. Multiply the cost of treatment for hypertension by the typical patient panel size and the cost of overtreatment could be as much as \$135,000 per each practicing physician. The impact of overmedication can also result in complications such as dizziness, fainting, exhaustion, or even patient falls. The Midmark 626 is the first exam chair of its kind designed to promote proper patient positioning for a more accurate blood-pressure measurement.

- Easily position patients for a better BP measurement
- Powered movement of the back section ensures the patient's back is supported
- A Barrier-Free® low chair-height allows the majority of patients, even female patients with heights in the 3rd percentile, to place their feet flat on the floor

#626-001, Base Model (133-0301).....	ea
#626-002, With Heated Upholstery (133-0312).....	ea
#626-003, With IQ Scale (133-0323).....	ea
#626-004, With Heated Upholstery and Scale (133-0334).....	ea



425 Air Lift Physician Stool

Designer Series

Designed for your comfort and support all day, every day. Choose the features you want. Start from the "bottom up" with the frame and seat, then add the back and armrests. Order everything you want at once, or add on later, as you prefer.

(151-4236).....	ea
-----------------	----



272 Exam Stool

Provides optimal comfort, support, and stability. Star-shaped base design allows for closer access to the patient, and stable, easy maneuverability.

- Black, plastic composite, 5-leg caster base and hood
- Pneumatic seat-height adjustment

Exam Stool.....	ea
Specify:	
#272-001-845, Lunar Gray	(132-2295)
#272-001-814, Stone	(132-2291)



Better care experiences for all ages and sizes.

MTI chairs, cabinets, and accessories are designed with specialty caregivers in mind.

- Designed to have lower entry heights smaller bases and backrests for easier access;
- Made with motorized lifts, bringing every patient to working height;
- Manufactured for long term dependability with 41% lower repair costs on average than our competitor.

Specialty equipment for Podiatry, ENT, Dermatology, Oral Surgery, Wound Care, Ophthalmology, Plastic Surgery

For more information, contact your local Henry Schein Representative.



25 MTI
CELEBRATING
25 YEARS



- 725-lb patient load capacity
- 4 programmable positions
- Programmable home and ADA-compliant home positions
- Power lift, back, tilt, and foot extension
- Hand and foot controls
- Oval articulating headrest
- 22" or 25" seat width
- 18¾"–39" seat height
- Leg rest features knee-break with 7" power foot extension
- Optional patient arms: removable, floating (compensating) arms, or seat rail or back rail mounted arms
- Removable upholstery: seamless or designer
- Tapered backrest for easier caregiver access
- Sturdy all-steel frame
- Multiple G2 accessory options

#829-001-80, Standard Base (22" seat width)
(135-8121).....ea

#829-002-80, Swivel Base (22" seat width)
(139-4043).....ea

#829-003-22, Mobile Base (22" seat width)
(140-6171).....ea

#830-001-80, Standard Base (25" seat width)
(133-8866).....ea

#830-002-80, Swivel Base (25" seat width)
(136-9044).....ea

#830-003-22, Mobile Base (25" seat width)
(133-9127).....ea



Specialized Equipment
for Specialty Healthcare
Medical Technology Industries

463 and 464 Exam/Procedure Chairs

The most feature-packed, ADA compliant* chair on the market. Standard features include ComfortSync™ coordinated footrest, EasyView™ chair-back mounted switches, removable slide-back floating armrests, push-button 333° swivel base and oval articulating headrest.

- 725-lb patient load capacity
- 4 programmable positions
- Programmable home and ADA-compliant home positions
- Power lift, recline, and tilt (no tilt on 463)
- 21" seat width
- 18¾"–37" seat height (2" lower than any other ENT chair on the market)
- Removable upholstery: seamless or designer
- Tapered backrest for easier caregiver access
- Sturdy, all-steel frame

#463-002-80, Swivel Base
(136-8741).....ea

#464-002-80, Swivel Base
(136-8744).....ea

*ADA compliant with transfer support accessory for easy patient ingress/egress




Specialized Equipment
for Specialty Healthcare
Medical Technology Industries

830 Series Procedure Chair

The highly versatile 830 Series chairs are ADA compliant in all base versions (standard, swivel, or mobile) using the transfer support, providing easy patient ingress/egress, and equal access for patients with disabilities and caregivers. All 830 chairs come equipped with SmartTech™ technology to enhance safety, reduce energy consumption, and provide smoother movement, additional memory options, and more user preferences.



An AED is only helpful in an emergency if it's visible and accessible to bystanders. Be prepared for the school year ahead.



ZOLL

ZOLL

ZOLL AED Plus®

Bringing Lifesaving Defibrillation to Everyone

- Provides a one-piece electrode for accurate placement
- Uses consumer lithium batteries that are easy to replace, long-lasting, and inexpensive

Semi-Automatic Unit

(866-7559).....ea

Contains: AED, electrodes, batteries, carry case, ECG, wall-mount bracket & operator's guide.

Fully Automatic Unit

(119-7910).....ea

Contains: AED, electrodes, batteries, carry case, ECG, wall-mount bracket & operator's guide.

CPR-D padz®—Adult Defibrillator Pad

(253-2611).....ea

Pedi padz® II—Pediatric Electrode Pad

(203-7054).....ea

Type 123 Lithium Batteries

(610-4769).....10/pkg

Soft Carrying Case

(106-1937).....ea

3D Wall Sign

(110-5384).....ea

Wall Cabinet

(108-7520).....ea

AED Plus Simulator

(499-7385).....ea



ZOLL

Powerheart® G5 AED

The Powerheart® G5 AED with Intellisense® CPR provides corrective feedback to ensure delivery of high-quality compressions. One-button dual-language functionality and comprehensive self-tests give rescuers the confidence to act.

- Intellisense® CPR Feedback provides real-time voice and text corrective prompts for CPR quality
- Rescue Ready® conducts regular, automatic self-tests on pads, battery, and internal circuitry
- RescueCoach® guides users through each critical step of the rescue with user-paced voice/text prompts
- Dual language change from primary to secondary language in the push of a button

#G5S80C-S, Semi-Automatic

(139-5101).....ea

#G5A80C-S, Fully Automatic

(139-5102).....ea



ZOLL

AED Pro® Automated External Defibrillator

This external defibrillator provides advanced capabilities for BLS and ALS users with innovative Real CPR Help® technology for rate and depth of chest compressions. The See-Thru CPR® capability filters CPR chest compressions, displaying the underlying cardiac rhythm. The AED Pro is unmatched in its ruggedness and durability and is ideal for any environment.

- Available as semiautomatic only or semiautomatic with manual override for ALS professionals
- Built-in internal memory and external USB memory technology allows expanded data storage capacity and maximum flexibility for data management and transfer
- Electrodes and battery compatibility with ZOLL M Series®, E Series®, and R Series® professional defibrillators
- ECG monitoring with 3-lead cable and high-quality display
- Unrivaled battery options (long-life disposable or rechargeable) that meet specific needs

(869-6270).....ea

IT'S THAT TIME OF YEAR

graham
medical



Back to School shopping
means
Back to School Appointments.

**Make sure your offices
are ready to go.**

Henry Schein Product #	Product	Material	Color	Size	Packaging
(162-8990)	MediShorts®	Nonwoven	Blue	Pediatric	50/CS
(104-9190)	Gown	Tissue/Poly	White	Pediatric	50/CS
(123-7369)	Gown	Tissue/Poly	Animal Parade	Pediatric	50/CS
(118-3930)	Gown	Tissue/Poly	Blue	30" x 42"	50/CS
(395-0890)	Tablepaper	Crepe	Animal Parade	18" x 125'	12RL/CS
(991-7184)	Tablepaper	Crepe	PediaPals®	18" x 125'	12RL/CS
(118-3932)	Drape	2-Ply	White	40" x 48"	100/CS
(395-7226)	Tablepaper/Scale liner	Poly-Perf	White	2" x 125'	12RL/CS



HENRY SCHEIN®

Premium Exam Drapes

2-Ply, All-Virgin Tissue

- Designed to blanket the patient in lightweight comfort

White, 36" x 48"	(101-3658)	100/case
White, 40" x 48"	(100-1797)	100/case
Blue, 40" x 48"	(100-3360)	100/case
Mauve, 40" x 48"	(100-6428)	100/case
White, 40" x 60"	(101-8521)	100/case



HENRY SCHEIN®

Premium Exam Drapes

3-Ply, All-Virgin Tissue

- Soft, opaque, and heavyweight
- Embossed tissue drape with dependable wet strength is ideal for short to moderate exams
- Conforms to the body and offers the ultimate in absorbency

White, 40" x 48"	(101-7084)	100/case
------------------	------------	----------



HENRY SCHEIN®

Dri-Gard® Towels

Teddy Bear Design

- Add color to the operatory and bring smiles to younger patients
- Made in the USA
- Not made with natural rubber latex

2-Ply Tissue/1-Ply Poly		
13" x 19"	(100-9712)	500/case



HENRY SCHEIN®

Dri-Gard® Towels

2-Ply Tissue/1-Ply Poly

- Provide positive moisture barrier and quick absorbency
- Unique horizontal embossing pattern protects patient's clothing and increases the strength of the towel
- Made in the USA
- Not made with natural rubber latex

13" x 19"		500/case
Specify:		
Aqua	(100-9276)	Lavender (101-0476)
Beige	(100-3501)	Peach (100-7384)
Blue	(100-7265)	Rose (100-6341)
Gray	(100-1313)	White (100-5715)
Green	(100-1350)	Yellow (100-9704)



HENRY SCHEIN®

Dri-Gard® Plus Towels

3-Ply Tissue/1-Ply Poly

- Provide a positive moisture barrier with maximum absorbency
- Unique horizontal embossing pattern protects patient's clothing and increases strength of the towel
- Made in the USA
- Not made with natural rubber latex

13" x 19"		500/case
Specify:		
Aqua	(100-9851)	Lavender (101-6810)
Beige	(100-5015)	Peach (100-7395)
Blue	(100-9825)	Rose (100-2024)
Gray	(100-9545)	White (100-7908)
Green	(100-4479)	Yellow (100-3395)



HENRY SCHEIN®

Standard Professional Towels

Specially treated to increase the wet and dry tensile strength of towel. Provides the perfect blend of softness and absorbency. Made in the USA. Not made with natural rubber latex.

13" x 19"		500/case
Specify:		
White	(100-3161)	Blue (100-0183)



HENRY SCHEIN®

Pillowcases

Tissue/Poly/Tissue

- Provide linen-like softness for patient comfort and an inner poly barrier for extra protection
- Strong and durable

White, 21" x 30"	(104-9206)	100/case
------------------	------------	----------



HENRY SCHEIN®

Premium FABRICEL® Pillowcase

1-Ply Tissue/Poly

- Disposable pillowcase is super soft and strong because it is made from Fabricel® material
- By bonding supple tissue to polyethylene film with a heat-lamination process, a unique material that is soft, strong, and tear resistant is created
- Deluxe material prevents soak-through

White, 21" x 30"	(101-2526)	100/case
------------------	------------	----------



HENRY SCHEIN®

Headrest Covers

Tissue/Poly

Soft tissue outside feels smooth and comfortable to the touch. The poly backing helps prevent moisture soak-through to headrest.

- Tissue outer surface provides greater comfort for your patients
- Poly lined to help protect your equipment

Small Size		
10" x 10"		500/case
Specify:		
White	(101-7104)	Mauve (101-3198)
Blue	(101-4595)	Green (102-8451)
Lavender	(102-5543)	
Regular Size (13" Open Side)		
10" x 13"		500/case
Specify:		
White	(101-8627)	Mauve (101-1051)
Blue	(101-8955)	Green (102-5542)
Lavender	(102-9374)	

HENRY SCHEIN®

Table Paper

- Quiet, absorbent barrier
- Designed and manufactured to provide strength, protection, and value
- White

Smooth, 225 ft

18"

(900-4978)..... 12 rolls/case

21"

(900-4979)..... 12 rolls/case



HENRY SCHEIN®

Premium Table Paper Rolls

- Designed and manufactured to provide strength, protection, and value
- White

Crepe, 125 ft

18"

(100-9865)..... 12/case

21"

(100-1657)..... 12/case

Smooth, 225 ft

14"

(104-4501)..... 12/case

18"

(100-9491)..... 12/case

21"

(100-0177)..... 12/case

Smooth, Heavy Duty, 225 ft

18"

(104-6092)..... 12/case

21"

(104-5273)..... 12/case



HENRY SCHEIN®

Poly-Perf® Exam Table Paper

- Absorbent top tissue layer with poly back for fluid barrier
- For exam table coverage; can also be used as scale and table liners
- Perforated every 18" or 21"
- 125-ft rolls
- White

18"

(570-0728)..... 9/case

21"

(570-0727)..... 9/case



HENRY SCHEIN®

Facial Tissues

- 2-ply: strong, and absorbent
- Soft and elegantly embossed

White, 9" x 8"

(900-4214)..... 100/box



HENRY SCHEIN®

Multifold Towels

- 1-ply: strong, and absorbent
- Interlocking technology
- Whitened without chemicals containing chlorine bleach
- Fits a wide range of universal dispensers
- 250 sheets/pack
- Made in the USA

White, 9½" x 9¼"

(570-0703)..... 4000/case

Contains: 250 sheets per pack, 16 packs per case.



HENRY SCHEIN®

C-Fold Towels

- 1-ply: strong, and absorbent
- Individually folded towels
- Whitened without chemicals containing chlorine bleach
- Fits a wide range of universal dispensers
- 150 sheets/pack
- Made in the USA

White, 12.87" x 10.5"

(570-0704)..... 2400/case

Contains: 150 sheets per pack, 16 packs per case.



HENRY SCHEIN®

Disposable Airlaid Washcloths

- Airlaid is made from synthetic staple fibers and natural cellulose pulp
- Offers absorbency, strength, and softness
- Ideal for incontinence cleanups, bathing, and other patient applications

White, 10" x 10"

(570-2403)..... 1000/case

Contains: 50 washcloths per bag, 20 bags per case.



HENRY SCHEIN®

Plastic Cups

Embossed Drinking Cups with Rolled Lips

- Unique double lip allows for consistent, single-cup dispensing
- Heavy, sturdy plastic: leakproof; won't crumple in hands
- Not made with natural rubber latex
- Made in the USA

5 oz..... 1000/case

Specify:

White	(100-3902)	Dusty Rose	(100-2104)
Blue	(100-4361)	Green	(101-2443)
Yellow	(100-1596)	Peach	(101-6812)
Beige	(100-3193)	Lavender	(101-2156)
Gray	(100-0473)		



HENRY SCHEIN®

Plastic Cups

Embossed Drinking Cups with Rolled Lips

- Heavy, sturdy plastic; leakproof and will not crumple in hands
- Releases easily from dispenser; easily accessible

Translucent, 5 oz

(104-9869)..... 100/pkg



HENRY SCHEIN®

Graduated Plastic Medicine Cups

Plastic medicine cups with easy-to-read graduations.

(112-5555)..... 100/bag



HENRY SCHEIN®

Emesis Basins

- Plastic
- Kidney shaped
- 16-oz capacity

Mauve, 8½"

(112-5809)..... 25/box



Cotton Rolls

Starched, 100% cotton, absorbent nasal roll.

- Starched
- 100% cotton
- Absorbent nasal roll

#2 Medium, Unwrapped, Non-sterile, 1½" x ¾"	(102-5455).....	2000/box
#2 Medium, Wrapped, Non-sterile, 1½" x ¾"	(101-6327).....	2000/box
#2 Medium, Wrapped, Sterile, 1½" x ¾"	(101-8764).....	2000/box



Premium Cotton Rolls

#2 Medium, 1½" x ¾"

Made from 100% cotton for greater absorbency. More absorbent than the leading brands. Longer fibers reduce linting. Longer bleaching process produces a whiter cotton roll. Can easily be cut into smaller pieces with scissors for effective use in tight oral-cavity areas. Nonlinting, non-sterile, and pliable. Maintains shape when wet. Will not adhere to mucous membranes.

Non-sterile	(896-9286).....	2000/box
Sterile	(104-8192).....	2000/box



Cotton Bulk Absorbent Rolls, USP

1-lb Bulk Rolls

Sterile	(100-0145).....	ea
Non-Sterile	(100-0147).....	ea



Cotton Balls

- 100% pure-white cotton
- Soft and absorbent

Sterile, Medium	(100-4752).....	500/box
Non-sterile, Medium	(100-0468).....	4000/case
Non-sterile, Large	(100-1491).....	2000/case



Tongue Depressors

- Highly polished
- Smooth
- Designed for oral use or topical application

Non-Sterile		
Junior, 5½"	(100-0758).....	500/box
Adult, 6"	(100-2416).....	500/box
Sterile		
Adult, 6"	(100-0238).....	100/box



Cotton-Tipped Applicators

- Applicators are highlighted by a soft cotton tip that is absorbent and will not lint

Non-sterile, 3"	(100-6015).....	1000/box
Non-sterile, 6"	(100-9175).....	1000/box
Sterile, 6", 100 Packs of 2	(100-9249).....	200/box



Wooden Applicator/Mixing Sticks

Plain Wood, 6"

- International peel symbols: ensure easy opening from proper end of package
- Long shelf life (more than 10 years): economical

(100-9485).....	1000/box
-----------------	----------



Protection Plus Disposable Underpads

- Made of waterproof, blue poly backing
- Absorbent core of tissue layers
- White, hypoallergenic spunbond facing sheet
- Not made with natural rubber latex

17" x 24"	(112-5633).....	300/case
23" x 24"	(112-5634).....	200/case
23" x 36"	(570-1182).....	150/case
30" x 30"	(570-1183).....	150/case
30" x 36"	(570-1181).....	100/case



Patient Belongings Bags

- Convenient bags for moving and storing patient belongings
- Plastic, opaque white
- Print: Patient Belongings, Name, Room
- 20" x 20" x 4"

Bags with Handles	(570-1179).....	250/case
Bags with Drawstrings	(570-1180).....	250/case

Exam Gloves



CRITERION

CRITERION® Nitrile N100

Formerly known as Criteron® Nitrile. Same Great Product, New Criterion® Packaging\nNitrile Powder-Free Exam

Durable nitrile construction for reliable protection and tactility. Chlorinated for reduced friction and easy donning. Textured fingers provide an enhanced grip. Chemotherapy tested to the ASTM D6978 standard. Beaded cuff. 3.15-mil thickness at the palm. Blue.

Gloves.....	100/box
Specify:	
X-small	(111-8533)
Small	(111-8535)
Medium	(111-8536)
Large	(111-8537)
X-large (90/box)	(111-8539)



CRITERION

CRITERION® Low Derma

Powder-Free Nitrile Exam, Non-Sterile
Made from a soft, durable nitrile formulation without the use of chemical accelerators or sulfur, preventing Type-IV allergic reactions. Reduces the potential for dermatitis for glove wearers, yet still provides softness and comfort and maintains superior protection. Chlorinated for reduced friction and easy donning. Textured fingers promote sure handling of instruments and materials. Chemotherapy tested to the ASTM D6978 standard. Beaded cuff. 2.0-mil thickness at the palm. Slate Blue.

Gloves.....	200/box
Specify:	
X-Small	(570-0630)
Small	(570-0631)
Medium	(570-0632)
Large	(570-0633)
X-Large (180/box)	(570-0634)



CRITERION

CRITERION® N300

Powder-Free Nitrile Exam, Non-Sterile
Ultrathin formulation for excellent tactile sensitivity and dexterity. Chlorinated for reduced friction and easy donning. Textured fingers promote sure handling of instruments and materials. Chemotherapy tested to the ASTM D6978 standard. Ecofriendly packaging—with 300 gloves in a box—reduces waste and save space. Beaded cuff. 2.0-mil thickness at the palm. Ice Blue.

Gloves.....	300/box
Specify:	
X-Small	(112-6829)
Small	(112-6830)
Medium	(112-6831)
Large	(112-6832)
X-Large (250/box)	(112-6833)



CRITERION

CRITERION® N200

Powder-Free Nitrile Exam, Non-Sterile
Slightly thinner than the N100, the durable nitrile construction provides reliable protection and tactility. Chlorinated for reduced friction and easy donning. Textured fingers provide an enhanced grip. Chemotherapy tested to the ASTM D6978 standard. Beaded cuff. 2.36-mil thickness at the palm. Available in both Periwinkle Blue and Sky Blue.

Periwinkle Blue.....	200/box
Specify:	
X-Small	(900-7437)
Small	(900-7438)
Medium	(900-7439)
Large	(900-7440)
X-large (180/box)	(900-7441)
Sky Blue.....	200/box
Specify:	
Small	(570-2098)
Medium	(570-2099)
Large	(570-2100)
X-Large (180/box)	(570-2101)



CRITERION

CRITERION® Aloe + Vitamin E

Powder-Free Nitrile Exam, Non-Sterile

- Extra-soft stretch nitrile with soothing aloe and vitamin E lining that provides an enhanced coating to rejuvenate and nourish hands
- Beyond the benefit of enhanced skin protection, this glove has the added benefit of easy donning and the strength of nitrile
- Not made with natural rubber latex
- Thin feel fits like latex and provides improved tactile sensitivity
- Textured fingers provide excellent grip in a wet/dry dental setting
- "Smart Dispense" technology allows a single-glove delivery with fewer wrinkles
- Radiant Razz color brightens up any procedure and patient experience
- 2.4-mil thick

Gloves.....	200/box
Specify:	
X-Small	(112-7143)
Small	(112-7144)
Medium	(112-7146)
Large	(112-7147)



Together, we can build a healthier world for all.

Practice Green is designed to empower the health care community to positively impact the future of the planet by reducing the ecological footprint and promoting sustainability—both in daily life and while at work.



A product included in the Practice Green portfolio has been carefully selected by Henry Schein as “green” based on the product information provided by the supplier partner that it is either reusable, contains recycled materials, can be recycled, is biodegradable, or is from a sustainable source.



Scan the QR code to learn more, or visit
www.henryschein.com/PracticeGreen

© 2024 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors. 24DS3152

**NEW!****Precision Eco Bio Face Masks,
Level 3**

Made from biopolymers that are produced from renewable crop resources.

- Hypoallergenic
- Absorbent, non-woven inner layer and hydrophobic outer layer
- Non-latex and graphene-free
- Bendable nose-clip wire for custom fit

Earloop Masks..... 50 per Box

Specify:
White (145-3980) Black (145-3981)

**BeeSure® Radiant Face Masks,
Level 3**

Designed to maintain a noncollapsible space for maximum breathing room. Forms an ergonomic seal around nose, chin, and cheeks to minimize air leakage.

- Fluid resistant
- Not made with natural rubber latex
- Lint-free and dye-free

Earloop Masks..... 40 per Box

Specify:
Aqua (143-3260) Lime (143-3261)
Orange (143-3262) Violet (143-3263)

**BeeSure® Floral Earloop Masks
Level 3**

- Quad-Fold® design offers 15% more breathing volume
- RoundEdge® corners and cushioned earloops for comfort
- Non-latex CoolSilk® and dye-free inner layer for a smooth, nonirritating feeling

Earloop Masks..... 50/box


Specify:
Hibiscus Pink (142-4158)
Plumeria Blue (142-4159)
Orchid Lavender (142-4163)
Fern Green (142-4165)
Daisy Yellow (142-4164)

**Ansell****Microflex® Neogard™ Touch
73-737**

Chloroprene Powder-Free Exam
Comfortable neoprene (polychloroprene) exam gloves that are certified by US Ergonomics for enhanced wet grip capabilities. Designed with ERGOFORM™ Technology to support musculoskeletal health and increase worker productivity. 3.5-mil fingertip thickness and 2.8-mil palm thickness. Green.

Gloves..... 200/box

Specify:
X-small (142-5299)
Small (142-5300)
Medium (142-5301)
Large (142-5302)
X-large (142-5303)

**Ansell****MICRO-TOUCH® DermaClean®**

Latex Powder-Free Exam

- Textured outer surface provides an excellent grip in wet and dry conditions
- Fully textured surface for excellent grip
- Superior elasticity provides great comfort and flexibility
- Natural rubber latex formulation ensures a comfortable, well-fitting glove
- Smooth chlorinated inner surface facilitates ease of donning
- 9½" glove length
- 5.5-mil palm thickness
- 6.3-mil finger thickness

Gloves..... 100/box

Specify:
X-Small (143-2488)
Small (143-2489)
Medium (143-2490)
Large (143-2491)
X-Large (80/box) (143-2492)

**Ansell****Soft White Nitrile**


Nitrile Powder-Free Exam

- Soft nitrile formulation provides a comfortable fit and feel
- The lightweight nitrile construction provides extra comfort and dexterity while working
- Protection from Type I latex allergy in HCW's or patients
- HYDRASOFT™ Skin Moisturizing Technology: Helps prevent dryness

Gloves..... 100/box

Specify:
X-Small (565-4510)
Small (565-8087)
Medium (565-7431)
Large (565-9481)
X-large (565-1575)


Ansell

LOW
DERMATITIS
POTENTIALNEW &
IMPROVED

MICROFLEX® XCEED® XC-310

YOUR NEW FAVORITE GLOVE IS NOW BETTER THAN EVER

SOFTER & STRONGER THAN LEADING COMPETITORS



Note: Data collected at Ansell Scientific Affairs labs in Reno, NV on March 2, 2022. ASTM test methodology used for tensile strength, elongation and modulus tests.
EN test methodology used for EN force at break test.

*Please see specially marked packaging

For more information, please reach out to your sales representative or visit [ansell.com](https://www.ansell.com).

Ansell, ® and ™ are trademarks owned by Ansell Limited or one of its affiliates. © 2024 Ansell Limited. All Rights Reserved.

bee sure
be good

**more gloves per box =
more patients served**



**less transportation cost =
less carbon footprint**



**overall decrease in
gloves cost to you**



The First and Only 350 Gloves/box UltraSlim Nitrile Exam Gloves

350 gloves/box mean more savings and being more green
XS (141-4267) S (141-4268) M (141-4269) *L (141-4270) *300/box



Super Slim Nitrile Exam Gloves

Superior second skin comfort, Super soft & stretchy for greater tactile sensitivity
XS (353-0052) S (353-0053) M (353-0054) L (353-0055)



Slim Nitrile Gloves

Immaculate tactile sensitivity + White Color
XS (108-3835) S (108-7932) M (108-4480) L (108-9150) XL (108-3559)



[beesure usa](https://www.beesureusa.com) [beesureusa](https://www.beesureusa.com)

For more information, contact your local Henry Schein Medical Representative or visit [henryschein.com/medical](https://www.henryschein.com/medical)

CRITERION

CRITERION™ Sensitive Skin Earloop Masks, Level 3

Same great features as the Criterion™ Sensitive Skin Earloop Masks Level 1, but with the highest barrier of protection.

- Fluid resistant
- Earloops and flexible nosepiece for a secure fit
- Hypoallergenic and free of chemicals, inks, and dyes
- Made in the USA
- Not made with natural rubber latex

Masks.....50/box

Specify:

Blue (570-3276) Pink (570-3277)



CRITERION

CRITERION™ Procedure Masks Level 3

- Light and comfortable for extended wear
- Filter media is surrounded by soft, nonirritating inner and outer facing
- Full-width flexible nosepiece provides fit and security
- Not made with natural rubber latex

Masks.....50/box

Specify:

Black (570-3273) Teal (570-2042)
Blue (570-2044) White (570-2043)
Lavender (570-2047) Yellow (570-2046)
Pink (570-2045)



CRITERION

CRITERION™ Earloop Masks Level 2

Formerly known as Henry Schein Procedure Masks Level 2. Same Great Product, New Criterion™ Packaging

- Made from spunbond outer/inner layers with filter media
- Flat earloops for greater comfort
- Meets ASTM standards as a low-barrier (fluid resistant) face mask
- Not made with natural rubber latex
- Made in the USA
- ASTM Level 2

Masks.....50/box

Specify:

Blue (104-7799) Pink (104-9200)
Lavender (104-8071)



CRITERION

CRITERION™ Surgical Tie-On Masks

4-ply surgical mask adds an extra layer of protection and comfort. Same great features as Criterion™ Surgical Tie-On Masks Level 1, but with the highest barrier of protection.

- Fluid resistant
- Flexible nosepiece
- Pleated design for enhanced breathability
- Not made with natural rubber latex

Level 3

(570-3294).....50/box



CRITERION

CRITERION™ Earloop Masks Level 1

Formerly known as Henry Schein Earloop Level 1 Procedure Masks. Same Great Product, New Criterion™ Packaging

- Light, comfortable masks for extended wear when necessary
- Ear loops are made from special stretch yarn: no irritating rubber or plastic touches the skin
- 3-ply, glass-free filter media
- Soft, nonirritating inner and outer facing
- Full-width flexible nosepiece provides fit and security
- Manufactured to exceed all current ASTM F2100-04 standards
- Not made with natural rubber latex
- ASTM Level 1

Masks.....50/box

Specify:

Blue (104-3809) Yellow (104-2849)
Lavender (104-8600) White (104-6611)
Pink (104-3730)



CRITERION

CRITERION™ Sensitive Skin Earloop Masks Level 1

Non-Latex, Blue

- Fluid resistant
- Earloops and a flexible nosepiece for a secure fit
- Hypoallergenic and free of chemicals, inks, and dyes
- Made in the USA
- Not made with natural rubber latex

Ultrasoft.....50/box

Specify:

Blue (104-8072) Pink (104-8073)



HENRY SCHEIN®

Earloop Face Masks with Eye Shields

- Optically clear shield is treated on both sides with antifog material
- Not made with natural rubber latex
- Teal
- ASTM Level 3

(104-7321).....25/box



CRITERION

CRITERION™ Anti-Fog Earloop Masks, Level 3

Lightweight and comfortable, these 3-ply procedural face masks offer ASTM Level 3 protection. Each mask features an integrated antifog strip to prevent eyewear from fogging and a soft, nonirritating inner and outer layer.

- Fluid resistant
- Pleated design for enhanced breathability
- Earloops and a flexible nosepiece for a secure fit
- Not made with natural rubber latex

Masks.....50/box

Specify:

Blue (570-3281) White (570-3283)
Lavender (570-3282) Pink (570-3284)



CRITERION

CRITERION™ Extra Breathing Room Earloop Masks

Formerly known as Maxi-Gard™ Air Procedure Masks With Air Space™ Technology Level 1. Same Great Product, New Criterion™ Packaging

With a contoured design that enhances airflow and breathability, these 3-ply procedural face masks offer ASTM Level 1 and Level 3 protection. Each mask features a soft, nonirritating inner and outer layer.

- Fluid resistant
- Earloops and a flexible nosepiece for a secure fit
- Not made with natural rubber latex

Masks.....50/box

Specify:

Blue (570-0722) Lavender (570-0724)
White (570-0723)



HENRY SCHEIN®

Pleated Surgical Face Masks

- Triple-layer, fluid-resistant, and high filtration rate
- Soft, lightweight, comfortable design
- Manufactured to exceed all current ASTM F2100-03a standards
- Cloth-like ties
- Shingle pleats: eliminates ledge for fluid pooling
- Full-width covered nosepiece: ensures proper closure
- Polypropylene outer facing, polypropylene inner: less chance of wire dislodging
- Not made with natural rubber latex
- ASTM Level 1

Tie-On Masks

(104-6809).....50/box





LEVEL UP! S3+

It's time to level up your game plan with our brand new **S3+**. This face mask offers **ASTM level 3** protection for remarkable fluid resistance, along with Quad-Fold for increased and smooth breathability. S3+ includes our Cool-Silk lint free and dye free inner layers to deliver a lighter feel and less irritation throughout your day.

S3+ Face Masks 50/box
White (138-8759), Blue (138-8758), Black (139-2761)

© CRANBERRY.USA f CRANBERRYUSA c CRANBERRYUSA.COM
For more information, contact your local Henry Schein Medical Representative or visit henryschein.com/medical

Cranberry

Trusted & Trustworthy

It's true. The PURELL® brand is trusted and preferred by healthcare workers.

But more importantly, we work tirelessly to earn that trust.

Our hand hygiene is gentle on hands, even with frequent use. And our powerful surface products can even be used near patients and on electronics.

The Difference Is Clear.



Contact your distributor or GOJO representative today.

Distributed by:



©2023 GOJO Industries, Inc. All rights reserved. | 34993 (05/2023)





HENRY SCHEIN® **Hand Sanitizer**

Waterless, alcohol-based gel hand sanitizer. Ideal for use where soap and water are not readily available. Rinse-free: no water required. Kills 99.9% of most common disease-causing germs spread by hands. Soothing emollients and a pleasant, fresh, clean scent. Leaves hands soft and clean. 66.5% ethyl alcohol. Made in the USA.

- Waterless, alcohol-based gel hand sanitizer
- Kills 99.9% of most common disease-causing germs spread by the hands
- Fresh clean scent and soothing emollients leave hands soft and clean

8-oz Pump Bottle	
(900-4443).....	ea
16-oz Pump Bottle	
(900-4442).....	ea



HENRY SCHEIN® **Maxiclens™ Antiseptic Skin Cleanser**

4.0% chlorhexidine gluconate formulation. Fast acting in just 15 seconds, while inhibiting bacterial regrowth for up to 6 hours. Skin emollients protect against irritation, keeping hands soft and soothed. Residue-free, nonsticky formulation with no added dyes. All-in-one solution is tough on bacteria, yet gentle on the skin. Pleasant floral scent.

1-qt Bottle	
(570-3268).....	ea
1-gal Refill	
(570-3269).....	ea



HENRY SCHEIN® **Foaming Antibacterial Hand Soap**

- Provides a rich, dense lather
- Rinses easily using minimal water
- Leaves hands feeling soft and smooth
- Pleasant Spring Water Clean fragrance
- Contains 0.6% chloroxylenol (PCMX) for antimicrobial action, and aloe vera to add moisture with each use

18.6-oz Pump Bottle	
(900-4441).....	ea



HENRY SCHEIN® **Antibacterial Liquid Hand Soap**

- Provides a luxurious feel, leaving hands smooth and soft
- Pleasant Cucumber and Aloe fragrance
- Contains 0.6% chloroxylenol (PCMX) for antimicrobial action, and aloe vera to add moisture with each use

8-oz Pump Bottle	
(900-4439).....	ea
16-oz Pump Bottle	
(900-4438).....	ea
1-gal Refill Bottle	
(900-4440).....	ea



HENRY SCHEIN® **Quick Wipes**

- Soft, strong towelettes: allow easy and gentle cleansing of hands and face
- Effective and gentle formulation that contains aloe: formulated for fast drying and gentleness
- Convenient packaging: towelettes easily dispensed through orifice; each towelette is saturated with cleaning liquid

Pop-Up Container, 9" x 5½"	
(104-4176).....	150/cont



HENRY SCHEIN® **Maxiclens® Instant Hand Sanitizer**

- Kills 99.9% of germs in 15 seconds without the use of soap and water
- Contains aloe vera to add moisture with every use
- Henry Schein's hand sanitizer helps meet OSHA, APIC, and CDC hand-washing recommendations
- Citrus scent

800-mL Bag for Manual Dispenser	
(112-5988).....	ea
800-mL Bag for Automatic Dispenser	
(112-5991).....	ea
Manual Dispenser	
(112-5984).....	ea
Automatic Dispenser, Hand Sanitizer	
(112-5986).....	ea
Stand for Dispenser	
(112-5989).....	ea



HENRY SCHEIN® **Maxiclens® Instant Antimicrobial Soap**

- Contains 0.6% PCMX for antimicrobial action and aloe vera to add moisture with each use
- Aids in reducing cross contamination
- Helps meet APIC and OSHA recommendations for hand washing

800-mL Bag for Manual Dispenser	
(112-5987).....	ea
800-mL Bag for Automatic Dispenser	
(112-5990).....	ea
Manual Dispenser	
(112-5983).....	ea
Automatic Dispenser, Soap	
(112-5985).....	ea
Stand for Dispenser	
(112-5989).....	ea

Infection Prevention



HENRY SCHEIN® MaxiSpray® Plus

Surface Disinfectant Cleaner

- Hospital-level, general-purpose disinfectant/cleaner proven to kill TB in 5 minutes, and HBV, HIV, and HCV in 3 minutes
- EPA-registered as a broad-spectrum disinfectant for both surface and immersion use
- Bactericidal, virucidal, tuberculocidal, and fungicidal

24-oz Trigger-Spray Bottle

(101-4114).....ea

1-gal Refill

(101-0285).....ea



HENRY SCHEIN®

MAXIWIPE® Germicidal Cloths

Disinfectant wipe that kills TB, RSV, HBV, HCV, HIV-1, rhinovirus, MRSA, VRE, influenza A (H1N1) virus, in 2 minutes.

- Bactericidal, tuberculocidal, and virucidal
- Compatible with a broad range of surfaces and equipment in healthcare
- Thick, absorbent, material helps remove soil from surfaces without tearing.
- Unique deep-well lid on canister seals securely to prevent moisture loss.
- Meets CDC and OSHA regulations
- Intermediate-level disinfectant

Large

(106-6794).....160/cont

X-Large

(112-5305).....65/cont



HENRY SCHEIN®

MAXIWIPE™ Alcohol-Free Germicidal Disposable Wipes

Formulated with care: Alcohol- and fragrance-free wipes for staff and patients with respiratory sensitivities. Efficient, one-step cleaner and disinfectant. Effective bactericidal, tuberculocidal, and virucidal, including SARS-CoV-2 (the virus that causes COVID-19), in 3 minutes. EPA Registration No. 9480-9, now qualified for List Q*. Intermediate level.

Large, 6" x 6¾"

(570-3056).....160/cont

X-Large, 7½" x 15"

(570-3057).....75/cont

*MaxiWIPE™ Alcohol-Free has demonstrated effectiveness against emerging viral pathogens, such as the monkeypox virus, on hard, nonporous surfaces.

The Premium Package for Pediatric/Back to School.



Surface Disinfection



Hand Hygiene



Screen Cleaning



Floorstand



#1
DISINFECTANT
WIPE IN
HEALTHCARE

Super Sani-Cloth® Wipes - (113-5423)
Sani-Cloth® Plus Wipes - (267-0721)
Sani-Cloth® AF3 Wipes - (267-0039)



Sani-Hands® Wipes - (610-8876)



Easy Screen® Wipes - (125-0248)



Hygiene Floorstand - (P014100)

Learn more by contacting your Henry Schein Rep

*Market leadership claims based on Q1 2024 Clarivate Dataa

©2023 PDI PD103249001



CaviWipes™

Easy as

1 2 3

CaviWipes disinfecting wipes

make surface disinfection easy with quick, effective protection giving you peace of mind.



CaviWipes™ HP

Hydrogen Peroxide Wipe

for alcohol-sensitive equipment and fast turnarounds

Part #	Item Description
(640-0122)	Large (6" x 6.75") 160 count canister
(142-5192)	Extra Large (9" x 12") 65 count canister



CaviWipes™ 2.0

Low Alcohol/Quat Wipe

for everyday surfaces and most medical devices

Part #	Item Description
(140-3081)	Large (6" x 6.75") 160 count canister
(140-3082)	Extra Large (9" x 12") 65 count canister



CaviWipes™ Bleach

Bleach Wipe

for high-risk areas and outbreak control

Part #	Item Description
(129-6430)	Standard (6" x 10.5") 90 count canister



Scan the QR Code
to Learn More

©2024 Metrex Research, LLC. All trademarks are property of Metrex Research, LLC. All Rights Reserved. MKT-24-0578 Rev 0

Metrex™
INFECTION PREVENTION

Visibility and Control



ClearLook®
Lighted Ear Curette™



MADE IN THE USA



InfantScoop®

FlexLoop®

AngleLoop®

MicroLoop®

VersaLoop®



Safer, more effective cerumen removal procedures

For more information, contact your local Henry Schein Medical Representative or visit henryschein.com/medical.



Summit Medical™ ENT Products **Back to School Essentials**



**RhinoGuard™ Disposable
Atomizer Tips**
(106-4157)



**Ventilation
Tubes**
(142-9477 | 142-9338)



**Apdyne Phenol
Applicator Kit**
(143-5289)

Speak to your Henry Schein representative for additional sizes.



Alcohol Prep Pads

- Sterile, 2-ply
- Soft, absorbent, nonwoven pads impregnated with 70% isopropyl alcohol USP


Medium
(570-2219)..... 200/box
Large
(570-1307) 200/box



BZK Towelettes

- Premoistened antiseptic towelettes used for general cleansing for patients and staff
- Alcohol-free
- Unscented
- Cleans and cools simultaneously
- Sting-free cleansing of all minor cuts, scrapes and burns
- BZK is nontoxic and nonirritating

(112-6735).....100/box



Povidone-Iodine Swabsticks and Prep Pads

Sterile
Topical antiseptic used for skin preparation helps to reduce bacteria that can cause infection. Saturated with 10% USP povidone-iodine.
Single Pack, 4"L
(570-1549).....25 pkg/box
Contains: 1 swabstick per pkg.
Triple Pack, 4"L
(570-1548).....25 pkg/box
Contains: 3 swabsticks per pkg.



Ear Curettes

- Disposable
- Color coded
- Individually wrapped
- Non-sterile
- Loop or spoon tip available


Curettes.....50/box
Specify:
Loop, Pink, Small (900-4669)
Loop, White, Large (900-4670)
Spoon, Gray, Small (900-4671)
Spoon, Peach, Large (900-4672)



Disposable Scalpels

- Stainless steel
- Plastic handle with millimeter scale and durable stainless steel blade
- Strong: reduced breakage
- Sharp: clean, true, consistent cut
- Packaging: superior protection and handling
- Individually peel-backed for convenience
- Testing: consistent high quality; assured quality to international standards

Scalpels.....10/box
Specify:
#10 (112-6196) #15 (112-6187)
#11 (112-6189) #15S (112-7296)
#12 (112-6194) #20 (112-6193)



Maxima® Safety Scalpels

- Sterile
- Single use
- Retractable
- Easily extend or retract the blades while avoiding contact with the sharp blade edges
- Color coded for easy identification of blade size

Safety Scalpels..... 10/box
Specify:
Size 10 (570-0243) Size 12 (570-0244)
Size 11 (570-0245) Size 15 (570-0246)



Non-Fenestrated Drapes

Sterile, Polylined, Tissue/Poly/Tissue

- Heat sealed along 2 outer edges to ensure a pinhole-free plastic layer that guarantees maximum protection against bacterial strike-through from below
- Features polyethylene lining between 2 layers of absorbent-strength tissue, virtually eliminating risk of bacterial migration and strike-through
- All sterile fields are available in convenient dispenser boxes of 50 each

White, 18" x 26"
(900-4686).....50/box



Fenestrated Drapes

Sterile, with 3" Round Fenestration
Tissue/Poly/Tissue

- Allows access to the operative site while it protects surrounding area
- 3" round fenestration option is die cut and heat sealed, leaving a smooth, even edge that prevents linting
- Not made with natural rubber latex

White, 18" x 26"
(900-4685).....50/box



Surgical Skin Markers

These surgical marking pens provide an accurate and reliable way to mark the incision site prior to surgery.

- Prep-resistant
- Nonabrasive tip
- Large ink reservoir

Non-Sterile, Fine Tip, Ruler, and Up to 9 Labels
(900-7018)..... 100/box
Sterile, Fine Tip, Ruler, and Up to 6 Labels
(900-7019).....50/box
Mini (WriteSite Jr.), Dark Blue
(900-7020)..... 100/box



Povidone-Iodine Solution

Sterile
10% PVP solution. Used for patient preoperative skin preparation.
1-oz Pouches
(141-5340).....250/case



Povidone-Iodine Swabsticks

Sterile, 1s
10% PVP solution. Used for patient preoperative skin preparation.
4" Swabsticks
(141-5338).....1000/case



Butterfly Wound Closure Strips

Sterile
Absorbent, nonadherent pads that protect and cushion wounds without sticking. Individually wrapped in easy-to-peel pouches.
 $\frac{3}{8}$ " x $1\frac{13}{16}$ "
(681-3633).....100/box

Sklar
A Legacy of Excellence

When Quality Matters Most

Sklar Econo™ Sterile single-use instruments offer both convenience and safety, with proven effectiveness in reducing Healthcare-Associated Infections (HAIs) and minimizing cross-contamination risks.

Please contact your local Henry Schein representative for samples and more product information.



Integra® Miltex® Products Limit Uncertainty



Exam Room Instruments



Scissors

[953-8705]	Operating Scissors, 5¾", straight, sharp/blunt points
[953-1865]	Iris Scissors, 4⅞", straight, fine, sharp points
[953-3771]	Iris Scissors, 4½", curved, fine, sharp points
[953-4908]	Universal Bandage Scissors, 7½", serrated blade, black handle
[953-2248]	Spenser Stitch Scissors, 3½", delicate
[127-5355]	Miltex Wire Cutting Scissors, 4¾", angled to side, one serrated blade

Forceps

[953-7280]	Adson Dressing Forceps, 4¾", delicate, serrated
[953-8252]	Adson Tissue Forceps, 4¾", delicate, 1x2 teeth
[953-2940]	Brown-Adson Tissue Forceps, 4¾", straight, 7x7 teeth
[953-2346]	Halsted Mosquito Forceps, 4⅞", straight
[953-3372]	Halsted Mosquito Forceps, 4¾", curved
[953-9662]	Swiss Cilia & Suture Forceps, 3¾", 2.8mm wide, slanted smooth jaws

Needle Holders

[953-2929]	Webster Needle Holder, 5¼", smooth jaws
[953-8903]	Halsey Needle Holder, 5¼", smooth jaws
[953-0667]	Halsey Needle Holder, 5¼", smooth jaws, Tungsten Carbide
[953-8985]	Olsen-Hegar Needle Holder combined w/ Suture Scissors, 4¾", serrated jaws, extra delicate, Tungsten Carbide
[953-5741]	Mayo-Hegar Needle Holder, 5¼", serrated jaws

Exam Room Specialties

[953-8998]	Taylor Percussion Hammer, 7½", standard size, solid handle, chrome
[953-2331]	Disposable Scalpel, safety sleeve, Size 10, 10/box
[100-0729]	Tuning Fork, aluminum alloy, C-256 vibrations
[100-2323]	Tuning Fork, aluminum alloy, C-512 vibrations
[953-5856]	Disposable Ear Curette, ring tip, pink, 50/box
[953-9486]	Nail Nipper, 5½", concave jaws, double spring, stainless steel
[100-4520]	Vantage Wartenberg Pinwheel
[365-0042]	Vantage Finger Ring Cutter, 6¼"

To place an order, please contact:

HENRY SCHEIN®
MEDICAL

www.henryschein.com/medical

Integra, the Integra logo and Miltex are registered trademarks of Integra LifeSciences Corporation or its subsidiaries in the United States and/or other countries. ©2022 Integra LifeSciences Corporation. All rights reserved. 0505707-1-EN 2022-01



HENRY SCHEIN®

Surgical Blades

- Sterile
- Individual foil packs
- Computer-controlled grinding ensures a consistently sharp edge

Stainless Steel.....100/box
Specify:

#10	(112-6199) #20	(112-6184)
#11	(112-6190) #21	(112-6185)
#12	(112-6177) #22	(112-6195)
#15	(112-6178) #23	(112-6200)
#15S	(112-7294) #25	(112-6201)

Carbon Steel.....10/pkg
Specify:

#10	(112-6183) #12	(112-6176)
#11	(112-6175) #15	(112-6181)



HENRY SCHEIN®

Maxima® Lister Bandage Scissors

German Stainless Steel, Extra Fine

4 1/2"	(104-8240).....ea
7 1/4"	(104-7364).....ea
8", Large Ring	(104-5841).....ea



HENRY SCHEIN®

Maxima® Iris Scissors

German Stainless Steel

3 1/2", Curved, Delicate Pattern	(102-3439).....ea
3 1/2", Straight Delicate Pattern	(102-5248).....ea
4", Curved	(101-8111).....ea
4", Straight	(101-6893).....ea
4 1/2", Straight	(100-8118).....ea



HENRY SCHEIN®

Maxima® Littauer Jr. Stitch Scissors

German Stainless Steel, Slender Pattern

4 1/2", Straight	(100-1420).....ea
------------------	-------------------



HENRY SCHEIN®

Maxima® Operating Room Scissors

German Stainless Steel, Straight, Blunt Points

4 3/4", Sharp	(102-3440).....ea
---------------	-------------------



HENRY SCHEIN®

Maxima® Spencer Stitch Scissors

German Stainless Steel

3 1/2"	(100-6024).....ea
--------	-------------------



HENRY SCHEIN®

Vienna Nasal Speculum

Standard Steel, Standard Pattern

Small	(100-0773).....ea
Medium	(100-1290).....ea
Large	(100-9171).....ea



HENRY SCHEIN®

Maxima® Adson Tissue Forceps

German Stainless Steel

4 3/4"	(101-6279).....ea
4 3/4", 1 x 2 Teeth	(101-4803).....ea




HENRY SCHEIN®

Hartman-Noyes Alligator Forceps

Standard Steel, Serrated Jaws

3 1/2"	(100-7328).....ea
5 1/2"	(100-0886).....ea
6 1/2"	(100-6937).....ea



HENRY SCHEIN®

Hartman-Noyes Ear Forceps

Standard Stainless Steel, Alligator Type, Serrated Jaws

3 1/4"	(104-5508).....ea
--------	-------------------



HENRY SCHEIN®

Disposable Mirror

High-quality plane-surface patient mirror made of white, nonautoclavable ABS plastic.

(100-6640).....72/box



HENRY SCHEIN®

Maxima® Taylor Percussion Hammer 7"

German Stainless Steel, Solid Handle, Pediatric Size

(104-8462).....ea



HENRY SCHEIN®

Plain Splinter Forceps 3 1/2"

Standard Stainless Steel

(100-3799).....ea

Actimove®

Delta-Cast®

AquaCast®
LINER

ORTHO-GLASS®

Your Trusted Partner In Orthopedic Solutions for Positive Care Outcomes

In Orthopedics, Essity provides a range of therapeutic solutions, clinical expertise and world class brands. This means one partner to support your needs.

Essity has a full line of:

- quality sports tapes and bandages
- braces and supports
- casting and splinting

Excited to now offer Coach®, Zonas® and, Elastikon® all by Actimove®



Reformable, Adjustable Bracing System



©2021 DJO, LLC

MKT1501-11263 Rev A

Formed for
Kids
upper extremity
bracing

Exos® Pediatric Short Arm Fracture Brace Open Thumb

HSI PART #	DJO® PART #	DESCRIPTION	SIZING
(118-5299)	311-11-1111	PED SAFB-OT-LEFT, 6" LENGTH, 5.75"-<6.5" CIRCUM.	3XS (PEDS)
(118-5301)	311-21-1111	PED SAFB-OT-LEFT, 7" LENGTH, 6.5" - <7.25" CIRCUM.	2XS (PEDS)
(126-3229)	311-31-1111	PED SAFB-OT-LEFT, 8" LENGTH, 6.5" - <7.25" CIRCUM.	1XS (PEDS)
(118-5300)	311-12-1111	PED SAFB-OT-RIGHT, 6" LENGTH, 5.75"-6.5" CIRCUM.	3XS (PEDS)
(118-5302)	311-22-1111	PED SAFB-OT-RIGHT, 7" LENGTH, 6.5"-<7.25" CIRCUM.	2XS (PEDS)
(126-3228)	311-32-1111	PED SAFB-OT-RIGHT, 8" LENGTH, 6.5"-<7.25" CIRCUM.	1XS (PEDS)

Black part numbers shown.
Other colors are available.

For more information or to place an order,
please contact your Henry Schein Representative
Phone 800-772-4346 / Fax 800-329-9109
or visit: www.henryschein.com/medical



Individual results may vary. Neither DJO, LLC, nor any of its subsidiaries dispense medical advice. The contents of this AD do not constitute medical, legal, or any other type of professional advice. Contact your physician.

*BOA® is a registered trademark of BOA Technology, Inc.

California customers please use website when ordering for Prop 65 information.

3M Medical Securement Solutions

Drive efficiency in your clinical practice and reduce tape SKUs.

A choice that's critical to care. At Solventum, formerly 3M Health Care, we've made it easy to choose the right securement for each patient and application, with four classes of securement that can work harder and smarter.



Non-adhesive



3M™ Coban™ NL Non-Latex Self-Adherent Wrap†

Non-adhesive securement for clinical applications such as blood draws, dressings, immobilization, securement for difficult to dress areas (head, fingers, toes) and support, and mild compression for soft tissue injuries (i.e. strains, sprains).*



Multi-purpose



3M™ Micropore™ S Surgical Tape†

Multi-purpose securement for clinical applications such as blood draws, light-weight dressings, IV lines and tubing (secondary securement), and non-critical tubes.



Flexible



3M™ Medipore™ H Soft Cloth Surgical Tape

Flexible securement for clinical application such as dressing and added pressure, IV lines and tubing (secondary securement), when swelling or movement is anticipated, central venous catheters (secondary securement), chest tubes, and surgical drain tubes.



High-adhesion



3M™ Durapore™ Surgical Tape

High-strength securement for clinical applications such as urinary catheters, nasogastric tubes, orogastric tubes, and patient positioning.



Contact your Henry Schein representative or visit [3M.com/MedicalSecurement](https://www.3M.com/MedicalSecurement) to learn more.

†Skin performance: Designates products that deliver the securement power you need while minimizing damage to skin.

*May use as a component in a compression wrap system support only under supervision of a wound care specialist.

© Solventum 2024. All rights reserved. Solventum and the S logo are trademarks of Solventum or its affiliates. 3M and the 3M logo are trademarks of 3M. Other trademarks are the property of their respective owners.



Steri-Strip™ Skin Closures

Sterile

Uses a pressure-sensitive, hypoallergenic adhesive. Adhesive style is breathable, made of microporous adhesive, and reinforced with filaments for added strength. It is appropriate for low-tension wounds. Contains iodophor: ideal for use when additional antimicrobial activity is desired.

- Less scarring, enhanced patient comfort, and less chance of infection than with sutures or staples
- Fast application
- Reduces likelihood of skin irritation
- Use to close lacerations, surgical incisions, and following early suture/staple removal
- Provides wound support and assists in increasing tensile strength of wound

#R1541, ¼" x 3"	
(777-6538).....	150/box
Contains: 3 strips per envelope, 50 envelopes per box.	
#R1542, ¼" x 1½"	
(777-0843).....	300/box
Contains: 6 strips per envelope, 50 envelopes per box.	
#R1546, Adhesive, ¼" x 4"	
(777-8250).....	500/box
Contains: 10 strips per envelope, 50 envelopes per box.	
#R1547, Adhesive, ½" x 4"	
(777-1358).....	300/box
Contains: 6 strips per envelope, 50 envelopes per box.	
#R1548, Adhesive, 1" x 5"	
(777-8945).....	100/box
Contains: 4 strips per envelope, 25 envelopes per box.	



Soft Flexible Fabric Bandages

- Highly absorbent, nonstick pad that protects the wound
- Prefolded wrapper tabs make it simple to open and apply
- Long-lasting adhesive gently sticks to the skin
- Sterile
- Not made with natural rubber latex

Strip, ¾" x 3"	
(127-1251).....	100/box
Strip, 1" x 3"	
(127-1254).....	100/box
Spot, ¾"	
(127-1247).....	100/box
Oval, ¾" x 1¼"	
(127-1252).....	100/box
Toe Shield, 1½" x 2¾"	
(127-1253).....	100/box
Knuckle, 1½" x 3"	
(127-1246).....	100/box
Small Fingertip, 1¾" x 2"	
(127-1243).....	100/box
Large Fingertip, 1¾" x 3"	
(127-1244).....	50/box
Patch, 1½" x 2"	
(127-1248).....	100/box
Patch, 2" x 3"	
(127-1357).....	50/box
4-Wing, 2¾" x 3"	
(127-1255).....	50/box



3M™ Micropore™ Surgical Tape

Gentle adhesion to skin holds well on damp skin for secure placement. Porous and breathable design for added comfort.

- Conformable
- Highly breathable
- Hypoallergenic: not made with natural rubber latex

Single-Use Rolls

#1530S-1, 1", 100/box	
(777-6928).....	box
Standard Rolls	
10 yd, White.....	box
Specify:	
#1530-0, ½", 24/box	(777-5044)
#1530-1, 1", 12/box	(777-6080)
#1530-2, 2", 6/box	(777-9382)
#1530-3, 3", 4/box	(777-9909)
10 yd, Tan.....	box
Specify:	
#1533-0, ½", 24/box	(777-4451)
#1533-1, 1", 12/box	(777-5133)
#1533-2, 2", 6/box	(777-5407)



Designer Adhesive Stat Strip Bandages

- Used to cover and protect minor cuts and wounds
- Sterile and not made with natural rubber latex
- Long-lasting and highly absorbent
- Features fun designs

Strips, ¾" x 3".....	100/box
Specify:	
Emoji	(127-5610)
Glitter Stars & Stripes	(127-1245)
Blue & Pink Camouflage	(127-1294)
Glitter	(127-1262)
Herbie the Dinosaur	(127-1336)
Pink	(127-1300)
Neon	(127-1338)



Stat Strip® Bandages

Breathable bandage with a highly absorbent, nonstick pad that protects the wound as it heals. Specifically designed to use with gloves.

Flexible Fabric Strips

¾" x 3"	
(127-1279).....	100/box
1" x 3"	
(127-1280).....	100/box
Sheer Strips	
¾" x 3"	
(127-1277).....	100/box
1" x 3"	
(127-1278).....	100/box



3M™ Transpore™ Surgical Tape

Transparent, perforated plastic tape that tears easily, works reliably.

- Easy to handle with gloves
- Good adhesion to skin and tubing
- Transparent
- Porous
- Water-resistant
- Hypoallergenic: not made with natural rubber latex
- 10 yd, Transparent

#1527-0, ½", 24/box	
(777-9735).....	box
#1527-1, 1" Width, 12/box	
(777-7305).....	box
#1527-2, 2", 6/box	
(777-9422).....	box
#1527-3, 3", 4/box	
(635-1724).....	box



Pediatric Character Adhesive Bandages

- To cover minor cuts, scrapes, and wounds
- Sterile and not made with natural rubber latex
- Long-lasting adhesive
- Nonstick wound pad
- Highly breathable

Bandages, ¾" x 3".....	100/box
Specify:	
Strawberry Shortcake	(127-1320)
Care Bears	(127-1322)
Transformers	(127-1316)
Candy Land	(127-1319)
Tonka	(127-1318)
My Little Pony	(127-1317)
Scooby Doo	(127-1257)



Super Hero Stat Strip Bandages, ¾" x 3"

Strip Bandages.....	100/box
Specify:	
Spider-Man	(127-1324)
Iron Man & Captain America	(127-1325)
Hulk & Thor	(127-1326)
Superman, Wonder Woman & Flash	(127-1311)
Batman, Green Lantern & Aquaman	(127-1312)



Herbie
(127-1336)




Neon
(127-1338)



Farm Friends
(145-8472)



Emoji
(127-5610)



Matey's
(146-0988)



Planet & Stars
(127-1282)



New modern designs with Dukal's trusted Stat Strip® Technology

Save time and decrease infection risk by reducing potential skin contact with Dukal's easy-to-open technology. Dukal's Stat Strip® Adhesive Bandages' specifically designed wrapper makes hygienic application as easy as grip, pull, apply – even with gloves.

› To learn more visit dukal.info/design-series





HENRY SCHEIN® Animal Adhesive Bandages

Sterile

- Exclusive design for Henry Schein
- Individually wrapped adhesive bandage is sterilized for maximum protection of minor cuts, punctures, and abrasions
- Nonstick, highly absorbent pad will not stick to wound
- Air vents let air in and keep moisture out
- Not made with natural rubber latex

Spots, 7/8"..... 100/box
(900-4506)..... 100/box
3/4" x 3"..... 100/box
(900-4507)..... 100/box

*Made in the USA



HENRY SCHEIN® Crayon Adhesive Bandages

Sterile

- Individually wrapped adhesive bandage is sterilized for maximum protection of minor cuts, punctures, and abrasions
- Pad will not stick to wounds
- Permits wound to drain
- Not made with natural rubber latex

3/4" x 3"..... 100/box
(112-6995)..... 100/box



HENRY SCHEIN® Flexible Fabric Elastic Strips

- Stretch
- Conform
- Permit free movement
- Will not stick to wound
- Not made with natural rubber latex

3/4" x 3"..... 100/box
(112-6132)..... 100/box
1" x 3"..... 100/box
(112-6133)..... 100/box
Knuckle, 1 1/2" x 3"..... 100/box
(112-6144)..... 100/box

HENRY SCHEIN® Sheer Plastic Bandages

Sterile

- Individually wrapped adhesive bandages for protecting minor cuts, punctures, and abrasions
- Highly absorbent pads will not stick to wound
- Air vents allow air in and keep moisture out
- Not made with natural rubber latex

3/4" x 3"..... 100/box
(112-6138)..... 100/box
1" x 3"..... 100/box
(112-6142)..... 100/box
7/8" Spots..... 100/box
(112-6143)..... 100/box
Professional Tray Pack, 3/4" x 3"..... 1500/case
(112-6134)..... 1500/case
Professional Tray Pack, 1" x 3"..... 1500/case
(112-6135)..... 1500/case



HENRY SCHEIN® Conforming Stretch Gauze Rolls

4.1-yd Rolls

- Economical knitted bandage in various widths
- Made from high-quality materials: nylon, polyester, and cotton
- Provides light compression as it cushions and protects
- Offers softness, breathability, and alleviation of irritation
- Low lint content reduces contamination to wound site
- Improved packaging is resealable, durable, easy to store, and remains clean
- Not made with natural rubber latex

Non-Sterile

1"..... 12/bag
(104-2704)..... 12/bag
2"..... 12/bag
(104-8825)..... 12/bag
3"..... 12/bag
(104-8997)..... 12/bag
4"..... 12/bag
(104-2728)..... 12/bag
6"..... 6/bag
(104-9424)..... 6/bag
Sterile
2"..... 12/bag
(104-2739)..... 12/bag
3"..... 12/bag
(104-9724)..... 12/bag
4"..... 12/bag
(104-6708)..... 12/bag
6"..... 6/bag
(104-4293)..... 6/bag



HENRY SCHEIN® Premium Hypoallergenic Cloth Surgical Tape

10-yd Rolls

- Reliable adhesion
- Highly porous
- Thin, flexible backing
- Gentle removal
- Not made with natural rubber latex
- Non-sterile

Tape..... box
Specify:
1/2", 24/box (101-3025) 2", 6/box (101-5979)
1", 12/box (101-2323)



HENRY SCHEIN® Premium Surgical Hypoallergenic Paper Tape

10-yd Rolls

- Excellent adhesion
- Pinked, nonfraying edges
- Strong woven backing
- Gentle removal
- Not made with natural rubber latex
- Non-sterile

Tape..... box
Specify:
1/2", 24/box (101-4305) 2", 6/box (101-1332)
1", 12/box (101-9673)



HENRY SCHEIN® Premium Hypoallergenic Transparent Surgical Tape

10-yd Rolls

- Excellent adhesion
- Porous, with flexible backing
- Tears easily in both directions
- Not made with natural rubber latex
- Non-sterile

Tape..... box
Specify:
1/2", 24/box (101-6655) 2", 6/box (101-9226)
1", 12/box (101-4645)



HENRY SCHEIN®

Self-Adherent Bandages

Non-Sterile, 5-yd Rolls

- Self-adherent: sticks to itself, not the skin
- Comfortable porous material: pain-free removal that will not damage skin

Latex

Tan, 1".....ea
(570-0113).....ea

Tan, 2".....ea
(900-4322).....ea

Tan, 3".....ea
(900-4323).....ea

Tan, 4".....ea
(900-4324).....ea

Tan, 6".....12 rolls/pkg
(900-4328).....12 rolls/pkg

Multi, 1".....30 rolls/pkg
(112-7366).....30 rolls/pkg

Multi, 2".....36 rolls/pkg
(900-4331).....36 rolls/pkg

Multi, 3".....24 rolls/pkg
(900-4332).....24 rolls/pkg

Non-Latex

Tan, 1".....ea
(112-7355).....ea

Tan, 2".....ea
(112-7357).....ea

Tan, 3".....ea
(112-7358).....ea

Tan, 4".....ea
(112-7359).....ea

Tan, 6".....ea
(112-7362).....ea

Multi, 1".....30 rolls/pkg
(112-7363).....30 rolls/pkg

Multi, 2".....30 rolls/pkg
(112-7364).....30 rolls/pkg

Multi, 3".....24 rolls/pkg
(112-7365).....24 rolls/pkg



HENRY SCHEIN®

Unna's Boot

10-yd Rolls

- Cotton gauze formulated with a zinc-oxide emulsion
- Used for treating venous leg ulcers, edema, and dermatitis
- Soothes inflammation and rashes
- Not made with natural rubber latex

3".....ea
(100-3444).....ea

4".....ea
(100-3256).....ea

3" with Calamine.....ea
(900-7455).....ea

4" with Calamine.....ea
(900-7456).....ea



HENRY SCHEIN®

Elastic Bandages with Velcro®

Non-sterile, 4½-yd Rolls

- Elastic bandage with Velcro® brand fasteners
- Provides even compression with strength, comfort, and breathability
- Premium honeycomb-knitted material
- Not made with natural rubber latex

2".....10/box
(900-4681).....10/box

3".....10/box
(900-4682).....10/box

4".....10/box
(900-4683).....10/box

6".....10/box
(900-4684).....10/box



HENRY SCHEIN®

Elastic Bandages

Non-Sterile, 4½-yd Rolls

- High cotton content to keep patients cool and comfortable
- Durable bandages can be laundered, sterilized, and reused
- Ideal for support and compression applications, such as treatment of varicosities, sprains, strains, and torn or injured ligaments
- Not made with natural rubber latex

2".....10 rolls/pkg
(900-4677).....10 rolls/pkg

3".....10 rolls/box
(900-4678).....10 rolls/box

4".....10 rolls/box
(900-4679).....10 rolls/box

6".....10 rolls/box
(900-4680).....10 rolls/box



HENRY SCHEIN®

Packing Strips

5-yd Rolls

- Sterile
- 100% cotton
- Nonfraying
- Nonlinting
- Selvage edge
- Absorbs fluid; less chance of leaving foreign matter in the wound
- Not made with natural rubber latex

¼".....ea
(104-7976).....ea

½".....ea
(104-6304).....ea

1".....ea
(104-8748).....ea

2".....ea
(104-6302).....ea

Packing Strip with Iodoform

¼".....ea
(112-5821).....ea

½".....ea
(112-5820).....ea



HENRY SCHEIN®

Gauze Pads

Highly absorbent. Can be used for general cleaning, debriding, and prepping. 100% woven cotton. Sterile 1's in peel-back packages. 12-ply.

- Highly absorbent
- 100% woven cotton
- Can be used for general cleaning, debriding, and prepping
- Box is perforated at bottom for easier gauze removal

2" x 2".....100/pkg
(570-2659).....100/pkg

3" x 3".....100/pkg
(570-2629).....100/pkg

4" x 4".....100/pkg
(570-2649).....100/pkg



HENRY SCHEIN® All-Gauze Sponges

Soft, pliable, and absorbent. Available in 3 sizes in both 8-ply and 12-ply styles. Non-sterile.

- Precision sizes and folding: consistency within functions
- Soft and pliable: nonabrasive touch
- Highly absorbent: soaks up unwanted fluid
- Can be used for general cleaning, debriding, and prepping

8 Ply, 2" x 2"	
(570-2657).....	200/pkg
12 Ply, 2" x 2"	
(570-2658).....	200/pkg
8 Ply, 3" x 3"	
(570-2650).....	200/pkg
12 Ply, 3" x 3"	
(570-2648).....	200/pkg
8 Ply, 4" x 4"	
(570-2656).....	200/pkg
12 Ply, 4" x 4"	
(570-2630).....	200/pkg



HENRY SCHEIN® Premium All-Gauze Sponges

Manufactured from 100% USP type-VII cotton gauze that is bleached and processed to established standards. Soft and pliable. Protects as a wound covering as it provides maximum absorbency.

Non-Sterile	
8 Ply, 2" x 2"	
(570-2647).....	200/pkg
12 Ply, 2" x 2"	
(570-2628).....	200/pkg
12 Ply, 3" x 3"	
(570-2640).....	200/pkg
8 Ply, 4" x 4"	
(570-2646).....	200/pkg
12 Ply, 4" x 4"	
(570-2645).....	200/pkg
Sterile	
8 Ply, 2" x 2"	
(570-2642).....	100/pkg
12 Ply, 3" x 3"	
(570-2631).....	80/pkg
8 Ply, 4" x 4"	
(570-2644).....	50/pkg
12 Ply, 4" x 4"	
(570-2637).....	50/pkg

HENRY SCHEIN® Premium Non-Woven Sponges

Non-sterile, 4 Ply

- General-use, non-woven rayon/polyester-blend material
- Quarter-fold design
- Less linting: minimizes lint contamination
- Soft and pliable: prevents irritation to mucosa
- Greater absorbency: reduces number of sponges required
- Closed weave
- Not made with natural rubber latex

2" x 2"	
(101-9278).....	200/pkg
3" x 3"	
(101-8748).....	200/pkg
4" x 4"	
(101-4336).....	200/pkg



HENRY SCHEIN® CRITERION™ Non-Woven Sponges

Non-sterile, 4 Ply

- High-quality, excellent construction, rayon/poly-blend sponge that is soft and gentle to touch and very comfortable against skin
- Thicker and more absorbent than standard sponges
- Multiuse sponge can be used for all applications of wound-care treatments
- Easy-to-open, peel-back packaging
- Not made with natural rubber latex

2" x 2"	
(104-2930).....	200/pkg
3" x 3"	
(104-2410).....	200/pkg
4" x 4"	
(104-2411).....	200/pkg

HENRY SCHEIN® CRITERION™ Non-Woven Sponges

Sterile 2's, 4 Ply

- High quality
- Excellently constructed
- Rayon/poly-blend that is soft and gentle to the touch
- High-quality, excellently constructed rayon/poly-blend sponge that is soft and gentle to the touch and very comfortable against skin
- Thicker and more absorbent than standard sponges
- Multiuse sponge that can be used for all applications of wound-care treatments
- Easy-to-open, peel-back packaging
- Not made with natural rubber latex

2" x 2"	
(104-7897).....	50/pkg
3" x 3"	
(104-4132).....	50/pkg
4" x 4"	
(104-8665).....	50/pkg



HENRY SCHEIN® Exodontia Sponges

100% cotton gauze sponges with maximal absorbency cotton inner liner. Superior and consistent quality. Available in 3 sizes. 8-ply.

Non-sterile—Filled

2" x 2"	
(570-2636).....	5000/case
3" x 3"	
(570-2602).....	4000/case
4" x 4"	
(570-2651).....	2000/case
Sterile—Filled	
2" x 2"	
(570-2653).....	200/pkg



HENRY SCHEIN® Premium Exodontia Sponges

Manufactured from all USP 100% cotton. Soft, absorbent, and pure white. Highest quality. 8-ply. Non-Sterile—Filled

2" x 2"	
(570-2635).....	5000/case
Sterile—Filled	
(570-2652).....	200/pkg



HENRY SCHEIN® Syngaue 50™

Non-sterile, 4 Ply

- General-use nonwoven material with a 70/30% blend of rayon and polyester. 47–50% more absorbent than 12-ply gauze
- Absorbs 10 times its weight: greater absorbency reduces the number of sponges required
- Soft and pliable: prevents irritation to mucosa; allows shaping for a better fit
- Sterilizable: lint-free mesh minimizes lint contamination
- Nonraveling: less mess
- Not made with natural rubber latex

2" x 2"	
(100-2125).....	200/pkg
3" x 3"	
(101-9668).....	200/pkg
4" x 4"	
(101-4583).....	200/pkg
(107-9502).....	2000/case



HENRY SCHEIN® Type VII Gauze

- Sterile
- Not made with natural rubber latex
- 4" x 4"

8 Ply	
(570-2655).....	10/pkg
12 Ply	
(570-2633).....	10/pkg
16 Ply	
(570-2638).....	10/pkg
12 Ply, X-ray Detectable	
(570-2632).....	10/pkg
16 Ply, X-ray Detectable	
(570-2634).....	10/pkg





SUPERIOR POINT

Needles are honed with a unique multiedged point to provide greater ease of penetration.

SUPERIOR ATTACHMENT

Needles are drilled for a seamless needle-suture juncture, providing a smoother transition with less tissue disruption.

SUPERIOR PACKAGING

Unique double-sterile barrier packaging ends costly waste.



NON-ABSORBABLE SUTURES 12/Box

MATERIAL	NEEDLE	LENGTH	SIZE				
			6-0	5-0	4-0	3-0	2-0
Silk Black Braided	C-9 24mm 1/2 Circle	18"/45cm				632G* (100-6320)	
	C-43 15mm 1/2 Circle	18"/45cm				(100-7762)	
	C-3 13mm 3/8 Circle	18"/45cm		640G* (101-4148)	641G* (101-2590)		
	C-6 18mm 3/8 Circle	18"/45cm		682G* (101-3207)	683G* (100-6830)	1679H* (100-7772)	
	C-7 24mm 3/8 Circle	18"/45cm				684G* (101-2049)	
	C-22 16mm 3/8 Circle	18"/45cm	(101-2636)				
Nylon-Black Monofilament		18"/45cm	697G* (900-7482)				
	C-3 13mm 3/8 Circle	18"/45cm	1698G* (900-7473)	698G* (900-7470)	699G* (900-7466)		
	C-6 18mm 3/8 Circle	18"/45cm		661G* (101-7137)	662G* (101-1209)		
		10"/25cm			(900-7464)		
	C-7 24mm 3/8 Circle	18"/45cm				663G* (101-3728)	
	C-8 26mm 3/8 Circle	18"/45cm					664G* (900-7476)
	C-17 12mm 3/8 Circle	10"/25cm	(900-7460)	(900-7459)	(900-7461)		
	C-22 16mm 3/8 Circle	18"/45cm	1665G* (101-6053)				
	PC-31 18mm 3/8 Circle	18"/45cm		1666G* (900-7472)	1667G* (900-7487)		
	DSM-19 18mm 3/8 Circle	18"/45cm				(570-0479)	
Polypropylene-Blue Monofilament	DSM-24 24mm 3/8 Circle	18"/45cm				(570-0480)	
	DSCM-16 16mm 3/8 Precision	18"/45cm		(570-0458)	(570-0459)		
	C-1 11mm 3/8 Circle	18"/45cm	8697G* (900-7483)				
	C-3 13mm 3/8 Circle	18"/45cm	8695G* (900-7479)	8698G* (900-7484)	8699G* (900-7486)		
	C-6 18mm 3/8 Circle	18"/45cm		8661G* (101-6409)	8683G* (101-7477)		
	C-7 24mm 3/8 Circle	18"/45cm				8684G* (101-7700)	
	C-22 16mm 3/8 Circle	18"/45cm	(101-4378)				
	DS-18 18mm 3/8 Circle	18"/45cm				8665G* (900-7481)	
	PC-31 18mm 3/8 Circle	18"/45cm			8682G* (900-7478)		
	DSM-16 16mm 3/8 Circle	18"/45cm		(570-0471)			
	DSM-19 18mm 3/8 Circle	18"/45cm		(570-0468)		(570-0467)	
	DSM-24 24mm 3/8 Circle	18"/45cm				(570-0470)	
Polyester- Green Braided	DSM-26 16mm 3/8 Circle	18"/45cm					(570-0473)
	DSCM-13 16mm 3/8 Precision	18"/45cm	(570-0472)				
	DSCM-16 16mm 3/8 Precision	18"/45cm		(570-0469)			
	D-14 17mm 1/2 Circle	18"/45cm		X951H* (100-7773)			

*Comparable Ethicon™ material

* - Ethicon Equivalent

ABSORBABLE SUTURES 12/Box

MATERIAL	NEEDLE	LENGTH	SIZE				
			6-0	5-0	4-0	3-0	2-0
Plain Gut	C-9 24mm 1/2 Circle	18"/45cm				612G* (100-6120)	
	C-3 13mm 3/8 Circle	18"/45cm		686G* (101-3192)			
	C-6 18mm 3/8 Circle	18"/45cm			(100-7771)		
	C-6 18mm 3/8 Circle	27"/70cm			H821H* (101-0668)	H822H* (101-2976)	
	DSCM-13 13mm 3/8 circle	18"/45cm	1916G* (900-7485)	1915G* (900-7471)			
Chronic Gut	C-3 13mm 3/8 Circle	18"/45cm		687G* (101-8824)			
	C-6 18mm 3/8 Circle	18"/45cm			1637G* (900-7637)	1638H* (900-7636)	
	C-6 18mm 3/8 Circle	27"/70cm			635H* (101-2787)	636H* (101-1582)	
	HR-26 26mm 1/2 Circle	30"/76cm				G122H* (570-0466)	
	DSM-13 13mm 3/8 Circle	18"/45cm			1654G* (570-0475)		
	DSM-19 19mm 3/8 Circle	18"/45cm			1637G* (570-0476)		
Polyglycolic Acid (PGA) Undyed-Braided	C-3 13mm 3/8 Circle	18"/45cm	J492G* (900-7468)	J493G* (900-4417)	J494G* (900-4418)		
	C-6 18mm 3/8 Circle	27"/70cm			J422H* (900-7463)	J423H* (900-7462)	
	PC-30 18mm 1/2 Circle	18"/45cm			J507G* (900-7477)		
	PC-31 18mm 3/8 Circle	18"/45cm			J496G* (900-4419)	J497G* (900-4421)	
(PGA) Undyed-Braided	HR-26 26mm 1/2 Circle	27"/45cm		(570-0481)			
Polyglycolic Acid (PGA) Undyed-Braided-Fast Absorbing	DSM-13 13mm 3/8 Circle	18"/45cm		VR493* (570-0477)	VR494* (570-0478)		
Polyglycolic Acid (PGA) Violet	C-6 18mm 3/8 Circle	18"/45cm			J392H* (900-4420)		
	DSM-19 19mm 3/8 Circle	18"/45cm				J393H* (570-0474)	
PGA PCL Undyed-Mono	DSM-13 13mm 3/8 Circle	18"/45cm		Y493G* (900-7467)	Y494G* (900-7465)		
	DSM-19 19mm 3/8 Circle	18"/45cm			Y496G* (900-7488)	Y497G* (900-7474)	
	DSM-13 13mm 3/8 Circle	18"/45cm	Y492G* (570-0465)				
	DSM-19 18mm 3/8 Circle	27"/45cm			(570-0460)	(570-0461)	
	DSM-24 18mm 3/8 Circle					(570-0432)	
PGA PCL Violet-Mono	DSM-13 13mm 3/8 Circle	18"/45cm		Y463G* (570-0462)	(570-0464)		


*Comparable Ethicon™ material

* - Ethicon Equivalent


NEEDLE CHART

REVERSE CUTTING NEEDLES

3/8 Circle




1/2 Circle



PRECISION REVERSE CUTTING NEEDLES

3/8 Circle



PRECISION CONVENTIONAL CUTTING NEEDLES

3/8 Circle



TAPER CUTTING NEEDLES

1/2 Circle





HENRY SCHEIN®

Topical Skin Adhesive

Sterile, 0.5-mL Single-Patient-Use Plastic Ampules
The proven choice of wound closure with over 40 years of clinical experience, this FDA-approved liquid adhesive consists of n-butyl-2 cyanoacrylate.

- Highest bond strength of all topical tissue adhesives on the market
- Will not clog or dry out when opened
- Provides a precise, controlled application
- Inhibits growth of bacteria
- Available in dyed and clear versions
- Flexible version provides more flexible coverage for longer incisions

Ampule Without Tip.....	ea
Specify:	
Clear	(112-6106)
Blue	(112-6111)
Flexible Ampule with Tip, Blue	
(570-0351).....	ea
Flexible, with Tip, Blue	
(112-7158).....	10/box



HENRY SCHEIN®

Disposable Skin Stapler

- Enables the user to work with higher efficiency and assists in avoiding cross infection effectively
- Specially designed for surgical skin suturing
- Novel structure and application security allow for quick speed when sewing up during operation, modest histological reaction, and good wound-matching and healing
- No pain during removal
- Scarring is minimal

Skin Stapler, 35 Wide	
(570-2956).....	ea
Skin Stapler, 35 Regular	
(570-2957).....	ea
Skin Staple Remover	
(900-7011).....	ea

HENRY SCHEIN®

Wound Closure Strips




- Adhesive wound closures with an exceptionally porous, spun-bond fiber backing include bilateral nylon filaments for high tensile strength
- To be used in addition to or in place of sutures
- Each size is in a color-coded package for quick and easy size identification
- Medical-grade adhesive properties: Affords excellent adhesion, resistance to moisture, superior longevity, and is gentle to the skin
- Hypoallergenic: Designed for patient comfort and longer-term exposure
- Not made with natural rubber latex

1/8" x 3"	
(900-4317).....	250/box
Contains: 5 strips per envelope, 50 envelopes per box.	
1/4" x 3"	
(900-4315).....	150/box
Contains: 3 strips per envelope, 50 envelopes per box.	
1/4" x 4"	
(900-4318).....	500/box
Contains: 10 strips per envelope, 50 envelopes per box.	
1/4" x 1 1/2"	
(900-4316).....	300/box
Contains: 6 strips per envelope, 50 envelopes per box.	
1/2" x 4"	
(900-4314).....	300/box
Contains: 6 strips per envelope, 50 envelopes per box.	
1" x 4"	
(900-7506).....	100/box
Contains: 4 strips per envelope, 25 envelopes per box.	




Blue = Common Sizes Used



30 X 45
Right Fit

23 Gallon Slim Jim Container



33 X 39
Right Fit

32 Gallon Container



40 X 46
Right Fit

44 Gallon Container



38 X 58
Right Fit

55 Gallon Container

"Right-sizing" Your Can Liner Adds Up!

Why **"Super-Size"** when you can **"Right-Size"** and save money & contribute to source reductions?



For more information, contact your local Henry Schein Medical Representative or visit henryschein.com/medical



AP medical
your safety choice

Sharps Container

- 4½"H x 12"W x 10¾"L
- Body: 1.8 mm; Lid: 1.53 mm
- Counterbalance lid
- Irreversible permanent closure
- Maximum filling level indicated
- Resistant to sharp and cutting objects
- Red, 5.4 qt

#242468
(139-8802).....ea

Phlebotomy Collector

One-step system eliminates recapping. Rigid, puncture-resistant, polypropylene container. May be autoclaved and compacted or incinerated according to local regulations.

#305487, 1.5 qt
(987-5812).....ea



Sharps Containers

- Durable and puncture-resistant
- Horizontal drop maximizes container volume
- Locks for final disposal
- Transparent lid helps reduce risk of overfilling
- Red

Slide Lid
#4623, 2 qt
(139-8161).....ea
Mailbox Lid
#4624, 5 qt
(139-8162).....ea



Sharpstainer® Disposable Containers

Red, autoclavable, and puncture-resistant. Tamper-proof locking cap. One-way safety insert.

Small, 0.7 qt, 3½" x 3½" x 5"
(153-3569).....ea
Tall Tray, 1.7 qt
(153-8085).....ea
10 qt
(120-4364).....ea

HENRY SCHEIN®
BRAND PRODUCTS

Rely on Us™

Sharps Containers

From a Name You Can Trust.

- For use by professional health-care givers
- Meets all NIOSH and OSHA regulations
- FDA approved
- Puncture-resistant
- Transparent body and/or top allows for easy viewing
- High-impact and shatter-resistant polypropylene



#0310-1500-HS, 1 qt (112-7149).....ea	#0354-150M-HS, 5 qt (112-7152).....ea	#0320-150R-HS, 2 gal, Rotary Lid (570-0377).....ea
#0310-150T-HS, 1.5 qt (570-0376).....ea	#0319-1500-HS, 1 gal, Slide Lid (570-0378).....ea	#0330-150R-HS, 3 gal (570-0212).....ea
#0322-150R-HS, 2.2 qt (570-0213).....ea	#0319-150R-HS, 1 gal, Rotary Lid (570-0379).....ea	#0380-150R-HS, 8 gal (570-0380).....ea
#0354-150B-HS, 5.4 qt (112-7151).....ea	#0320-1500-HS, 2 gal, Slide Lid (112-7150).....ea	



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!



Ready, set, vaccinate

with BD hypodermic solutions

It's simple: the more immunity you can give your patients, the better. That means having enough of the right vaccines on hand *and* the high-quality needles and syringes you need to deliver them. BD hypodermic solutions are designed to deliver consistent quality and reliable performance.

BD SafetyGlide® Needles are designed to help keep clinicians safe from accidental needlestick injuries and feature Activation-Assist® Technology—making safe injections easier and more efficient.



To purchase your
BD vaccine devices,
scan this QR code.



BD, Franklin Lakes, NJ, 07417, U.S.

bd.com

BD, the BD Logo, Activation-Assist and SafetyGlide are trademarks of Becton, Dickinson and Company or its affiliates.
© 2024 BD. All rights reserved. BD-120838 (04/24)





Low-Dead-Space Syringes with Needles

22 ga x 1 1/2"

Low-dead-space syringe eliminates the hub, since needle is attached directly to syringe barrel. This configuration eliminates pharmaceutical wastage because virtually all of the medication is delivered. 1 1/2" needle allows all pharmaceuticals to be accessed without inverting the vial. Large, permanent markings are easy to read. Nontoxic and latex-free.

1cc

(606-1100)..... 100/box



Hypodermic Disposable Needles

SUS-304 Stainless Steel

Sterile, Nontoxic, Pyrogen-Free, Non-Latex
Needles and boxes are color coded for easy identification. Each needle is individually packaged, sterilized, and ready for use.

Needles..... 100/box

Specify:

#26403, Dark Blue, 25 ga x 5/8"	(337-9855)
#26405, Dark Blue, 25 ga x 1"	(337-5760)
#26406, Dark Blue, 25 ga x 1 1/2"	(337-2640)
#26407, Light Blue, 23 ga x 3/4"	(337-1604)
#26408, Light Blue, 23 ga x 1"	(337-2432)
#26411, Black, 22 ga x 1"	(337-1881)
#26412, Black, 22 ga x 1 1/2"	(337-3130)



VanishPoint® Safety Needles and Syringes

- Needle retracts into the syringe when the plunger is fully depressed or after the injection
- Easy, one-handed activation
- Lubricated needle
- Requires less disposal space in sharps containers
- Regular wall needle with tribeveled tip

U-100 Insulin Syringes

0.5 mL

#15291, 30 ga x 3/16"	
(137-1779).....	100/box
#10161, 25 ga x 1"	
(701-5933).....	100/box
#10331, 22 ga x 1"	
(701-9859).....	100/box
#10351, 21 ga x 1"	
(125-4798).....	100/box



The Sharp Choice

MYCO Sterile Sharps Supplies

Reli® Wealy® Retractable Safety Syringe

- OSHA-compliant safety-engineered sharp
- Passive activation of the needle retraction safety feature eliminates capping step
- Easy, single-handed operation
- Low dead space for 1 mL/cc safety syringes
- Sterile

#SS1ML25G101

(140-6632)..... 100/box



Building a Healthier Tomorrow


SOL-CARE™ Luer Lock Safety Syringe with Exchangeable Needle

- Highly transparent barrel with bold, precise scale for perfect visualization and accurate dosage
- Integrated plunger stopper prevents plunger from being pulled out, avoiding exposure to hazardous fluids

3-mL Syringe

#100078IM, 25 ga x 5/8", Orange	
(261-0258).....	100/box
#100079IM, 25 ga x 1", Orange	
(139-6144).....	100/box
#100077IM, 23 ga x 1", Blue	
(261-0260).....	100/box

“NOT ALL SAFETY NEEDLES ARE CREATED EQUAL”



USE **REAL SAFETY™**
SYRINGES THAT
SAVE MONEY
TIME & LIVES

FROM



RETRACTABLE
TECHNOLOGIES, INC.

www.retractable.com



Doing Well by Doing Good



Visit www.HenrySchein.com/HSCares
to find out more about the Calendar of Caring Programs.



Thank you for helping us raise more than \$3 million
to help health happen!



Adrenalin® Injection

1 mg/mL, 1-mL SDV
(248-0883).....ea
1 mg/mL, 1-mL SDV
(890-0257).....25/pkg
1 mg/mL, 30-mL MDV
(129-4847).....ea



Albuterol Inhalation Solution,

0.042%, 3 mL
(116-4177).....30/box



Albuterol Sulfate Inhaler HFA

Compare to ProAir®
90 mcg
(140-9653).....ea



Amoxicillin Suspension b.i.d.

Suspension b.i.d., 400/5 mL
(277-0717).....100 mL/btl



Aplisol® Injectables


Tuberculin PPD Stabilized Solution, 5 TU/0.1 mL
1-mL Vial
(110-5199).....10/pkg
5-mL Vial
(110-5099).....50/pkg
CPT® Code: 86580*
Item stored under refrigeration. May be shipped separately.



kaléo

Auvi-Q Epinephrine Auto-Injectors

0.1 mg
(142-4090).....2/pkg
0.15 mg
(142-4091).....2/pkg
0.3 mg
(142-4092).....2/pkg



Bicillin® C-R Injectable* Prefilled Syringes

Penicillin G Benzathine and Penicillin G Procaine,
IM Suspension
Pediatric, 21 ga x 1"
900/300 U/mL, 2 mL
(321-6738).....10/pkg
Pediatric, 21 ga x 1"
1,200,000 U/mL, 2 mL
(321-2437).....10/pkg
Adult, 21 ga x 1½"
1,200,000 U/mL, 2 mL
(321-7671).....10/pkg
*Tubex® Injector is not required.



Bicillin® L-A Injectable* Prefilled Syringes

Penicillin G Benzathine IM Suspension
Adult, 21 ga x 1½"
1,200,000 U/mL, 2 mL
(321-7445).....10/pkg
Adult, 18 ga x 1½"
2,400,000 U/mL, 4 mL
(321-1663).....10/pkg
*Tubex injector is not required.



Budesonide Inhalation Suspension

Compare to Pulmicort®
0.5 mg/2 mL
(143-5369).....ea



Ceftriaxone Sodium Injectable

Compare to Rocephin
500-mg SDV
(135-9645).....10/box
1-gram Vial
(133-6436).....10/box
2-gram Vial
(135-7800).....10/box



Ceftriaxone Sodium Injectable

Compare to Rocephin®
250-mg SDV
(133-6543).....ea
250-mg SDV
(131-4702).....10/box
500-mg SDV
(133-6542).....ea



Celestone® Soluspan® Injectable Suspension

6 mg/mL, 5-mL MDV

(140-7856).....ea



Cyanocobalamin (Vitamin B-12) Injectable

Compare to Rubramin PC®

1000 mcg/mL, 1-mL MDV

(244-1902).....25/pkg



Cyanocobalamin Injectable (B-12)

1000 mcg/mL

1-mL SDV

(141-8752).....10/box

10-mL MDV

(129-3812).....25/box

30-mL MDV

(129-3810).....25/box



Depo-Medrol® Injectables

Methylprednisolone Acetate

20 mg/mL, 5-mL MDV

(908-3069).....ea

40 mg/mL, 1-mL SDV

(908-5362).....ea

40 mg/mL, 1-mL SDV

(231-1203).....25/box

40 mg/mL, 5-mL MDV

(908-9783).....ea

40 mg/mL, 5-mL MDV

(908-8955).....25/box

40 mg/mL, 10-mL MDV

(908-1188).....ea

40 mg/mL, 10-mL MDV

(908-1177).....25/box

80 mg/mL, 1-mL SDV

(908-7936).....ea

80 mg/mL, 1-mL SDV

(908-8623).....25/box

80 mg/mL, 5-mL MDV

(908-3787).....ea

80 mg/mL, 5-mL MDV

(908-7913).....25/box



Dexamethasone Sodium Phosphate Injection

Compare to Decadron®

4 mg/mL, 30-mL MDV

(248-0894).....ea

4 mg/mL, 30-mL MDV

(140-1494).....25/box



Dexamethasone Sodium Phosphate Injection

4 mg/mL, 5-mL MDV

(375-1975).....25/pkg




Dexamethasone Sodium Phosphate Injectable

Compare to Decadron®

Preservative-free

10 mg/mL, 1-mL SDV

(406-7616).....25/pkg



Dexamethasone Sodium Phosphate

10 mg/mL, 10-mL MDV

(139-6050).....10/box



Dexamethasone Sodium Phosphate Injection

10 mg/mL, 10-mL MDV

(248-0875).....ea



Dexamethasone Sodium Phosphate Injectable

Compare to Decadron®

10 mg/mL, 1-mL SDV

(138-6758).....25/pkg





Dexamethasone Tablets

Compare to Decadron®

0.5 mg	
(277-0039).....	100/btl
0.75 mg	
(277-0525).....	100/btl
1.5 mg	
(277-0049).....	100/btl
4 mg	
(277-0050).....	100/btl



Epinephrine Auto-Injector

Junior, 0.15 mg	
(138-1212).....	2/pkg
Adult, 0.3 mg	
(143-6151).....	2/pkg



Epinephrine Auto Injectors

Authorized generic of Adrenaclick®.

0.15 mg	
(127-6483).....	2/pkg
0.3 mg	
(127-9954).....	2/pkg



Epinephrine Injectable Ampules

Compare to Adrenalin®

1:1000, 1 mL	
(248-0899).....	ea
1:1000, 1 mL	
(128-9991).....	10/box



Kenalog®-40 Injectables

Triamcinolone Acetonide

40 mg/mL	
1-mL SDV	
(196-8300).....	ea
5-mL MDV	
(196-1737).....	ea
10-mL MDV	
(196-9429).....	ea



Kenalog®-80 Injection

80 mg/mL	
1-mL Vial	
(136-4157).....	ea
5-mL Vial	
(136-4158).....	ea



FOSUN PHARMA

Ketorolac Injection

15 mg/mL, 1-mL SDV	
(146-1679).....	25/box
30 mg/mL, 1-mL SDV	
(258-0457).....	ea
30 mg/mL, 1-mL SDV	
(143-6324).....	25/box
30 mg/mL, 2-mL SDV	
(258-0456).....	ea
30 mg/mL, 2-mL SDV	
(143-6293).....	25/box



Levalbuterol Inhalation Solution

0.31 mg	
(138-1259).....	30/box
0.63 mg	
(138-1150).....	30/box
1.25 mg	
(133-9591).....	30/box



Medroxyprogesterone

Compare to Depo-Provera®


150 mg/mL, 1-mL SDV	
(146-6392).....	ea
(146-6382).....	25/box



Methylprednisolone Acetate Injection

Compare to Depo-Medrol®

40 mg/mL, 5-mL MDV	
(145-7924).....	ea
40 mg/mL, 10-mL MDV	
(143-8564).....	ea
80 mg/mL, 5-mL MDV	
(142-2720).....	ea



Methylprednisolone Acetate Injection

Compare to Depo-Medrol®

1-mL SDV	
40 mg/mL	
(140-9021).....	ea
40 mg/mL	
(140-9026).....	25/box
80 mg/mL	
(140-9048).....	25/box



Silver Nitrate Applicators, 6" Wood

- Applications: skin or mucous membrane cauterization; removing warts and granulated tissue
 - 75% silver nitrate; 25% potassium nitrate
 - 50 mg active ingredients/applicator
 - Silver nitrate tips are soluble in all secretions
 - Control strength by dipping in sterile, distilled water
 - Rigid wooden applicator for precision placement
 - Convenient tube vial
- (112-6994)..... 100/vial



Solu-Cortef® Injectables

Preservative-Free

100 mg, 2-mL SDV	
(908-0011).....	ea
100 mg, 2-mL SDV	
(908-0012).....	25/pkg
250 mg, 2-mL SDV	
(908-0013).....	ea
250 mg, 2-mL SDV	
(114-8910).....	25/pkg
500 mg, 4-mL SDV	
(908-0014).....	ea
1000 mg, 8-mL SDV	
(908-0015).....	ea
100-mg, 2-mL SDV	
(835-5169).....	ea



sanofi

Tubersol®

Tuberculin PPD (Mantoux), 5 TU/0.1 mL

1-mL Vial	
(546-9191).....	10/pkg

CPT® Code: 86580†



Triamcinolone Acetonide Injection

Compare to Kenalog®-40

40 mg/mL	
1-mL SDV	
(142-2167).....	ea
(142-3667).....	25/box
5-mL MDV	
(142-3668).....	ea
10-mL MDV	
(142-2227).....	ea



Ventolin HFA Inhaler

90 mcg	
8 g	
(131-5535).....	ea
18 g	
(228-3204).....	ea



Solu-Medrol® Injectables

Act-O-Vial®, Preservative-Free

40 mg, 1-mL Vial	
(382-0014).....	25/pkg
125 mg, 2-mL Vial	
(382-0015).....	25/pkg
500 mg, 4-mL Vial	
(382-0016).....	ea
1 g, 8-mL Vial	
(382-0017).....	ea



Pain Ease® Mist Spray

Topical Anesthetic Skin Refrigerant

3.9-oz Mist Spray (791-0634).....	ea
1-oz Medium-Stream Spray (132-9467).....	ea
3.9-oz Medium-Stream Spray (761-0029).....	ea

Can be used in facilities that restrict use of flammable components.



Ethyl Chloride Spray

Topical Anesthetic Skin Refrigerant

• Fast: Swab the target area with antiseptic. Spray 4–7 seconds immediately prior to injections, IV start, or venipuncture	
• Low cost: pennies per application	
Fine-Stream Spray, 3.9 oz (138-9019).....	ea
Fine Spray–3.9 oz. (138-9016).....	4/pkg
Medium-Stream Spray, 3.9 oz (138-5777).....	ea
Medium-Stream Spray, 3.9 oz (138-5778).....	4/pkg



Ethyl Chloride Spray

Topical Anesthetic Skin Refrigerant

• Fast: Swab the target area with antiseptic. Spray 4–7 seconds immediately prior to injections, IV start, or venipuncture	
• Low cost: pennies per application	
Mist Spray Can, 3.9 oz (138-6406).....	ea
(138-6360).....	4/pkg
(138-6330).....	12/case

Gebauer's
PainEase®
The INSTANT topical anesthetic



Flu season is getting cool this
year with Pain Ease

Visit www.Gebauer.com/HSPE to learn more



4444 East 153rd. Street • Cleveland, Ohio 44128

sanofi

ActHIB®

Haemophilus b Ped Conjugate (Tetanus Toxoid Conjugate) Vaccine, Thimerosal-Free

0.5-mL SDV with Diluent

(546-0024)..... 5/pkg

CPT® Code: 90648†n



sanofi

Adacel®

Tetanus Toxoid, Reduced Diphtheria Toxoid, and Acellular Pertussis Vaccine Adsorbed

0.5-mL SDV

(546-7789)..... 10/pkg

0.5-mL Prefilled Syringe Without Needle

Not made with natural rubber latex.

(132-5841)..... 5/pkg

CPT® Code: 90715†



GSK


Bexsero®

Meningococcal Group B Vaccine

0.5-mL Prefilled Syringe

(127-2151)..... 10/pkg

CPT® Code: 90620†



GSK

Boostrix®

Tetanus Toxoid, Reduced Diphtheria, and Acellular Pertussis Adsorbed Vaccine

0.5-mL SDV

(254-6475)..... 10/pkg

0.5-mL Tip-Lok® Prefilled Syringe

(254-0031)..... 10/pkg

CPT® Code: 90715†



sanofi

DAPTACEL®

Diphtheria and Tetanus Toxoids and Acellular Pertussis (DTaP) Vaccine (Adsorbed), Thimerosal-Free

0.5-mL SDV, Pediatric

(546-1136)..... 10/pkg

CPT® Code: 90700†



GSK

Engerix-B®—Pediatric/Adolescent

Hepatitis B Recombinant Vaccine

10 µg, 0.5-mL Syringe

(254-0026)..... 10/pkg

CPT® Code: 90744†



MERCK

Gardasil® 9

Human Papillomavirus 9-Valent Vaccine, Recombinant

0.5-mL Prefilled Luer-lock Syringe

(558-0044)..... 10/pkg

CPT® Code: 90651†



GSK

HAVRIX®—Pediatric/Adolescent

Hepatitis A (Inactivated) Vaccine, Thimerosal-Free, 720 ELU/0.5 mL

(254-0023)..... 10/pkg

CPT® Code: 90633†





GSK

Hiberix®

Haemophilus b Conjugate Vaccine (Tetanus Toxoid Conjugate) for injection

With Diluent, 0.5-mL SDV

(146-1811).....10/box

CPT® Code: 90648†



sanofi

IPOL®

Poliovirus (Salk) (Inactivated) Vaccine

5-mL, 10-Dose Vial

(546-7189).....ea

CPT® Code: 90713†

No credit will be given for return of this product.



GSK

Kinrix™

Diphtheria and Tetanus Toxoids and Acellular Pertussis Adsorbed and Inactivated Poliovirus Vaccine

0.5-mL Prefilled TipLok Syringe

(254-0025).....5/pkg

CPT Code: 90696†



GSK

MENVEO®

Meningococcal (Groups A, C, Y, and W-135) Oligosaccharide Diphtheria CRM197 Conjugate Vaccine Solution for Intramuscular Injection

0.5-mL SDV

(127-9616).....5/pkg

CPT Code: 90734†



sanofi

MenQuadfi™


Meningococcal (Groups A, C, Y, W) Conjugate Vaccine Solution for Intramuscular Injection

0.5-mL SDV

(140-3989).....5/carton

CPT® Code: 90619†

Item stored under refrigeration. May be shipped separately.



MERCK

M-M-R® II

Measles, Mumps, and Rubella Vaccine

0.5-mL SDV with Diluent

(558-0110).....10/pkg

CPT® Code: 90707†

M-M-R® II is a registered trademark of Merck & Co., Inc.



GSK

INFANRIX®

Diphtheria and Tetanus Toxoids and Acellular Pertussis (DTaP) Vaccine (Adsorbed), Thimerosal-Free

0.5-mL Prefilled Syringe

(124-0010).....10/pkg

CPT® Code: 90700†



GSK

Menveo One-Vial

Meningococcal (Groups A, C, Y, and W-135) Oligosaccharide Diphtheria CRM197 Conjugate Vaccine

0.5-mL SDV

(145-4145).....10/pkg

CPT® Code: 90734†



GSK

Pediarix®

Diphtheria and Tetanus Toxoids and Acellular Pertussis (DTaP) (Adsorbed), Hepatitis B (Recombinant), and Inactivated Poliovirus Vaccine (Combined)

0.5-mL Tip-Lok® Prefilled Syringe

(254-0028).....10/pkg

CPT® Code: 90723†



PedvaxHIB®

Haemophilus b Conjugate (Meningococcal Protein Conjugate) Liquid Vaccine

0.5-mL SDV
(558-3778).....10/pkg

CPT® Code: 90647†

PedvaxHIB® is a registered trademark of Merck & Co., Inc.



Pneumovax® 23

Pneumococcal Vaccine (Polyvalent)

0.5-mL Prefilled Syringe
(122-2087).....10/pkg

CPT® Code: 90732†

Pneumovax® 23 is a registered trademark of Merck & Co., Inc.




sanofi

Pentacel®

Diphtheria and Tetanus Toxoids and Acellular Pertussis Adsorbed, Inactivated Poliovirus and Haemophilus b Conjugate (Tetanus Toxoid Conjugate) Vaccine

0.5-mL SDV, Thimerosal-Free
(140-7444).....5/pkg

CPT® Code: 90698†




Prenar 20™

Pneumococcal 20-Valent Conjugate Vaccine

0.5-mL Prefilled Syringe
(140-8153).....ea

0.5-mL Prefilled Syringe
(140-8152).....10/box

CPT® Code: 90677†



Priorix

Measles, Mumps, and Rubella Vaccine, Live

0.5-mL SDV
(142-3773).....10/box

CPT® Code: 90707†



ProQuad®

Measles, Mumps, Rubella, and Varicella Virus Vaccine Live

0.5-mL SDV
(558-0053).....10/pkg

CPT Code 90701†

ProQuad® is a Registered Trademark of Merck & Co., Inc.



sanofi


Quadracel™

Diphtheria and Tetanus Toxoids and Acellular Pertussis Adsorbed and Inactivated Poliovirus Vaccine

0.5-mL SDV
(141-7656).....10/pkg

0.5-mL Prefilled Syringe
(141-7655).....10/pkg

CPT® Code: 90696†



RECOMBIVAX HB® – Pediatric/Adolescent

Hepatitis B Vaccine, Thimerosal-Free

5 µg/mL, 0.5-mL SDV
(558-6254).....10/pkg

0.5-mL Prefilled Luer-Lock Syringe without Safety Needle
(122-2084).....10/pkg

CPT® Code: 90744*

RECOMBIVAX HB® is a registered trademark of Merck & Co., Inc.





GSK

Rotarix®
Rotavirus Vaccine, Live, Oral
1.5-mL Pediatric Prefilled Syringe
(143-4023)..... 10/pkg
CPT® Code: 90681†




MERCK

RotaTaq®
Rotavirus Vaccine (live), Oral, Pentavalent,
2-mL Single-Dose Tube
(558-8763)..... 10/pkg
(558-0009)..... 25/pkg
CPT® Code: 90680†
RotaTaq® is a registered trademark of Merck & Co., Inc.



GRIFOLS

TDVAX™
Tetanus and Diphtheria Toxoids, Adsorbed
7+ years. Not made with natural rubber latex.
Preservative-free.
0.5-mL SDV
(125-9839)..... 10/pkg
CPT® Code: 90714†




sanofi

Tenivac™
(Tetanus and Diphtheria Toxoids Adsorbed)
0.5-mL SDV
(546-0016)..... 10/pkg
0.5-mL Prefilled Syringe
(546-0017)..... 10/pkg
CPT® Code: 90714†
Includes Federal Excise Tax (FET).



Pfizer

Trumenba®
Meningococcal Group B Vaccine
0.5-mL Prefilled Syringe
(248-0649)..... 5/pkg
(248-0648)..... 10/pkg
CPT® Code: 90621†n



MERCK

Vaqta® Pediatric/Adolescent
Hepatitis A (Inactivated) Vaccine, Thimerosal-Free,
25 U
0.5-mL Prefilled Luer-lock Syringe without Safety
Needle
(122-2080)..... 10/pkg
CPT® Code: 90633
VAQTA® is a registered trademark of Merck & Co., Inc.



MERCK

VARIVAX®
Varicella (Chicken Pox) Virus (Live) Vaccine,
Thimerosal-Free, Lyophilized
0.5-mL SDV
(558-1592)..... 10/box
CPT® Code: 90716†
VARIVAX® is a registered trademark of Merck & Co., Inc.



sanofi

Vaxelis™
Diphtheria and Tetanus Toxoids and Acellular
Pertussis, Inactivated Poliovirus, Haemophilus b
Conjugate and Hepatitis B Vaccine
0.5-mL SDV
(140-3160)..... 10/pkg
0.5-mL Prefilled Syringe
(140-3161)..... 10/pkg
CPT® Code: 90697†



MERCK

Vaxneuvance™
Pneumococcal 15-Valent Conjugate Vaccine
0.5-mL Prefilled Syringe
(141-1001)..... ea
(141-1002)..... 10/box
CPT® Code: 90671†
Vaxneuvance™ is a trademark of Merck & Co., Inc.

HENRY SCHEIN®
**A&D Ointment,
Unit Dose**

Helps treat chafed skin associated with diaper rash and helps seal out wetness. Temporarily protects minor cuts, scrapes, and burns, and temporarily protects and relieves chapped and cracked skin and lips. Also helps prevent and protect from the drying effects of wind and cold weather.

5 gm
(570-3126).....144/box
Contains: 144–5-gm packets per box.



HENRY SCHEIN®
**APAP Extra Strength Tablets–
Unit Dose**

Compare to Tylenol® Extra Strength

- Manufactured and packaged in the USA
- Highest quality-control standards under the FDA
- Packaged in 250 sealed-in pouches of 2 tablets each
- Convenient dispenser box
- Color-coded inner and outer packaging

500 mg
(900-4690).....500/box
Contains: 250 pouches of 2 tablets.



HENRY SCHEIN®
**APAP Tablets,
Unit Dose**

Compare to Tylenol®
Nonaspirin pain reliever and fever reducer.

- 2 tablets per pouch
- Convenient dispenser box
- Made in the USA

325 mg
(900-4689).....500/box
Contains: 250 pouches of 2 tablets.



HENRY SCHEIN®
Bacitracin Zinc Ointment, USP

1-oz Tube
(570-3129).....ea
0.9-gm Unit-Dose Packs
(570-3125).....144/box



HENRY SCHEIN®
Ear Wash System

A safe and effective option for the removal of earwax (cerumen).



HENRY SCHEIN®
Hydrocortisone Cream 1%, USP

Compare to Active Ingredient in Cortaid®
Itch relief for minor skin irritations, inflammation, and rashes. Active ingredient: hydrocortisone, USP, 10 mg.

1-oz Tube
(570-3124).....ea



HENRY SCHEIN®
**Hydrogen Peroxide Solution,
USP 3%**

8-oz Bottle
(112-7012).....ea
16-oz Bottle
(112-7069).....ea
1-gal Bottle
(112-7070).....ea

Features and benefits:

- Flexible hose and rigid tube options include the bottle of earwash fluid, trigger mechanism, and 3 tips
- Accessories such as an ear wash basin, additional tips, and splash shields for the rigid tube unit are sold separately*
- Ear wash basin is ergonomically molded to fit under the ear
- Simply fill the bottle with warm water, secure the trigger mechanism onto the bottle, twist a disposable tip into the end of the nozzle, and press the trigger to pump the water through the nozzle into the ear

Ear Wash Kit with Flexible Hose and 3 Tips
(570-1884).....ea

Ear Wash Kit with Rigid Tube and 3 Tips
(570-1885).....ea

Ear Wash Tips
(570-1883).....20/bag

Plastic Ear Wash Basin
(570-1903).....ea

Ear Wash Splash Shields for Rigid Tube
(570-1917).....10/bag


*There is no option for Ear Wash Splash Shields for the Flexible Hose unit.



HENRY SCHEIN®
Isopropyl Alcohol ♦

70% Isopropyl
16-oz Bottle
(112-7067).....ea
32-oz Bottle
(112-7011).....ea
1-gal Bottle
(112-7068).....ea

♦ Hazardous material



HENRY SCHEIN®
Triple Antibiotic Ointment, USP

Compare to Active Ingredients in Neosporin®

1-oz Tube
(570-3130).....ea
0.9-gm Unit-Dose Packs
(570-3128).....144/box



Castile® Soap Towelettes

- Ideal for midstream clean-catch urine specimen procedures
- Gentle cleansing agent, pH balanced
- Contains 2% coconut oil-based soap solution

2" x 2.25"..... 100/box
(610-0148)..... 100/box



MEDI-FIRST®

Cherry Cough Drops UD

Compare to Halls®

For the temporary relief of coughing and dry throats. Individually wrapped.

(839-3443)..... 125/box

Contains: 125 packets of 1.

Active ingredient: menthol, 6.1 mg and eucalyptus.



Debrox® Earwax Removal Kit

Provides a safe, nonirritating method of removing earwax through the power of microfoam action. Active ingredient: carbomide peroxide 6.5%.

Kit..... ea
(131-5909)..... ea

Contains: soft rubber bulb ear syringe & ½-oz bottle of Debrox.

½-oz Bottle Refill

(135-5624)..... ea



Little Remedies Cough Syrup

Honey, 4-oz Bottle

(133-6107)..... ea



Little Remedies® Saline Spray & Drops

1-oz Bottle

(133-6108)..... ea



Medique®

Medi-First® Honey Lemon Cough Drops

Menthol/eucalyptus vapor-action cough drops that clear nasal passages and provide temporary relief of minor throat irritation and cough due to colds or inhaled irritants. Individually wrapped.

#84025

(136-3464)..... 125/box



Nice'n Clean Baby Wipes, Unscented

(113-4964)..... 40/box
(602-5044)..... 80/box



Pedia-Lax® Enema for Children

Gives children the same fast-acting relief as a regular Fleet® enema in an easy, child-size dose.

2.25 fl oz
(248-0570)..... ea

GSK

Robitussin® Honey Cough & Chest Congestion DM, Maximum Strength

4-oz Bottle

(132-6215)..... ea



GSK

Robitussin® Cough+Chest Congestion DM, Non-Drowsy

Temporarily relieves coughs due to minor throat and bronchial irritations that may occur with a cold. Helps loosen phlegm (mucus) and thin bronchial secretions to drain bronchial tubes. Active ingredients: 20 mg of dextromethorphan HBr cough suppressant and 200 mg of guaifenesin expectorant. For ages 12+.

8-oz Bottle

(132-5718)..... ea

Gentle skincare for baby



Aquaphor® Baby Healing Ointment

- Skin Protectant
- One essential solution for Baby's many skincare needs.

Aquaphor® Baby Wash & Shampoo

- Gently cleanses baby skin and hair without drying.

Aquaphor® Baby Diaper Rash Cream

- 15% Zinc Oxide Skin Protectant
- Helps treat and prevent diaper rash

Beiersdorf Description	Size	Format	UPC	HSI Item Code
Aquaphor Baby Healing Ointment	7 oz.	Tube	0-72140-01946-4	124-7620
Aquaphor Baby Wash & Shampoo	16.9 fl. oz.	Pump Bottle	0-72140-02109-2	128-5313
Aquaphor Baby 3 in 1 Diaper Rash Cream	3.5 oz.	Tube in Carton	0-72140-00250-3	138-1290

Aquaphor®
Baby



There's before, and then there's **Aquaphor**

© 2024 Beiersdorf, Inc.

Trusted By
Parents.
Feared
by Lice



(143-7060)

(143-7055)

Nix®

Heading back
to school?
So are Lice.

#1

Pediatrician
Recommended Brand

Be Prepared! Treat &
NOW Repel with Nix®



Aquaphor
HEALING OINTMENT

Aquaphor® Baby 3 in 1 Diaper Rash Cream

3½-oz Tube
(138-1290).....ea



PrestigeConsumer
HEALTHCARE

Boudreaux's Butt Paste®

4-oz Tube
(131-8864).....ea



PrestigeConsumer
HEALTHCARE

Nix Ultra® Lice Treatment

- Kills lice and their eggs
- Kills super lice
- Fast and effective
- Pesticide-free and non-toxic
- Easy application
- Tough on lice; gentle on skin

Lice Shampoo
(135-5184)..... 4 oz/box
Lice Treatment Liquid
(135-5618)..... 3.4 oz/box



Aquaphor
HEALING OINTMENT

Aquaphor® Healing Ointment

Advanced Therapy
• Dermatologist recommended
• Advanced therapy for dry, cracked, or irritated skin
• Clinically proven to restore smooth, healthy skin
#11047, 0.35-oz Tubes
(709-3153).....2/pkg



PrestigeConsumer
HEALTHCARE

Nix Ultra® Lice Removal Kit

(136-9432).....ea
Contains: 1 of each: Nix Ultra® Lice Removal Comb, 3.4-oz Nix Ultra® solution & 5-oz Nix® Lice Killing Spray.



PrestigeConsumer
HEALTHCARE

Nix® Lice Prevention Daily Leave-In Spray

6-oz Bottle
(143-7055).....ea



Aquaphor
HEALING OINTMENT

Aquaphor® Healing Ointment, Unit Dose

(911-0011).....



PrestigeConsumer
HEALTHCARE

Nix® Lice Crème Rinse

2-oz Bottle
(131-8942).....ea
(131-8943).....2/pkg



PrestigeConsumer
HEALTHCARE

Nix Ultra® Lice Treatment & Prevention Kit

(143-7060).....ea



ACCLEAN®

Acclean® 5% Sodium Fluoride Varnish

5% sodium-fluoride varnish is simple to apply. Premeasured unit-dose cups fit easily into enclosed finger ring, and varnish applies effortlessly with included Microbrush® applicator. This easy-to-use, pearlescent white varnish sets on contact with saliva and dries to a natural tooth color. Formulated with only natural sweeteners and xylitol. Delivers the highest fluoride concentration available (compared to nonvarnish fluorides.)

Unit Dose—0.5 mL.....50/box

Specify:
Bubblegum (112-4869) Mint (112-7169)
Cherry (112-4868) Raspberry (112-5963)

Contains: 50 unit-dose cups, 50 Microbrush® applicators & 1 finger ring.

Bulk Package—0.5 mL.....200/box

Specify:
Bubblegum (112-5999) Mint (112-7170)
Cherry (112-5998) Raspberry (112-6001)



Colgate
Oral Pharmaceuticals, Inc.

Colgate® PreviDent® Varnish 5% Sodium Fluoride (22,600 ppm)

Delivers a smooth, even, and quick application. Convenient unit-dose packages each contain 0.4 mL of varnish and a brush. Dries to a natural enamel color and contains xylitol. Available in great-tasting Caramel Sundae, Raspberry, Mint, and Grape Bubblegum flavors.

Unit Dose—0.4 mL.....50/box

Specify:
Grape (543-0148)
Mint (543-0120)
Raspberry (543-0121)
Caramel Sundae (129-8852)

Bulk Pack—0.4 mL.....150/pkg

Specify:
Mint (543-4811)
Raspberry (543-2267)

*Advance Care Program Pricing Available. Annual commitment: minimum 5 boxes. For full prescribing information, go to www.colgateprofessional.com.



Medicom

Duraflor® 5% Sodium Fluoride Varnish

Highest concentration of fluoride available. Sets rapidly when it comes to contact with saliva. Clinically proven to provide relief of hypersensitivity. May be used a cavity liner. Sweetened with xylitol. Does not contain tree nuts, peanuts, sesame, corn, shellfish, egg, soy, dairy products, or gluten.

Bubblegum, 10-mL Tube
(107-9307).....ea

Unit Doses
Bubblegum, 0.25 mL
(107-8603).....32/box
(107-8931).....200/box



premier

Enamel Pro® Varnish with ACP

The only varnish formulated to deliver amorphous calcium phosphate (ACP). Desensitizes dentin—it has been shown to diminish hydraulic conductance by 73%.¹ The fluoride and ACP provide a semipermanent occlusion with fluorapatite.² The ACP technology delivers 4X more fluoride uptake into the enamel than the leading brand of transparent varnish.³ Easy-to-use, single-unit doses with an applicator brush.

Unit Dose, 0.25 mL, Bubblegum
(378-0166).....35/box

Contains: 35 single-unit-dose application sticks, instructions & 35 post-treatment patient instructions.

Unit Dose—0.4 mL, Standard Pack35/box

Specify:
Bubblegum (378-0162)
Strawberry Cream (378-0156)
Vanilla Mint (378-0372)

Contains: 35 single-unit-dose application sticks, instructions & 35 post-treatment patient instructions.

Unit Dose—0.4 mL, Bulk Pack200/box

Specify:
Bubblegum (378-0171)
Strawberry Cream (378-0169)
Vanilla Mint (378-0371)

Contains: 200 single-unit-dose application sticks, instructions & 35 post-treatment patient instructions.

¹⁻³ Data on file.

3M

Vanish™ 5% Sodium Fluoride White Varnish with Tri-Calcium Phosphate

White varnish, enhanced with tricalcium phosphate, becomes virtually invisible after application. Saliva-activated formula adheres to dry or moist teeth and will spread after application due to the patented formula, which allows migration to tooth surfaces. Indicated to treat hypersensitive teeth, exposed dentin, and root-surface sensitivity.

Trial Pack—0.5 mL.....50/box

Specify:
Cherry (737-0016) Mint (737-0017)

Melon (737-0018)
Clinic Pack—0.5 mL.....100/box

Specify:
Cherry (737-0019) Mint (737-0020)

Melon (737-0021)
Bulk Pack—0.5 mL.....1000/box

Specify:
Cherry (737-0022)
Melon (737-0024)
Mint (737-0023)
Assorted (500 Cherry & 500 Melon) (737-0025)
Assorted (400 Melon, 400 Cherry & 200 Mint) (737-0026)



PREVENTECH

Vella® 5% Sodium Fluoride Varnish with Xylitol

Features a smooth, thin-layer formula and upgraded applicator brush with more bristles to reduce clumping during application. Provides relief from hypersensitivity and sets quickly when it comes in contact with saliva. Dries clear. Made in the USA.

Unit Dose—0.5 mL.....35 per Box

Specify:
Caramel (124-7651) Spearmint (113-6486)
Bubblegum (896-0377) Strawberry (113-6488)

Melon (115-5242)
Bulk Pack—0.5 mL.....100 per Box

Specify:
Caramel (124-7649) Spearmint (113-6487)
Bubblegum (896-0378) Strawberry (113-6489)
Melon (115-5241)



Space-Saving Cold Storage for Pediatric Care

Trusted by pediatricians nationwide, ABS delivers peace of mind, reliability, and compliance with stringent industry standards

American Biotech Supply (ABS) ensures uncompromised integrity of pediatric supplies with state-of-the-art cold storage solutions. Certified to the NSF/ANSI 456 Standard for vaccine storage, refrigerators and freezers protect pharmaceuticals at optimal temperatures.

Our cutting-edge technology guarantees precise temperature control, prioritizing patient safety, preventing waste, allowing for peak delivery of vaccines.



Unrivalled Precision

Ensures medication efficacy, safeguarding patient safety and reducing wastage.



Space-Saving Design

Optimize facility space, while ensuring reliable storage for vaccines.



Energy Efficient

Largest selection of models with natural refrigerants and Energy Star certification.



Fastest Lead Times

On-hand inventory for immediate fulfillment needs with logistics professionals to manage nationwide freight and installation.



**For more information, contact your
Henry Schein Sales Representative.**

Are you prepared for the back-to-school vaccine demand wave?

Back-to-school routine immunizations require additional inventory of vaccines. Safe and effective vaccine storage is critical in maintaining the efficacy of vaccines and making sure every dose counts.

- GX Solutions refrigerators and freezers feature OptiCool™ cooling technology which offers optimized temperature uniformity, stability, and recovery, essential in ensuring that vaccines are stored in an optimal storage environment.
- NSF/ANSI 456 Vaccine Storage certified - Ensures your equipment includes the features and temperature performance to safeguard critical vaccine inventory.

For more information, contact your Henry Schein Sales Representative.



Helmer
SCIENTIFIC



i.Series® iPR120-GX Pharmacy Refrigerator

20.2-cu-ft Capacity, Glass Door, 2 to 10°C
GX Solutions and the i.C3® Information Center set i.Series® models apart by providing a professional medical-grade refrigerator with advanced monitoring and control. This pharmaceutical-grade refrigerator is designed for storing medications, vaccines, and other pharmaceuticals. It is NSF/ANSI 456 vaccine storage-certified and provides superior stability, uniformity, and recovery. It also has the added benefits of quiet operation, low heat output, and energy savings, all thanks to the advanced OptiCool™ technology.

- Provides continuous temperature monitoring, event logs, temperature graphs, and more
- Designed specifically for the central pharmacy, medication rooms, clean rooms, satellite locations, and outpatient clinics
- 50–65% more efficient than traditional technology
- Quiet operation (47 dB) to create fewer distractions and better working environment
- Audible and visual alarms include power failure, high- and low-temperature alarms, door ajar, condenser temperature, and low battery
- Includes Certificate of Calibration generated with an ISO 17025 reference thermometer
- Energy Star® certified

#5115120-1
(141-8672).....ea




i.Series® iLF105-GX Laboratory Freezer

5.3-cu-ft Capacity, Solid Door, -15 to -30°C
A GX Solutions professional medical-grade, -30°C freezer. Designed for the unique needs of critical health care applications, such as storing samples, reagents, controls, and frozen vaccines. It delivers optimized control in 3 important areas: temperature, noise, and energy management.

- Powered by OptiCool™ technology, this freezer provides improved performance and user experience through energy efficiency, low noise levels, and stable storage environments to safeguard valuable patient samples, frozen vaccines, and more
- 30% more efficient than traditional technology
- Digital microprocessor temperature controller, with built-in alarms, includes minimum/maximum display and reset
- Energy Star® certified
- Powerful i.C3® Information Center is exclusive to i.Series® models and provides constant temperature monitoring, multiple information logs, and security features
- GX Solutions is the first professional medical-grade line of freezers for clinical laboratories, pharmacies, and other health care facilities
- NSF/ANSI 456 Vaccine Storage Standard certified

#5222105-1
(141-9439).....ea



High-Performance Undercounter Medical Freezer

Designed to maintain temperatures between 5°F and -13°F (-15°C and -25°C) with a tight uniformity of $\pm 5^\circ\text{C}$ for storing critical vaccines and biologics, this model is engineered to meet CDC guidelines in a small footprint. Precise temperature control and uniformity protect vaccine and biologic viability whether the unit is loaded or empty.

- Capacity: 5 cubic ft
- 3 coated steel wire shelves: 19.2" x 14.75"
- Temperature control range: 5°F to -13°F (-15°C to -25°C) in 0.1°C increments
- Alarms: high/low temperature
- Refrigerant: 6-600a
- Door swing: right hinge
- Exterior dimensions: 34.1"H x 23.8"W x 23.7"D
- Interior dimensions: 25.6"H x 19.9"W x 17.1"D

#PF-L5181W-PA

(141-5626).....ea




Pharmaceutical Refrigerator/Freezer Combo

Designed specifically for the storage of vaccines and other pharmaceuticals, this purpose-built refrigerator/freezer combo is designed for temperature uniformity and stability that is not commonly found with household or commercial cabinets. Separate refrigerator and freezer compartments provide a wider range of temperature storage points to accommodate a variety of vaccine product cold-chain requirements.

- Shelves: Refrigerator: 3 adjustable tempered glass; Freezer: 2 PE coated wire
- Refrigerant: Refrigerator: SNAP approved HC R-600A; Freezer: HC R-290
- Temperature control: Refrigerator: 36°F to 57°F (2°C to 14°C); Freezer: -31°F to 5°F (-35°C to -15°C)
- Alarms: High/low temperature
- Refrigerator 11.5 cu ft (326 L); Freezer 4.8 cu ft (136 L)

#MPR-N450FH-PA

(136-5777).....ea



Small Medical Refrigerator

This laboratory undercounter refrigerator combines high-performance temperature control with environmentally friendly refrigerants and energy efficiency. Precision temperature control and uniformity prevent the interior from cooling below freezing and risking vaccine and biologic viability, whether the unit is loaded or empty. Engineered to meet the strict regulatory guidelines of the CDC, ANSI, and NSF.

- Exterior dimensions: 34.1"H x 23.5"W x 23.7"D
- Interior dimensions: 25.6"H x 19.9"W x 17.1"D
- Capacity: 5 cubic ft
- Refrigerant: R-600a
- 36° to 46°F (2°C to 8°C)
- High/low temperature alarm

#PR-L5181W-PA

(140-3362).....ea



Life Science Innovator
Since 1966

Precise Temperature Control to Safeguard Pediatric Patients


Pharmaceutical Refrigerators & Freezers

Our refrigerators and freezers are designed for precise temperature control and uniformity to store critical pharmaceuticals, vaccines, enzymes, and other medical products.

Contact your Henry Schein sales representative for more details.



MPR-S300H-PA | 1374201



PR-L5181W-PA | 1403362



PF-L5181W-PA | 1415626



MPR-N250FH-PA | 1405115

PHC Corporation of North America





**PRECISE, STABLE, EFFICIENT
AND WHISPER QUIET**

Accucold is proud to lead the industry with the largest line of cold storage product sizes certified to NSF/ANSI 456 Vaccine Storage standards

— Proven stability!



accucold[®]
Performance Series
Purpose built for medical, pharmacy, and laboratory
Best in class delivery - Largest inventory

OPTIGEN[®] Allergen-Specific IgE Assay



**Test up to 36 allergens with a small
volume of patient blood sample**

Introducing New Panels: Environmental 36 and Food 29

- ✓ Comprehensive pre-configured panels
- ✓ Streamlined inventory management
- ✓ Cost savings

AVAILABLE REAGENT PANELS

#85067 OPTIGEN Food 29 (142-6007)	20/kit
#85068 OPTIGEN Environmental 36 (142-6010)	20/kit
#80109 Control Reagents for OPTIGEN (131-1156)	1/kit

CPT[®] CODE: 86003

© 2024 Minaris Medical America, Inc. All rights reserved.
DOC.NO.1230 Rev.00





OneStep™ + Ultra Flu Test Kit

From A Name You Can Trust.

Simple to use, results in minutes.

Intended for the qualitative detection of influenza type A and type B nucleoprotein antigens directly from nasal and nasopharyngeal swab specimens from patients with signs and symptoms of respiratory infection.

(570-1886)..... 25/box

Contains: 25 test sticks, 25 sterile nasal swabs, 27 extraction buffer vials, 1 influenza A control swab & 1 influenza B control swab.

CPT® Codes: Influenza A, 87804QW*/Influenza B, 87804QW-59†



Features & Benefits:

- CLIA waived
- Results in 10 minutes
- Sensitivity (A/B): 90.3%/88%
- Specificity (A/B): 96.7%/99.2%
- For use with nasal and nasopharyngeal swab specimens
- Made in the USA



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

OneStep+ Strep A Test Strips

Convenient, compact packaging with everything needed to run 25 tests including positive and negative controls.

- CLIA waived
- Sensitivity: 97%; specificity: 95%
- Room-temperature storage
- Canister packaging for ease of multiple testing
- Walk-away procedure
- 25 tests including external negative and positive controls



(107-4050)..... 25/box

(570-0122)..... 50/box

CPT® Code: 87880QW†



OneStep+ Ultra Mono Test

- 3-step process
- 5-minute time to result
- Dipstick test
- Whole blood, serum, or plasma sample type
- CLIA-waived for whole blood
- Moderate complexity for serum and plasma
- Performance (reconciled with EBV serology): 100% sensitivity; 95.9% specificity

(570-1728)..... 25/box

Contains: 25 each of: test sticks in a container, test tubes, transfer pipettes, 50-µL capillary pipettes; and 1 each of: diluent, mono positive control, mono negative control, workstation & package insert.

CPT® Code: 86308QW† (whole blood); 86308† (serum and plasma)

OneStep Pro+ Strep A Cassette Kit

- CLIA-waived
- Easy-to-use flip design; all-in-one extraction, incubation, and test
- Results in 5–10 minutes
- Sensitivity of 96.2% and specificity of 98.7%
- External negative and positive controls included
- Made in the USA



(900-7570)..... 25 tests/kit

CPT® Code: 87880QW†



No instrument, no problem

QuickVue® Rapid Lateral-Flow Tests produce accurate results in minutes, allowing for faster diagnosis and treatment



Easy to interpret

Visually read test that incorporates internal procedural controls and a two-color result to help ensure accuracy. Everything you need is included in the kit.



Easy to implement

Test kits do not require an instrument, which reduces implementation set-up and counter-top space, and removes recurring maintenance and software updates.



QuickVue®

Broad testing menu

The broad testing menu allows you to customize your test-to-treat program based on your organization's needs.

SARS Antigen | Strep A | H. pylori
Influenza A+B | RSV | iFOB



Learn more about QuickVue tests by scanning the QR code, visiting quidelortho.com, calling us at **800.874.1517**, or reaching out to your Henry Schein representative.

New QuidelOrtho branding may not be available in all markets, subject to country-specific regulatory approval. Please confirm with your local commercial team.

© 2024 QuidelOrtho Corporation PR-105741-NA-EN-US-v1



QuickVue® Influenza Tests

An immunoassay for detecting influenza types A and B. Three simple steps with results in 10 minutes or less. Color-coded procedure card helps demonstrate each step of test. Kit contains positive and negative controls, and stores at room temperature.

- Detects both A and B types of influenza
- CLIA waived
- Sensitivity (A/B): 94%/70%
- Specificity (A/B): 90%/97%
- Color-coded procedure card helps demonstrate each step of the test
- No instrumentation needed
- Results in 10 minutes or less, right in the exam room
- Room-temperature storage with up to a 24-month shelf life from date of manufacture
- Made in the USA

#20183, A+B Test, Differentiates Type A or Type B (227-5611)..... 25/box

#20171, Transport Tubes* with 25 Sterile Foam Swabs Inside Plastic Transport Tubes (227-3550)..... 25/box

CPT® Codes: Influenza A, 87804QW*/

Influenza B, 87804-59†



QuickVue In-Line® One-Step Strep A Tests

True one-step test with a unique in-line extraction and rapid procedure that allows testing to occur right at the point of care. Two-color results give you clear readability, as well as proven accuracy and dependable results. Kit stores at room temperature and includes positive and negative external controls.



- Built-in controls for CLIA waived
- Easy to use and easy to read
- 2-color, rapid results in 5 minutes
- Sensitivity: 92%
- Specificity: 99%
- Room-temperature storage
- Made in the USA

#0343

(227-4207)..... 25/box

CPT® Code: 87880QW†



QuickVue+® Strep A Tests

- Built-in controls for CLIA compliance
- Easy test procedure
- ± endpoint
- Built-in workstation
- Sensitivity: 95%
- Specificity: 98%
- Overall agreement: 98%
- Room-temperature storage
- Moderately complex
- Made in the USA

#20122

(489-0075)..... 25/box

CPT® Code: 87880†



QuickVue® Dipstick Strep A Tests

- CLIA waived
- Results in 5 minutes or less
- 2-color endpoint: results are easy to read and interpret
- Sensitivity: 92%
- Specificity: 98%
- Overall agreement: 95%
- Room-temperature storage
- Built-in positive and negative controls
- Made in the USA

#20108

(227-0611)..... 50/box

CPT® Code: 87880QW†





Test for COVID-19*, Flu A and Flu B with one sample

Quickly detect COVID-19* and Influenza A or Influenza B with a single test with clear, reliable, digital results in 15 minutes.

See what's new at BDVeritor.com

* In the USA, this product has not been FDA cleared or approved but has been authorized by the FDA under an Emergency Use Authorization for use by authorized laboratories; use by laboratories certified under the CLIA, 42 U.S.C. §263a, that meet requirements to perform moderate, high, or waived complexity tests. The product is authorized for use at the Point of Care (POC), i.e., in patient care settings operating under a CLIA Certificate of Waiver, Certificate of Compliance, or Certificate of Accreditation.

This product has been authorized only for the detection of proteins from SARS-CoV-2, influenza A and influenza B, not for any other viruses or pathogens; and, in the USA, the emergency use of this product is only authorized for the duration of the declaration that circumstances exist justifying the authorization of emergency use of in vitro diagnostics for detection and/or diagnosis of COVID 19 under Section 564(b)(1) of the Federal Food, Drug and Cosmetic Act, 21 U.S.C. §360bbb-3(b)(1), unless the declaration is terminated or authorization is revoked sooner.

BD Life Sciences, 7 Loveton Circle, Sparks, MD 21152-0999, U.S.
800-638-8663


BD, the BD Logo and BD Veritor are trademarks of Becton, Dickinson and Company or its affiliates.
© 2024 BD. All rights reserved. (BD-88448 3307-US-0424 April 2024)



BD Veritor™ Plus System

- Faster access to test results: Help inform appropriate and targeted treatment decisions
- Improved workflow: The BD Veritor™ Plus Analyzer improves workflow with an Analyze Now or Walk Away mode, providing test results (with a print-enabled function) and aiding in decision-making while the patient is on site
- Improved traceability: The BD Veritor™ InfoScan module can scan 3 different identifiers (operator ID, specimen ID, and test kit/lot) that can be locally recorded, printed, or downloaded, and reduces hands-on time through automated documentation
- Modularity: Provides a unique solution to meet each testing site's preference
- Connectivity: BD Veritor™ InfoSync module provides the same benefits as BD Veritor™ InfoScan. Connects and documents tests results across multiple sites and can integrate with a hospital electronic medical record or laboratory information system
- Advanced particle technology: Enhances sensitivity by using a proprietary process to produce highly stable, modified colloidal metal particles.
- Adaptive read technology: Compensates for background. Nonspecific binding helps improve specificity to reduce false-positive results

#256066, Veritor Plus Analyzer
(126-5527).....ea



BD Veritor™ System RSV CLIA-Waived Kit

The first CLIA-waived RSV test referenced against a higher sensitivity standard than culture. Test device is ready to insert into reader 10 minutes after sample is added. The same sample processed for RSV can also be used for flu A+B. Digital results displayed in 10 seconds.

- Sample type: NP swab sample
- Positive percent agreement (PPA): 81.6%
- Negative percent agreement (NPA): 99.1%
- 6 years and younger

#256038, RSV Test
(121-7553).....30/box
CPT* Code: 87807QW†



BD Veritor™ System Flu A+B Test Kit

The first CLIA-waived flu A+B test referenced against PCR, a higher sensitivity standard than culture. Device is ready to insert into reader 10 minutes after sample is added.

- Positive percent agreement (A/B): 82.1%/74.6%
- Negative percent agreement (A/B): 98.6%/99.6%
- Positive agreement vs. PCR for current visual read rapid tests ranges from 10% to 70%
- Streamlined workflow: Provides a digital result in <11 minutes, with <50 seconds of hands-on time

#256045, Flu A+B Test Kit
(992-0002).....30/box
CPT* Codes: Influenza A: 87804QW† / Influenza B: 87804-59†; Using the Objective and Differentiated Reader System: 87999†



BD Veritor™ System Group A Strep Test

The first CLIA-waived digital immunoassay (DIA) for the rapid detection of Group A strep (GAS) with an instrumented result. Test device is ready to insert into reader 5 minutes after sample is added. Once the test device is inserted in the reader, an objective, digital test result is displayed in 10 seconds.


- Sample type: Throat swab sample
- Sensitivity*: 95.4%
- Specificity*: 95.7%

#256040, Strep A Test
(992-0006).....30/box
CPT* Code: 87880QW†
*Reference method: Bacterial culture; data from package insert.



Meet Our Awesome Family

BV • STREP • FLU A&B • TRICHOMONAS • H. PYLORI • MONO • RSV • HCG • IFOB



Not all products are available in all countries.

© 2023 SEKISUI Diagnostics, LLC. All rights reserved. OSOM® is a registered U.S. trademark of SEKISUI Diagnostics, LLC. Because every result matters™ is a trademark of SEKISUI Diagnostics, LLC.

Today's healthcare environment is challenging. Let our team help you **minimize variability**, **provide standardization options**, and **reduce overall costs** with our rapid diagnostics products and services.

Answers for Healthcare. Awesome...yes, OSOM®.

Contact your Henry Schein Sales Rep for more information or visit [HenrySchein.com/Medical](https://www.henryschein.com/Medical).



SEKISUI
DIAGNOSTICS
Because every result matters™



OSOM® Ultra Strep A Tests

Faster results for faster treatment—statistically equivalent to single-swab culture with results in just 5 minutes.

- Simple to use—Get easy-to-read results in 4 easy steps and no need for additional equipment
- Greater sensitivity—Antibody-labeled color particles coated at 2 separate locations allows for greater sensitivity than conventional single-label technology without sacrificing specificity
- Easy-to-read 2-color results
- Simplify reagent addition—single reagent delivery system in 25-test kit eliminates the need to handle 2 separate reagents
- QC Inside®—2 extra tests included in the 50-test kit at no additional charge to run external quality-control testing for more confidence in results
- Made in the USA

#147

(392-7389)..... 25/box

#149 (Includes 2 additional tests for external QC testing.)

(102-7617)..... 52/box

CPT® Code: 87880QW*

Ask about RGS. Genzyme compares to rigorous Gold Standard.



OSOM® Strep A Test Kit

Achieve fast results for better treatment. Aids in the diagnosis of strep A, preventing the unnecessary use of antibiotics.

- Simple to use: get easy-to-read results in 4 easy steps, with no need for additional equipment
- Trust the results: high sensitivity and specificity with results in just 5 minutes
- Easy-to-read, 2-color results: blue test line and red control line
- Collect and detect: qualitative detection of group A streptococci antigens from throat swabs; can also confirm presumptive group A streptococcal colonies recovered from cultures less than 72 hours old
- QC Inside®: 2 extra tests included at no additional charge to run external quality control testing for more confidence in the results
- Made in the USA

#141

(224-7320)..... 52/box

CPT® Code: 87880QW*



Sofia® 2 Platform

The Sofia 2 analyzer represents the next generation of rapid diagnostic testing

Proven lateral-flow technology, proprietary advanced fluorescence chemistry and assay development techniques are integrated into two small and portable benchtop analyzers. Our high-performance analyzers deliver point-of-care testing for better patient outcomes.

QuidelOrtho™



Simplify your workflow

Our Sofia benchtop systems and fluorescent immunoassays deliver highly accurate, objective and automated results in three to 15 minutes.



Efficient data storage and retrieval

View results remotely from your computer with Virena® System wireless connectivity. Patient results are automatically stored and results never expire.



Learn more about Sofia by scanning the QR code, visiting **quidelortho.com**, calling us at **800.874.1517** or reaching out to your Henry Schein representative.



Convenient and easy to use

The Sofia systems are small in form, have intuitive interfaces and only require about one minute of total hands-on time.

New QuidelOrtho branding may not be available in all markets, subject to country-specific regulatory approval. Please confirm with your local commercial team.

© 2024 QuidelOrtho Corporation PR-10572-NA-EN-US-v1



Sofia® Rapid Strep A Test

Easy-to-use test with enhanced assay sensitivity and results in 5 minutes. With results available in minutes, you can rapidly diagnose and treat, all in the same visit.

- Accurate detection with throat swab samples
- Sensitivity: 90.6 %
- Specificity: 96.1%
- Room temperature storage
- 24-month shelf life from date of manufacture
- For use with Sofia Analyzer (sold separately)

#20274, Sofia Strep A+
(122-7865)..... 25/box

CPT® Code: 87880*



Sofia® Influenza A+B FIA Test

This CLIA-waived test uses advanced immunofluorescence-based lateral-flow technology to detect influenza A and B viral nucleoprotein antigens in nasal and nasopharyngeal swab and nasopharyngeal aspirate/wash specimens taken directly from symptomatic patients.

- Enhanced assay sensitivity
- Results in 15 minutes
- Accurate detection with nasal swab, NP swab and NPA/wash samples
- Nasal swab sensitivity (A/B): 90%/89%; specificity (A/B): 95%/96%
- NP swab sensitivity (A/B): 97%/90%; specificity (A/B): 95%/97%
- NPA/wash sensitivity (A/B): 99%/88%; specificity (A/B): 96%/96%
- Room temperature storage
- 21-month shelf life from date of manufacture
- For use with Sofia Analyzer (sold separately)
- CLIA-waived for nasal and nasopharyngeal swabs

#20299, Sofia 2 Analyzer
(128-5292).....ea

#20218, Sofia Influenza A+B Test
(122-4092)..... 25/box

CPT® Codes: Influenza A: 87804QW*; Influenza B: 87804-59*





ID NOW™ Molecular Instrument

ID NOW is a rapid, instrument-based, isothermal system for the qualitative detection of infectious diseases. The unique ID NOW isothermal nucleic acid amplification technology provides molecular results in just minutes, allowing you to make effective clinical decisions sooner.

- Small footprint enables ease of use at point of care
- Simple operation via visual touchscreen
- Requires minimal training
- Connectivity can be provided to your network

#NAT-024, ID NOW Molecular Instrument (607-0052).....ea
#IDNOWPRINT, ID NOW Universal Printer (133-7015).....ea



ID NOW™ Strep A 2

Now the fastest Strep A test without the need for culture confirmation, the CLIA-waived ID NOW Strep A 2 test provides molecular results in 6 minutes or less on the unique ID NOW platform. This allows you to confidently test, and with a positive result, treat, during the patient's visit.

- Early detection: positive results in as few as 2 minutes, with 99% detected within 3 minutes
- Negative results in 6 minutes
- No culture confirmation of negative results
- Nearly 90% fewer false negatives than RADTs
- No incubation time or thermocycling required
- CLIA waived
- Allows for the accurate detection of group A strep
- Provides the confidence to make clinical decisions sooner
- Facilitates effective patient management
- Aids targeted antiviral therapy and antimicrobial stewardship
- Room temperature storage

#734000, ID NOW Strep A 2 (132-5149).....24/box
#734-080 ID NOW Control Swab Kit (12 positive, 12 negative) (131-9966).....ea



ID NOW™ RSV

Molecular. In minutes.
ID NOW RSV is the first molecular rapid test to detect RSV in 13 minutes or less, using our unique ID NOW platform. ID NOW allows you to make real-time clinical decisions that impact patient care.

- CLIA waived
- Room temperature storage

#435000, ID NOW RSV (126-8988).....24/box
#435080, ID NOW RSV Control Swab Kit (12 positive, 12 negative) (126-8992).....ea



ID NOW™ COVID-19 2.0

Provides reliable, real-time molecular test results in minutes for timely clinical management and infection control decisions—during the patient encounter. Detects SARS-CoV-2 in 6–12 minutes with the option to add on the ID NOW™ Influenza A & B 2 test without collecting another sample.

- Highly accurate, rapid molecular COVID-19 test
- Easily add on flu A & B testing based on clinical necessity, supporting diagnostic stewardship
- CLIA-waived intuitive procedure allows for easy standardization across care settings

#192-000, Test Kit (145-6334).....24/box
#192-080, Control Swab Kit (142-3602).....24/box
#190-010, Swab Transport Tube Accessory Pack (140-3061).....24/box
CPT® Code: 87635†



ID NOW™ Influenza A & B 2

Delivers molecular flu results in less than 5 minutes on our unique ID NOW platform, making it significantly faster than other molecular methods and more accurate than conventional rapid tests.

- Small footprint enables ease of use at point of care
- Simple operation via visual touchscreen
- CLIA waived
- Positive results detected in as few as 5 minutes, with 95% detected within 7 minutes, negative results in 13 minutes
- Detects up to 42% more true positives than RADT
- Specimen type: nasal or nasopharyngeal swabs, direct or eluted in viral transport media
- Room temperature storage

#427-000, ID NOW Influenza A & B 2 Test Kit (129-8422).....24/box
#425-080, ID NOW Influenza A & B Control Swab Kit (12 positive, 12 negative) (607-0055).....ea





ANSWER THE DOUBT.

SARS-CoV-2, Flu A, Flu B, RSV?
Detect and differentiate with **MORE** clarity.



Get **MORE** than results
with the **PCRplus** advantage.
From Cepheid.

Learn more at cepheid.com

US-IVD. In Vitro Diagnostic Medical Device. May not be available in all countries.



Xpert® Xpress Flu and Xpert® Xpress Flu/RSV

- CLIA-waived, real-time RT-PCR molecular test results
- One cartridge, one sample, three results
- Single-use, self-contained cartridge reduces waste and risk of contamination
- Room temperature (36°F–82°F, 2°C–28°C) reagent storage
- Sample types: nasal swab, nasopharyngeal swab
- Integrated QC in every cartridge
- Detects viral RNA, enabling better detection of seasonal mutations of the flu virus

#XPRSFLU-10, Xpress Flu

(127-8873).....10/box



Xpert® Xpress Strep A

- CLIA-waived, real-time RT-PCR molecular test results
- No confirmation required for negative results*
- Sample type: throat swab (ESwab™ in liquid Amies)
- Integrated QC in every cartridge
- Single-use, self-contained cartridge reduces waste and risk of contamination
- Room temperature (36°F–82°F, 2°C–28°C) reagent storage

#XPRSTREPA-10, Xpress Strep A

(129-6941).....10/box

*Follow-up testing by culture is required if the Xpert Xpress Strep A test result is negative and clinical symptoms persist, or there is an outbreak of acute rheumatic fever (ARF).



GeneXpert® Xpress

Unlock the Power of PCR at Point of Care
The Cepheid GeneXpert Xpress brings the power of CLIA-waived PCR testing to the physician office lab and beyond.

- CLIA-waived PCR testing: Flu, FLU/RSV, and Strep A
- Scalable design: 2-module and 4-module options available for customized throughput needs
- Random access: True on-demand, walkaway testing system that runs up to 4 tests at one time
- Operator and QC lockout: Customizable features to meet the needs of every practice
- Built-in redundancy: Each module operates as its own independent analyzer
- No specialized training or daily maintenance—simple to operate with minimal setup needs
- Bidirectional interface: Interface to LIS (HL7, ASTM) available

#GXIV-2-CLIA, 2-Module Configuration

(137-5675).....ea

#GXIV-4-CLIA, 4-Module Configuration

(137-5676).....ea

HELP KIDS GROW INTO HEALTHY ADULTS

Screen early with Cholestech LDX™

Early detection of familial hypercholesterolemia (FH) is key to CVD prevention, yet **only 4-5% of children 9-11 years old receive pediatric lipid screening** per NIH guidelines^{1,2}.

Benefits of screening with Cholestech LDX™



Small, 40mL
fingerstick
sample



Convenient
for kids,
parents, and
HCPs



Accurate,
reliable
results in
5 minutes



Contact your
Henry Schein representative
today.



REFERENCES: 1. Singh A, et al. *J Am Coll Cardiol*. 2019;73(19):2439-2450. 2. Knowles JW, et al. *JAMA*. 2017;318:381-382.
©2023 Abbott. All rights reserved. COL-19910



Cholestech LDX™ System

Complete Lipid Profile, Glucose

Lab-accurate lipid profile and glucose results in about 5 minutes from a simple fingerstick. Allows provider to get results while patient is still in office—provider can counsel and intervene face-to-face. This improves patient compliance and increases control and management of heart disease and diabetes. Reimbursable by Medicare and most private payers.

- CLIA waived
- Procedure saves time and enhances productivity
- Small, tough, and portable
- Eliminates anxiety of venipuncture; requires only a fingerstick
- On-the-spot adjustment of therapy and improved patient compliance
- Accuracy comparable to laboratory tests
- Saves time: fewer follow-up visits required
- Enhances compliance through teachable moments between healthcare providers and patients



#10-959, Analyzer
(738-2917)..... ea
Contains: 1 of each: analyzer, printer, power supply, starter pack, optics check cassette, printer cable, user's manual, procedure manual & training video. Starter pack contains 200 lancets, 200 capillary tubes, 200 capillary plungers, 1 MiniPet® Pipette, 50 pipette tips, 1 accessory tray & 1 fingerstick tip pad.

#97986, Test for TC
(143-3847).....10/box

#97988, Test for TC and Glucose
(143-3860).....10/box

#97987, Test for TC and HDL Panel
(143-3849)..... 10/box
#97990, Test for TC, HDL, and Glucose Panel
(143-3862)..... 10/box
#97989, Test for Lipid Profile
(143-3866)..... 10/box
#97991, Test for Lipid Profile and Glucose
(143-3872)..... 10/box
#88770, Calibration Verification Material (Lipids/
Glucose), 2 mL
(116-7301).....ea
Contains: 1 vial of each level (1–4).
#88769, Multi-Analyte Controls, 2 mL
(116-4232).....ea
Contains: 1 vial of each level (1–2).
#88773, Multi-Analyte Control, 0.25 mL
(825-0041).....ea
Contains: 1 vial of each level (1–2).
#11-010, Disposable Pipette Tips
For use with 35-µL, 40-µL, and 50-µL MiniPet Pi-
pettes.
(267-4596)..... 50/box
#10-311, Capillary Plungers
For use with all types of Cholestech capillary tubes.
(960-1587)..... 50/box
#11-787, Printer Cable, Connects Cholestech LDX to
Thermal Label Printer
(803-3213).....ea
CPT® Codes for lipid profile: #80061QW†; for lipid profile/ glucose:
80061QW†/839471QW†

THE MARKET LEADER
FOR HbA1c TESTING

Afinion™ 2 Simply More Efficient



The Afinion 2 Analyzer enables fast and easy quantitative determinations of hemoglobin A1c (HbA1c) and albumin-creatinine ratio (ACR). Its compact size and rapid result time make it ideal for managing patients with diabetes.

Benefits of HbA1c testing with Afinion 2



Small, 1.5 µL
fingerstick
sample



Easy to
use in just
3 steps



Highly
accurate
results in
as little as
3 minutes



No
maintenance
required



Contact your Henry Schein
representative today.



©2024 Abbott. All rights reserved. COL-23834



Afinion™ 2 Analyzer

With its compact size and panel of tests, the Afinion 2 system is ideal for point-of-care testing in physician offices, clinics, community health centers, and hospital outpatient clinics. From just a small urine or fingerstick whole-blood sample, lab-accurate results for ACR and HbA1c are made available during the consultation. With built-in connectivity, the Afinion 2 Analyzer can communicate with LIS and HIS systems to share and manage data more effectively, contributing to better patient management everywhere.

- Size: 7.2" x 7.5" x 13"
- Weight: 7½ lb
- Patient/operator ID and patient/control results stored in memory
- No maintenance required
- No operator calibration
- Extensive and comprehensive failsafe systems prevent the possibility of erroneous results
- QC and operator lockout
- Afinion tests can be used with both the Afinion 2 Analyzer and the Alere Afinion AS100 Analyzer

#1116985, Afinion 2 Analyzer	
(136-7869).....	ea
#1116784, Analyzer Cleaning Kit	
(136-4626).....	50/box
#14-716AFI, Universal Printer with Power Supply	
(119-8807).....	ea
#26333, Universal Printer Labels	
(119-8813).....	ea
#1115349, Barcode Scanner	
(119-8889).....	ea



Afinion™ ACR Test

Test for quantitative determination of albumin, creatinine, and albumin/creatinine ratio (ACR) in human urine, used for early identification of renal disease in patients with diabetes and/or hypertension.

- Lab-quality results in 5 minutes
- 3.5 µL of random/spot urine sample
- Quantitative results
- Increased accuracy by measuring creatinine

#1116972, ACR Test	
(136-4619).....	15/box
#1116973, ACR Controls	
(136-4620).....	ea



Afinion™ HbA1c Test

Test for quantitative determination of glycated hemoglobin (HbA1c) in human whole blood, used to monitor metabolic control in patients with diabetes.

- 3-minute test time
- Excellent precision (<3% CV)
- No interference from Hb variants
- NGSP and IFCC certified
- CLIA-waived

#1116974, HbA1c Tests	
(136-4596).....	15/box
#1116975, HbA1c Control	
(136-4618).....	ea



Childhood is Full of Firsts

Make sure kids get their first cholesterol screening between the ages of 9 and 11.¹



Make it easy with CardioChek® Plus.

Results in as little as 90 seconds from a fingerstick.



HS Item # 570-1523

Lipid Panel Starter Bundle

¹ Sleeth, C. M., Peterson, A. L., & Hotanson, J. (2020, July 1). Pediatric dyslipidemia screening by pediatricians and family medicine physicians: Current practices and future directions. American Academy of Pediatrics. https://publications.aap.org/pediatrics/article/146/1_MeetingAbstract/636/4732/Pediatric-dyslipidemia-screening-by-pediatricians?autologincheck=redirected

©2024 Polymer Technology Systems, Inc. MKG 003137/Rev. 4 04/2024



HENRY SCHEIN®

CardioChek® Plus Analyzer

Fast, Accurate, Efficient Lipid and Glucose Results at the Point of Care

Recognizing the strong connection between heart attack, stroke, and diabetes, the CardioChek® Plus Analyzer provides rapid, on-site blood-chemistry values and measurements to help identify at-risk individuals.

- Fast: On-site lipid and glucose results in as little as 90 seconds; simple, quality-control testing
- Portable: compact and lightweight, easy to transport and store; handheld and battery-operated; and room-temperature test strips
- Patient-friendly: fingerstick with sample size of 1.1 µL to 40 µL
- Wireless communication: easy data transfer
- Measures total and HDL cholesterol, triglycerides, and glucose
- Calculates LDL and non-HDL cholesterol as well as TC/HDL and LDL/HDL ratios
- Certified accuracy
- CLIA waived
- Meets NCEP guidelines for accuracy and precision
- FCC- and ISO-certified wireless communications
- CRMLN-certified for CHOL and HDL tests
- FDA cleared

#PTS-2700-HS

(570-0411).....ea

HENRY SCHEIN®

PTS Panels® Starter Kit

#PTS-4060-HS

(570-1036).....ea

Contains: Analyzer, 60 lipid/glucose tests, multichemistry controls & 40-µL cap tubes.



HENRY SCHEIN®

PTS Panels® Lipid Test Strips

#PTS-1710-HS

(570-0414).....15/box

Contains: 15 test strips.



HENRY SCHEIN®

PTS Panels® eGlucose Test Strips

#PTS-2713-HS

(570-0415).....50/box

Contains: 50 test strips.



HENRY SCHEIN®

PTS Panels® Lipid/Glucose Smart Bundle 2-Pack

#PTS-2729-HS

(570-0412).....2/pkg

Contains: 1 of each: 15-count lipid test-strip vials & 15-count vial of glucose test strips.



HENRY SCHEIN®

PTS Panels® 40-µL Capillary Tubes

#PTS-2866-HS

(570-0413).....16/bag

LAB-ACCURATE RESULTS—ON SITE. IN MINUTES.

POINT OF CARE

i-STAT® SYSTEM + PICCOLO XPRESS CHEMISTRY ANALYZER

FAST

Results in minutes*, accelerating patient care decisions.

ACCURATE

Lab-accurate results for a wide range of tests.*

VERSATILE

Supports testing in a variety of care settings.

IMPACTFUL

Improve patient satisfaction with convenient, on-site testing.

*Data on file Abbott Point of Care Inc, Instructions for use, test comparison tables.



i-STAT SYSTEM


- Data easily integrated into EMR and/or LIS
- Expert support for successful implementation, including 24/7 technical support

CLIA WAIVED TESTS:

Creatinine and Glucose

i-STAT

i-STAT Part No: 04J60-20 Waived Kit, 04J48-50 Moderately Complex Kit, 06F20-20 Waived Analyzer, 04P75-03 Moderately Complex Analyzer



PICCOLO XPRESS

- Comprehensive CLIA-waived menu available
- Intuitive function to increase operational efficiency

CLIA WAIVED TESTS:

CMP, BMP, Lipid Panel, Lipid Panel Plus, General Chemistry, Kidney Check, Renal Panel, MetLyte 8 and Electrolytes

piccolo xpress

Piccolo Xpress Part No: 07P05-50 Piccolo Xpress with 3YR Serv, 07P05-51 Piccolo Xpress with 5YR Serv



TO LEARN MORE

contact your Abbott Point of Care representative, your distribution representative, or visit www.globalpointofcare.abbott

©Abbott Point of Care Inc. | i-STAT is a trademark of Abbott. | For in vitro diagnostic use only. This material is intended for a U.S. audience only. | 400 College Road East, Princeton, NJ 08540 USA | (609) 454-9000 | (609) 419-9370 (fax) | www.globalpointofcare.abbott | Piccolo Xpress is a registered trademark of Abaxis, Inc. and distributed by Abbott Point of Care. | Abaxis 888-3355 Rev A Henry Schein Ad — Abbott POC i-STAT/Piccolo | 4930.Rev1.APOC.EN-US 07/23



Piccolo Xpress® Chemistry Analyzer

This fully automated chemistry analyzer provides the accuracy of larger laboratory analyzers in a compact and portable system.

- Only 0.1 cc of whole blood, serum or plasma needed to perform the test
- Small footprint occupies very little counter space
- Easy to use: load sample, insert disc and read results
- Self-contained reagent discs deliver complete results in 12 minutes
- 11 panels are CLIA waived—largest menu on a single platform
- 5 bidirectional USB ports for easy connectivity to EMR/LIS or peripherals
- Barcoded discs self-calibrate with every run and minimize data input errors
- iQC™ system provides continuous quality control checks on every run to eliminate the potential for erroneous results

#07P05-01, Piccolo Xpress® Chemistry Analyzer (106-6368).....ea

Contains: 1–0.1cc Minipette, 1–96-pack of disposable Minipette tips, 1 operator's manual, training CD, 1 power cord, 1 power supply, 10 lithium-heparin tubes, 1 screwdriver & 2–12-packs of results cards.

CLIA-Waived Test Panels

#07P02-08, Comprehensive Metabolic Panel: ALB, ALP, ALT, AST, BUN, Ca, Cl-, CRE, GLU, K+, Na+, TBIL, tCO₂, TP (174-6978).....10/box

#07P02-04, Basic Metabolic Panel: BUN, Ca, Cl-, CRE, GLU, K+, Na+, tCO₂, eGFR* (137-3987).....10/box

#07P02-05, Lipid Panel: CHOL, CHOL/HDL*, nHDLc*, HDL, LDL*, TRIG, VLDL* (* Calculated) (146-0189).....10/box

#07P02-12, Lipid Panel Plus: ALT, AST, CHOL, CHOL/HDL*, nHDLc*, GLU, HDL, LDL*, TRIG, VLDL* (* Calculated) (104-6269).....10/box

#07P02-10, Liver Panel Plus: ALB, ALP, ALT, AMY, AST, GGT, TBIL, TP (104-3998).....10/box

#07P02-01, General Chemistry 6: ALT, AST, BUN, CRE, GGT, GLU, eGFR* (108-5450).....10/box

#07P02-11, General Chemistry 13: ALB, ALP, ALT, AMY, AST, BUN, Ca, CRE, GGT, GLU, TBIL, TP, UA, eGFR* (982-0407).....10/box

#07P02-02, Electrolyte Panel: Cl-, K+, Na+, tCO₂ (106-1669).....10/box

#07P02-09, Kidney Check: BUN, CRE, eGFR* (109-7255).....10/box

#07P02-07, Renal Function Panel: ALB, BUN, Ca, Cl-, CRE, GLU, K+, Na+, PHOS, tCO₂, eGFR* (166-0886).....10/box

#07P02-03, MetLyte 8 Panel: BUN, CK, Cl-, CRE, GLU, K+, Na+, tCO₂, eGFR* (114-1284).....10/box

Moderately Complex Test Panels

#07P02-06, Hepatic Function Panel: ALB, ALP, ALT, AST, DBIL, TBIL, TP (195-7347).....10/box



TRUE METRIX® PRO Meter and Strips




HENRY SCHEIN® True Metrix® Pro Meter

- No coding
- Tiny, 0.5-µL sample size
- Stores 500 results
- As fast as 4 seconds
- BG range: 20–600
- Data management
- Download capabilities
- Glucose control detection
- Strip-release button

#RE4099P-43

(570-3418)..... ea

Contains: 1 True Metrix® Pro meter, installed 3-volt battery, instructions for use, logbook & compact carrying case.



HENRY SCHEIN®

True Metrix® Pro Strips

Advanced technologies—the meter, a complex algorithm, chemistry, and electrodes—on the test strip work together as part of the True Metrix® Pro system to produce accurate results. Featuring Triple Sense Technology®, the system provides proven accuracy and confidence in results.

- No-coding technology eliminates need for coding of meter
- Code is assigned during manufacturing process
- Code is printed onto contacts on end of test strip
- Meter reads code on insertion

#R3099P-450, Strips

(570-0327)..... 50/box

Controls

#R5H01-1, Level 1

(125-4034)..... ea

#R5H01-2, Level 2

(125-4035)..... ea

#R5H01-3, Level 3

(125-4036)..... ea



OneStep™ + Pro Hb Analyzer

The Henry Schein OneStep™ + Pro Hb Analyzer provides an accurate hemoglobin and hematocrit measurement.

1
test

2
results

3
seconds

This point-of-care hemoglobin analyzer is designed to provide quantitative, lab-quality hemoglobin results in only 3 seconds.

- Reagent-free microcuvettes
- Microcuvette shelf life of up to 2 years, even after the canister has been opened
- Cleared for use on patients 6 months or older
- Measurement range 5.0 to 25.6 g/dL
- Small 15 µL sample volume
- 3.25" Color LCD display
- CLIA Waived

Starter Kit

(570-1314)ea

Contains: 1 analyzer, 1 box of 100 microcuvettes, 1 optical system check, 4 AA batteries & 1 USB cable.

Analyzer

(570-1311)ea

Contains: 1 analyzer, 4 AA batteries & 1 USB cable.

Microcuvettes

(570-1312) 200/box

Contains: 4 canisters of 50 cuvettes.

Controls Set

(570-1313) ea

Contains: 3 levels of controls: Low, Mid & High.

Optical System Check

(570-1322) ea

Contains: ID chip & microcuvette.



RELY ON US FOR QUALITY, SELECTION, PERFORMANCE, AND VALUE!

More Healthcare Professionals Choose HemoCue¹

HemoCue systems set the industry standard for point-of-care hemoglobin testing. Our systems are proven through rigorous clinical study. Nearly 4,000 publications include HemoCue hemoglobin testing systems.²



Your Partner in Improving Care

Our field-based team of experts will support you every step of the way, from every day questions to detailed hands-on training.

Choose the HemoCue Point of Care Hemoglobin System to Meet Your Needs

Contact your HemoCue or Henry Schein representative today
Learn more at hemocue.us

HemoCue America | California, USA | customerservice@hemocue.com
©2024 HemoCue America. All rights reserved. MMUS-1065



1. Sanchis-Gomar, F., Cortell-Ballester, J., Pareja-Galeano, H., Banfi, G., & Lippi, G. (2013). Hemoglobin point-of-care testing: the HemoCue system. *Journal of laboratory automation*, 18(3), 198-205.
2. U.S. National Library of Medicine. (n.d.). Search: Hemocue Hemoglobin - NLM. National Center for Biotechnology Information. Retrieved February 1, 2024, from <https://www.ncbi.nlm.nih.gov/search/all/?term=hemocue+hemoglobin>



Hb 801 Analyzer

The HemoCue[®] Hb 801 system seamlessly integrates point-of-care testing into your workflow.

- Wireless connectivity through middleware
- Results in under 1 second
- Reagent-free cuvettes, improving open-vial stability
- Large display with intuitive on-screen symbols
- Internal quality control self-test
- USB adapter (included) or 3 AA batteries (not included) power supply

Hb 801 Analyzer
(135-3349)..... ea

Hb 801 Hemoglobin Microcuvettes
(135-3352)..... 200/box

Contains: 25 per bottle, 4 bottles per box.



Hb 201+ Handheld Hemoglobin Analyzer

Offers the same important benefits as its predecessor, the B-Hemoglobin Photometer. Smaller size gives more flexibility and mobility to use the system anytime and anywhere. Performs an internal self-test. No control cuvette required.

#121721, Analyzer
(345-2344)..... ea

Contains: 1 of each: analyzer, AC adapter, manual & training CD-ROM.

#111715, #3201 Microcuvettes
(110-9167)..... 100/box

Contains: 4 packs of 25.

#111716, Microcuvettes
(345-8733)..... 200/box

Contains: 4 packs of 50.


- Produces results in less than 1 minute with lab-quality accuracy
- 3 easy steps
- Precalibrated, portable analyzer that requires minimal maintenance
- No calibration or instrument manipulation (between cuvette/reagent batches) is needed
- Safe and easy to use
- Reimbursable
- CLIA '88 waived
- 3-year warranty

POC-22-NAM-3567

Meet the Quadruple Aim in Diabetes Care with In-office HbA1c and uACR

Better outcomes. Lower costs. Better patient experience. Better clinician experience.

siemens-healthineers.us/chronicdisease



SIEMENS
Healthineers



SIEMENS

DCA™ Vantage Analyzer

Point-of-care immunoassay analyzer for diabetes management. Provides precise and accurate quantitative test results for HbA1c in whole blood during a patient's visit to help improve decision making. With a single urine sample, the analyzer also provides quantitative test results for low levels of albumin, creatinine, and A:C ratio in the physician's office or clinic.

#5075US, Analyzer.....ea
 (110-0743).....ea
 #5035C, HbA1c Reagent Kit
 (147-8045).....10/pkg

#5068A, HbA1c Control Kit
 (147-6067).....ea
 Contains: 2 vials of each control: high & low.



SIEMENS

Microalbumin/Creatinine Reagent Kits

Detects and monitors albuminuria, an early indication of kidney disease. Achieve test results in just 7 minutes from a urine specimen, using the DCA™ Vantage Analyzer.

• Made in the USA

#6011A
 (147-8917).....10/pkg
 CPT® Codes: 83036QW†/82043†/82570*

SIEMENS

DCA Microalbumin/Creatinine Control Kit

• Made in the USA

#6012A
 (147-8477).....ea
 Contains: 2 vials of ea control (high & low) & reconstitution fluid.

Resilience + Passion

We make
diagnostics that

matter



We understand it is crucial to immediately diagnose and treat vulvovaginal disorders which are linked to a broad range of associated health complications, from preterm birth to an increase in Sexually Transmitted Infections.¹

The OSOM® BVBLUE® and OSOM® Trichomonas Tests are easy-to-use, CLIA-waived and produce accurate and objective results for two of the most common causes of vulvovaginal disorders so you can get the answers fast and your patients back to doing what they love.

Like you, we understand there is a patient behind every answer—and that's what matters most.

For more information, contact your local Henry Schein Medical Representative or visit henryschein.com/medical.



The Perfect Pair

OSOM® TRICHOMONAS RAPID TEST

OSOM® BVBLUE® TEST

SEKISUI
DIAGNOSTICS
Because every result matters™

¹ Brown, H. Improving the Diagnosis of Vulvovaginitis. Population Health Management. Vol. 23, suppl 1, 2020

© 2024 SEKISUI Diagnostics, LLC. All rights reserved. OSOM® is a registered trademark of SEKISUI Diagnostics, LLC. Because every result matters™ is a trademark of SEKISUI Diagnostics, LLC. BVBLUE® is a registered trademark of Gryphus Diagnostics, LLC.



osom

OSOM® Trichomonas Rapid Test

The only CLIA-waived lateral-flow rapid test for the detection of trichomoniasis while providing 95% agreement with culture and wet mount combined.

- Simple-to-use—Get easy-to-read results in just 10 minutes and 4 easy steps
- Eliminate the need for wet mount—Antigen-based rapid test does not require visualization of live organisms
- CDC-approved method of testing for the detection of *T. vaginalis*—unaffected by other common infections
- Trust your results—95% agreement against Composite Reference Standard (culture and wet mount)
- Easy-to-use—Requires minimal training and no additional equipment
- QC Inside®—2 extra tests included at no additional charge to run external quality-control testing for more confidence in results
- Made in the USA

#181, Test
(406-7631)..... 25/box
#182, Control
(611-1157)..... ea

CPT® Code: 87808QW†



osom

OSOM® BVBLUE® Test

For *Gardnerella vaginalis*

Diagnose bacterial vaginosis (BV) in just 10 minutes for faster treatment and education in the same office visit.

- Test for the 4 major pathogens that cause BV in just one test—reacts with sialidase produced by *Gardnerella vaginalis*, *Bacteroides* spp., *Prevotella*, and *Mobiluncus* spp.
- Simple to use—Get easy-to-read results in 4 easy steps and no need for additional equipment
- Decrease the chance of human error with minimal hands-on requirements
- Quicker results—Read results in just 10 minutes
- Easy to read—Instant color change for accurate results
- Get more patient tests per kit—Exclusive QC Inside® provides 2 extra tests at no additional charge to run external quality control testing
- Made in the USA

#183, Test
(414-4170)..... 25/box
#184, Control
(327-4117)..... ea

CPT® Code: 87905QW†



HENRY SCHEIN® Urispec™ Plus Urinalysis Reader & Urispec™ 11-Way Urine Strip Promotion

Provides dependable urine status results to detect early stages of many diseases, such as diabetes, kidney disease, and urinary tract infections. Instrument-read results have long proven to be advantageous for both busy health care professionals and patients. Comprehensive interface options and the optimized printouts minimize risk associated with manual transcriptions. Reliable results can be obtained virtually immediately at the point of care. Easy switch between different test strips: Urispec Plus is for use with the Urispec 11-Way, 5-Way Plus, 4-Way, and GP+A urine test strips for up to 11 parameters for the semiquantitative detection of occult blood, urobilinogen, bilirubin, protein, nitrate, ascorbic acid, ketones, glucose, pH, specific gravity, and leukocytes. Use of these test strips on the Urispec Plus analyzer is CLIA-waived.

- CLIA-waived
- Readings eliminate subjectivity of visual color interpretation, which makes urinalysis simpler and more reliable
- Features no-touch operation: After strip placement, measurement starts automatically, which makes urinalysis faster and more hygienic
- Easy switch between different test strips: Urispec Plus is for use with Urispec 11-Way, 5-Way Plus, 4-Way, and GP+A urine test strips

(570-1309).....ea
Contains: 1 Urispec™ Plus Urinalysis Reader & 12 bottles of Urispec™ 11-Way Urine Strips.



HENRY SCHEIN® Urispec™ 5-Way Plus Urine Reagent Test Strips

- Tests for glucose, blood, protein, nitrite, leukocytes, and urobilinogen
- Strips for visual read or for use with the Henry Schein Urispec Plus Urinalysis Reader only

(900-7508).....100/btl
CPT® Codes: 81000†/81001†/81002†/81003†

HENRY SCHEIN® Urispec™ 4-Way Urine Reagent Strips

- Protein, nitrite, glucose, and leukocytes
- Ideal for OB/GYN
- Infection, diabetes, and kidney function
- Dip and read in 60 seconds
- Strips for visual read or for use with Urispec Plus Urinalysis Reader

(107-7157).....100/btl
CPT® Codes: 81000†/81001†/81002†/81003†



HENRY SCHEIN® Urispec™ GP+A Urine Reagent Test Strips

Test for glucose and protein, with an ascorbic acid test for accuracy.

- Protein and glucose
- Ideal for metabolic, kidney, and diabetic monitoring
- Dip and read in 60 seconds
- Strips for visual read or for use with Urispec Plus Urinalysis Reader

(100-6567).....100/btl
CPT® Codes: 81000†/81001†/81002†/81003†



HENRY SCHEIN® Urispec™ 11-Way Urine Reagent Test Strips

Compare to Hydratrend Urine Reagent Test Strips

For urinalysis by the dip-and-read technique. Similar to all Multistix® 10 SG Plus Ascorbic Acid. Tests for glucose, bilirubin, ketone, specific gravity, blood, pH, protein, urobilinogen, nitrite, and leukocytes. Also includes additional ascorbic acid analyte pad. The importance of the addition of an ascorbic acid analyte pad is that the glucose test must be repeated if the ascorbic acid reaction is positive. An ascorbic acid content of as little as 5 mg/dL (0.3 mmol) can interfere with a glucose assay in low concentrations. As with glucose detection, blood detection is also disturbed by low concentrations of ascorbic acid, whereas high ascorbic acid concentrations can interfere with the nitrate and bilirubin test zones.

- Unique "read all" feature at 60 seconds
- Similar to Multistix® 10 SG Plus Ascorbic Acid Test and all brands of 10 SG products
- Glucose, bilirubin, ketone, specific gravity, blood, pH, protein, urobilinogen, nitrite, and leukocytes
- Reader available
- Ascorbic acid pad to ensure uncompromised results
- Room-temperature storage
- Strips for visual read or for use with Urispec Plus Urinalysis Reader

(102-2285).....100/btl
CPT® Codes: 81000†/81001†/81002†/81003†



HENRY SCHEIN® Urispec™ Pro 2 GP Urine Strips

Strips for Visual Read
Tests for glucose and protein.

- Made in the USA

(900-4996).....100/btl
CPT® Codes: 81000†/81001†/81002†/81003†



HENRY SCHEIN® Urispec™ 5 OB Urine Strips

Strips for Visual Read
Tests for glucose, protein, blood, nitrites, and leukocytes.

- Made in the USA

(900-4995).....100/btl
CPT® Codes: 81000†/81001†/81002†/81003†



HENRY SCHEIN® Urispec™ 10SG Urine Strips

Strips for Visual Read
Tests for glucose, bilirubin, ketone, urobilinogen, specific gravity, blood, pH, protein, nitrites, and leukocytes.

- Made in the USA

(900-4994).....100/btl
CPT® Codes: 81000†/81001†/81002†/81003†



HENRY SCHEIN® Urispec™ Urinalysis Controls

For Most Brands of Routine Urine Dip Sticks and Pregnancy Tests

A bilevel product packaged in plastic tubes and dropper bottles. Controls are liquid, ready-to-use, and are designed to validate urinalysis reagent strips. Packaged in plastic dip tubes and dropper bottles, they are convenient and easy to use. Long open- and closed-vial stability at 2°C–8°C for 2 years unopened, or 90 days once the Dip Control is opened, and 180 days once the Drop Control is opened. This allows the user ample time to utilize the product before expiration, resulting in cost savings.

- Values for all pregnancy test brands
- Assays for common brands of urine dip for routine urinalysis (2 parameter-10SG)
- Includes values for microalbumin, creatinine, and hCG
- Two distinct levels (normal and abnormal)
- Closed-vial dating up to 2 years with refrigerator storage
- 30-day stability at room temperature

3 Bilevel Sets, 15 mL/Tube (6 Vials)
Dip stick in tube, with open-vial stability of 90 days at 2°C–8°C.

(112-6841).....ea
2 Bilevel Sets, 25 mL/Dropper Bottle (4 Bottles)
Open-vial stability of 180 days at 2°C–8°C.
(112-6842).....ea





HENRY SCHEIN®

Specimen Containers

These 4-oz specimen containers have blue, tamper-evident, leak-resistant tops and are suitable for routine urine or specimen collection and transport. Container includes patient information label.

- Sterile, with Tab, Unwrapped, 4 oz
(900-7433).....100/case
- Sterile, Individually Wrapped, 4 oz
(900-7431).....100/case
- Non-Sterile, Bulk, 4 oz
(900-7432).....500/case



globe
SCIENTIFIC INC.

Uri-Pak Urine Collection Kits with Tubes

The Uri-Pak® Kit offers all the components necessary for specimen collection, handling and testing. Each kit consists of 5 plastic bags. Each bag contains:

- 100–12 mL PS centrifuge tubes
- 100–PE snap-caps
- 100–3.5 oz collection cups
- 100–adhesive patient ID labels
- Optional 10-place disposable tube rack

- #112018R, Kit with 12 mL Tubes, Blue Flat-Top Caps, Cups, Labels, and Cardboard Rack
(424-0042)..... 500/case
Contains: 100 tests per bag.
- #112019R, Kit with 12 mL Tubes, Yellow Caps with Sanitary Grips, Cups, Labels, and Cardboard Rack
(395-6061)..... 500/case
Contains: 100 tests per bag.
- #112019R, Kit with 12 mL Tubes, Yellow Caps with Sanitary Grips, Cups, and Labels
(120-4241)..... 500/case
Contains: 100 tests per bag.
- #112192, Patient ID Labels, 25/Sheet
(781-6446).....2000/case



globe
SCIENTIFIC INC.

Specimen Transport Bag

These specimen bags are available with a zipper closure or tamper-evident glue seal. The bags are printed with a biohazard symbol and feature an attached document pouch that keeps important paperwork safe.

- #4918, 6" x 9" Bag with Score Line & Document Pouch
(106-5227).....1000/case



globe
SCIENTIFIC INC.

Urine Collection Containers and Cups

The 4-oz (120-mL) containers are available in both sterile and non-sterile versions. Non-sterile containers have a frosted writing area for patient information. Sterile containers are individually wrapped and have a trilingual patient ID label. They also have easy-to-turn, leak-resistant screw caps and molded graduation marks. The 6.5-oz (190-mL) cups come in both polypropylene and polystyrene versions and are available with or without paper lids.

- #5912, Trilingual Patient ID Label, Sterile, Individually Wrapped
(606-6028).....100/case
- #5913, 4-oz Container, Sterile, Individually Wrapped
(868-0112).....100/case
- #5915, 4-oz Container, Non-sterile, Etched Writing Space, Bulk
(642-6596).....500/case
- #5916, 4-oz Container, Non-sterile, with Affixed Thermometer Strip
(413-4578).....500/case
- #5920, Polypropylene (PP) Beakers
(868-8877).....500/case



globe
SCIENTIFIC INC.

TransferTop™ Urine Collection and Transfer Containers

- Requires no open transferring of samples for reduced spills, improved lab hygiene, and limited exposure to dangerous samples
 - Reduces risk of contamination for improved test results
 - Directly transfers samples to a tube suitable for analysis and transport
 - Works with most vacuum urine tubes
 - Provides consistent sample transfer volumes when used with like vacuum tubes
 - Improves efficiencies and sample throughput when compared to traditional open systems
- #3856L, 120 mL
(129-2541).....300/case



HENRY SCHEIN®

Top Hat Specimen Collector

- 800 cc
(112-5810).....25/box



To Order: **1-800-772-4346** 8:00am–8:00pm (et)

henryschein.com/medical

*California customers please use website
when ordering for Prop 65 information.*

Due to manufacturer price increases, prices are subject to change without notification. Offer valid on our full array of products. Upon any request by Medicare, Medicaid, Tricare or other payer, you must report the value of any benefit received from a discount program (e.g., points redemptions, gifts or other special awards). ©2024 Henry Schein, Inc. No copying without permission.

Not responsible for typographical errors. No adjustments from prior sales. Offer not to be combined with other promotions and/or special contract pricing. Void where prohibited. Local regulatory requirements may apply to use or installation of certain products. It is your responsibility to understand and comply with any such requirements prior to purchase, use or installation of products.

[†]For informational purposes only.

Source: ReimbursementCodes.com and/or www.cms.gov.

Customer is responsible for verification of billing/coding in accordance with applicable specific circumstances.



RELY ON US™ FOR YOUR 2024–2025
FLU VACCINE NEEDS!

A Flu Vaccine Program you can trust from **A COMPANY YOU CAN RELY ON**

This is our 18th year of offering the Worry-Free Guarantee Flu Vaccine Program,
more than 90 years of expertise in pharmaceuticals!



Market-Leading Brands

- Seqirus Flucelvax® Trivalent, Afluria® Trivalent, and Flud® Trivalent
- Sanofi Fluzone® Trivalent, Fluzone® HD Trivalent, and Flublok® Trivalent
- GSK Fluarix® Trivalent and Flulaval® Trivalent



Product Choice to Meet All Your Needs

- Adult and Pediatric
- Multi-dose Vials and Prefilled Syringes
- Cell Culture and Egg-based Options Available



sanofi

GSK

WORRY-FREE PROGRAM COMPONENTS

ORDER TODAY!

Ask Your Sales Consultant for Details!

- Guaranteed Ship By Date
- Extended Payment Terms to 12/19/2024[†]
- Return Privilege^{*}
- Option to increase quantity at time of order or later in the season^{**}
- Discounts on ancillary and companion items during flu season
- Flu Practice Marketing Kit^{***}

Visit www.henryschein.com/flu

^{*}Offer good on minimum purchase of 7 vials and/or packs of Flucelvax, Afluria, Flud, Fluzone, Flublok, Fluarix and Flulaval. Returns must be in complete unopened vials or packs of syringes and allowable return percentage is specific to brand. Ask your sales consultant for more information. Returns will be in a form of a credit toward next year's flu purchases. Allowable flu vaccine returns must be received back to Henry Schein, Inc. between March 1, 2025 and April 30, 2025 for full returns credit.

^{**}Subject to availability. ^{***}Not available in California. [†]Subject to credit approval.

©2024 Henry Schein, Inc. No copying without permission. Not responsible for typographical errors.