

# Web Tutorial: Enabling the Barcode Scanner

1

## Access Keychain Barcode

**Note:** This feature only works with Internet Explorer browser.

- A. Login
- B. Click **My Account**
- C. Click **Barcode Scanner**

A

B

C

**Barcode Scanner**

For more detailed information about the Coronavirus please visit [www.henryscheindental.com/coronavirus](http://www.henryscheindental.com/coronavirus). THIS WEBSITE PROCESSES U.S. CUSTOMER ORDERS ONLY. [CLICK HERE FOR IMPORTANT SHIPPING INFORMATION.](#)

**My Account** / My Account / Barcode Scanner

Welcome to our new and easy to use Barcode Scanner-just Point, Click and Scan!

Here is what you will need:

**System Requirements**

- Browser: Internet Explorer 6.0+ [Download Here](#)
- Operating System: Windows XP, Vista, 7, 8
- USB Drivers
- ActiveX Installation

**Hardware Requirements**

- Barcode Scanner: Symbol CS1504 - OR - Opticon OPN2003/2004
- USB Port: 1
- Cables: USB or Serial

**Activation Process**

- Step 1. Download and Install the USB Driver:
  - If you are using the Symbol CS1504 Barcode Reader, download the Symbol USB Driver [Download Here](#)
  - If you are using the Opticon OPN2003 or OPN2004, download the Opticon USB Driver [Download Here](#)
- Step 2. Install the Scanner's barcode reader Activex control for XP [Download Here](#)
- Step 3. For any issues using our ActiveX plugin in Internet Explorer on your computer, please install the following Windows application to upload your items to the Henry Schein website, for ordering and inventory management purposes. [Download Here](#)
- Step 4. Click Enable Scanner

**Need to order the Barcode Scanner?**

To order the Barcode Scanner or for more information on how this new device can make your ordering process simpler, please contact your Henry Schein Dental Sales Consultant.

[Print a list of Item Quantity Barcodes](#)  
[Print Barcode from a Shopping List](#)  
[Print Barcodes from your items purchased](#)

**Already Own the Barcode Scanner?**

If you have already purchased the Barcode Scanner, simply click the 'Enable scanner' button below. Enabling your scanner will provide you with the one-time installation of the scanner's barcode reader called ActiveX. If the installation is not automatic, you may click the 'ActiveX Installation' link available within System Requirements or refer to the Setup Instructions.

**ENABLE SCANNER**

2

## Enable Keychain Barcode

**Note:** You may need to contact ECS support to load this software.

- A. Review and follow **Activation Process** steps
- B. Click **Enable Scanner** button to complete installation

A

B