



# ETCO<sub>2</sub> Sampling Cannulas



MedSource ETCO<sub>2</sub> Sampling Cannulas allow for a breath-by-breath sampling of ETCO<sub>2</sub> while simultaneously delivering oxygen.



Sampling of ETCO<sub>2</sub> in exhaled respiratory gases allows for healthcare workers to assess adequacy of patient perfusion, ventilation and oxygenation in patients receiving critical care. These cannulas pair with Capnograph Sidestream Monitoring Systems and come in a variety of configurations to meet hospital, acute care and prehospital ETCO<sub>2</sub> monitoring needs.

## SPECIFICATION

- Dip molded PVC nasal cannula
- Male luer lock
- Not made with natural rubber latex
- PVC material not made with DEHP
- Kink resistant tubing
- Pair with any sidestream monitor (low flow 50ml/min)



[www.medsourcelabs.com](http://www.medsourcelabs.com)

Item #	Description	UOM
MS-35400A7	Adult Nasal CO <sub>2</sub> Sampling Cannula W/O <sub>2</sub> delivery Luer Conn. 7ft.	25 EA/BX, 200 EA/CS
MS-35400A10	Adult Nasal CO <sub>2</sub> Sampling Cannula W/O <sub>2</sub> delivery Luer Conn.10ft.	25 EA/BX, 150 EA/CS
MS-35400A14	Adult Nasal CO <sub>2</sub> Sampling Cannula W/O <sub>2</sub> delivery Luer Conn.14ft.	25 EA/BX, 100 EA/CS

# ETCO<sub>2</sub> Sampling Cannulas



## Button Filter

Minimizes entrance of moisture from sampling ports to monitor

MedSource ETCO<sub>2</sub> Sampling Cannulas allow for a breath-by-breath sampling of ETCO<sub>2</sub> while simultaneously delivering oxygen.

## Universal Connector

Pairs with all oxygen delivery systems

**Kink Resistant Tubing**

## Softed Nares


Ultra soft PVC material for patient comfort

## Reflective Orange Connector

Pairs with sidestream low flow capnography monitors

## Softed Tear Shaped Scoop

Lifts up and away from patient's lip for extended comfort



Sampling of ETCO<sub>2</sub> in exhaled respiratory gases allows for healthcare workers to assess adequacy of patient perfusion, ventilation and oxygenation in patients receiving critical care. These cannulas pair with Capnograph Sidestream Monitoring Systems and come in a variety of configurations to meet hospital, acute care and prehospital ETCO<sub>2</sub> monitoring needs.

## SPECIFICATION

- Dip molded PVC nasal cannula and oral scoop
- Reflective orange connector
- Not made with natural latex rubber
- PVC material not made with DEHP
- 7' of kink resistant tubing
- Pair with any sidestream monitor (low flow 50ml/min)



[www.medsourcelabs.com](http://www.medsourcelabs.com)

Item #	Description	UOM
MS-35601A	Adult Nasal/Oral Dual CO <sub>2</sub> Sampling Cannula W/ O <sub>2</sub> Delivery, Orange Luer W/Reflective Conn., Sampling Scoop, Filter, 7ft.	25 EA/BX, 200 EA/CS
MS-35601P	Pediatric Nasal/Oral Dual CO <sub>2</sub> Sampling Cannula W/ O <sub>2</sub> Delivery, Orange Luer W/Reflective Conn., Sampling Scoop, Filter, 7ft.	25 EA/BX, 200 EA/CS
MS-35700	Intubated Sampling Line W/T-Connector, Orange Luer W/Reflective Conn., Sampling Scoop, Filter, 7ft.	25 EA/BX, 100 EA/CS